

# BENTLEY LT



The *lightweight Bentley LT* offers up to 20° of tilt from the front of the seat. This allows the user to keep their *feet on the ground* while the chair is fully tilted. The rigid frame provides *superior stability* and maneuverability.



Hand Propulsion



Foot Propulsion



Self-Tilt



Desk/Table Accessible



Transit Tie Downs

**Bentley LT** 20° knee pivot tilt chair with a lightweight aluminum frame

## Bentley exclusive features

- ✓ Foot and hand propel in tilt position with no knee-rise
- ✓ Short wheelbase for maneuverability
- ✓ Adjustable rear axle position to move center of gravity



## Available options

- ✓ Transit tie downs WC19 & ISO7176-19
- ✓ 30° manual dynamic recline
- ✓ Power tilt
- ✓ Easy-Out front rigging



DYNAMIC TILT  
Up to 20°



WEIGHT CAPACITY  
250 lbs



SEAT WIDTH  
14" to 22"



SEAT DEPTH  
16" to 22"



SEAT HEIGHT  
13" to 20"



BACK SEAT ANGLE  
90° to 122°



OVERALL WIDTH  
+9.5" to 10.5"



CHAIR WEIGHT  
33 lbs

visit [pdgmobility.com/blt](https://pdgmobility.com/blt)  
for more information