FUZE T20



The *Fuze T20* offers up to 20° of tilt with ease of propelling for occupants and caregivers. The chair has a sleek aluminum frame and is *easy to set up and adjust*. Its forward pivot location provides *low knee-rise* during tilt.











Fuze T20 A sleek 20° knee pivot tilt chair

visit pdgmobility.com/t20

for more information

Available options See order form

✓ Adjustable width frame



DYNAMIC TILT Up to 20°



WEIGHT CAPACITY 250 lbs



SEAT WIDTH 14" to 22"



SEAT DEPTH 15" to 22"



SEAT HEIGHT 13" to 20"



BACK SEAT ANGLE 84° to 120°



0VERALL WIDTH +9.5" to 10.5"



CHAIR WEIGHT 42 lbs

