

# POWER TILT




The *Power Tilt* is an optional upgrade for the PDG Fuze T50, Bentley, and select Stellar tilt chairs. It allows users and/or attendants to **operate the tilt feature with ease**. With two included battery packs, the power tilt is **always ready**.

## Exclusive features

- ✓ Factory installed by PDG - no third-party shipping, expenses, or delays involved
- ✓ System can also be installed on chairs in the field with a kit available from PDG
- ✓ All Power Tilts include 2 battery packs and an offboard charger to charge one battery while the other one is in use on the chair
- ✓ Toggle, button, and egg switches with auto-reversing function available
- ✓ Covered by HCPCS Code E1002 in the United States


## Available on these PDG tilt chairs

visit [pdgmobility.com/power-tilt](https://pdgmobility.com/power-tilt)  
for more information



### *Fuze T50 and T50n*

- ✓ Up to 350 lbs capacity
- ✓ Fixed and adjustable frame



### *Bentley and Bentley LT*

- ✓ Up to 250 lbs capacity



### *Stellar GL, GLT, and LEAP*

- ✓ Up to 250 lbs capacity