

STELLAR GLT



The *lightweight Stellar GLT* is a versatile 45° low knee-rise tilt chair. It promotes self propulsion and *participation in mobility related activities of daily living*. The rigid frame provides superior stability and maneuverability.



Hand Propulsion



Foot Propulsion



Self-Tilt



Desk/Table Accessible



Transit Tie Downs

Stellar GLT 45° tilt chair with a lightweight aluminum frame

Stellar exclusive features

- ✓ Foot and hand propel in tilt position with low knee-rise
- ✓ Short wheelbase for maneuverability
- ✓ Adjustable rear axle position to move center of gravity

Available options

- ✓ Transit tie downs WC19 & ISO7176-19
- ✓ 30° manual dynamic recline
- ✓ Power tilt
- ✓ Easy-Out front rigging



DYNAMIC TILT
Up to 45°



WEIGHT CAPACITY
250 lbs



SEAT WIDTH
14" to 22"



SEAT DEPTH
16" to 22"



SEAT HEIGHT
13" to 20"



BACK SEAT ANGLE
90° to 122°



OVERALL WIDTH
+9.5" to 10.5"



CHAIR WEIGHT
34 lbs

visit pdgmobility.com/glt
for more information