

# STELLAR IMPACT



The *Stellar Impact* was designed for the most demanding chair users. It offers up to 45° tilt, and its *ultra-tough steel construction* provides *maximum durability* for high tone, highly agitated, and/or continuously rocking occupants.



Hand Propulsion



Foot Propulsion



Self-Tilt



Desk/Table Accessible



Transit Tie Downs

**Stellar Impact** *The toughest tilt chair in the world*

*Available options* See order form

- ✓ Transit tie downs WC-19 & ISO-7176-19
- ✓ 30° manual dynamic recline
- ✓ Dynamic Rocker Back

visit [pdgmobility.com/impact](https://pdgmobility.com/impact) for more information



DYNAMIC TILT  
Up to 45°



WEIGHT CAPACITY  
250 lbs



SEAT WIDTH  
14" to 22"



SEAT DEPTH  
17" to 22"



SEAT HEIGHT  
14" to 21"



BACK SEAT ANGLE  
90° to 122°



OVERALL WIDTH  
+9.5" to 10.5"



CHAIR WEIGHT  
64 lbs