



# BENTLEY LT-R

**Order Form - Canada**  
Prices in Canadian Dollars. January 1, 2024

## Customer and Order Information

Purchase Order *	Tag	v.0101
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in
		<i>*required</i>

## Base Chair

**Bentley LT-R 20° Manual Tilt Wheelchair** Reverse Wheel Layout enhances user access for self-propulsion. [pdgmobility.com/blt-r](http://pdgmobility.com/blt-r)

<input checked="" type="checkbox"/> 2701000	Bentley 20° Manual Tilt Wheelchair	\$ 4,300
<input checked="" type="checkbox"/> 2701009	Lightweight Aluminum Frame and Reverse Wheel Layout Upgrade	\$ 300

**i** Choose the standard **Bentley** with a steel frame for agitated or heavier users. See separate order form.

## Weight Capacity

2702000 **250 lbs Chair Weight Capacity**

**User Weight Range** Required for setup of tilt uplift force.

<input type="checkbox"/> 2702430	Up to 75 lbs	<input type="checkbox"/> 2702431	75 lbs to 100 lbs	<input type="checkbox"/> 2702432	100 lbs to 115 lbs
<input type="checkbox"/> 2702433	115 lbs to 150 lbs	<input type="checkbox"/> 2702434	150 lbs to 200 lbs	<input type="checkbox"/> 2702435	200 lbs to 250 lbs

## Frame Options

2702300 **Transit Tie Down System (TTDS)**

TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details.

## Seat Width

Overall chair width is +9.5".

**Standard Seat Widths**

<input type="checkbox"/> 2703014	14" Seat Width	<input type="checkbox"/> 2703015	15" Seat Width	<input type="checkbox"/> 2703016	16" Seat Width
<input type="checkbox"/> 2703017	17" Seat Width	<input type="checkbox"/> 2703018	18" Seat Width	<input type="checkbox"/> 2703019	19" Seat Width
<input type="checkbox"/> 2703020	20" Seat Width	<input type="checkbox"/> 2703021	21" Seat Width	<input type="checkbox"/> 2703022	22" Seat Width

## Seat Depth

**Standard Seat Depths (Adjustable)**

<input type="checkbox"/> 2704016	16" Seat Depth	<input type="checkbox"/> 2704017	17" Seat Depth	<input type="checkbox"/> 2704018	18" Seat Depth
<input type="checkbox"/> 2704019	19" Seat Depth	<input type="checkbox"/> 2704020	20" Seat Depth		

**Custom Seat Depths**

<input type="checkbox"/> 2704021	21" Seat Depth	<input type="checkbox"/> 2704022	22" Seat Depth	\$ 400
----------------------------------	----------------	----------------------------------	----------------	--------

## Frame Color

**PDG Premium Colors** View color swatches on PDG's website. [pdgmobility.com/colors](http://pdgmobility.com/colors)

<input type="checkbox"/> 2705001	Black Shimmer	<input type="checkbox"/> 2705015	Sky Blue	<input type="checkbox"/> 2705016	Poppy Red
<input type="checkbox"/> 2705003	Fine Texture Black	<input type="checkbox"/> 2705010	Lollypop Blue	<input type="checkbox"/> 2705011	Lollypop Red
<input type="checkbox"/> 2705017	Arctic Grey	<input type="checkbox"/> 2705018	Tropical Teal	<input type="checkbox"/> 2705019	Tangerine Orange
<input type="checkbox"/> 2705020	Olive Green	<input type="checkbox"/> 2705021	Lush Lilac	<input type="checkbox"/> 2705012	Bubble Gum Pink
<input type="checkbox"/> 2705022	Leaf Green	<input type="checkbox"/> 2705023	Boysenberry Purple	<input type="checkbox"/> 2705014	Daffodil Yellow

2705200 **Custom Color** Contact PDG for specifications. \$ 300

### Armrests

2706000 **No Armrests**

**Height Adjustable Armrests (Single Post)** \$ 295

2706111 8"-11" Range       2706114 10"-14" Range       2706116 12"-16" Range

2706100 **Fixed Height Armrests (12")** Full Length Armpad only.

**Armrest Height** If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.  
Seat Pan to top of Armrest: \_\_\_\_\_ Does not include height of Armpad

**Armrest Position** Indicate preference for Armrest position. If left blank, the default setting will be used.  
In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

2706401 Rear       2706402 Mid       2706403 Front

### Armpads

**Regular Armpads**

2707100 Desk Length - 2" x 10"       2707101 Full Length - 2" x 14"

**Waterfall Armpads** \$ 115

2707105 Desk Length - 2" x 10"       2707106 Full Length - 2" x 14"

**Wide Armpads** May interfere with wheels in some tilt configurations. \$ 200

2707110 Desk Length - 3" x 10"       2707111 Full Length - 3" x 12"

**Gel Armpads** May interfere with wheels in some tilt configurations. \$ 135

2707120 Full Length - 2" x 12"       2707121 Full Length - 3.5" x 12"

**Contour Arm Trough with Hand Pad** May interfere with wheels in some tilt configurations. \$ 125

2707130 Right (4.5" x 18")       2707131 Left (4.5" x 18")

### Armrest Options

**Armrest Wideners** Move Armrests outward by up to 1" on each side. Seat size does not change. May cause tilt interference.

2707210 1" Total Additional Width \$ 95

2707211 2" Total Additional Width \$ 195

2707310 **Clothing Guards** Not compatible with 8"-11" Height range. \$ 110

**Armrest Extensions** Not compatible with Desk Length Armpads. \$ 75

2707400 3" Extension - Forward       2707401 3" Extension - Rearward

### Seat Pan

2708000 **No Seat Pan**

2708001 **Solid Aluminum Seat Pan** \$ 135

### Seat Height

**Don't forget** **Select desired Seat-to-Floor Height** STF Height is adjustable within the range of the chosen Subframe Plate (Low/High).

6" Casters	Low						High						
	RW	12"	13"	14"	15"	16"	RW	14"	15"	16"	17"	18"	19"
20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

8" Casters	Low						High						
	RW	12"	13"	14"	15"	16"	RW	14"	15"	16"	17"	18"	19"
24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

### Casters

**Urethane Casters** 8" Casters are compatible only with 24" Large Wheels.

<input type="checkbox"/>	2711061	6" x 1.25" Caster size	<input type="checkbox"/>	2711081	8" x 1.25" Caster size	
<input type="checkbox"/>	2711063	6" x 1.50" Caster size				\$ 85
<input type="checkbox"/>	2711084	8" x 2.00" Caster size				\$ 135

**Pneumatic Casters** 8" Casters are compatible only with 24" Large Wheels.

<input type="checkbox"/>	2711161	6" x 1.25" Caster size				\$ 110
<input type="checkbox"/>	2711181	8" x 1.25" Caster size				\$ 120
<input type="checkbox"/>	2711184	8" x 2.00" Caster size				\$ 150

<input type="checkbox"/>	2711260	<b>Frog Legs 6" x 1.50" Aluminum Soft-Roll Casters</b>				\$ 275
--------------------------	---------	--	--	--	--	--------

<input type="checkbox"/>	2711400	<b>Frog Legs Suspension Forks</b>	Only compatible with 24" Large Wheels. Requires 6" Soft-Roll Caster.			\$ 595
--------------------------	---------	-----------------------------------	--	--	--	--------

### Large Wheels

**Composite Wheels**

<input type="checkbox"/>	2712020	20" Wheel size	<input type="checkbox"/>	2712022	22" Wheel size	<input type="checkbox"/>	2712024	24" Wheel size
--------------------------	---------	----------------	--------------------------	---------	----------------	--------------------------	---------	----------------

**Spoke Wheels**

<input type="checkbox"/>	2712120	20" Wheel size	<input type="checkbox"/>	2712122	22" Wheel size	<input type="checkbox"/>	2712124	24" Wheel size	\$ 75
--------------------------	---------	----------------	--------------------------	---------	----------------	--------------------------	---------	----------------	-------



**Large Wheel Position** Indicate Large Wheel positioning preference. If left blank, the default setting will be used.

In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

<input type="checkbox"/>	2712401	Rear	<input type="checkbox"/>	2712402	Mid	<input type="checkbox"/>	2712403	Front	Not compatible with 20"-22" Seat Depth
--------------------------	---------	------	--------------------------	---------	-----	--------------------------	---------	-------	--

### Axle Type

<input type="checkbox"/>	2713000	<b>Quick Release Axles</b>
--------------------------	---------	----------------------------

<input type="checkbox"/>	2713010	<b>Fixed Axles</b>
--------------------------	---------	--------------------

### Large Tires

<input type="checkbox"/>	2714000	<b>Urethane Tires</b>		
--------------------------	---------	-----------------------	--	--

<input type="checkbox"/>	2714010	<b>Pneumatic Tires</b>	May result in a finished Seat Height 1/2" greater than stated.	\$ 60
--------------------------	---------	------------------------	--	-------

<input type="checkbox"/>	2714020	<b>Pneumatic Tires with Flat-Free Inserts</b>	May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels.	\$ 95
--------------------------	---------	---	---	-------

### Handrims

<input type="checkbox"/>	2715000	<b>No Handrims</b>
--------------------------	---------	--------------------

<input type="checkbox"/>	2715010	<b>Standard Aluminum Handrims</b>
--------------------------	---------	-----------------------------------

<input type="checkbox"/>	2715020	<b>Plastic Coated Handrims</b>	\$ 150
--------------------------	---------	--------------------------------	--------

<input type="checkbox"/>	2715030	<b>Projection Handrims</b>	Composite Wheels only. 10 Projections, 30°, 1-1.8" Length.	\$ 495
--------------------------	---------	----------------------------	--	--------

### Wheel Locks & Brakes

**User Wheel Locks**

<input type="checkbox"/>	2716010	Push to Lock	<input type="checkbox"/>	2716020	Pull to Lock
--------------------------	---------	--------------	--------------------------	---------	--------------

**User-activated Unilateral Wheel Lock** Pull to Lock.

<input type="checkbox"/>	1016310	Left Activation	<input type="checkbox"/>	1016311	Right Activation	\$ 795
--------------------------	---------	-----------------	--------------------------	---------	------------------	--------

### Wheel Lock Options

<input type="checkbox"/>	2716100	<b>Wheel Lock Extensions</b>	\$ 60
--------------------------	---------	------------------------------	-------

<input type="checkbox"/>	2716200	<b>Hill Holders</b>	Requires Pneumatic Large Tires with Flat-Free Inserts.	\$ 450
--------------------------	---------	---------------------	--	--------



# BENTLEY LT-R

**Order Form - Canada**  
Prices in Canadian Dollars. January 1, 2024

## Anti-Tippers

### Front Anti-Tippers Required Safety Option. Must choose one type.

- 2717010 Standard \$ 75
- 2717015 Spring-loaded Dampens the forward-tip movement of the chair during transfers and allows for improved ground clearance. Not removable. \$ 200

## Back Posts

### Back Post Style

- 2718010 Stroller Handle Back

### Back Post Height Measured from the Seat Pan to the Handle.

- 2718425 Short (22")  2718431 Tall (25")

### Back Post Folding Fold-down Back is not compatible with Reclining Backrest.

- 2719220 Fold-down

### Seat to Back Angle Seat Angles above 98° may cause interference at Seat Depths above 19".

- 2719090 Set at 90°  2719098 Set at 98°  2719106 Set at 106°
- 2719114 Set at 114°  2719122 Set at 122°

## Back Options

### Angle/Height Adjustable Stroller Handle Not compatible with Headrest Mounting Plate.

- 2719010 Button on Outside  2719011 Button on Inside \$ 150

### Endomorph Back Backrest will be narrower than Seat by total of 2x the Offset selected below. Minimum width of Backrest is 12".

- 2718102 1" Offset per side  2718103 1.5" Offset per side  2718104 2" Offset per side \$ 395

### 10" Backrest Extension Not compatible with Angle/Height Adjustable Stroller Handle.

- 2719210 \$ 150

### Headrest Mounting Plate

Not compatible with Backrest Extension, Angle/Height Adjustable Stroller Handle. Not required if mounting to a solid Back.

- 2719212 \$ 95

## Back Upholstery

- 2720000 **No Back Upholstery**

### Standard Upholstery

- 2720418 Short (18" Height)  2720421 Regular (21" Height)  2720424 Tall (24" Height) \$ 195

### Tension Adjustable Upholstery Height of Upholstery is 16".

- 2720450 \$ 295

### Position of Upholstery Select the desired position of the top of the Upholstery relative to the Seat Pan.

		22" Back Post									
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"		
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>					
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>					
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>					

		25" Back Post									
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"		
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Tall							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

## Tilt Activation

### Attendant Tilt only Select trigger mounting location below.

- 2723011 Right Back Post  2723012 Left Back Post

### User Tilt only Full Length Arm Pads only. Select trigger mounting location below.

- 2723021 Right Armrest  2723022 Left Armrest \$ 100

### Attendant & User Tilt Full Length Arm Pads only. Select trigger mounting locations below.

- 2723031 Right Back Post  2723032 Left Back Post \$ 195
- 2723041 Right Armrest  2723042 Left Armrest

- 2723110 **Tilt Trigger Extension** \$ 50

See next page for more options >

<b>Power Tilt</b> 250 lbs Weight Capacity. 20° tilt range. Includes 2 batteries for off-chair charging. <a href="http://pdgmobility.com/power-tilt">pdgmobility.com/power-tilt</a>		
Minimum Seat Width is 15". Power Tilt is compatible with either O2 Holder or IV Pole on the same chair, not both. Not available for 12" STF Height.		
<input type="checkbox"/>	2723400 Standard Power Tilt	\$ 1,795
<input type="checkbox"/>	2723410 Power Tilt with Manual Tilt Override Allows for instant manual release of tilt in case of power tilt failure.	\$ 2,295
Select Control Switch	Toggle Switch <input type="checkbox"/> 2723421 Installed on Right side <input type="checkbox"/> 2723422 Installed on Left side	
	Button Switch <input type="checkbox"/> 2723431 Installed on Right side <input type="checkbox"/> 2723432 Installed on Left side	
	Stealth Egg Switch <input type="checkbox"/> 2723440 Not installed	
	Optional for Egg Switch: <input type="checkbox"/> 2723445 Swing Away Mounting Hardware Not compatible with all Headrest makes and models.	\$ 625
Specify Headrest to be used: _____		

Recline	
<input type="checkbox"/>	2723310 <b>Reclining Backrest (Manual)</b> \$ 895
30° maximum recline. Not compatible with Backrest Extension, Fold-down Back, Tension Adjustable Upholstery. Requires longer Frame. Maximum Seat Depth is 19". Finished Back Post Height is 3" taller than selected Back Post Height.	

Front Rigging	
<input type="checkbox"/>	2725000 <b>No Front Rigging</b>
<b>Easy-Out™ Front Rigging</b> Effortless insertion and removal. No Hanger Bracket remains on frame. <a href="http://pdgmobility.com/easy-out">pdgmobility.com/easy-out</a>	
<input type="checkbox"/>	2725410 Fixed 70° Not available for 12" STF Height. \$ 95
<input type="checkbox"/>	2725420 Universal Angle Adjustable (0°-110°) Set at angle: _____ (use increments of 10°) 12" STF Height is only compatible with 0°-60°. \$ 295
Adjust angle by loosening one bolt. Tool-free leg-length adjustment.	
<input type="checkbox"/>	2725110 <b>Standard Plug In 70° Front Rigging</b> Legacy Design. Not compatible with Easy-Out.

Footplates	
<input type="checkbox"/>	2726010 <b>Composite Footplates</b>
<input type="checkbox"/>	2726020 <b>Aluminum Angle Adjustable Footplates</b> \$ 165
<input type="checkbox"/>	2726030 <b>One Piece Fixed Footplate</b> Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable. \$ 350
<input type="checkbox"/>	2726040 <b>One Piece Flip-up Footplate</b> Requires selection of Composite or Aluminum Angle Adjustable. \$ 395
<input type="checkbox"/>	2726050 <b>Interlocking Angle Adjustable Footplates</b> Available for Seat Widths up to 21" only. \$ 195

Footrest/Legrest Accessories	
<input type="checkbox"/>	2727012 <b>Calf Pads</b> Includes reinforced premium Calf Pad mount. \$ 200
<input type="checkbox"/>	2727032 <b>Gel Knee Pads</b> Compatible with Universal Angle Adjustable Front Rigging only. \$ 200
<input type="checkbox"/>	2727040 <b>Length Adjustable Heel Loops</b> \$ 50
<input type="checkbox"/>	2727050 <b>Bodypoint Padded Calf Strap</b> \$ 195

Front Rigging Extension Tubes	
<b>Regular Extension Tubes</b>	
<input type="checkbox"/>	2728010 Short 14"-18" for Easy-Out Front Rigging 13"-17" for Standard Plug In Front Rigging
<input type="checkbox"/>	2728020 Long 18"-24" for Easy-Out Front Rigging 15"-19" for Standard Plug In Front Rigging \$ 75
	<b>Lower Leg Length</b> If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.
Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.	

Miscellaneous Options	
<b>Spoke Guards</b> Made from durable and lightweight aluminum. \$ 195	
<input type="checkbox"/>	2730020 Black <input type="checkbox"/> 2730021 Frame Color

See next page for more options >




# BENTLEY LT-R

**Order Form - Canada**  
Prices in Canadian Dollars. January 1, 2024

<input type="checkbox"/>	2730050	<b>Protective Lower Padding</b>				\$	795
		<b>Stealth Flip Down Abductor</b>	Requires Solid Seat Pan.			\$	650
<input type="checkbox"/>	2730060	Small Pad	3"x4"x3"x2" size	<input type="checkbox"/>	2730062	Large Pad	5"x5"x4"x3" size
<input type="checkbox"/>	2730080	<b>Gel Hip Guides</b>	Reduces available Seat Width by 2".			\$	495
<input type="checkbox"/>	2730110	<b>IV Pole</b>				\$	295
<input type="checkbox"/>	2730120	<b>O2 Holder</b>				\$	395
<input type="checkbox"/>	2730200	<b>Q'Strain Lap Belt</b>	Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" Belt Length.			\$	245
<input type="checkbox"/>	2730222	<b>Wheelchair Backpack</b>	Black with PDG Logo. 17"x14" size.			\$	145
<input type="checkbox"/>	2730170	<b>Adhesive Back Hook Strip, 2" x 30"</b>	For Cushion mounting.			\$	10

**On-Chair Seating & Positioning**

 Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. [pdgmobility.com/seating](https://pdgmobility.com/seating)

**Order Notes**

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.