



# BENTLEY LT

**Order Form - Canada**  
Prices in Canadian Dollars. January 1, 2024

## Customer and Order Information

Purchase Order *	Tag	v.0401
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in
		<i>*required</i>

## Base Chair

**Bentley LT 20° Manual Tilt Wheelchair**

<input checked="" type="checkbox"/> 2401000	Bentley 20° Manual Tilt Wheelchair	\$ 4,300
<input checked="" type="checkbox"/> 2401009	Lightweight Aluminum Frame Upgrade Aluminum Frame is not recommended for agitated users.	\$ 300
<b>i</b> Choose the standard <b>Bentley</b> with a steel frame for agitated or heavier users. See separate order form.		

## Weight Capacity

2402000 **250 lbs Chair Weight Capacity**

**User Weight Range** Required for setup of tilt uplift force.

<input type="checkbox"/> 2402430	Up to 75 lbs	<input type="checkbox"/> 2402431	75 lbs to 100 lbs	<input type="checkbox"/> 2402432	100 lbs to 115 lbs
<input type="checkbox"/> 2402433	115 lbs to 150 lbs	<input type="checkbox"/> 2402434	150 lbs to 200 lbs	<input type="checkbox"/> 2402435	200 lbs to 250 lbs

## Frame Options

2402300 **Transit Tie Down System (TTDS)**

TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details.

## Seat Width

Overall chair width is +9.5" with 12/16" Rear Wheels and +10.5" with 20/22/24" Rear Wheels.

**Standard Seat Widths**

<input type="checkbox"/> 2403014	14" Seat Width	<input type="checkbox"/> 2403015	15" Seat Width	<input type="checkbox"/> 2403016	16" Seat Width
<input type="checkbox"/> 2403017	17" Seat Width	<input type="checkbox"/> 2403018	18" Seat Width	<input type="checkbox"/> 2403019	19" Seat Width
<input type="checkbox"/> 2403020	20" Seat Width	<input type="checkbox"/> 2403021	21" Seat Width	<input type="checkbox"/> 2403022	22" Seat Width

## Seat Depth

**Standard Seat Depths (Adjustable)**

<input type="checkbox"/> 2404016	16" Seat Depth	<input type="checkbox"/> 2404017	17" Seat Depth	<input type="checkbox"/> 2404018	18" Seat Depth
<input type="checkbox"/> 2404019	19" Seat Depth	<input type="checkbox"/> 2404020	20" Seat Depth		
<b>Custom Seat Depths</b>					
<input type="checkbox"/> 2404021	21" Seat Depth	<input type="checkbox"/> 2404022	22" Seat Depth	\$ 400	

## Frame Color

**PDG Premium Colors** View color swatches on PDG's website. [pdgmobility.com/colors](http://pdgmobility.com/colors)

<input type="checkbox"/> 2405001	Black Shimmer	<input type="checkbox"/> 2405015	Sky Blue	<input type="checkbox"/> 2405016	Poppy Red
<input type="checkbox"/> 2405003	Fine Texture Black	<input type="checkbox"/> 2405010	Lollypop Blue	<input type="checkbox"/> 2405011	Lollypop Red
<input type="checkbox"/> 2405017	Arctic Grey	<input type="checkbox"/> 2405018	Tropical Teal	<input type="checkbox"/> 2405019	Tangerine Orange
<input type="checkbox"/> 2405020	Olive Green	<input type="checkbox"/> 2405021	Lush Lilac	<input type="checkbox"/> 2405012	Bubble Gum Pink
<input type="checkbox"/> 2405022	Leaf Green	<input type="checkbox"/> 2405023	Boysenberry Purple	<input type="checkbox"/> 2405014	Daffodil Yellow
<input type="checkbox"/> 2405200	<b>Custom Color</b> Contact PDG for specifications.				\$ 300

### Armrests

2406000 **No Armrests**

**Height Adjustable Armrests (Single Post)** \$ 295

2406111 8"-11" Range       2406114 10"-14" Range       2406116 12"-16" Range

2406100 **Fixed Height Armrests (12")** Full Length Armpad only.

**Armrest Height** If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.  
Seat Pan to top of Armrest: \_\_\_\_\_ Does not include height of Armpad

**Armrest Position** Indicate preference for Armrest position. If left blank, the default setting will be used.  
In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

2406401 Rear       2406402 Mid       2406403 Front

### Armpads

**Regular Armpads**

2407100 Desk Length - 2" x 10"       2407101 Full Length - 2" x 14"

**Waterfall Armpads** \$ 115

2407105 Desk Length - 2" x 10"       2407106 Full Length - 2" x 14"

**Wide Armpads** May interfere with wheels in some tilt configurations. \$ 200

2407110 Desk Length - 3" x 10"       2407111 Full Length - 3" x 12"

**Gel Armpads** May interfere with wheels in some tilt configurations. \$ 135

2407120 Full Length - 2" x 12"       2407121 Full Length - 3.5" x 12"

**Contour Arm Trough with Hand Pad** May interfere with wheels in some tilt configurations. \$ 125

2407130 Right (4.5" x 18")       2407131 Left (4.5" x 18")

### Armrest Options

**Armrest Wideners** Move Armrests outward by up to 1" on each side. Seat size does not change. May cause tilt interference.

2407210 1" Total Additional Width \$ 95

2407211 2" Total Additional Width \$ 195

2407310 **Clothing Guards** Not compatible with 8"-11" Height range. \$ 110

**Armrest Extensions** Not compatible with Desk Length Armpads. \$ 75

2407400 3" Extension - Forward       2407401 3" Extension - Rearward

### Seat Pan

2408000 **No Seat Pan**

2408001 **Solid Aluminum Seat Pan** \$ 135

### Seat Height

**Don't forget** **Select desired Seat-to-Floor Height** STF Height is adjustable within the range of the chosen Subframe Plate (Low/High).

Casters	Low						High						
	RW	13"	14"	15"	16"	17"	RW	15"	16"	17"	18"	19"	20"
6" Casters	12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
8" Casters	RW	13"	14"	15"	16"	17"	RW	15"	16"	17"	18"	19"	20"
	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

## Casters

### Urethane Casters

8" Casters are compatible only with 16" and 24" Rear Wheels.

<input type="checkbox"/>	2411061	6" x 1.25" Caster size	<input type="checkbox"/>	2411081	8" x 1.25" Caster size	
<input type="checkbox"/>	2411063	6" x 1.50" Caster size				\$ 85
<input type="checkbox"/>	2411084	8" x 2.00" Caster size				\$ 135

### Pneumatic Casters

8" Casters are compatible only with 16" and 24" Rear Wheels.

<input type="checkbox"/>	2411161	6" x 1.25" Caster size				\$ 110
<input type="checkbox"/>	2411181	8" x 1.25" Caster size				\$ 120
<input type="checkbox"/>	2411184	8" x 2.00" Caster size				\$ 150

<input type="checkbox"/>	2411260	<b>Frog Legs 6" x 1.50" Aluminum Soft-Roll Casters</b>				\$ 275
--------------------------	---------	--	--	--	--	--------

<input type="checkbox"/>	2411400	<b>Frog Legs Suspension Forks</b>	Only compatible with 16 and 24" Rear Wheels. Requires 6" Soft-Roll Caster.			\$ 595
--------------------------	---------	-----------------------------------	--	--	--	--------

## Rear Wheels

### Composite Wheels

<input type="checkbox"/>	2412012	12" Wheel size	<input type="checkbox"/>	2412016	16" Wheel size	<input type="checkbox"/>	2412020	20" Wheel size
<input type="checkbox"/>	2412022	22" Wheel size	<input type="checkbox"/>	2412024	24" Wheel size			

### Spoke Wheels

<input type="checkbox"/>	2412120	20" Wheel size	<input type="checkbox"/>	2412122	22" Wheel size	<input type="checkbox"/>	2412124	24" Wheel size	\$ 75
--------------------------	---------	----------------	--------------------------	---------	----------------	--------------------------	---------	----------------	-------



### Rear Wheel Position

Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.

In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

<input type="checkbox"/>	2412401	Rear	<input type="checkbox"/>	2412402	Mid	<input type="checkbox"/>	2412403	Front
--------------------------	---------	------	--------------------------	---------	-----	--------------------------	---------	-------

## Axle Type

<input type="checkbox"/>	2413000	<b>Quick Release Rear Axles</b>
<input type="checkbox"/>	2413010	<b>Fixed Rear Axles</b>

## Rear Tires

<input type="checkbox"/>	2414000	<b>Urethane Tires</b>		
<input type="checkbox"/>	2414010	<b>Pneumatic Tires</b>	May result in a finished Seat Height 1/2" greater than stated.	\$ 60
<input type="checkbox"/>	2414020	<b>Pneumatic Tires with Flat-Free Inserts</b>	May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels.	\$ 95

## Handrims

Handrims are only available with 20", 22", and 24" Rear Wheels.

<input type="checkbox"/>	2415000	<b>No Handrims</b>		
<input type="checkbox"/>	2415010	<b>Standard Aluminum Handrims</b>		
<input type="checkbox"/>	2415020	<b>Plastic Coated Handrims</b>	Not compatible with Drum Brakes.	\$ 150
<input type="checkbox"/>	2415030	<b>Projection Handrims</b>	Composite Wheels only. 10 Projections, 30°, 1-1.8" Length.	\$ 495

## Wheel Locks & Brakes

### User Wheel Locks

<input type="checkbox"/>	2416010	Push to Lock	<input type="checkbox"/>	2416020	Pull to Lock	
<input type="checkbox"/>	2416030	<b>Attendant Wheel Locks (Push to Lock)</b>				\$ 85
<input type="checkbox"/>	2416040	<b>User &amp; Attendant Wheel Locks</b>	Not compatible with Drum Brakes.			\$ 125
<input type="checkbox"/>	2416300	<b>Drum Brakes</b>	Available with 12" Composite or 20" Spoke Rear Wheels. Urethane Tires only.			\$ 795
		<b>User-activated Unilateral Wheel Lock</b>	Pull to Lock.			\$ 795
<input type="checkbox"/>	2416310	Left Activation	<input type="checkbox"/>	2416311	Right Activation	



Wheel Lock Options			
<input type="checkbox"/>	2416100	<b>Wheel Lock Extensions</b> Not compatible with Attendant Wheel Lock.	\$ 60
<input type="checkbox"/>	2416200	<b>Hill Holders</b> Requires Pneumatic Rear Tires with Flat-Free Inserts.	\$ 450

Anti-Tippers			
<input checked="" type="checkbox"/>	2417010	<b>Rear Anti-Tippers</b> Required Safety Option.	\$ 75

Back Posts			
<b>Back Post Style</b>			
<input checked="" type="checkbox"/>	2418010	Stroller Handle Back	
<b>Back Post Height</b> Measured from the Seat Pan to the Handle.			
<input type="checkbox"/>	2418425	Short (22")	<input type="checkbox"/> 2418431 Tall (25")
<b>Back Post Folding</b> Fold-down Back is not compatible with Reclining Backrest.			
<input checked="" type="checkbox"/>	2419220	Fold-down	
<b>Seat to Back Angle</b>			
<input type="checkbox"/>	2419090	Set at 90°	<input type="checkbox"/> 2419098 Set at 98° <input type="checkbox"/> 2419106 Set at 106°
<input type="checkbox"/>	2419114	Set at 114°	<input type="checkbox"/> 2419122 Set at 122°

Back Options			
<b>Angle/Height Adjustable Stroller Handle</b> Not compatible with Headrest Mounting Plate.			\$ 150
<input type="checkbox"/>	2419010	Button on Outside	<input type="checkbox"/> 2419011 Button on Inside
<b>Endomorph Back</b> Backrest will be narrower than Seat by total of 2x the Offset selected below. Minimum width of Backrest is 12".			\$ 395
<input type="checkbox"/>	2418102	1" Offset per side	<input type="checkbox"/> 2418103 1.5" Offset per side <input type="checkbox"/> 2418104 2" Offset per side
<input type="checkbox"/>	2419210	<b>10" Backrest Extension</b> Not compatible with Angle/Height Adjustable Stroller Handle.	\$ 150
<input type="checkbox"/>	2419212	<b>Headrest Mounting Plate</b>	\$ 95
Not compatible with Backrest Extension, Angle/Height Adjustable Stroller Handle. Not required if mounting to a solid Back.			

Back Upholstery																																																																									
<input type="checkbox"/>	2420000	<b>No Back Upholstery</b>																																																																							
<b>Standard Upholstery</b>			\$ 195																																																																						
<input type="checkbox"/>	2420418	Short (18" Height)	<input type="checkbox"/> 2420421 Regular (21" Height) <input type="checkbox"/> 2420424 Tall (24" Height)																																																																						
<input type="checkbox"/>	2420450	<b>Tension Adjustable Upholstery</b> Height of Upholstery is 16".	\$ 295																																																																						
<b>Position of Upholstery</b> Select the desired position of the top of the Upholstery relative to the Seat Pan.																																																																									
<table border="1"> <thead> <tr> <th rowspan="2">Style</th> <th colspan="10">22" Back Post</th> </tr> <tr> <th>17"</th> <th>18"</th> <th>19"</th> <th>20"</th> <th>21"</th> <th>22"</th> <th>23"</th> <th>24"</th> <th>25"</th> </tr> </thead> <tbody> <tr> <td>Short</td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Regular</td> <td></td> <td></td> <td></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Tension Adj</td> <td></td> <td></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td></td> <td></td> <td></td> <td></td> </tr> </tbody> </table>										Style	22" Back Post										17"	18"	19"	20"	21"	22"	23"	24"	25"	Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>					Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>					Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>															
Style	22" Back Post																																																																								
	17"	18"	19"	20"	21"	22"	23"	24"	25"																																																																
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																																																			
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																																																			
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																																																			
<table border="1"> <thead> <tr> <th rowspan="2">Style</th> <th colspan="10">25" Back Post</th> </tr> <tr> <th>17"</th> <th>18"</th> <th>19"</th> <th>20"</th> <th>21"</th> <th>22"</th> <th>23"</th> <th>24"</th> <th>25"</th> </tr> </thead> <tbody> <tr> <td>Short</td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td>Regular</td> <td></td> <td></td> <td></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td>Tall</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td>Tension Adj</td> <td></td> <td></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </tbody> </table>										Style	25" Back Post										17"	18"	19"	20"	21"	22"	23"	24"	25"	Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Tall							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Style	25" Back Post																																																																								
	17"	18"	19"	20"	21"	22"	23"	24"	25"																																																																
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																																															
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																																															
Tall							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																																															
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																																															

Tilt Activation			
<b>Attendant Tilt only</b> Select trigger mounting location below.			
<input type="checkbox"/>	2423011	Right Back Post	<input type="checkbox"/> 2423012 Left Back Post
<b>User Tilt only</b> Full Length Arm Pads only. Select trigger mounting location below.			\$ 100
<input type="checkbox"/>	2423021	Right Armrest	<input type="checkbox"/> 2423022 Left Armrest
<b>Attendant &amp; User Tilt</b> Full Length Arm Pads only. Select trigger mounting locations below.			\$ 195
<input type="checkbox"/>	2423031	Right Back Post	<input type="checkbox"/> 2423032 Left Back Post
<input type="checkbox"/>	2423041	Right Armrest	<input type="checkbox"/> 2423042 Left Armrest

See next page for more options >



<input type="checkbox"/>	2423110	<b>Tilt Trigger Extension</b>			\$ 50
<p><b>Power Tilt</b> 250 lbs Weight Capacity. 20° tilt range. Includes 2 batteries for off-chair charging. <a href="http://pdgmobility.com/power-tilt">pdgmobility.com/power-tilt</a>          Minimum Seat Width is 15". Power Tilt is compatible with either O2 Holder or IV Pole on the same chair, not both.</p>					
<input type="checkbox"/>	2423400	Standard Power Tilt			\$ 1,795
<input type="checkbox"/>	2423410	Power Tilt with Manual Tilt Override	Allows for instant manual release of tilt in case of power tilt failure.		\$ 2,295
		Toggle Switch	<input type="checkbox"/> 2423421 Installed on Right side	<input type="checkbox"/> 2423422 Installed on Left side	
		Button Switch	<input type="checkbox"/> 2423431 Installed on Right side	<input type="checkbox"/> 2423432 Installed on Left side	
		Stealth Egg Switch	<input type="checkbox"/> 2423440 Not installed		
		Optional for Egg Switch:	<input type="checkbox"/> 2423445 Swing Away Mounting Hardware	Not compatible with all Headrest makes and models.	\$ 625
Specify Headrest to be used: _____					

Recline					
<input type="checkbox"/>	2423310	<b>Reclining Backrest (Manual)</b>			\$ 895
30° maximum recline. Not compatible with Backrest Extension, Fold-down Back, Tension Adjustable Upholstery, Dynamic Rocker Back. Requires longer Frame. Maximum Seat Depth is 19". Finished Back Post Height is 3" taller than selected Back Post Height.					


Front Rigging					
<input type="checkbox"/>	2425000	<b>No Front Rigging</b>		<input type="checkbox"/> 2425001 <b>Sector Blocks only</b>	
<p><b>Easy-Out™ Front Rigging</b> Effortless insertion and removal. No Hanger Bracket remains on frame. <a href="http://pdgmobility.com/easy-out">pdgmobility.com/easy-out</a></p>					
<input type="checkbox"/>	2425410	Fixed 70°			\$ 95
<input type="checkbox"/>	2425420	Universal Angle Adjustable (0°-110°)	Set at angle: _____ (use increments of 10°)		\$ 295
Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings.					
<b>Swing Away Front Rigging</b> Legacy Design. Not compatible with Easy-Out.					
<input type="checkbox"/>	2425130	Swing Away 60°	<input type="checkbox"/> 2425140 Swing Away 70°		
<input type="checkbox"/>	2425200	Swing Away Elevating Legrests	Includes Calf Pads.		\$ 375
<input type="checkbox"/>	2425300	1" Offset Hangers for Swing Away			\$ 75
<input type="checkbox"/>	2425110	<b>Standard Plug In 70° Front Rigging</b>	Legacy Design. Not compatible with Easy-Out.		

Footplates					
<input type="checkbox"/>	2426010	<b>Composite Footplates</b>			
<input type="checkbox"/>	2426020	<b>Aluminum Angle Adjustable Footplates</b>			\$ 165
<input type="checkbox"/>	2426030	<b>One Piece Fixed Footplate</b>	Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.		\$ 350
<input type="checkbox"/>	2426040	<b>One Piece Flip-up Footplate</b>	Requires selection of Composite or Aluminum Angle Adjustable.		\$ 395
<input type="checkbox"/>	2426050	<b>Interlocking Angle Adjustable Footplates</b>	Available for Seat Widths up to 21" only.		\$ 195

Footrest/Legrest Accessories					
<b>Calf Pads</b> Not compatible with Gel Padding.					
<input type="checkbox"/>	2427012	For Easy-Out™ Front Rigging	Includes reinforced premium Calf Pad mount.		\$ 200
<input type="checkbox"/>	2427010	For Swing Away Front Rigging			\$ 135
<input type="checkbox"/>	2427020	<b>Gel Padding</b>	Not compatible with Calf Pad or Easy-Out Front Rigging.		\$ 125
<b>Gel Knee Pads</b>					
<input type="checkbox"/>	2427032	For Easy-Out™ Angle Adjustable Front Rigging			\$ 200
<input type="checkbox"/>	2427030	For Swing Away ELR Front Rigging			\$ 100
<input type="checkbox"/>	2427040	<b>Length Adjustable Heel Loops</b>			\$ 50
<input type="checkbox"/>	2427050	<b>Bodypoint Padded Calf Strap</b>			\$ 195

Front Rigging Extension Tubes					
<b>Regular Extension Tubes</b>					
<input type="checkbox"/>	2428010	Short	14"-18" for Easy-Out Front Rigging	13"-17" for Swing Away and Standard Plug In Front Rigging	
<input type="checkbox"/>	2428020	Long	18"-24" for Easy-Out Front Rigging	15"-19" for Swing Away and Standard Plug In Front Rigging	\$ 75

See next page for more options >

	<p><b>Lower Leg Length</b> If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.</p> <p>Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.</p>
--	---

## Miscellaneous Options

<p><b>One-Arm Drive</b> Available with 22" and 24" Composite Rear Wheels. Minimum 17" STF and 18" STF, respectively. Limits Tilt. Not compatible with User &amp; Attendant Wheel Locks, Drum Brakes, Quick Release Rear Axles, Pneumatic Tires. <span style="float: right;">\$ 3,500</span></p> <p><input type="checkbox"/> 2430010 Installed on Right side      <input type="checkbox"/> 2430011 Installed on Left side</p> <p><input type="checkbox"/> 2430012 Large Hand Setup 24" Rear Wheels with 22" Driving Handrim.</p>	
<p><b>Spoke Guards</b> Available for 20"-24" Composite Rear Wheels. Made from durable and lightweight aluminum. <span style="float: right;">\$ 195</span></p> <p><input type="checkbox"/> 2430020 Black      <input type="checkbox"/> 2430021 Frame Color</p>	
<p><input type="checkbox"/> 2430050 <b>Protective Lower Padding</b> <span style="float: right;">\$ 795</span></p>	
<p><b>Stealth Flip Down Abductor</b> Requires Solid Seat Pan. <span style="float: right;">\$ 650</span></p> <p><input type="checkbox"/> 2430060 Small Pad 3"x4"x3"x2" size      <input type="checkbox"/> 2430062 Large Pad 5"x5"x4"x3" size</p>	
<p><input type="checkbox"/> 2430080 <b>Gel Hip Guides</b> Reduces available Seat Width by 2". <span style="float: right;">\$ 495</span></p>	
<p><input type="checkbox"/> 2430110 <b>IV Pole</b> Not compatible with Attendant Wheel Locks. Requires frame modification and cannot be added to chair later. <span style="float: right;">\$ 295</span></p>	
<p><input type="checkbox"/> 2430120 <b>O2 Holder</b> Not compatible with Attendant Wheel Locks. Requires frame modification and cannot be added to chair later. <span style="float: right;">\$ 395</span></p>	
<p><input type="checkbox"/> 2430200 <b>Q'Straint Lap Belt</b> Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" Belt Length. <span style="float: right;">\$ 245</span></p>	
<p><input type="checkbox"/> 2430222 <b>Wheelchair Backpack</b> Black with PDG Logo. 17"x14" size. <span style="float: right;">\$ 145</span></p>	
<p><input type="checkbox"/> 2430170 <b>Adhesive Back Hook Strip, 2" x 30"</b> For Cushion mounting. <span style="float: right;">\$ 10</span></p>	

## On-Chair Seating & Positioning

 Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. [pdgmobility.com/seating](http://pdgmobility.com/seating)

## Order Notes

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.