



### Customer and Order Information

Purchase Order *	Tag	v.0401
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in

### Base Chair

<input checked="" type="checkbox"/> 1001000	<b>Bentley 20° Manual Tilt Wheelchair</b>	\$ 4,300
<input type="checkbox"/> ⓘ	Also available: the lightweight <b>Bentley LT</b> with aluminum frame. See separate order form. <a href="http://pdgmobility.com/blt">pdgmobility.com/blt</a>	

### Weight Capacity

<b>Chair Weight Capacity</b>		
<input type="checkbox"/> 1002000	250 lbs Capacity Available for 14" to 22" Seat Width.	
<input type="checkbox"/> 1002050	300 lbs Capacity Available for 14" to 22" Seat Width.	\$ 200
<input type="checkbox"/> 1005006	450 lbs Capacity Available for 20" to 32" Seat Width.	\$ 500
<b>User Weight Range</b> Required for setup of tilt uplift force.		
<input type="checkbox"/> 1002430	Up to 75 lbs	
<input type="checkbox"/> 1002431	75 lbs to 100 lbs	
<input type="checkbox"/> 1002432	100 lbs to 115 lbs	
<input type="checkbox"/> 1002433	115 lbs to 150 lbs	
<input type="checkbox"/> 1002434	150 lbs to 200 lbs	
<input type="checkbox"/> 1002435	200 lbs to 250 lbs	
<input type="checkbox"/> 1002436	250 lbs to 300 lbs	
<input type="checkbox"/> 1002437	300 lbs to 375 lbs	
<input type="checkbox"/> 1002438	375 lbs to 450 lbs	

### Frame Options

<input type="checkbox"/> 1002200	<b>Extended Base with Rear Floor Jacks</b>	\$ 900
Lifts Rear Wheels off ground for increased stability. Limited ground clearance. Not compatible with Attendant Wheel Lock, HD Anti-Tipper, IV Pole, O2 Holder.		
<input type="checkbox"/> 1002300	<b>Transit Tie Down System (TTDS)</b>	
TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Not compatible with Seat Widths greater than 22". Tested to max. Weight Capacity of 300 lbs.		

### Seat Width

Overall chair width is +9.5" with 12/16" Rear Wheels and +10.5" with 20/22/24" Rear Wheels.

<b>Standard Seat Widths</b>		
<input type="checkbox"/> 1003014	14" Seat Width	
<input type="checkbox"/> 1003015	15" Seat Width	
<input type="checkbox"/> 1003016	16" Seat Width	
<input type="checkbox"/> 1003017	17" Seat Width	
<input type="checkbox"/> 1003018	18" Seat Width	
<input type="checkbox"/> 1003019	19" Seat Width	
<input type="checkbox"/> 1003020	20" Seat Width	
<input type="checkbox"/> 1003021	21" Seat Width	
<input type="checkbox"/> 1003022	22" Seat Width	
<b>Custom Seat Widths</b> Requires 450 lbs Weight Capacity Frame. Extended delivery lead times may apply.		
<input type="checkbox"/> 1003023	23" Seat Width	\$ 375
<input type="checkbox"/> 1003024	24" Seat Width	\$ 375
<input type="checkbox"/> 1003025	25" Seat Width	\$ 700
<input type="checkbox"/> 1003026	26" Seat Width	\$ 700
<input type="checkbox"/> 1003027	27" Seat Width	
<input type="checkbox"/> 1003028	28" Seat Width	
<input type="checkbox"/> 1003029	29" Seat Width	
<input type="checkbox"/> 1003030	30" Seat Width	
<input type="checkbox"/> 1003031	31" Seat Width	
<input type="checkbox"/> 1003032	32" Seat Width	

### Seat Depth

<b>Standard Seat Depths (Adjustable)</b>		
<input type="checkbox"/> 1004016	16" Seat Depth	
<input type="checkbox"/> 1004017	17" Seat Depth	
<input type="checkbox"/> 1004018	18" Seat Depth	
<input type="checkbox"/> 1004019	19" Seat Depth	
<input type="checkbox"/> 1004020	20" Seat Depth	
<b>Custom Seat Depths</b>		
<input type="checkbox"/> 1004021	21" Seat Depth	
<input type="checkbox"/> 1004022	22" Seat Depth	\$ 400

Frame Color			
<b>PDG Premium Colors</b> View color swatches on PDG's website. <a href="http://pdgmobility.com/colors">pdgmobility.com/colors</a>			
<input type="checkbox"/> 1005001	Black Shimmer	<input type="checkbox"/> 1005015	Sky Blue
<input type="checkbox"/> 1005003	Fine Texture Black	<input type="checkbox"/> 1005010	Lollypop Blue
<input type="checkbox"/> 1005017	Arctic Grey	<input type="checkbox"/> 1005018	Tropical Teal
<input type="checkbox"/> 1005020	Olive Green	<input type="checkbox"/> 1005021	Lush Lilac
<input type="checkbox"/> 1005022	Leaf Green	<input type="checkbox"/> 1005023	Boysenberry Purple
<input type="checkbox"/> 1005016	Poppy Red	<input type="checkbox"/> 1005011	Lollypop Red
		<input type="checkbox"/> 1005019	Tangerine Orange
		<input type="checkbox"/> 1005012	Bubble Gum Pink
		<input type="checkbox"/> 1005014	Daffodil Yellow
<input type="checkbox"/> 1005200	<b>Custom Color</b>	Contact PDG for specifications.	\$ 300

Armrests			
<input type="checkbox"/> 1006000	<b>No Armrests</b>		
<b>Single Post Height Adjustable Armrests</b>			\$ 295
<input type="checkbox"/> 1006111	8"-11" Range	<input type="checkbox"/> 1006114	10"-14" Range
		<input type="checkbox"/> 1006116	12"-16" Range
<b>Dual Post Height Adjustable Armrests</b>			\$ 295
<input type="checkbox"/> 1006211	8"-11" Range	<input type="checkbox"/> 1006214	10"-14" Range
		<input type="checkbox"/> 1006216	12"-16" Range
<b>Fixed Height Armrests (12")</b> Full Length Armpad only.			
<input type="checkbox"/> 1006100	Single Post	<input type="checkbox"/> 1006200	Dual Post
	<b>Armrest Height</b>	If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.	
	Seat Pan to top of Armrest: _____	Does not include height of Armpad	
	<b>Armrest Position</b>	Indicate preference for Armrest position. If left blank, the default setting will be used.	
In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.			
<input type="checkbox"/> 1006401	Rear	<input type="checkbox"/> 1006402	Mid
		<input type="checkbox"/> 1006403	Front
<b>Locking Flip-up Armrests</b> Not recommended for users with excessive strength or tone.			\$ 295
Not compatible with Fold-down Back, Clothing Guards, Back Upholstery. Available with Full Length and Desk Length Armpads.			
<input type="checkbox"/> 1006310	Set at 10" Height	<input type="checkbox"/> 1006311	Set at 11" Height
		<input type="checkbox"/> 1006312	Set at 12" Height
<input type="checkbox"/> 1006313	Set at 13" Height	<input type="checkbox"/> 1006314	Set at 14" Height

Armpads			
<b>Regular Armpads</b>			
<input type="checkbox"/> 1007100	Desk Length - 2" x 10"	<input type="checkbox"/> 1007101	Full Length - 2" x 14"
<b>Waterfall Armpads</b> Not compatible with Dual Post Armrests			\$ 115
<input type="checkbox"/> 1007105	Desk Length - 2" x 10"	<input type="checkbox"/> 1007106	Full Length - 2" x 14"
<b>Wide Armpads</b> May interfere with wheels in some tilt configurations.			\$ 200
<input type="checkbox"/> 1007110	Desk Length - 3" x 10"	<input type="checkbox"/> 1007111	Full Length - 3" x 12"
<b>Gel Armpads</b> May interfere with wheels in some tilt configurations.			\$ 135
<input type="checkbox"/> 1007120	Full Length - 2" x 12"	<input type="checkbox"/> 1007121	Full Length - 3.5" x 12"
<b>Contour Arm Trough with Hand Pad</b> May interfere with wheels in some tilt configurations.			\$ 125
<input type="checkbox"/> 1007130	Right (4.5" x 18")	<input type="checkbox"/> 1007131	Left (4.5" x 18")

Armrest Options			
<b>Armrest Wideners</b> Move Armrests outward by up to 1" on each side. Seat size does not change. May cause tilt interference.			
<input type="checkbox"/> 1007210	1" Total Additional Width		\$ 95
<input type="checkbox"/> 1007211	2" Total Additional Width		\$ 195
<b>Clothing Guards</b> Not compatible with Flip-up Armrests.			
<input type="checkbox"/> 1007310	For Single Post Armrests	Not compatible with 8"-11" Height range.	\$ 110
<input type="checkbox"/> 1007320	For Dual Post Armrests		\$ 165
<b>Armrest Extensions</b> Not compatible with Flip-up Armrests, Desk Length Armpads.			\$ 75
<input type="checkbox"/> 1007400	3" Extension - Forward	<input type="checkbox"/> 1007401	3" Extension - Rearward

### Seat Pan

- 1008000 **No Seat Pan**
- 1008001 **Solid Aluminum Seat Pan** Required for 22" Seat Depth. \$ 135

### Seat Height

**Don't forget** **Select desired Seat-to-Floor Height** STF Height is adjustable within the range of the chosen Subframe Plate (Low/High).

	Low						High						
	RW	13"	14"	15"	16"	17"	RW	15"	16"	17"	18"	19"	20"
6" Casters	12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
8" Casters	RW	13"	14"	15"	16"	17"	RW	15"	16"	17"	18"	19"	20"
	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

### Casters

- Urethane Casters** 8" Casters are compatible only with 16" and 24" Rear Wheels.
- 1011061 6" x 1.25" Caster size  1011081 8" x 1.25" Caster size
  - 1011063 6" x 1.50" Caster size \$ 85
  - 1011084 8" x 2.00" Caster size \$ 135
- Pneumatic Casters** 8" Casters are compatible only with 16" and 24" Rear Wheels.
- 1011161 6" x 1.25" Caster size \$ 110
  - 1011181 8" x 1.25" Caster size \$ 120
  - 1011184 8" x 2.00" Caster size \$ 150
  - 1011260 **Frog Legs 6" x 1.50" Aluminum Soft-Roll Casters** \$ 275
  - 1011400 **Frog Legs Suspension Forks** Only compatible with 16 and 24" Rear Wheels. Requires 6" Soft-Roll Caster. \$ 595

### Rear Wheels

- Composite Wheels**
- 1012012 12" Wheel size  1012016 16" Wheel size  1012020 20" Wheel size
  - 1012022 22" Wheel size  1012024 24" Wheel size
- Spoke Wheels** \$ 75
- 1012120 20" Wheel size  1012122 22" Wheel size  1012124 24" Wheel size
  - 1012224 **24" Heavy Duty Spoke Wheels** \$ 350
- Rear Wheel Position** Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.  
In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.
- 1012401 Rear  1012402 Mid  1012403 Front

### Axle Type

- 1013000 **Quick Release Rear Axles**
- 1013010 **Fixed Rear Axles**

### Rear Tires

- 1014000 **Urethane Tires**
- 1014010 **Pneumatic Tires** May result in a finished Seat Height 1/2" greater than stated. \$ 60
- 1014020 **Pneumatic Tires with Flat-Free Inserts** May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels. \$ 95



Handrims		Handrims are only available with 20", 22", and 24" Rear Wheels.
<input type="checkbox"/>	1015000 <b>No Handrims</b>	
<input type="checkbox"/>	1015010 <b>Standard Aluminum Handrims</b>	
<input type="checkbox"/>	1015020 <b>Plastic Coated Handrims</b> Not compatible with Drum Brakes.	\$ 150
<input type="checkbox"/>	1015030 <b>Projection Handrims</b> Composite Wheels only. 10 Projections, 30°, 1-1.8" Length.	\$ 495

Wheel Locks & Brakes			
<b>User Wheel Locks</b>			
<input type="checkbox"/>	1016010 Push to Lock	<input type="checkbox"/>	1016020 Pull to Lock
<input type="checkbox"/>	1016030 <b>Attendant Wheel Locks (Push to Lock)</b>		\$ 85
<input type="checkbox"/>	1016040 <b>User &amp; Attendant Wheel Locks</b> Not compatible with Drum Brakes.		\$ 125
<input type="checkbox"/>	1016300 <b>Drum Brakes</b> Available with 12" Composite or 20" Spoke Rear Wheels. Urethane Tires only.		\$ 795
<b>User-activated Unilateral Wheel Lock</b> Pull to Lock. Available up to 22" Seat Width only.			\$ 795
<input type="checkbox"/>	1016310 Left Activation	<input type="checkbox"/>	1016311 Right Activation

Wheel Lock Options		
<input type="checkbox"/>	1016100 <b>Wheel Lock Extensions</b> Not compatible with Attendant Wheel Lock.	\$ 60
<input type="checkbox"/>	1016200 <b>Hill Holders</b> Requires Pneumatic Rear Tires with Flat-Free Inserts.	\$ 450

Anti-Tippers		
<input checked="" type="checkbox"/>	1017010 <b>Rear Anti-Tippers</b> Required Safety Option.	\$ 75
<input type="checkbox"/>	1017020 <b>Front Anti-Tippers</b>	\$ 795
<input type="checkbox"/>	1017030 <b>Heavy Duty Rear Anti-Tippers</b> Provide additional rearward stability for users with excessive movement. Not compatible with Attendant Wheel Locks, Rear Floor Jacks, IV Pole, O2 Holder.	\$ 895

Back Posts					
<b>Back Post Style</b>					
<input type="checkbox"/>	1018010 Stroller Handle Back	<input type="checkbox"/>	1018020 Cane Handle Back	Seat Widths of up to 22" only. May reduce tilt range.	
<b>Back Post Height</b> Measured from the Seat Pan to the Handle.					
<input type="checkbox"/>	1018425 Short (22")	<input type="checkbox"/>	1018431 Tall (25")		
<b>Back Post Folding</b> Fold-down Back is not compatible with Dynamic Rocker Back Interface, Reclining Backrest.					
<input type="checkbox"/>	1019220 Fold-down	<input type="checkbox"/>	1019221 Non-folding		
<b>Seat to Back Angle</b>					
<input type="checkbox"/>	1019090 Set at 90°	<input type="checkbox"/>	1019098 Set at 98°	<input type="checkbox"/>	1019106 Set at 106°
<input type="checkbox"/>	1019114 Set at 114°	<input type="checkbox"/>	1019122 Set at 122°		

Back Options					
<b>Angle/Height Adjustable Stroller Handle</b> Not compatible with Headrest Mounting Plate or Cane Handle Back.		\$ 150			
<input type="checkbox"/>	1019010 Button on Outside	<input type="checkbox"/>	1019011 Button on Inside		
<b>Endomorph Back</b> Backrest will be narrower than Seat by total of 2x the Offset selected below. Minimum width of Backrest is 12".		\$ 395			
<input type="checkbox"/>	1018102 1" Offset per side	<input type="checkbox"/>	1018103 1.5" Offset per side	<input type="checkbox"/>	1018104 2" Offset per side
<input type="checkbox"/>	1019210 <b>10" Backrest Extension</b> Not compatible with Angle/Height Adjustable Stroller Handle.		\$ 150		
<input type="checkbox"/>	1019212 <b>Headrest Mounting Plate</b> Not compatible with Backrest Extension, Angle/Height Adjustable Stroller Handle, Cane Handle Back. Not required if mounting to a solid Back.		\$ 95		
<input checked="" type="checkbox"/>	1019230 <b>Spring Back Suspension</b> Not compatible with Reclining Back.				
<input type="checkbox"/>	1019240 <b>Seating Dynamics Dynamic Rocker Back Interface</b> Not compatible with Fold-down Back, Spring Back Suspension, 30° Recline, Cane Handle, Standard Upholstery. Finished Back Post Height is 4" taller than selected Back Post Height.		\$ 2,995		

## Back Upholstery

1020000 **No Back Upholstery**

**Standard Upholstery** Not available with Cane Handle Back. \$ 195

1020418 Short (18" Height)       1020421 Regular (21" Height)       1020424 Tall (24" Height)

**Tension Adjustable Upholstery** Height of Upholstery is 16". Not available with Cane Handle Back.

1020450 For 14" to 22" Seat Width \$ 295

1020451 For 24" to 32" Seat Width \$ 375

**Position of Upholstery** Select the desired position of the top of the Upholstery relative to the Seat Pan.

22" Back Post									
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

25" Back Post									
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tall							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

## Tilt Activation

**Attendant Tilt only** Select trigger mounting location below. Chairs with 450 lbs Weight Capacity require two tilt triggers (L and R) for tilt activation.

1023011 Right Back Post       1023012 Left Back Post

**User Tilt only** Not available for 450 lbs Weight Capacity. Full Length Arm Pads only. Select trigger mounting location below. \$ 100

1023021 Right Armrest       1023022 Left Armrest

**Attendant & User Tilt** Not available for 450 lbs Weight Capacity. Full Length Arm Pads only. Select trigger mounting locations below. \$ 195

1023031 Right Back Post       1023032 Left Back Post

1023041 Right Armrest       1023042 Left Armrest

1023110 **Tilt Trigger Extension** \$ 50

**Power Tilt** 250 lbs Weight Capacity. 20° tilt range. Includes 2 batteries for off-chair charging. [pdgmobility.com/power-tilt](http://pdgmobility.com/power-tilt)  
Minimum Seat Width is 15". Power Tilt is compatible with either O2 Holder or IV Pole on the same chair, not both.

1023400 Standard Power Tilt \$ 1,795

1023410 Power Tilt with Manual Tilt Override Allows for instant manual release of tilt in case of power tilt failure. \$ 2,295

1023421 Installed on Right side       1023422 Installed on Left side  
 1023431 Installed on Right side       1023432 Installed on Left side  
 1023440 Not installed  
 1023445 Swing Away Mounting Hardware Not compatible with all Headrest makes and models. \$ 625

Optional for Egg Switch: \_\_\_\_\_

Specify Headrest to be used: \_\_\_\_\_

## Recline

1023310 **Reclining Backrest (Manual)** \$ 895

30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold-down Back, Tension Adjustable Upholstery, Dynamic Rocker Back. Requires longer Frame. Maximum Seat Depth is 19". Finished Back Post Height is 3" taller than selected Back Post Height.

## Front Rigging

1025000 **No Front Rigging**       1025001 **Sector Blocks only**

**Easy-Out™ Front Rigging** Effortless insertion and removal. No Hanger Bracket remains on frame. [pdgmobility.com/easy-out](http://pdgmobility.com/easy-out)

1025410 Fixed 70° \$ 95

1025420 Universal Angle Adjustable (0°-110°) Set at angle: \_\_\_\_\_ (use increments of 10°) \$ 295

Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings.

**Swing Away Front Rigging** Legacy Design. Not compatible with Easy-Out.

1025130 Swing Away 60°       1025140 Swing Away 70°


1025200 Swing Away Elevating Legrests Includes Calf Pads. \$ 375

1025300 1" Offset Hangers for Swing Away \$ 75

1025110 **Standard Plug In 70° Front Rigging** Legacy Design. Not compatible with Easy-Out.

Footplates			
<input type="checkbox"/>	1026010	<b>Composite Footplates</b>	
<input type="checkbox"/>	1026020	<b>Aluminum Angle Adjustable Footplates</b>	\$ 165
<input type="checkbox"/>	1026030	<b>One Piece Fixed Footplate</b> Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.	\$ 350
<input type="checkbox"/>	1026040	<b>One Piece Flip-up Footplate</b> Requires selection of Composite or Aluminum Angle Adjustable.	\$ 395
<input type="checkbox"/>	1026050	<b>Interlocking Angle Adjustable Footplates</b> Available for Seat Widths up to 21" only.	\$ 195

Footrest/Legrest Accessories			
<b>Calf Pads</b> Not compatible with Gel Padding.			
<input type="checkbox"/>	1027012	For Easy-Out™ Front Rigging Includes reinforced premium Calf Pad mount.	\$ 200
<input type="checkbox"/>	1027010	For Swing Away Front Rigging	\$ 135
<input type="checkbox"/>	1027020	<b>Gel Padding</b> Not compatible with Calf Pad or Easy-Out Front Rigging.	\$ 125
<b>Gel Knee Pads</b>			
<input type="checkbox"/>	1027032	For Easy-Out™ Angle Adjustable Front Rigging	\$ 200
<input type="checkbox"/>	1027030	For Swing Away ELR Front Rigging	\$ 100
<input type="checkbox"/>	1027040	<b>Length Adjustable Heel Loops</b>	\$ 50
<input type="checkbox"/>	1027050	<b>Bodypoint Padded Calf Strap</b> Sizes greater than 22" are not padded.	\$ 195
<b>Miller's Dynamic Footrest Mounts</b> Not compatible with Easy-Out Front Rigging.			
<input type="checkbox"/>	1027060	For Standard Plug In Front Rigging	\$ 445
<input type="checkbox"/>	1027065	For Swing Away 70° Front Rigging Not available on 13" STF Height.	\$ 575

Front Rigging Extension Tubes			
<b>Regular Extension Tubes</b>			
<input type="checkbox"/>	1028010	Short 14"-18" for Easy-Out Front Rigging 13"-17" for Swing Away and Standard Plug In Front Rigging	
<input type="checkbox"/>	1028020	Long 18"-24" for Easy-Out Front Rigging 15"-19" for Swing Away and Standard Plug In Front Rigging	\$ 75
<b>Extra Strength Extension Tubes (Solid)</b> Available on Plug In Front Rigging only. Not compatible with Swing Away or Easy-Out Front Rigging.			\$ 100
<input type="checkbox"/>	1028110	Standard (13"-17")	<input type="checkbox"/> 1028120 Long (15"-19")
	<b>Lower Leg Length</b> If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.		
	Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.		

Miscellaneous Options			
<b>One-Arm Drive</b> Available with 22" and 24" Composite Rear Wheels. Minimum 17" STF and 18" STF, respectively. Limits Tilt. Not compatible with User & Attendant Wheel Locks, Drum Brakes, Quick Release Rear Axles, Pneumatic Tires.			\$ 3,500
<input type="checkbox"/>	1030010	Installed on Right side	<input type="checkbox"/> 1030011 Installed on Left side
<input type="checkbox"/>	1030012	Large Hand Setup 24" Rear Wheels with 22" Driving Handrim.	
<b>Spoke Guards</b> Available for 20"-24" Composite Rear Wheels. Made from durable and lightweight aluminum.			\$ 195
<input type="checkbox"/>	1030020	Black	<input type="checkbox"/> 1030021 Frame Color
<input type="checkbox"/>	1030050	<b>Protective Lower Padding</b>	\$ 795
<b>Stealth Flip Down Abductor</b> Requires Solid Seat Pan.			\$ 650
<input type="checkbox"/>	1030060	Small Pad 3"x4"x3"x2" size	<input type="checkbox"/> 1030062 Large Pad 5"x5"x4"x3" size
<input type="checkbox"/>	1030080	<b>Gel Hip Guides</b> Reduces available Seat Width by 2".	\$ 495
<input type="checkbox"/>	1030110	<b>IV Pole</b> Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.	\$ 295
<input type="checkbox"/>	1030120	<b>O2 Holder</b> Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.	\$ 395
<input type="checkbox"/>	1030200	<b>Q'Strain Lap Belt</b> Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" Belt Length.	\$ 245
<input type="checkbox"/>	1030222	<b>Wheelchair Backpack</b> Black with PDG Logo. 17"x14" size.	\$ 145
<input type="checkbox"/>	1030170	<b>Adhesive Back Hook Strip, 2" x 30"</b> For Cushion mounting.	\$ 10



# BENTLEY

**Order Form - Canada**  
Prices in Canadian Dollars. January 1, 2024

## On-Chair Seating & Positioning



Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. [pdgmobility.com/seating](http://pdgmobility.com/seating)

## Order Notes

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.