



Customer and Order Information

Purchase Order *	Tag	v.0401
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in
		<i>*required</i>

Base Chair

Stellar GLT 45° Manual Tilt Wheelchair

<input checked="" type="checkbox"/> 2501000	Stellar GL 45° Manual Tilt Wheelchair	\$ 4,300
<input checked="" type="checkbox"/> 2501009	Lightweight Aluminum Frame Upgrade Aluminum Frame is not recommended for agitated users.	\$ 300
i Choose the standard Stellar GL with a steel frame for agitated or heavier users. See separate order form.		

Weight Capacity

<input checked="" type="checkbox"/> 2502000	250 lbs Chair Weight Capacity				
User Weight Range Required for setup of tilt uplift force.					
<input type="checkbox"/> 1902430	Up to 75 lbs	<input type="checkbox"/> 2502431	75 lbs to 100 lbs	<input type="checkbox"/> 2502432	100 lbs to 115 lbs
<input type="checkbox"/> 1902433	115 lbs to 150 lbs	<input type="checkbox"/> 2502434	150 lbs to 200 lbs	<input type="checkbox"/> 2502435	200 lbs to 250 lbs

Frame Options

<input type="checkbox"/> 2502300	Transit Tie Down System (TTDS)
TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details.	

Seat Width

Overall chair width is +9.5" with 12/16" Rear Wheels and +10.5" with 20/22/24" Rear Wheels.

Standard Seat Widths

<input type="checkbox"/> 2503014	14" Seat Width	<input type="checkbox"/> 2503015	15" Seat Width	<input type="checkbox"/> 2503016	16" Seat Width
<input type="checkbox"/> 2503017	17" Seat Width	<input type="checkbox"/> 2503018	18" Seat Width	<input type="checkbox"/> 2503019	19" Seat Width
<input type="checkbox"/> 2503020	20" Seat Width	<input type="checkbox"/> 2503021	21" Seat Width	<input type="checkbox"/> 2503022	22" Seat Width

Seat Depth

Standard Seat Depths (Adjustable)

<input type="checkbox"/> 2504016	16" Seat Depth	<input type="checkbox"/> 2504017	17" Seat Depth	<input type="checkbox"/> 2504018	18" Seat Depth
<input type="checkbox"/> 2504019	19" Seat Depth	<input type="checkbox"/> 2504020	20" Seat Depth		
Custom Seat Depths					
<input type="checkbox"/> 2504021	21" Seat Depth	<input type="checkbox"/> 2504022	22" Seat Depth	\$ 400	

Frame Color

PDG Premium Colors View color swatches on PDG's website. pdgmobility.com/colors

<input type="checkbox"/> 2505001	Black Shimmer	<input type="checkbox"/> 2505015	Sky Blue	<input type="checkbox"/> 2505016	Poppy Red
<input type="checkbox"/> 2505003	Fine Texture Black	<input type="checkbox"/> 2505010	Lollypop Blue	<input type="checkbox"/> 2505011	Lollypop Red
<input type="checkbox"/> 2505017	Arctic Grey	<input type="checkbox"/> 2505018	Tropical Teal	<input type="checkbox"/> 2505019	Tangerine Orange
<input type="checkbox"/> 2505020	Olive Green	<input type="checkbox"/> 2505021	Lush Lilac	<input type="checkbox"/> 2505012	Bubble Gum Pink
<input type="checkbox"/> 2505022	Leaf Green	<input type="checkbox"/> 2505023	Boysenberry Purple	<input type="checkbox"/> 2505014	Daffodil Yellow
<input type="checkbox"/> 2505200	Custom Color Contact PDG for specifications.				\$ 300

Armrests

2506000 **No Armrests**

Height Adjustable Armrests (Single Post) \$ 295

2506111 8"-11" Range 2506114 10"-14" Range 2506116 12"-16" Range

2506100 **Fixed Height Armrests (12")** Full Length Armpad only.

Armrest Height If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.
Seat Pan to top of Armrest: _____ Does not include height of Armpad

Armrest Position Indicate preference for Armrest position. If left blank, the default setting will be used.
In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

2506401 Rear 2506402 Mid 2506403 Front

Armpads

Regular Armpads

2507100 Desk Length - 2" x 10" 2507101 Full Length - 2" x 14"

Waterfall Armpads \$ 115

2507105 Desk Length - 2" x 10" 2507106 Full Length - 2" x 14"

Wide Armpads May interfere with wheels in some tilt configurations. \$ 200

2507110 Desk Length - 3" x 10" 2507111 Full Length - 3" x 12"

Gel Armpads May interfere with wheels in some tilt configurations. \$ 135

2507120 Full Length - 2" x 12" 2507121 Full Length - 3.5" x 12"

Contour Arm Trough with Hand Pad May interfere with wheels in some tilt configurations. \$ 125

2507130 Right (4.5" x 18") 2507131 Left (4.5" x 18")

Armrest Options

Armrest Wideners Move Armrests outward by up to 1" on each side. Seat size does not change. May cause tilt interference.

2507210 1" Total Additional Width \$ 95

2507211 2" Total Additional Width \$ 195

2507310 **Clothing Guards** Not compatible with 8"-11" Height range. \$ 110

Armrest Extensions Not compatible with Desk Length Armpads. \$ 75

2507400 3" Extension - Forward 2507401 3" Extension - Rearward

Seat Pan

2508000 **No Seat Pan**

2508001 **Solid Aluminum Seat Pan** \$ 135

Seat Height

Don't forget **Select desired Seat-to-Floor Height** STF Height is adjustable within the range of the chosen Subframe Plate (Low/High).

6" Casters	Low							High						
	RW	13"	14"	15"	16"	17"	18"	RW	15"	16"	17"	18"	19"	20"
12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

8" Casters	Low							High						
	RW	13"	14"	15"	16"	17"	18"	RW	15"	16"	17"	18"	19"	20"
16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Casters

Urethane Casters 8" Casters are compatible only with 16" and 24" Rear Wheels.

<input type="checkbox"/> 2511061	6" x 1.25" Caster size	<input type="checkbox"/> 2511081	8" x 1.25" Caster size	
<input type="checkbox"/> 2511063	6" x 1.50" Caster size			\$ 85
<input type="checkbox"/> 2511084	8" x 2.00" Caster size			\$ 135

Pneumatic Casters 8" Casters are compatible only with 16" and 24" Rear Wheels.

<input type="checkbox"/> 2511161	6" x 1.25" Caster size			\$ 110
<input type="checkbox"/> 2511181	8" x 1.25" Caster size			\$ 120
<input type="checkbox"/> 2511184	8" x 2.00" Caster size			\$ 150

<input type="checkbox"/> 2511260	Frog Legs 6" x 1.50" Aluminum Soft-Roll Casters			\$ 275
----------------------------------	--	--	--	--------

<input type="checkbox"/> 2511400	Frog Legs Suspension Forks	Only compatible with 16 and 24" Rear Wheels. Requires 6" Soft-Roll Caster.		\$ 595
----------------------------------	-----------------------------------	--	--	--------

Rear Wheels

Composite Wheels

<input type="checkbox"/> 2512012	12" Wheel size	<input type="checkbox"/> 2512016	16" Wheel size	<input type="checkbox"/> 2512020	20" Wheel size
<input type="checkbox"/> 2512022	22" Wheel size	<input type="checkbox"/> 2512024	24" Wheel size		

Spoke Wheels

<input type="checkbox"/> 2512120	20" Wheel size	<input type="checkbox"/> 2512122	22" Wheel size	<input type="checkbox"/> 2512124	24" Wheel size	\$ 75
----------------------------------	----------------	----------------------------------	----------------	----------------------------------	----------------	-------



Rear Wheel Position Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.

In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

<input type="checkbox"/> 2512401	Rear	<input type="checkbox"/> 2512402	Mid	<input type="checkbox"/> 2512403	Front
----------------------------------	------	----------------------------------	-----	----------------------------------	-------

Axle Type

<input type="checkbox"/> 2513000	Quick Release Rear Axles
<input type="checkbox"/> 2513010	Fixed Rear Axles

Rear Tires

<input type="checkbox"/> 2514000	Urethane Tires		
<input type="checkbox"/> 2514010	Pneumatic Tires	May result in a finished Seat Height 1/2" greater than stated.	\$ 60
<input type="checkbox"/> 2514020	Pneumatic Tires with Flat-Free Inserts	May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels.	\$ 95

Handrims

Handrims are only available with 20", 22", and 24" Rear Wheels.

<input type="checkbox"/> 2515000	No Handrims		
<input type="checkbox"/> 2515010	Standard Aluminum Handrims		
<input type="checkbox"/> 2515020	Plastic Coated Handrims	Not compatible with Drum Brakes.	\$ 150
<input type="checkbox"/> 2515030	Projection Handrims	Composite Wheels only. 10 Projections, 30°, 1-1.8" Length.	\$ 495

Wheel Locks & Brakes

User Wheel Locks

<input type="checkbox"/> 2516010	Push to Lock	<input type="checkbox"/> 2516020	Pull to Lock	
<input type="checkbox"/> 2516030	Attendant Wheel Locks (Push to Lock)			\$ 85
<input type="checkbox"/> 2516040	User & Attendant Wheel Locks	Not compatible with Drum Brakes.		\$ 125
<input type="checkbox"/> 2516300	Drum Brakes	Available with 12" Composite or 20" Spoke Rear Wheels. Urethane Tires only.		\$ 795
	User-activated Unilateral Wheel Lock	Pull to Lock.		\$ 795
<input type="checkbox"/> 2516310	Left Activation	<input type="checkbox"/> 2516311	Right Activation	



Wheel Lock Options		
<input type="checkbox"/>	2516100 Wheel Lock Extensions	Not compatible with Attendant Wheel Lock. \$ 60
<input type="checkbox"/>	2516200 Hill Holders	Requires Pneumatic Rear Tires with Flat-Free Inserts. \$ 450

Anti-Tippers		
<input checked="" type="checkbox"/>	2517010 Rear Anti-Tippers	Required Safety Option. \$ 75

Back Posts		
Back Post Style		
<input checked="" type="checkbox"/>	2518010	Stroller Handle Back
Back Post Height Measured from the Seat Pan to the Handle.		
<input type="checkbox"/>	2518425	Short (22")
<input type="checkbox"/>	2518431	Tall (25")
Back Post Folding Fold-down Back is not compatible with Reclining Backrest.		
<input checked="" type="checkbox"/>	2519220	Fold-down
Seat to Back Angle		
<input type="checkbox"/>	2519090	Set at 90°
<input type="checkbox"/>	2519098	Set at 98°
<input type="checkbox"/>	2519106	Set at 106°
<input type="checkbox"/>	2519114	Set at 114°
<input type="checkbox"/>	2519122	Set at 122°

Back Options		
Angle/Height Adjustable Stroller Handle		Not compatible with Headrest Mounting Plate. \$ 150
<input type="checkbox"/>	2519010	Button on Outside
<input type="checkbox"/>	2519011	Button on Inside
Endomorph Back		Backrest will be narrower than Seat by total of 2x the Offset selected below. Minimum width of Backrest is 12". \$ 395
<input type="checkbox"/>	2518102	1" Offset per side
<input type="checkbox"/>	2518103	1.5" Offset per side
<input type="checkbox"/>	2518104	2" Offset per side
<input type="checkbox"/>	2519210	10" Backrest Extension Not compatible with Angle/Height Adjustable Stroller Handle. \$ 150
<input type="checkbox"/>	2519212	Headrest Mounting Plate \$ 95
Not compatible with Backrest Extension, Angle/Height Adjustable Stroller Handle. Not required if mounting to a solid Back.		

Back Upholstery																																																																						
<input type="checkbox"/>	2520000	No Back Upholstery																																																																				
Standard Upholstery		\$ 195																																																																				
<input type="checkbox"/>	2520418	Short (18" Height)																																																																				
<input type="checkbox"/>	2520421	Regular (21" Height)																																																																				
<input type="checkbox"/>	2520424	Tall (24" Height)																																																																				
<input type="checkbox"/>	2520450	Tension Adjustable Upholstery Height of Upholstery is 16". \$ 295																																																																				
Position of Upholstery Select the desired position of the top of the Upholstery relative to the Seat Pan.																																																																						
<table border="1"> <thead> <tr> <th colspan="2"></th> <th colspan="9">22" Back Post</th> </tr> <tr> <th>Style</th> <th>17"</th> <th>18"</th> <th>19"</th> <th>20"</th> <th>21"</th> <th>22"</th> <th>23"</th> <th>24"</th> <th>25"</th> </tr> </thead> <tbody> <tr> <td>Short</td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Regular</td> <td></td> <td></td> <td></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Tension Adj</td> <td></td> <td></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td></td> <td></td> <td></td> </tr> </tbody> </table>												22" Back Post									Style	17"	18"	19"	20"	21"	22"	23"	24"	25"	Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>													
		22" Back Post																																																																				
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"																																																													
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																																																
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																																																
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																																																
<table border="1"> <thead> <tr> <th colspan="2"></th> <th colspan="9">25" Back Post</th> </tr> <tr> <th>Style</th> <th>17"</th> <th>18"</th> <th>19"</th> <th>20"</th> <th>21"</th> <th>22"</th> <th>23"</th> <th>24"</th> <th>25"</th> </tr> </thead> <tbody> <tr> <td>Short</td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td>Regular</td> <td></td> <td></td> <td></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td>Tall</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td>Tension Adj</td> <td></td> <td></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </tbody> </table>												25" Back Post									Style	17"	18"	19"	20"	21"	22"	23"	24"	25"	Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Tall							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		25" Back Post																																																																				
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"																																																													
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																																													
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																																													
Tall							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																																													
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																																													

Tilt Activation		
Attendant Tilt only Select trigger mounting location below.		
<input type="checkbox"/>	2523011	Right Back Post
<input type="checkbox"/>	2523012	Left Back Post
User Tilt only Not available on widths greater than 22". Full Length Arm Pads only. Select trigger mounting location below.		\$ 100
<input type="checkbox"/>	2523021	Right Armrest
<input type="checkbox"/>	2523022	Left Armrest
Attendant & User Tilt Not available on widths greater than 22". Full Length Arm Pads only. Select trigger mounting locations below.		\$ 195
<input type="checkbox"/>	2523031	Right Back Post
<input type="checkbox"/>	2523032	Left Back Post
<input type="checkbox"/>	2523041	Right Armrest
<input type="checkbox"/>	2523042	Left Armrest

See next page for more options >



<input type="checkbox"/>	2523110	Tilt Trigger Extension				\$ 50		
Tilt Limit Optional. If no tilt limit is selected, the full tilt range will be available.								
<input type="checkbox"/>	2523220	Set at 20°	<input type="checkbox"/>	2523230	Set at 30°	<input type="checkbox"/>	2523240	Set at 40°
Power Tilt 250 lbs Weight Capacity. 45° tilt range. Includes 2 batteries for off-chair charging. pdgmobility.com/power-tilt Minimum Seat Width is 16". Power Tilt is compatible with either O2 Holder or IV Pole on the same chair, not both.								
<input type="checkbox"/>	2523400	Standard Power Tilt				\$ 1,795		
<input type="checkbox"/>	2523410	Power Tilt with Manual Tilt Override	Allows for instant manual release of tilt in case of power tilt failure.			\$ 2,295		
➤ Select Control Switch		Toggle Switch	<input type="checkbox"/>	2523421	Installed on Right side	<input type="checkbox"/>	2523422	Installed on Left side
		Button Switch	<input type="checkbox"/>	2523431	Installed on Right side	<input type="checkbox"/>	2523432	Installed on Left side
		Stealth Egg Switch	<input type="checkbox"/>	2523440	Not installed			
		Optional for Egg Switch:	<input type="checkbox"/>	2523445	Swing Away Mounting Hardware	Not compatible with all Headrest makes and models.		\$ 625
<i>Specify Headrest to be used: _____</i>								

Recline						
<input type="checkbox"/>	2523310	Reclining Backrest (Manual)				\$ 895
30° maximum recline. Not compatible with Backrest Extension, Fold-down Back, Tension Adjustable Upholstery. Finished Back Post Height is 3" taller than selected Back Post Height. Limits Tilt by up to 20°.						

Front Rigging						
<input type="checkbox"/>	2525000	No Front Rigging	<input type="checkbox"/>	2525001	Sector Blocks only	
Easy-Out™ Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame. pdgmobility.com/easy-out						
<input type="checkbox"/>	2525410	Fixed 70°				\$ 95
<input type="checkbox"/>	2525420	Universal Angle Adjustable (0°-110°)	Set at angle: _____ (use increments of 10°)			\$ 295
Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings.						
Swing Away Front Rigging Legacy Design. Not compatible with Easy-Out.						
<input type="checkbox"/>	2525130	Swing Away 60°	<input type="checkbox"/>	2525140	Swing Away 70°	
<input type="checkbox"/>	2525200	Swing Away Elevating Legrests	Includes Calf Pads.			\$ 375
<input type="checkbox"/>	2525300	1" Offset Hangers for Swing Away				\$ 75
<input type="checkbox"/>	2525110	Standard Plug In 70° Front Rigging	Legacy Design. Not compatible with Easy-Out.			

Footplates						
<input type="checkbox"/>	2526010	Composite Footplates				
<input type="checkbox"/>	2526020	Aluminum Angle Adjustable Footplates				\$ 165
<input type="checkbox"/>	2526030	One Piece Fixed Footplate	Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.			\$ 350
<input type="checkbox"/>	2526040	One Piece Flip-up Footplate	Requires selection of Composite or Aluminum Angle Adjustable.			\$ 395
<input type="checkbox"/>	2526050	Interlocking Angle Adjustable Footplates	Available for Seat Widths up to 21" only.			\$ 195

Footrest/Legrest Accessories						
Calf Pads Not compatible with Gel Padding.						
<input type="checkbox"/>	2527012	For Easy-Out™ Front Rigging	Includes reinforced premium Calf Pad mount.			\$ 200
<input type="checkbox"/>	2527010	For Swing Away Front Rigging				\$ 135
<input type="checkbox"/>	2527020	Gel Padding	Not compatible with Calf Pad or Easy-Out Front Rigging.			\$ 125
Gel Knee Pads						
<input type="checkbox"/>	2527032	For Easy-Out™ Angle Adjustable Front Rigging				\$ 200
<input type="checkbox"/>	2527030	For Swing Away ELR Front Rigging				\$ 100
<input type="checkbox"/>	2527040	Length Adjustable Heel Loops				\$ 50
<input type="checkbox"/>	2527050	Bodypoint Padded Calf Strap				\$ 195

Front Rigging Extension Tubes

Regular Extension Tubes

- | | | | | | | | |
|--------------------------|---------|-------|---------|----------------------------|---------|---|-------|
| <input type="checkbox"/> | 2528010 | Short | 14"-18" | for Easy-Out Front Rigging | 13"-17" | for Swing Away and Standard Plug In Front Rigging | |
| <input type="checkbox"/> | 2528020 | Long | 18"-24" | for Easy-Out Front Rigging | 15"-19" | for Swing Away and Standard Plug In Front Rigging | \$ 75 |



Lower Leg Length If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.

Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.

Miscellaneous Options

Spoke Guards Available for 20"-24" Composite Rear Wheels. Made from durable and lightweight aluminum.

- | | | | | | | |
|--------------------------|---------|-------|--------------------------|---------|-------------|--------|
| <input type="checkbox"/> | 2530020 | Black | <input type="checkbox"/> | 2530021 | Frame Color | \$ 195 |
|--------------------------|---------|-------|--------------------------|---------|-------------|--------|

Stealth Flip Down Abductor Requires Solid Seat Pan.

- | | | | | | | | | |
|--------------------------|---------|-----------|------------------|--------------------------|---------|-----------|------------------|--------|
| <input type="checkbox"/> | 2530060 | Small Pad | 3"x4"x3"x2" size | <input type="checkbox"/> | 2530062 | Large Pad | 5"x5"x4"x3" size | \$ 650 |
|--------------------------|---------|-----------|------------------|--------------------------|---------|-----------|------------------|--------|

- | | | | | | | | |
|--------------------------|---------|-----------------------|-------------------------------------|--|--|--|--------|
| <input type="checkbox"/> | 2530080 | Gel Hip Guides | Reduces available Seat Width by 2". | | | | \$ 495 |
|--------------------------|---------|-----------------------|-------------------------------------|--|--|--|--------|

- | | | | | | | | |
|--------------------------|---------|----------------|--|--|--|--|--------|
| <input type="checkbox"/> | 2530110 | IV Pole | Not compatible with Attendant Wheel Locks. Requires frame modification and cannot be added to chair later. | | | | \$ 295 |
|--------------------------|---------|----------------|--|--|--|--|--------|

- | | | | | | | | |
|--------------------------|---------|------------------|--|--|--|--|--------|
| <input type="checkbox"/> | 2530120 | O2 Holder | Not compatible with Attendant Wheel Locks. Requires frame modification and cannot be added to chair later. | | | | \$ 395 |
|--------------------------|---------|------------------|--|--|--|--|--------|

- | | | | | | | | |
|--------------------------|---------|---------------------------|--|--|--|--|--------|
| <input type="checkbox"/> | 2530200 | Q'Straint Lap Belt | Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" Belt Length. | | | | \$ 245 |
|--------------------------|---------|---------------------------|--|--|--|--|--------|

- | | | | | | | | |
|--------------------------|---------|----------------------------|------------------------------------|--|--|--|--------|
| <input type="checkbox"/> | 2530222 | Wheelchair Backpack | Black with PDG Logo. 17"x14" size. | | | | \$ 145 |
|--------------------------|---------|----------------------------|------------------------------------|--|--|--|--------|

- | | | | | | | | |
|--------------------------|---------|---|-----------------------|--|--|--|-------|
| <input type="checkbox"/> | 2530170 | Adhesive Back Hook Strip, 2" x 30" | For Cushion mounting. | | | | \$ 10 |
|--------------------------|---------|---|-----------------------|--|--|--|-------|

On-Chair Seating & Positioning

 Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. pdgmobility.com/seating

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.