



### Customer and Order Information

Purchase Order *	Tag	v.0101
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in

### Base Chair

**Stellar HD 45° Manual Tilt Wheelchair** Extended delivery lead times may apply.

<input checked="" type="checkbox"/> 2101000	Stellar 45° Manual Tilt Wheelchair	\$ 4,300
<input checked="" type="checkbox"/> 2101201	Reinforced HD Frame Package Heavy-gauge steel reinforcement of frame for bariatric users. Required selection.	\$ 1,395

**i** Choose the **Stellar Impact** for highly agitated users up to 250 lbs. See separate order form.

### Weight Capacity

**Chair Weight Capacity**

<input type="checkbox"/> 2102045	450 lbs Capacity	
<input type="checkbox"/> 2102060	600 lbs Capacity Tilt is limited to maximum 25°.	\$ 300

**User Weight Range** Required for setup of tilt uplift force.

<input type="checkbox"/> 2102212	Up to 250 lbs	<input type="checkbox"/> 2102308	250 lbs to 350 lbs
<input type="checkbox"/> 2102310	350 lbs to 425 lbs	<input type="checkbox"/> 2102312	425 lbs to 600 lbs

### Seat Width

Overall chair width is +11"

**Standard Seat Widths**

<input type="checkbox"/> 2103022	22" Seat Width	<input type="checkbox"/> 2103024	24" Seat Width	<input type="checkbox"/> 2103026	26" Seat Width	\$ 375
<input type="checkbox"/> 2103028	28" Seat Width	<input type="checkbox"/> 2103030	30" Seat Width	<input type="checkbox"/> 2103032	32" Seat Width	\$ 800

### Seat Depth

**Standard Seat Depths (Adjustable)**

<input type="checkbox"/> 2104018	18" Seat Depth	<input type="checkbox"/> 2104019	19" Seat Depth	<input type="checkbox"/> 2104020	20" Seat Depth
<input type="checkbox"/> 2104021	21" Seat Depth	<input type="checkbox"/> 2104022	22" Seat Depth		\$ 300

### Frame Color

**PDG Premium Colors** View color swatches on PDG's website. [pdgmobility.com/colors](http://pdgmobility.com/colors)

<input type="checkbox"/> 2105001	Black Shimmer	<input type="checkbox"/> 2105015	Sky Blue	<input type="checkbox"/> 2105016	Poppy Red
<input type="checkbox"/> 2105003	Fine Texture Black	<input type="checkbox"/> 2105010	Lollypop Blue	<input type="checkbox"/> 2105011	Lollypop Red
<input type="checkbox"/> 2105017	Arctic Grey	<input type="checkbox"/> 2105018	Tropical Teal	<input type="checkbox"/> 2105019	Tangerine Orange
<input type="checkbox"/> 2105020	Olive Green	<input type="checkbox"/> 2105021	Lush Lilac	<input type="checkbox"/> 2105012	Bubble Gum Pink
<input type="checkbox"/> 2105022	Leaf Green	<input type="checkbox"/> 2105023	Boysenberry Purple	<input type="checkbox"/> 2105014	Daffodil Yellow

2105200 **Custom Color** Contact PDG for specifications. \$ 300

### Armrests

<input type="checkbox"/> 2106000	<b>No Armrests</b>				
<b>Single Post Height Adjustable Armrests</b>		\$ 295			
<input type="checkbox"/> 2106111	8"-11" Range	<input type="checkbox"/> 2106114	10"-14" Range	<input type="checkbox"/> 2106116	12"-16" Range

See next page for more options >

<b>Dual Post Height Adjustable Armrests</b>			\$ 295
<input type="checkbox"/> 2106211	8"-11" Range	<input type="checkbox"/> 2106214	10"-14" Range
		<input type="checkbox"/> 2106216	12"-16" Range
<b>Fixed Height Armrests (12")</b> Full Length Armpad only.			
<input type="checkbox"/> 2106100	Single Post	<input type="checkbox"/> 2106200	Dual Post
<b>Armrest Height</b> If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.			
Seat Pan to top of Armrest: _____ Does not include height of Armpad			
<b>Armrest Position</b> Indicate preference for Armrest position. If left blank, the default setting will be used.			
In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.			
<input type="checkbox"/> 2106401	Rear	<input type="checkbox"/> 2106402	Mid
		<input type="checkbox"/> 2106403	Front

<b>Armpads</b>			
<b>Regular Armpads</b>			
<input type="checkbox"/> 2107100	Desk Length - 2" x 10"	<input type="checkbox"/> 2107101	Full Length - 2" x 14"
<b>Waterfall Armpads</b> Not compatible with Dual Post Armrests			
<input type="checkbox"/> 2107105	Desk Length - 2" x 10"	<input type="checkbox"/> 2107106	Full Length - 2" x 14"
<b>Wide Armpads</b> May interfere with wheels in some tilt configurations.			
<input type="checkbox"/> 2107110	Desk Length - 3" x 10"	<input type="checkbox"/> 2107111	Full Length - 3" x 12"
<b>Gel Armpads</b> May interfere with wheels in some tilt configurations.			
<input type="checkbox"/> 2107120	Full Length - 2" x 12"	<input type="checkbox"/> 2107121	Full Length - 3.5" x 12"
<b>Contour Arm Trough with Hand Pad</b> May interfere with wheels in some tilt configurations.			
<input type="checkbox"/> 2107130	Right (4.5" x 18")	<input type="checkbox"/> 2107131	Left (4.5" x 18")

<b>Armrest Options</b>			
<b>Armrest Wideners</b> Move Armrests outward by up to 1" on each side. Seat size does not change. May cause tilt interference.			
<input type="checkbox"/> 2107210	1" Total Additional Width		\$ 95
<input type="checkbox"/> 2107211	2" Total Additional Width		\$ 195
<b>Clothing Guards</b>			
<input type="checkbox"/> 2107310	For Single Post Armrests	Not compatible with 8"-11" Height range.	\$ 110
<input type="checkbox"/> 2107320	For Dual Post Armrests		\$ 165
<b>Armrest Extensions</b> Not compatible with Desk Length Armpads.			
<input type="checkbox"/> 2107400	3" Extension - Forward	<input type="checkbox"/> 2107401	3" Extension - Rearward

<b>Seat Pan</b>			
<input checked="" type="checkbox"/> 2108001	<b>Solid Aluminum Seat Pan</b>		\$ 175

<b>Seat Height</b>			
<b>Seat to Floor Height</b>			
<input type="checkbox"/> 2110016	Set at 16" STF	<input type="checkbox"/> 2110017	Set at 17" STF
		<input type="checkbox"/> 2110018	Set at 18" STF

<b>Casters</b>			
<input checked="" type="checkbox"/> 2111050	<b>5" Casters with Heavy Duty Steel Fork</b>		

<b>Rear Wheels</b>			
<input checked="" type="checkbox"/> 2112224	<b>24" Heavy Duty Spoke Wheels</b>		
<b>Rear Wheel Position</b> Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.			
In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.			
<input type="checkbox"/> 2112401	Rear	<input type="checkbox"/> 2112402	Mid
		<input type="checkbox"/> 2112403	Front

<b>Axle Type</b>			
<input checked="" type="checkbox"/> 2113010	<b>Fixed Rear Axles</b>		



## Rear Tires

2114000 **Urethane Tires**

## Handrims

2115000 **No Handrims**

2115010 **Standard Aluminum Handrims**

2115020 **Plastic Coated Handrims** \$ 250

## Wheel Locks

**User Wheel Locks**

2116010 Push to Lock  2116020 Pull to Lock

## Wheel Lock Options

2116100 **Wheel Lock Extensions** \$ 60

## Anti-Tippers

2117010 **Rear Anti-Tippers** Required Safety Option. \$ 75

## Back Post

2118010 **Stroller Handle Back**

**Back Post Height** Measured from the Seat Pan to the Handle.

2118425 Short (22")  2118431 Tall (25")

2119090 **Seat to Back Angle set at 90°**

## Back Options

**Endomorph Back** Backrest will be narrower than Seat by total of 2x the Offset selected below. Minimum width of Backrest is 12". \$ 395

2118102 1" Offset per side  2118103 1.5" Offset per side  2118104 2" Offset per side

2119212 **Headrest Mounting Plate** Not required if mounting to a solid Back. \$ 95

2119240 **Seating Dynamics Dynamic Rocker Back Interface** \$ 2,995

Not compatible with 30° Recline, 21" & 22" Seat Depths. Finished Back Post Height is 4" taller than selected Back Post Height.

## Back Upholstery

2120000 **No Back Upholstery**

**Tension Adjustable Upholstery** \$ 375

2120421 Regular (18" Height)  2120424 Tall (22" Height)

**Position of Upholstery** Select the desired position of the top of the Upholstery relative to the Seat Pan.

		22" Back Post									
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"		
Regular	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>					
Tall				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>					

		25" Back Post									
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"		
Regular	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Tall					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

## Tilt Activation

2123010 **Attendant Tilt only** Two tilt triggers (L and R) are required for tilt activation.

**Tilt Limit** Optional. If no tilt limit is selected, the full tilt range will be available. Maximum Tilt Limit for 600 lbs Weight Capacity is 25°.

2123220 Set at 20°  2123230 Set at 30°  2123240 Set at 40°

## Recline

2123310 **Reclining Backrest (Manual)** \$ 895

30° maximum recline. Not compatible with 600 lbs Weight Capacity, Tension Adjustable Upholstery. Finished Back Post Height is 3" taller than selected Back Post

### Front Rigging

<input type="checkbox"/>	2125000	<b>No Front Rigging</b>	<input type="checkbox"/>	2125001	<b>Sector Blocks only</b>	
<b>Easy-Out™ Heavy Duty Front Rigging</b> Effortless insertion and removal. No Hanger Bracket remains on frame. <a href="http://pdgmobility.com/easy-out">pdgmobility.com/easy-out</a>						
<input type="checkbox"/>	2125410	Fixed 70°				\$ 95
<input type="checkbox"/>	2125420	Universal Angle Adjustable (0°-110°)	Set at angle: _____ (use increments of 10°)			\$ 295
Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings.						
<b>Swing Away Front Rigging</b> Legacy Design. Not compatible with Easy-Out.						
<input type="checkbox"/>	2125130	Swing Away 60°	<input type="checkbox"/>	2125140	Swing Away 70°	
<input type="checkbox"/>	2125200	Swing Away Elevating Legrests	Includes Calf Pads.			\$ 375
<input type="checkbox"/>	2125300	1" Offset Hangers for Swing Away				\$ 75
<input type="checkbox"/>	2125110	<b>Standard Plug In 70°</b>	Legacy Design. Not compatible with Easy-Out.			


### Footplates

<input checked="" type="checkbox"/>	2126020	<b>Aluminum Angle Adjustable Footplates</b>				\$ 175
-------------------------------------	---------	---	--	--	--	--------

### Footrest/Legrest Accessories

<b>Calf Pads</b> Not compatible with Gel Padding.						
<input type="checkbox"/>	2127012	For Easy-Out™ Front Rigging	Includes reinforced premium Calf Pad mount.			\$ 200
<input type="checkbox"/>	2127010	For Swing Away Front Rigging				\$ 135
<input type="checkbox"/>	2127020	<b>Gel Padding</b>	Not compatible with Calf Pad or Easy-Out Front Rigging.			\$ 125
<b>Gel Knee Pads</b>						
<input type="checkbox"/>	2127032	For Easy-Out™ Angle Adjustable Front Rigging				\$ 200
<input type="checkbox"/>	2127030	For Swing Away ELR Front Rigging				\$ 100
<input type="checkbox"/>	2127040	<b>Length Adjustable Heel Loops</b>				\$ 50
<input type="checkbox"/>	2127051	<b>Bodypoint Calf Strap</b>				\$ 195

### Front Rigging Extension Tubes

<b>Regular Extension Tubes</b>						
<input type="checkbox"/>	2128010	Short	14"-18" for Easy-Out Front Rigging	13"-17" for Swing Away and Standard Plug In Front Rigging		
<input type="checkbox"/>	2128020	Long	18"-24" for Easy-Out Front Rigging	15"-19" for Swing Away and Standard Plug In Front Rigging		\$ 75
<b>Extra Strength Extension Tubes (Solid)</b> Available on Plug In Front Rigging only. Not compatible with Easy-Out Front Rigging.						
<input type="checkbox"/>	2128110	Standard (13"-17")	<input type="checkbox"/>	2128120	Long (15"-19")	\$ 100
	<b>Lower Leg Length</b> If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.					
	Seat Pan to Foot Plates: _____		Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.			

### Miscellaneous Options

<input type="checkbox"/>	2130080	<b>Gel Hip Guides</b>	Reduces available Seat Width by 2".			\$ 495
<input type="checkbox"/>	2130110	<b>IV Pole</b>				\$ 295
<input type="checkbox"/>	2130120	<b>O2 Holder</b>				\$ 395
<input type="checkbox"/>	2130222	<b>Wheelchair Backpack</b>	Black with PDG Logo. 17"x14" size.			\$ 145
<input type="checkbox"/>	2130170	<b>Adhesive Back Hook Strip, 2" x 30"</b>	For Cushion mounting.			\$ 10

### On-Chair Seating & Positioning

	Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. <a href="http://pdgmobility.com/seating">pdgmobility.com/seating</a>					
---	---	--	--	--	--	--

### Order Notes

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.