



| Customer and Order Information | | |
|--|------------------|---|
| Purchase Order * | Tag | v.0401 |
| Dealer/Provider Name * | | |
| Purchaser Name and Contact Information * | | |
| Shipping Address * | | |
| User Information (if available) | Weight _____ lbs | Height _____ ft/in *required |

| Base Chair | | |
|--|---|----------|
| Stellar GLT 45° Manual Tilt Wheelchair | | |
| <input checked="" type="checkbox"/> 2501000 | Stellar GL 45° Manual Tilt Wheelchair | \$ 4,400 |
| <input checked="" type="checkbox"/> 2501009 | Lightweight Aluminum Frame Upgrade <small>Aluminum Frame is not recommended for agitated users.</small> | \$ 300 |
| ❗ Choose the standard Stellar GL with a steel frame for agitated or heavier users. <small>See separate order form.</small> | | |

| Weight Capacity | | |
|--|---|---|
| <input checked="" type="checkbox"/> 2502000 | 250 lbs Chair Weight Capacity | |
| User Weight Range <small>Required for setup of tilt uplift force.</small> | | |
| <input type="checkbox"/> 1902430 | Up to 75 lbs | <input type="checkbox"/> 2502431 75 lbs to 100 lbs |
| <input type="checkbox"/> 1902433 | 115 lbs to 150 lbs | <input type="checkbox"/> 2502432 100 lbs to 115 lbs |
| | <input type="checkbox"/> 2502434 150 lbs to 200 lbs | <input type="checkbox"/> 2502435 200 lbs to 250 lbs |

| Frame Options | | |
|--|---------------------------------------|--|
| <input type="checkbox"/> 2502300 | Transit Tie Down System (TTDS) | |
| <small>TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details.</small> | | |

| Seat Width | | |
|--|---|---|
| <small>Overall chair width is +9.5" with 12/16" Rear Wheels and +10.5" with 20/22/24" Rear Wheels.</small> | | |
| Standard Seat Widths | | |
| <input type="checkbox"/> 2503014 | 14" Seat Width | <input type="checkbox"/> 2503015 15" Seat Width |
| <input type="checkbox"/> 2503017 | 17" Seat Width | <input type="checkbox"/> 2503016 16" Seat Width |
| <input type="checkbox"/> 2503020 | 20" Seat Width | <input type="checkbox"/> 2503019 19" Seat Width |
| | <input type="checkbox"/> 2503021 21" Seat Width | <input type="checkbox"/> 2503022 22" Seat Width |

| Seat Depth | | |
|--|---|---|
| Standard Seat Depths (Adjustable) | | |
| <input type="checkbox"/> 2504016 | 16" Seat Depth | <input type="checkbox"/> 2504017 17" Seat Depth |
| <input type="checkbox"/> 2504019 | 19" Seat Depth | <input type="checkbox"/> 2504018 18" Seat Depth |
| | <input type="checkbox"/> 2504020 20" Seat Depth | |
| Custom Seat Depths | | |
| <input type="checkbox"/> 2504021 | 21" Seat Depth | <input type="checkbox"/> 2504022 22" Seat Depth |
| | | \$ 400 |

| Frame Color | | |
|---|--|---|
| PDG Premium Colors <small>View color swatches on PDG's website. pdgmobility.com/colors</small> | | |
| <input type="checkbox"/> 2505001 | Black Shimmer | <input type="checkbox"/> 2505015 Sky Blue |
| <input type="checkbox"/> 2505003 | Fine Texture Black | <input type="checkbox"/> 2505016 Poppy Red |
| <input type="checkbox"/> 2505017 | Arctic Grey | <input type="checkbox"/> 2505010 Lollypop Blue |
| <input type="checkbox"/> 2505020 | Olive Green | <input type="checkbox"/> 2505018 Tropical Teal |
| <input type="checkbox"/> 2505022 | Leaf Green | <input type="checkbox"/> 2505019 Tangerine Orange |
| | <input type="checkbox"/> 2505021 Lush Lilac | <input type="checkbox"/> 2505011 Lollypop Red |
| | <input type="checkbox"/> 2505023 Boysenberry Purple | <input type="checkbox"/> 2505012 Bubble Gum Pink |
| <input type="checkbox"/> 2505200 | Custom Color <small>Contact PDG for specifications.</small> | <input type="checkbox"/> 2505014 Daffodil Yellow |
| | | \$ 300 |

Armrests

2506000 **No Armrests**

Height Adjustable Armrests (Single Post) \$ 295

2506111 8"-11" Range 2506114 10"-14" Range 2506116 12"-16" Range

2506100 **Fixed Height Armrests (12")** Full Length Armpad only.

Armrest Height If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.
 Seat Pan to top of Armrest: _____ Does not include height of Armpad

Armrest Position Indicate preference for Armrest position. If left blank, the default setting will be used.
 In the interest of maintaining chair stably and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

2506401 Rear 2506402 Mid 2506403 Front

Armpads

Regular Armpads

2507100 Desk Length - 2" x 10" 2507101 Full Length - 2" x 14"

Waterfall Armpads \$ 115

2507105 Desk Length - 2" x 10" 2507106 Full Length - 2" x 14"

Wide Armpads May interfere with wheels in some tilt configurations. \$ 200

2507110 Desk Length - 3" x 10" 2507111 Full Length - 3" x 12"

Gel Armpads May interfere with wheels in some tilt configurations. \$ 135

2507120 Full Length - 2" x 12" 2507121 Full Length - 3.5" x 12"

Contour Arm Trough with Hand Pad May interfere with wheels in some tilt configurations. \$ 125

2507130 Right (4.5" x 18") 2507131 Left (4.5" x 18")

Armrest Options

Armrest Wideners Move Armrests outward by up to 1" on each side. Seat size does not change. May cause tilt interference.

2507210 1" Total Additional Width \$ 95

2507211 2" Total Additional Width \$ 195

2507310 **Clothing Guards** Not compatible with 8"-11" Height range. \$ 110

Armrest Extensions Not compatible with Desk Length Armpads. \$ 75

2507400 3" Extension - Forward 2507401 3" Extension - Rearward

Seat Pan

2508000 **No Seat Pan**

2508001 **Solid Aluminum Seat Pan** \$ 135

Seat Height

Don't forget **Select desired Seat-to-Floor Height** STF Height is adjustable within the range of the chosen Subframe Plate (Low/High).

| 6" Casters | Low | | | | | | | High | | | | | | |
|------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|
| | RW | 13" | 14" | 15" | 16" | 17" | 18" | RW | 15" | 16" | 17" | 18" | 19" | 20" |
| 12" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | 12" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 16" | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | 16" | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 20" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | 20" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 22" | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | 22" | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 24" | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | 24" | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |

| 8" Casters | Low | | | | | | | High | | | | | | |
|------------|-----|-----|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|------|-----|-----|--------------------------|--------------------------|--------------------------|--------------------------|
| | RW | 13" | 14" | 15" | 16" | 17" | 18" | RW | 15" | 16" | 17" | 18" | 19" | 20" |
| 16" | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | 16" | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 24" | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | 24" | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |

Casters

Urethane Casters 8" Casters are compatible only with 16" and 24" Rear Wheels.

| | | | | |
|----------------------------------|------------------------|----------------------------------|------------------------|--------|
| <input type="checkbox"/> 2511061 | 6" x 1.25" Caster size | <input type="checkbox"/> 2511081 | 8" x 1.25" Caster size | |
| <input type="checkbox"/> 2511063 | 6" x 1.50" Caster size | | | \$ 85 |
| <input type="checkbox"/> 2511084 | 8" x 2.00" Caster size | | | \$ 135 |

Pneumatic Casters 8" Casters are compatible only with 16" and 24" Rear Wheels.

| | | | | |
|----------------------------------|------------------------|--|--|--------|
| <input type="checkbox"/> 2511161 | 6" x 1.25" Caster size | | | \$ 110 |
| <input type="checkbox"/> 2511181 | 8" x 1.25" Caster size | | | \$ 120 |
| <input type="checkbox"/> 2511184 | 8" x 2.00" Caster size | | | \$ 150 |

| | | | | |
|----------------------------------|--|--|--|--------|
| <input type="checkbox"/> 2511260 | Frog Legs 6" x 1.50" Aluminum Soft-Roll Casters | | | \$ 275 |
|----------------------------------|--|--|--|--------|

| | | | | |
|----------------------------------|-----------------------------------|--|--|--------|
| <input type="checkbox"/> 2511400 | Frog Legs Suspension Forks | Only compatible with 16 and 24" Rear Wheels. Requires 6" Soft-Roll Caster. | | \$ 695 |
|----------------------------------|-----------------------------------|--|--|--------|

Rear Wheels

Composite Wheels

| | | | | | |
|----------------------------------|----------------|----------------------------------|----------------|----------------------------------|----------------|
| <input type="checkbox"/> 2512012 | 12" Wheel size | <input type="checkbox"/> 2512016 | 16" Wheel size | <input type="checkbox"/> 2512020 | 20" Wheel size |
| <input type="checkbox"/> 2512022 | 22" Wheel size | <input type="checkbox"/> 2512024 | 24" Wheel size | | |

Spoke Wheels

| | | | | | | |
|----------------------------------|----------------|----------------------------------|----------------|----------------------------------|----------------|-------|
| <input type="checkbox"/> 2512120 | 20" Wheel size | <input type="checkbox"/> 2512122 | 22" Wheel size | <input type="checkbox"/> 2512124 | 24" Wheel size | \$ 75 |
|----------------------------------|----------------|----------------------------------|----------------|----------------------------------|----------------|-------|



Rear Wheel Position Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.

In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

| | | | | | |
|----------------------------------|------|----------------------------------|-----|----------------------------------|-------|
| <input type="checkbox"/> 2512401 | Rear | <input type="checkbox"/> 2512402 | Mid | <input type="checkbox"/> 2512403 | Front |
|----------------------------------|------|----------------------------------|-----|----------------------------------|-------|

Axle Type

| | |
|----------------------------------|---------------------------------|
| <input type="checkbox"/> 2513000 | Quick Release Rear Axles |
| <input type="checkbox"/> 2513010 | Fixed Rear Axles |

Rear Tires

| | | | |
|----------------------------------|---|---|-------|
| <input type="checkbox"/> 2514000 | Urethane Tires | | |
| <input type="checkbox"/> 2514010 | Pneumatic Tires | May result in a finished Seat Height 1/2" greater than stated. | \$ 60 |
| <input type="checkbox"/> 2514020 | Pneumatic Tires with Flat-Free Inserts | May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels. | \$ 95 |

Handrims

Handrims are only available with 20", 22", and 24" Rear Wheels.

| | | | |
|----------------------------------|-----------------------------------|--|--------|
| <input type="checkbox"/> 2515000 | No Handrims | | |
| <input type="checkbox"/> 2515010 | Standard Aluminum Handrims | | |
| <input type="checkbox"/> 2515020 | Plastic Coated Handrims | Not compatible with Drum Brakes. | \$ 150 |
| <input type="checkbox"/> 2515030 | Projection Handrims | Composite Wheels only. 10 Projections, 30°, 1-1.8" Length. | \$ 595 |

Wheel Locks & Brakes

User Wheel Locks

| | | | | |
|----------------------------------|---|---|------------------|--------|
| <input type="checkbox"/> 2516010 | Push to Lock | <input type="checkbox"/> 2516020 | Pull to Lock | |
| <input type="checkbox"/> 2516030 | Attendant Wheel Locks (Push to Lock) | | | \$ 85 |
| <input type="checkbox"/> 2516040 | User & Attendant Wheel Locks | Not compatible with Drum Brakes. | | \$ 125 |
| <input type="checkbox"/> 2516300 | Drum Brakes | Available with 12" Composite or 20" Spoke Rear Wheels. Urethane Tires only. | | \$ 995 |
| | User-activated Unilateral Wheel Lock | Pull to Lock. | | \$ 995 |
| <input type="checkbox"/> 2516310 | Left Activation | <input type="checkbox"/> 2516311 | Right Activation | |



Wheel Lock Options

- 2516100 **Wheel Lock Extensions** Not compatible with Attendant Wheel Lock. \$ 60
- 2516200 **Hill Holders** Requires Pneumatic Rear Tires with Flat-Free Inserts. \$ 450

Anti-Tippers

- 2517010 **Rear Anti-Tippers** Required Safety Option. \$ 75

Back Posts

Back Post Style

- 2518010 Stroller Handle Back

Back Post Height

 Measured from the Seat Pan to the Handle.

- 2518425 Short (22") 2518431 Tall (25")

Back Post Folding

 Fold-down Back is not compatible with Reclining Backrest.

- 2519220 Fold-down

Seat to Back Angle

- 2519090 Set at 90° 2519098 Set at 98° 2519106 Set at 106°
- 2519114 Set at 114° 2519122 Set at 122°

Back Options

Angle/Height Adjustable Stroller Handle

 Not compatible with Headrest Mounting Plate. \$ 150

- 2519010 Button on Outside 2519011 Button on Inside

Endomorph Back

 Backrest will be narrower than Seat by total of 2x the Offset selected below. Minimum width of Backrest is 12". \$ 395

- 2518102 1" Offset per side 2518103 1.5" Offset per side 2518104 2" Offset per side

- 2519210 **10" Backrest Extension** Not compatible with Angle/Height Adjustable Stroller Handle. \$ 195

- 2519212 **Headrest Mounting Plate** \$ 95

Not compatible with Backrest Extension, Angle/Height Adjustable Stroller Handle. Not required if mounting to a solid Back.

Back Upholstery

- 2520000 **No Back Upholstery** \$ 265

Standard Upholstery

- 2520418 Short (18" Height) 2520421 Regular (21" Height) 2520424 Tall (24" Height)

- 2520450 **Tension Adjustable Upholstery** Height of Upholstery is 16". \$ 365

Position of Upholstery

 Select the desired position of the top of the Upholstery relative to the Seat Pan.

| Style | 22" Back Post | | | | | | | | |
|-------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|-----|-----|-----|
| | 17" | 18" | 19" | 20" | 21" | 22" | 23" | 24" | 25" |
| Short | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | |
| Regular | | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | |
| Tension Adj | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | |

| Style | 25" Back Post | | | | | | | | |
|-------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|
| | 17" | 18" | 19" | 20" | 21" | 22" | 23" | 24" | 25" |
| Short | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| Regular | | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| Tall | | | | | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| Tension Adj | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |

Tilt Activation

Attendant Tilt only

 Select trigger mounting location below.

- 2523011 Right Back Post 2523012 Left Back Post

User Tilt only

 Not available on widths greater than 22". Full Length Arm Pads only. Select trigger mounting location below. \$ 100

- 2523021 Right Armrest 2523022 Left Armrest

Attendant & User Tilt

 Not available on widths greater than 22". Full Length Arm Pads only. Select trigger mounting locations below. \$ 195

- 2523031 Right Back Post 2523032 Left Back Post
- 2523041 Right Armrest 2523042 Left Armrest

See next page for more options >



| | | | | | | | | |
|---|---------|--------------------------------------|--|---------|------------------------------|--|---------|------------------------|
| <input type="checkbox"/> | 2523110 | Tilt Trigger Extension | | | | \$ 50 | | |
| Tilt Limit Optional. If no tilt limit is selected, the full tilt range will be available. | | | | | | | | |
| <input type="checkbox"/> | 2523220 | Set at 20° | <input type="checkbox"/> | 2523230 | Set at 30° | <input type="checkbox"/> | 2523240 | Set at 40° |
| Power Tilt 250 lbs Weight Capacity. 45° tilt range. Includes 2 batteries for off-chair charging. pdgmobility.com/power-tilt Minimum Seat Width is 16". Power Tilt is compatible with either O2 Holder or IV Pole on the same chair, not both. | | | | | | | | |
| <input type="checkbox"/> | 2523400 | Standard Power Tilt | | | | \$ 2,095 | | |
| <input type="checkbox"/> | 2523410 | Power Tilt with Manual Tilt Override | Allows for instant manual release of tilt in case of power tilt failure. | | | \$ 2,695 | | |
| ➤ Select Control Switch | | Toggle Switch | <input type="checkbox"/> | 2523421 | Installed on Right side | <input type="checkbox"/> | 2523422 | Installed on Left side |
| | | Button Switch | <input type="checkbox"/> | 2523431 | Installed on Right side | <input type="checkbox"/> | 2523432 | Installed on Left side |
| | | Stealth Egg Switch | <input type="checkbox"/> | 2523440 | Not installed | | | |
| | | Optional for Egg Switch: | <input type="checkbox"/> | 2523445 | Swing Away Mounting Hardware | Not compatible with all Headrest makes and models. | | \$ 625 |
| <i>Specify Headrest to be used:</i> _____ | | | | | | | | |

| Recline | | | | | | |
|--|---------|------------------------------------|--|--|--|--------|
| <input type="checkbox"/> | 2523310 | Reclining Backrest (Manual) | | | | \$ 895 |
| 30° maximum recline. Not compatible with Backrest Extension, Fold-down Back, Tension Adjustable Upholstery. Finished Back Post Height is 3" taller than selected Back Post Height. Limits Tilt by up to 20°. | | | | | | |

| Front Rigging | | | | | | |
|---|---------|---|--|-------------------------|---------------------------|--------|
| <input type="checkbox"/> | 2525000 | No Front Rigging | <input type="checkbox"/> | 2525001 | Sector Blocks only | |
| Easy-Out™ Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame. pdgmobility.com/easy-out | | | | | | |
| <input type="checkbox"/> | 2525410 | Fixed 70° | | | | \$ 95 |
| <input type="checkbox"/> | 2525420 | Universal Angle Adjustable (0°-110°) | Set at angle: _____ | (use increments of 10°) | | \$ 295 |
| Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings. | | | | | | |
| Swing Away Front Rigging Legacy Design. Not compatible with Easy-Out. | | | | | | |
| <input type="checkbox"/> | 2525130 | Swing Away 60° | <input type="checkbox"/> | 2525140 | Swing Away 70° | |
| <input type="checkbox"/> | 2525200 | Swing Away Elevating Legrests | Includes Calf Pads. | | | \$ 375 |
| <input type="checkbox"/> | 2525300 | 1" Offset Hangers for Swing Away | | | | \$ 75 |
| <input type="checkbox"/> | 2525110 | Standard Plug In 70° Front Rigging | Legacy Design. Not compatible with Easy-Out. | | | |

| Footplates | | | | | | |
|--------------------------|---------|---|--|--|--|--------|
| <input type="checkbox"/> | 2526010 | Composite Footplates | | | | |
| <input type="checkbox"/> | 2526020 | Aluminum Angle Adjustable Footplates | | | | \$ 165 |
| <input type="checkbox"/> | 2526030 | One Piece Fixed Footplate | Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable. | | | \$ 350 |
| <input type="checkbox"/> | 2526040 | One Piece Flip-up Footplate | Requires selection of Composite or Aluminum Angle Adjustable. | | | \$ 395 |
| <input type="checkbox"/> | 2526050 | Interlocking Angle Adjustable Footplates | Available for Seat Widths up to 21" only. | | | \$ 195 |

| Footrest/Legrest Accessories | | | | | | |
|---|---------|--|---|--|--|--------|
| Calf Pads Not compatible with Gel Padding. | | | | | | |
| <input type="checkbox"/> | 2527012 | For Easy-Out™ Front Rigging | Includes reinforced premium Calf Pad mount. | | | \$ 200 |
| <input type="checkbox"/> | 2527010 | For Swing Away Front Rigging | | | | \$ 135 |
| <input type="checkbox"/> | 2527020 | Gel Padding | Not compatible with Calf Pad or Easy-Out Front Rigging. | | | \$ 125 |
| Gel Knee Pads | | | | | | |
| <input type="checkbox"/> | 2527032 | For Easy-Out™ Angle Adjustable Front Rigging | | | | \$ 200 |
| <input type="checkbox"/> | 2527030 | For Swing Away ELR Front Rigging | | | | \$ 100 |
| <input type="checkbox"/> | 2527040 | Length Adjustable Heel Loops | | | | \$ 50 |
| <input type="checkbox"/> | 2527050 | Bodypoint Padded Calf Strap | | | | \$ 195 |

Front Rigging Extension Tubes

Regular Extension Tubes

- | | | | | | | | |
|--------------------------|---------|-------|---------|----------------------------|---------|---|-------|
| <input type="checkbox"/> | 2528010 | Short | 14"-18" | for Easy-Out Front Rigging | 13"-17" | for Swing Away and Standard Plug In Front Rigging | |
| <input type="checkbox"/> | 2528020 | Long | 18"-24" | for Easy-Out Front Rigging | 15"-19" | for Swing Away and Standard Plug In Front Rigging | \$ 75 |



Lower Leg Length If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.

Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.

Miscellaneous Options

Spoke Guards Available for 20"-24" Composite Rear Wheels. Made from durable and lightweight aluminum. \$ 235

- | | | | | | |
|--------------------------|---------|-------|--------------------------|---------|-------------|
| <input type="checkbox"/> | 2530020 | Black | <input type="checkbox"/> | 2530021 | Frame Color |
|--------------------------|---------|-------|--------------------------|---------|-------------|

Stealth Flip Down Abductor Requires Solid Seat Pan. \$ 825

- | | | | | | | | |
|--------------------------|---------|-----------|------------------|--------------------------|---------|-----------|------------------|
| <input type="checkbox"/> | 2530060 | Small Pad | 3"x4"x3"x2" size | <input type="checkbox"/> | 2530062 | Large Pad | 5"x5"x4"x3" size |
|--------------------------|---------|-----------|------------------|--------------------------|---------|-----------|------------------|

- | | | | | | | |
|--------------------------|---------|-----------------------|-------------------------------------|--|--|--------|
| <input type="checkbox"/> | 2530080 | Gel Hip Guides | Reduces available Seat Width by 2". | | | \$ 595 |
|--------------------------|---------|-----------------------|-------------------------------------|--|--|--------|

- | | | | | | | |
|--------------------------|---------|----------------|--|--|--|--------|
| <input type="checkbox"/> | 2530110 | IV Pole | Not compatible with Attendant Wheel Locks. Requires frame modification and cannot be added to chair later. | | | \$ 345 |
|--------------------------|---------|----------------|--|--|--|--------|

- | | | | | | | |
|--------------------------|---------|------------------|--|--|--|--------|
| <input type="checkbox"/> | 2530120 | O2 Holder | Not compatible with Attendant Wheel Locks. Requires frame modification and cannot be added to chair later. | | | \$ 495 |
|--------------------------|---------|------------------|--|--|--|--------|

- | | | | | | | |
|--------------------------|---------|---------------------------|--|--|--|--------|
| <input type="checkbox"/> | 2530200 | Q'Straint Lap Belt | Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" Belt Length. | | | \$ 295 |
|--------------------------|---------|---------------------------|--|--|--|--------|

- | | | | | | | |
|--------------------------|---------|----------------------------|------------------------------------|--|--|--------|
| <input type="checkbox"/> | 2530222 | Wheelchair Backpack | Black with PDG Logo. 17"x14" size. | | | \$ 165 |
|--------------------------|---------|----------------------------|------------------------------------|--|--|--------|

- | | | | | | | |
|--------------------------|---------|---|-----------------------|--|--|-------|
| <input type="checkbox"/> | 2530170 | Adhesive Back Hook Strip, 2" x 30" | For Cushion mounting. | | | \$ 10 |
|--------------------------|---------|---|-----------------------|--|--|-------|

On-Chair Seating & Positioning

 Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. pdgmobility.com/seating

Order Notes

Please send your completed order form by email to info@pdgmobility.com.