



STELLAR HD

Order Form - Distributor

Prices in Canadian Dollars. January 1, 2024

Customer and Order Information

Purchase Order *	Tag	v.0101
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in

Base Chair

Stellar HD 45° Manual Tilt Wheelchair Extended delivery lead times may apply.

- 2101000 Stellar 45° Manual Tilt Wheelchair \$ 4,400
- 2101201 Reinforced HD Frame Package Heavy-gauge steel reinforcement of frame for bariatric users. Required selection. \$ 1,395
- Choose the **Stellar Impact** for highly agitated users up to 250 lbs. See separate order form.

Weight Capacity

Chair Weight Capacity

- 2102045 450 lbs Capacity
- 2102060 600 lbs Capacity Tilt is limited to maximum 25°. \$ 300

User Weight Range Required for setup of tilt uplift force.

- 2102212 Up to 250 lbs
- 2102308 250 lbs to 350 lbs
- 2102310 350 lbs to 425 lbs
- 2102312 425 lbs to 600 lbs

Seat Width

Overall chair width is +11"

Standard Seat Widths

- 2103022 22" Seat Width
- 2103024 24" Seat Width
- 2103026 26" Seat Width
- 2103028 28" Seat Width
- 2103030 30" Seat Width
- 2103032 32" Seat Width

Seat Depth

Standard Seat Depths (Adjustable)

- 2104018 18" Seat Depth
- 2104019 19" Seat Depth
- 2104020 20" Seat Depth
- 2104021 21" Seat Depth
- 2104022 22" Seat Depth

Frame Color



PDG Premium Colors View color swatches on PDG's website. pdgmobility.com/colors

- 2105001 Black Shimmer
- 2105003 Fine Texture Black
- 2105017 Arctic Grey
- 2105020 Olive Green
- 2105022 Leaf Green
- 2105015 Sky Blue
- 2105010 Lollypop Blue
- 2105018 Tropical Teal
- 2105021 Lush Lilac
- 2105023 Boysenberry Purple
- 2105016 Poppy Red
- 2105011 Lollypop Red
- 2105019 Tangerine Orange
- 2105012 Bubble Gum Pink
- 2105014 Daffodil Yellow
- 2105200 **Custom Color** Contact PDG for specifications. \$ 300

Armrests

- 2106000 **No Armrests**
- Single Post Height Adjustable Armrests** \$ 295
- 2106111 8"-11" Range
- 2106114 10"-14" Range
- 2106116 12"-16" Range

See next page for more options >

Dual Post Height Adjustable Armrests			\$ 295
<input type="checkbox"/> 2106211 8"-11" Range	<input type="checkbox"/> 2106214 10"-14" Range	<input type="checkbox"/> 2106216 12"-16" Range	
Fixed Height Armrests (12") Full Length Armpad only.			
<input type="checkbox"/> 2106100 Single Post	<input type="checkbox"/> 2106200 Dual Post		
 Armrest Height If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used. Seat Pan to top of Armrest: _____ Does not include height of Armpad			
 Armrest Position Indicate preference for Armrest position. If left blank, the default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.			
<input type="checkbox"/> 2106401 Rear	<input type="checkbox"/> 2106402 Mid	<input type="checkbox"/> 2106403 Front	


Armpads			
Regular Armpads			
<input type="checkbox"/> 2107100 Desk Length - 2" x 10"	<input type="checkbox"/> 2107101 Full Length - 2" x 14"		
Waterfall Armpads Not compatible with Dual Post Armrests			
<input type="checkbox"/> 2107105 Desk Length - 2" x 10"	<input type="checkbox"/> 2107106 Full Length - 2" x 14"		
Wide Armpads May interfere with wheels in some tilt configurations.			
<input type="checkbox"/> 2107110 Desk Length - 3" x 10"	<input type="checkbox"/> 2107111 Full Length - 3" x 12"		
Gel Armpads May interfere with wheels in some tilt configurations.			
<input type="checkbox"/> 2107120 Full Length - 2" x 12"	<input type="checkbox"/> 2107121 Full Length - 3.5" x 12"		
Contour Arm Trough with Hand Pad May interfere with wheels in some tilt configurations.			
<input type="checkbox"/> 2107130 Right (4.5" x 18")	<input type="checkbox"/> 2107131 Left (4.5" x 18")		

Armrest Options			
Armrest Wideners Move Armrests outward by up to 1" on each side. Seat size does not change. May cause tilt interference.			
<input type="checkbox"/> 2107210 1" Total Additional Width			
<input type="checkbox"/> 2107211 2" Total Additional Width			
Clothing Guards			
<input type="checkbox"/> 2107310 For Single Post Armrests	Not compatible with 8"-11" Height range.		
<input type="checkbox"/> 2107320 For Dual Post Armrests			
Armrest Extensions Not compatible with Desk Length Armpads.			
<input type="checkbox"/> 2107400 3" Extension - Forward	<input type="checkbox"/> 2107401 3" Extension - Rearward		

Seat Pan			
<input checked="" type="checkbox"/> 2108001 Solid Aluminum Seat Pan			
	\$ 175		

Seat Height			
Seat to Floor Height			
<input type="checkbox"/> 2110016 Set at 16" STF	<input type="checkbox"/> 2110017 Set at 17" STF	<input type="checkbox"/> 2110018 Set at 18" STF	

Casters			
<input checked="" type="checkbox"/> 2111050 5" Casters with Heavy Duty Steel Fork			

Rear Wheels			
<input checked="" type="checkbox"/> 2112224 24" Heavy Duty Spoke Wheels			
 Rear Wheel Position Indicate Rear Wheel positioning preference. If left blank, the default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.			
<input type="checkbox"/> 2112401 Rear	<input type="checkbox"/> 2112402 Mid	<input type="checkbox"/> 2112403 Front	

Axle Type			
<input checked="" type="checkbox"/> 2113010 Fixed Rear Axles			



Rear Tires

2114000 **Urethane Tires**

Handrims

2115000 **No Handrims**

2115010 **Standard Aluminum Handrims**

2115020 **Plastic Coated Handrims** \$ 250

Wheel Locks

User Wheel Locks

2116010 Push to Lock 2116020 Pull to Lock

Wheel Lock Options

2116100 **Wheel Lock Extensions** \$ 60

Anti-Tippers

2117010 **Rear Anti-Tippers** Required Safety Option. \$ 75

Back Post

2118010 **Stroller Handle Back**

Back Post Height Measured from the Seat Pan to the Handle.

2118425 Short (22") 2118431 Tall (25")

2119090 **Seat to Back Angle set at 90°**

Back Options

Endomorph Back Backrest will be narrower than Seat by total of 2x the Offset selected below. Minimum width of Backrest is 12". \$ 395

2118102 1" Offset per side 2118103 1.5" Offset per side 2118104 2" Offset per side

2119212 **Headrest Mounting Plate** Not required if mounting to a solid Back. \$ 95

2119240 **Seating Dynamics Dynamic Rocker Back Interface** \$ 2,995

Not compatible with 30° Recline, 21" & 22" Seat Depths. Finished Back Post Height is 4" taller than selected Back Post Height.

Back Upholstery

2120000 **No Back Upholstery**

Tension Adjustable Upholstery \$ 525

2120421 Regular (18" Height) 2120424 Tall (22" Height)

Position of Upholstery Select the desired position of the top of the Upholstery relative to the Seat Pan.

		22" Back Post								
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"	
Regular	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				
Tall				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				

		25" Back Post								
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"	
Regular	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Tall					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Tilt Activation

2123010 **Attendant Tilt only** Two tilt triggers (L and R) are required for tilt activation.

Tilt Limit Optional. If no tilt limit is selected, the full tilt range will be available. Maximum Tilt Limit for 600 lbs Weight Capacity is 25°.

2123220 Set at 20° 2123230 Set at 30° 2123240 Set at 40°

Recline

2123310 **Reclining Backrest (Manual)** \$ 895

30° maximum recline. Not compatible with 600 lbs Weight Capacity, Tension Adjustable Upholstery. Finished Back Post Height is 3" taller than selected Back Post

Front Rigging

<input type="checkbox"/>	2125000	No Front Rigging	<input type="checkbox"/>	2125001	Sector Blocks only	
Easy-Out™ Heavy Duty Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame. pdgmobility.com/easy-out						
<input type="checkbox"/>	2125410	Fixed 70°				\$ 95
<input type="checkbox"/>	2125420	Universal Angle Adjustable (0°-110°)	Set at angle: _____ (use increments of 10°)			\$ 295
Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings.						
Swing Away Front Rigging Legacy Design. Not compatible with Easy-Out.						
<input type="checkbox"/>	2125130	Swing Away 60°	<input type="checkbox"/>	2125140	Swing Away 70°	
<input type="checkbox"/>	2125200	Swing Away Elevating Legrests	Includes Calf Pads.			\$ 375
<input type="checkbox"/>	2125300	1" Offset Hangers for Swing Away				\$ 75
<input type="checkbox"/>	2125110	Standard Plug In 70°	Legacy Design. Not compatible with Easy-Out.			


Footplates

<input checked="" type="checkbox"/>	2126020	Aluminum Angle Adjustable Footplates				\$ 175
-------------------------------------	---------	---	--	--	--	--------

Footrest/Legrest Accessories

Calf Pads Not compatible with Gel Padding.						
<input type="checkbox"/>	2127012	For Easy-Out™ Front Rigging	Includes reinforced premium Calf Pad mount.			\$ 200
<input type="checkbox"/>	2127010	For Swing Away Front Rigging				\$ 135
<input type="checkbox"/>	2127020	Gel Padding	Not compatible with Calf Pad or Easy-Out Front Rigging.			\$ 125
Gel Knee Pads						
<input type="checkbox"/>	2127032	For Easy-Out™ Angle Adjustable Front Rigging				\$ 200
<input type="checkbox"/>	2127030	For Swing Away ELR Front Rigging				\$ 100
<input type="checkbox"/>	2127040	Length Adjustable Heel Loops				\$ 50
<input type="checkbox"/>	2127051	Bodypoint Calf Strap				\$ 195

Front Rigging Extension Tubes

Regular Extension Tubes						
<input type="checkbox"/>	2128010	Short	14"-18" for Easy-Out Front Rigging	13"-17" for Swing Away and Standard Plug In Front Rigging		
<input type="checkbox"/>	2128020	Long	18"-24" for Easy-Out Front Rigging	15"-19" for Swing Away and Standard Plug In Front Rigging		\$ 75
Extra Strength Extension Tubes (Solid) Available on Plug In Front Rigging only. Not compatible with Easy-Out Front Rigging.						
<input type="checkbox"/>	2128110	Standard (13"-17")	<input type="checkbox"/>	2128120	Long (15"-19")	\$ 100
	Lower Leg Length If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.					
	Seat Pan to Foot Plates: _____		Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.			

Miscellaneous Options

<input type="checkbox"/>	2130080	Gel Hip Guides	Reduces available Seat Width by 2".			\$ 595
<input type="checkbox"/>	2130110	IV Pole				\$ 345
<input type="checkbox"/>	2130120	O2 Holder				\$ 495
<input type="checkbox"/>	2130222	Wheelchair Backpack	Black with PDG Logo. 17"x14" size.			\$ 165
<input type="checkbox"/>	2130170	Adhesive Back Hook Strip, 2" x 30"	For Cushion mounting.			\$ 10

On-Chair Seating & Positioning

	Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. pdgmobility.com/seating					
---	---	--	--	--	--	--

Order Notes

Please send your completed order form by email to info@pdgmobility.com.