



| Customer and Order Information | | |
|--|------------------|---|
| Purchase Order * | Tag | v.0401 |
| Dealer/Provider Name * | | |
| Purchaser Name and Contact Information * | | |
| Shipping Address * | | |
| User Information (if available) | Weight _____ lbs | Height _____ ft/in *required |

| Base Chair | | |
|--|---|----------|
| Stellar Impact 45° Manual Tilt Wheelchair Extended delivery lead times may apply. | | |
| <input checked="" type="checkbox"/> 2301000 | Stellar 45° Manual Tilt Wheelchair | \$ 4,400 |
| <input checked="" type="checkbox"/> 2301201 | Reinforced Impact Frame Package Ultra-durable steel reinforcement of frame for highly agitated users. Required selection. | \$ 1,095 |
| i Choose the Stellar GL or the Stellar HD for users over 250 lbs. Custom reinforcement available. See separate order forms. | | |



| Weight Capacity | | |
|---|--------------------------------------|--|
| <input checked="" type="checkbox"/> 2302000 | 250 lbs Chair Weight Capacity | |

| Frame Options | | |
|----------------------------------|---|--------|
| <input type="checkbox"/> 2302200 | Extended Base with Rear Floor Jacks Lifts Rear Wheels off ground for increased stability. Limited ground clearance. Not compatible with Attendant Wheel Lock, HD Anti-Tipper, IV Pole, O2 Holder. | \$ 900 |
| <input type="checkbox"/> 2302300 | Transit Tie Down System (TTDS) TTDS includes 4 factory-installed transport brackets. See Owners Manual for complete details. | \$ 250 |

| Seat Width | | Overall chair width is +9.5" for 12/16" RW and +10.5" for 20/22/24" RW |
|----------------------------------|----------------|--|
| Standard Seat Widths | | |
| <input type="checkbox"/> 2303014 | 14" Seat Width | <input type="checkbox"/> 2303015 15" Seat Width |
| <input type="checkbox"/> 2303017 | 17" Seat Width | <input type="checkbox"/> 2303018 18" Seat Width |
| <input type="checkbox"/> 2303020 | 20" Seat Width | <input type="checkbox"/> 2303021 21" Seat Width |
| | | <input type="checkbox"/> 2303016 16" Seat Width |
| | | <input type="checkbox"/> 2303019 19" Seat Width |
| | | <input type="checkbox"/> 2303022 22" Seat Width |

| Seat Depth | | |
|--|----------------|---|
| Standard Seat Depths (Adjustable) | | |
| <input type="checkbox"/> 2304017 | 17" Seat Depth | <input type="checkbox"/> 2304018 18" Seat Depth |
| <input type="checkbox"/> 2304020 | 20" Seat Depth | <input type="checkbox"/> 2304019 19" Seat Depth |
| Custom Seat Depths | | |
| <input type="checkbox"/> 2304021 | 21" Seat Depth | <input type="checkbox"/> 2304022 22" Seat Depth |
| | | \$ 400 |

| Frame Color | | |
|--|---|---|
| PDG Premium Colors View color swatches on PDG's website. pdgmobility.com/colors | | |
| <input type="checkbox"/> 2305001 | Black Shimmer | <input type="checkbox"/> 2305015 Sky Blue |
| <input type="checkbox"/> 2305003 | Fine Texture Black | <input type="checkbox"/> 2305010 Lollypop Blue |
| <input type="checkbox"/> 2305017 | Arctic Grey | <input type="checkbox"/> 2305018 Tropical Teal |
| <input type="checkbox"/> 2305020 | Olive Green | <input type="checkbox"/> 2305021 Lush Lilac |
| <input type="checkbox"/> 2305022 | Leaf Green | <input type="checkbox"/> 2305023 Boysenberry Purple |
| | | <input type="checkbox"/> 2305016 Poppy Red |
| | | <input type="checkbox"/> 2305011 Lollypop Red |
| | | <input type="checkbox"/> 2305019 Tangerine Orange |
| | | <input type="checkbox"/> 2305012 Bubble Gum Pink |
| | | <input type="checkbox"/> 2305014 Daffodil Yellow |
| <input type="checkbox"/> 2305200 | Custom Color Contact PDG for specifications. | \$ 300 |

| Armrests | | | | |
|--|---|--------------------|-----------------------------------|-----------------------|
| <input type="checkbox"/> | 2306000 | No Armrests | | |
| Single Post Height Adjustable Armrests | | | \$ | 295 |
| <input type="checkbox"/> | 2306111 | 8"-11" Range | <input type="checkbox"/> | 2306114 10"-14" Range |
| <input type="checkbox"/> | 2306116 | 12"-16" Range | | |
| Dual Post Height Adjustable Armrests | | | \$ | 295 |
| <input type="checkbox"/> | 2306211 | 8"-11" Range | <input type="checkbox"/> | 2306214 10"-14" Range |
| <input type="checkbox"/> | 2306216 | 12"-16" Range | | |
| Fixed Height Armrests (12") Full Length Armpad only. | | | | |
| <input type="checkbox"/> | 2306100 | Single Post | <input type="checkbox"/> | 2306200 Dual Post |
|  | Armrest Height If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used. | | | |
| | Seat Pan to top of Armrest: _____ | | Does not include height of Armpad | |
|  | Armrest Position Indicate preference for Armrest position. If left blank, the default setting will be used. | | | |
| In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference. | | | | |
| <input type="checkbox"/> | 2306401 | Rear | <input type="checkbox"/> | 2306402 Mid |
| <input type="checkbox"/> | 2306403 | Front | | |

| Armpads | | | | |
|--|---------|------------------------|--------------------------|----------------------------------|
| Regular Armpads | | | | |
| <input type="checkbox"/> | 2307100 | Desk Length - 2" x 10" | <input type="checkbox"/> | 2307101 Full Length - 2" x 14" |
| Waterfall Armpads Not compatible with Dual Post Armrests | | | \$ | 115 |
| <input type="checkbox"/> | 2307105 | Desk Length - 2" x 10" | <input type="checkbox"/> | 2307106 Full Length - 2" x 14" |
| Wide Armpads May interfere with wheels in some tilt configurations. | | | \$ | 200 |
| <input type="checkbox"/> | 2307110 | Desk Length - 3" x 10" | <input type="checkbox"/> | 2307111 Full Length - 3" x 12" |
| Gel Armpads May interfere with wheels in some tilt configurations. | | | \$ | 135 |
| <input type="checkbox"/> | 2307120 | Full Length - 2" x 12" | <input type="checkbox"/> | 2307121 Full Length - 3.5" x 12" |
| Contour Arm Trough with Hand Pad May interfere with wheels in some tilt configurations. | | | \$ | 125 |
| <input type="checkbox"/> | 2307130 | Right (4.5" x 18") | <input type="checkbox"/> | 2307131 Left (4.5" x 18") |

| Armrest Options | | | | | | |
|---|---------|---------------------------|--|--------------------------|---------------------------------|-----|
| Armrest Wideners Move Armrests outward by up to 1" on each side. Seat size does not change. May cause tilt interference. | | | | | | |
| <input type="checkbox"/> | 2307210 | 1" Total Additional Width | | \$ | 95 | |
| <input type="checkbox"/> | 2307211 | 2" Total Additional Width | | \$ | 195 | |
| Clothing Guards | | | | | | |
| <input type="checkbox"/> | 2307310 | For Single Post Armrests | Not compatible with 8"-11" Height range. | | \$ | 110 |
| <input type="checkbox"/> | 2307320 | For Dual Post Armrests | | \$ | 165 | |
| Armrest Extensions Not compatible with Desk Length Armpads. | | | \$ | 75 | | |
| <input type="checkbox"/> | 2307400 | 3" Extension - Forward | | <input type="checkbox"/> | 2307401 3" Extension - Rearward | |

| Seat Pan | | | | | |
|--------------------------|---------|--------------------------------|--|----|-----|
| <input type="checkbox"/> | 2308000 | No Seat Pan | | | |
| <input type="checkbox"/> | 2308001 | Solid Aluminum Seat Pan | | \$ | 135 |

Seat Height



Select desired Seat-to-Floor Height STF Height is adjustable within the range of the chosen Lower Frame type (Low/High).

| 6" Casters | Low | | | | | | High | | | | | |
|------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|
| | RW | 14" | 15" | 16" | 17" | 18" | RW | 17" | 18" | 19" | 20" | 21" |
| | 12" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | 12" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 16" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | 16" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 20" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | 20" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 22" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | 22" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 24" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | 24" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |

| 8" Casters | Low | | | | | | High | | | | | |
|------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|
| | RW | 14" | 15" | 16" | 17" | 18" | RW | 17" | 18" | 19" | 20" | 21" |
| | 16" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | 16" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 24" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | 24" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |

Casters

Urethane Casters 8" Casters are compatible only with 16" and 24" Rear Wheels.

- 2311061 6" x 1.25" Caster size 2311081 8" x 1.25" Caster size
- 2311063 6" x 1.50" Caster size \$ 85
- 2311084 8" x 2.00" Caster size \$ 135

Pneumatic Casters 8" Casters are compatible only with 16" and 24" Rear Wheels.

- 2311161 6" x 1.25" Caster size \$ 110
- 2311181 8" x 1.25" Caster size \$ 120
- 2311184 8" x 2.00" Caster size \$ 150

- 2311260 **Frog Legs 6" x 1.50" Aluminum Soft-Roll Casters** \$ 275

- 2311400 **Frog Legs Suspension Forks** Only compatible with 16 and 24" Rear Wheels. Requires 6" Soft-Roll Caster. \$ 695

Rear Wheels

Composite Wheels

- 2312012 12" Wheel size 2312016 16" Wheel size 2312020 20" Wheel size
- 2312022 22" Wheel size 2312024 24" Wheel size

Spoke Wheels \$ 75

- 2312120 20" Wheel size 2312122 22" Wheel size 2312124 24" Wheel size

- 2312224 **24" Heavy Duty Spoke Wheels** \$ 350



Rear Wheel Position Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.

In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

- 2312401 Rear 2312402 Mid 2312403 Front

Axle Type

- 2313000 **Quick Release Rear Axles**

- 2313010 **Fixed Rear Axles**

Rear Tires

- 2314000 **Urethane Tires**

- 2314010 **Pneumatic Tires** May result in a finished Seat Height 1/2" greater than stated. \$ 60

- 2314020 **Pneumatic Tires with Flat-Free Inserts** May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels. \$ 95



| Handrims | | Handrims are available only with 20", 22", and 24" Rear Wheels. |
|--------------------------|---|---|
| <input type="checkbox"/> | 2315000 No Handrims | |
| <input type="checkbox"/> | 2315010 Standard Aluminum Handrims | |
| <input type="checkbox"/> | 2315020 Plastic Coated Handrims Not compatible with Drum Brakes. | \$ 150 |
| <input type="checkbox"/> | 2315030 Projection Handrims Composite Wheels only. 10 Projections, 30°, 1-1.8" Length. | \$ 595 |

| Wheel Locks and Brakes | | | | |
|---|--|--------------------------|--------------------------|--------|
| User Wheel Locks | | | | |
| <input type="checkbox"/> | 2316010 Push to Lock | <input type="checkbox"/> | 2316020 Pull to Lock | |
| <input type="checkbox"/> | 2316030 Attendant Wheel Locks (Push to Lock) | | | \$ 85 |
| <input type="checkbox"/> | 2316040 User & Attendant Wheel Locks Not compatible with Drum Brakes. | | | \$ 125 |
| <input type="checkbox"/> | 2316300 Drum Brakes Available with 12" Composite or 20" Spoke Rear Wheels. Urethane Tires only. | | | \$ 995 |
| User-activated Unilateral Wheel Lock Pull to Lock. | | | | \$ 995 |
| <input type="checkbox"/> | 2316310 Left Activation | <input type="checkbox"/> | 2316311 Right Activation | |

| Wheel Lock Options | | |
|--------------------------|---|--------|
| <input type="checkbox"/> | 2316100 Wheel Lock Extensions Not compatible with Attendant Wheel Lock. | \$ 60 |
| <input type="checkbox"/> | 2316200 Hill Holders Requires Pneumatic Rear Tires with Flat-Free Inserts. | \$ 450 |

| Anti-Tippers | | |
|-------------------------------------|--|--------|
| <input checked="" type="checkbox"/> | 2317010 Rear Anti-Tippers Required Safety Option. | \$ 75 |
| <input type="checkbox"/> | 2317020 Front Anti-Tippers | \$ 795 |
| <input type="checkbox"/> | 2317030 Heavy Duty Rear Anti-Tippers Provide additional rearward stability for users with excessive movement. Not compatible with Attendant Wheel Locks, Rear Floor Jacks, IV Pole, O2 Holder. | \$ 895 |

| Back Post | | | | | |
|---|------------------------------|--------------------------|---|--------------------------|---------------------|
| Back Post Style | | | | | |
| <input type="checkbox"/> | 2318010 Stroller Handle Back | <input type="checkbox"/> | 2318020 Cane Handle Back May reduce tilt range. | | |
| Back Post Height Measured from the Seat Pan to the Handle. | | | | | |
| <input type="checkbox"/> | 2318425 Short (22") | <input type="checkbox"/> | 2318431 Tall (25") | | |
| Back Post Folding Folding Back is not compatible with Dynamic Rocker Back Interface, Reclining Backrest. | | | | | |
| <input type="checkbox"/> | 2319220 Fold-down | <input type="checkbox"/> | 2319221 Non-folding | | |
| Seat to Back Angle | | | | | |
| <input type="checkbox"/> | 2319090 Set at 90° | <input type="checkbox"/> | 2319098 Set at 98° | <input type="checkbox"/> | 2319106 Set at 106° |
| <input type="checkbox"/> | 2319114 Set at 114° | <input type="checkbox"/> | 2319122 Set at 122° | | |

| Back Options | | | | | |
|---|--|--------------------------|------------------------------|--------------------------|----------------------------|
| Angle/Height Adjustable Stroller Handle Not compatible with Headrest Mounting Plate or Cane Handle Back. | | \$ 150 | | | |
| <input type="checkbox"/> | 2319010 Button on Outside | <input type="checkbox"/> | 2319011 Button on Inside | | |
| Endomorph Back Backrest will be narrower than Seat by total of 2x the Offset selected below. Minimum width of Backrest is 12". | | \$ 395 | | | |
| <input type="checkbox"/> | 2318102 1" Offset per side | <input type="checkbox"/> | 2318103 1.5" Offset per side | <input type="checkbox"/> | 2318104 2" Offset per side |
| <input type="checkbox"/> | 2319210 10" Backrest Extension Not compatible with Angle/Height Adjustable Stroller Handle. | | | \$ 195 | |
| <input type="checkbox"/> | 23219212 Headrest Mounting Plate Not compatible with Backrest Extension, Angle/Height Adjustable Stroller Handle, Cane Handle Back. Not required if mounting to a solid Back. | | | \$ 95 | |
| <input checked="" type="checkbox"/> | 2319230 Spring Back Suspension Not compatible with Reclining Back. | | | | |
| <input type="checkbox"/> | 2319240 Seating Dynamics Dynamic Rocker Back Interface Not compatible with Fold-down Back, 21" & 22" Seat Depths, Spring Back Suspension, 30° Recline, Cane Handle, Standard Upholstery. Finished Back Post Height is 4" taller than selected Back Post Height. | | | \$ 2,995 | |



Back Upholstery

2320000 **No Back Upholstery**

Standard Upholstery Not available with Cane Handle Back. \$ 265

2320418 Short (18" Height) 2320421 Regular (21" Height) 2320424 Tall (24" Height)

2320450 **Tension Adjustable Upholstery** Height of Upholstery is 16". Not available with Cane Handle Back. \$ 365

Position of Upholstery Select the desired position of the top of the Upholstery relative to the Seat Pan.

| 22" Back Post | | | | | | | | | |
|---------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|-----|-----|-----|
| Style | 17" | 18" | 19" | 20" | 21" | 22" | 23" | 24" | 25" |
| Short | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | |
| Regular | | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | |
| Tension Adj | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | |

| 25" Back Post | | | | | | | | | |
|---------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|
| Style | 17" | 18" | 19" | 20" | 21" | 22" | 23" | 24" | 25" |
| Short | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| Regular | | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| Tall | | | | | | | | <input type="checkbox"/> | <input type="checkbox"/> |
| Tension Adj | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |

Tilt Activation

Attendant Tilt only Select trigger mounting location below.

2323011 Right Back Post 2323012 Left Back Post

User Tilt only Full Length Arm Pads only. Select trigger mounting location below. \$ 100

2323021 Right Armrest 2323022 Left Armrest

Attendant & User Tilt Full Length Arm Pads only. Select trigger mounting locations below. \$ 195

2323031 Right Back Post 2323032 Left Back Post

2323041 Right Armrest 2323042 Left Armrest

2323110 **Trigger Extension** \$ 50

Tilt Limit Optional. If no tilt limit is selected, the full tilt range will be available.

2323220 Set at 20° 2323230 Set at 30° 2323240 Set at 40°

Recline

2323310 **Reclining Backrest (Manual)** \$ 895

30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold-down Back, Tension Adjustable Upholstery. Finished Back Post Height is 3" taller than selected Back Post Height.

Front Rigging

2325000 **No Front Rigging** 2325001 **Sector Blocks only**

Easy-Out™ Heavy Duty Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame. pdgmobility.com/easy-out

2325410 Fixed 70° \$ 95

2325420 Universal Angle Adjustable (0°-110°) Set at angle: _____ (use increments of 10°) \$ 295
Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings.

Swing Away Front Rigging Legacy Design. Not compatible with Easy-Out.

2325130 Swing Away 60° 2325140 Swing Away 70°

2325200 Swing Away Elevating Legrests Includes Calf Pads. \$ 375

2325300 1" Offset Hangers for Swing Away \$ 75

2325110 **Standard Plug In 70°** Legacy Design. Not compatible with Easy-Out.

Footplates

2326010 **Composite Footplates**

2326020 **Aluminum Angle Adjustable Footplates** \$ 165

2326030 **One Piece Fixed Footplate** Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable. \$ 350


2326040 **One Piece Flip-up Footplate** Requires selection of Composite or Aluminum Angle Adjustable. \$ 395

2326050 **Interlocking Angle Adjustable Footplates** Available for Seat Widths up to 21" only. \$ 195

Footrest/Legrest Accessories

| | | | |
|---|---------|--|--------|
| Calf Pads Not compatible with Gel Padding. | | | |
| <input type="checkbox"/> | 2327012 | For Easy-Out™ Front Rigging Includes reinforced premium Calf Pad mount. | \$ 200 |
| <input type="checkbox"/> | 2327010 | For Swing Away Front Rigging | \$ 135 |
| <input type="checkbox"/> | 2327020 | Gel Padding Not compatible with Calf Pad or Easy-Out Front Rigging. | \$ 125 |
| Gel Knee Pads | | | |
| <input type="checkbox"/> | 2327032 | For Easy-Out™ Angle Adjustable Front Rigging | \$ 200 |
| <input type="checkbox"/> | 2327030 | For Swing Away ELR Front Rigging | \$ 100 |
| <input type="checkbox"/> | 2327040 | Length Adjustable Heel Loops | \$ 50 |
| <input type="checkbox"/> | 2327050 | Bodypoint Padded Calf Strap Sizes greater than 22" are not padded. | \$ 195 |
| Miller's Dynamic Footrest Mounts Not compatible with Easy-Out Front Rigging. | | | |
| <input type="checkbox"/> | 2327060 | For Standard Plug In Front Rigging | \$ 600 |
| <input type="checkbox"/> | 2327065 | For Swing Away 70° Front Rigging | \$ 720 |

Front Rigging Extension Tubes

| | | | |
|--|--|--|--------------------------|
| Regular Extension Tubes | | | |
| <input type="checkbox"/> | 2328010 | Short 14"-18" for Easy-Out Front Rigging 13"-17" for Swing Away and Standard Plug In Front Rigging | |
| <input type="checkbox"/> | 2328020 | Long 18"-24" for Easy-Out Front Rigging 15"-19" for Swing Away and Standard Plug In Front Rigging | \$ 75 |
| Extra Strength Extension Tubes (Solid) Available on Plug In Front Rigging only. Not compatible with Easy-Out Front Rigging. | | | \$ 100 |
| <input type="checkbox"/> | 2328110 | Standard (13"-17") | <input type="checkbox"/> |
| | | 2328120 | Long (15"-19") |
|  | Lower Leg Length If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used. | | |
| | Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion. | | |

Miscellaneous Options

| | | | |
|--|---------|--|--------------------------|
| <input checked="" type="checkbox"/> | 2330410 | Strut Extension Damper Required seat damper for improved rocking dynamics. | \$ 255 |
| Spoke Guards Available for 20"-24" Composite Rear Wheels. Made from durable and lightweight aluminum. | | | \$ 235 |
| <input type="checkbox"/> | 2330020 | Black | <input type="checkbox"/> |
| | | 2330021 | Frame Color |
| <input type="checkbox"/> | 2330050 | Protective Lower Padding | \$ 995 |
| <input type="checkbox"/> | 2330080 | Gel Hip Guides Reduces available Seat Width by 2". | \$ 595 |
| <input type="checkbox"/> | 2330110 | IV Pole Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks. | \$ 345 |
| <input type="checkbox"/> | 2330120 | O2 Holder Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks. | \$ 495 |
| <input type="checkbox"/> | 2330200 | Q'Straint Lap Belt Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" Belt Length. | \$ 295 |
| <input type="checkbox"/> | 2330222 | Wheelchair Backpack Black with PDG Logo. 17"x14" size. | \$ 165 |
| <input type="checkbox"/> | 2330170 | Adhesive Back Hook Strip, 2" x 30" For Cushion mounting. | \$ 10 |

On-Chair Seating & Positioning

| | | | |
|---|--|--|--|
|  | Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. pdgmobility.com/seating | | |
|---|--|--|--|

Order Notes

Please send your completed order form by email to info@pdgmobility.com.