



BENTLEY LT-R

Order Form - USA

Prices in U.S. Dollars. January 1, 2024

Customer and Order Information

Purchase Order *	Tag	v.0101
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in
		*required

Base Chair

Bentley LT-R 20° Manual Tilt Wheelchair Reverse Wheel Layout enhances user access for self-propulsion. pdgmobility.com/blt-r

<input checked="" type="checkbox"/>	2701000 Bentley 20° Manual Tilt Wheelchair	E1161	\$ 3,150
<input checked="" type="checkbox"/>	2701009 Lightweight Aluminum Frame and Reverse Wheel Layout Upgrade	K0108	\$ 300
<input type="checkbox"/>	Choose the standard Bentley with a steel frame for agitated or heavier users. See separate order form.		

Weight Capacity

2702000 **250 lbs Chair Weight Capacity**

User Weight Range Required for setup of tilt uplift force.

<input type="checkbox"/> 2702430 Up to 75 lbs	<input type="checkbox"/> 2702431 75 lbs to 100 lbs	<input type="checkbox"/> 2702432 100 lbs to 115 lbs
<input type="checkbox"/> 2702433 115 lbs to 150 lbs	<input type="checkbox"/> 2702434 150 lbs to 200 lbs	<input type="checkbox"/> 2702435 200 lbs to 250 lbs

Frame Options

2702300 **Transit Tie Down System (TTDS)**

TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details.

Seat Width

Overall chair width is +9.5".

Standard Seat Widths

<input type="checkbox"/> 2703014 14" Seat Width	<input type="checkbox"/> 2703015 15" Seat Width	<input type="checkbox"/> 2703016 16" Seat Width
<input type="checkbox"/> 2703017 17" Seat Width	<input type="checkbox"/> 2703018 18" Seat Width	<input type="checkbox"/> 2703019 19" Seat Width
<input type="checkbox"/> 2703020 20" Seat Width	<input type="checkbox"/> 2703021 21" Seat Width	<input type="checkbox"/> 2703022 22" Seat Width
		E2201 \$ 375

Seat Depth

Standard Seat Depths (Adjustable)

<input type="checkbox"/> 2704016 16" Seat Depth	<input type="checkbox"/> 2704017 17" Seat Depth
<input type="checkbox"/> 2704018 18" Seat Depth	<input type="checkbox"/> 2704019 19" Seat Depth
<input type="checkbox"/> 2704020 20" Seat Depth	E2203 \$ 500
Custom Seat Depths	
<input type="checkbox"/> 2704021 21" Seat Depth	E2203 \$ 525
<input type="checkbox"/> 2704022 22" Seat Depth	E2204 \$ 800

Frame Color

PDG Premium Colors View color swatches on PDG's website. pdgmobility.com/colors

<input type="checkbox"/> 2705001 Black Shimmer	<input type="checkbox"/> 2705015 Sky Blue	<input type="checkbox"/> 2705016 Poppy Red
<input type="checkbox"/> 2705003 Fine Texture Black	<input type="checkbox"/> 2705010 Lollypop Blue	<input type="checkbox"/> 2705011 Lollypop Red
<input type="checkbox"/> 2705017 Arctic Grey	<input type="checkbox"/> 2705018 Tropical Teal	<input type="checkbox"/> 2705019 Tangerine Orange
<input type="checkbox"/> 2705020 Olive Green	<input type="checkbox"/> 2705021 Lush Lilac	<input type="checkbox"/> 2705012 Bubble Gum Pink
<input type="checkbox"/> 2705022 Leaf Green	<input type="checkbox"/> 2705023 Boysenberry Purple	<input type="checkbox"/> 2705014 Daffodil Yellow
<input type="checkbox"/> 2705200 Custom Color Contact PDG for specifications.	\$ 300	

Armrests

2706000 **No Armrests**

Height Adjustable Armrests (Single Post) E0973 \$ 345

2706111 8"-11" Range 2706114 10"-14" Range 2706116 12"-16" Range

2706100 **Fixed Height Armrests (12")** Full Length Armpad only.

Armrest Height If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.
 Seat Pan to top of Armrest: _____ Does not include height of Armpad

Armrest Position Indicate preference for Armrest position. If left blank, the default setting will be used.
 In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

2706401 Rear 2706402 Mid 2706403 Front

Armpads

Regular Armpads

2707100 Desk Length - 2" x 10" 2707101 Full Length - 2" x 14"

Waterfall Armpads \$ 115

2707105 Desk Length - 2" x 10" 2707106 Full Length - 2" x 14"

Wide Armpads May interfere with wheels in some tilt configurations. \$ 200

2707110 Desk Length - 3" x 10" 2707111 Full Length - 3" x 12"

Gel Armpads May interfere with wheels in some tilt configurations. \$ 135

2707120 Full Length - 2" x 12" 2707121 Full Length - 3.5" x 12"

Contour Arm Trough with Hand Pad May interfere with wheels in some tilt configurations. E2209 \$ 125

2707130 Right (4.5" x 18") 2707131 Left (4.5" x 18")

Armrest Options

Armrest Wideners Move Armrests outward by up to 1" on each side. Seat size does not change. May cause tilt interference.

2707210 1" Total Additional Width \$ 95

2707211 2" Total Additional Width \$ 195

2707310 **Clothing Guards** Not compatible with 8"-11" Height range. \$ 110

Armrest Extensions Not compatible with Desk Length Armpads. \$ 75

2707400 3" Extension - Forward 2707401 3" Extension - Rearward

Seat Pan

2708000 **No Seat Pan**

2708001 **Solid Aluminum Seat Pan** E2231 \$ 165

Seat Height

Don't forget **Select desired Seat-to-Floor Height** STF Height is adjustable within the range of the chosen Subframe Plate (Low/High).

6" Casters	Low						High						
	RW	12"	13"	14"	15"	16"	RW	14"	15"	16"	17"	18"	19"
20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

8" Casters	Low						High						
	RW	12"	13"	14"	15"	16"	RW	14"	15"	16"	17"	18"	19"
24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Casters

Urethane Casters

8" Casters are compatible only with 24" Large Wheels.

<input type="checkbox"/>	2711061	6" x 1.25" Caster size	<input type="checkbox"/>	2711081	8" x 1.25" Caster size		
<input type="checkbox"/>	2711063	6" x 1.50" Caster size					\$ 85
<input type="checkbox"/>	2711084	8" x 2.00" Caster size					\$ 135

Pneumatic Casters

8" Casters are compatible only with 24" Large Wheels.

<input type="checkbox"/>	2711161	6" x 1.25" Caster size				E2214 & E2215	\$ 110
<input type="checkbox"/>	2711181	8" x 1.25" Caster size				E2214 & E2215	\$ 120
<input type="checkbox"/>	2711184	8" x 2.00" Caster size				E2214 & E2215	\$ 150

<input type="checkbox"/>	2711260	Frog Legs 6" x 1.50" Aluminum Soft-Roll Casters					\$ 275
--------------------------	---------	--	--	--	--	--	--------

<input type="checkbox"/>	2711400	Frog Legs Suspension Forks	Only compatible with 24" Large Wheels. Requires 6" Soft-Roll Caster.			E1015	\$ 595
--------------------------	---------	-----------------------------------	--	--	--	-------	--------

Large Wheels

Composite Wheels

<input type="checkbox"/>	2712020	20" Wheel size	<input type="checkbox"/>	2712022	22" Wheel size	<input type="checkbox"/>	2712024	24" Wheel size
--------------------------	---------	----------------	--------------------------	---------	----------------	--------------------------	---------	----------------

Spoke Wheels

<input type="checkbox"/>	2712120	20" Wheel size	<input type="checkbox"/>	2712122	22" Wheel size	<input type="checkbox"/>	2712124	24" Wheel size	\$ 75
--------------------------	---------	----------------	--------------------------	---------	----------------	--------------------------	---------	----------------	-------



Large Wheel Position

Indicate Large Wheel positioning preference. If left blank, the default setting will be used.

In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

<input type="checkbox"/>	2712401	Rear	<input type="checkbox"/>	2712402	Mid	<input type="checkbox"/>	2712403	Front	Not compatible with 20"-22" Seat Depth
--------------------------	---------	------	--------------------------	---------	-----	--------------------------	---------	-------	--

Axle Type

<input type="checkbox"/>	2713000	Quick Release Axles
--------------------------	---------	----------------------------

<input type="checkbox"/>	2713010	Fixed Axles
--------------------------	---------	--------------------

Large Tires

<input type="checkbox"/>	2714000	Urethane Tires					
--------------------------	---------	-----------------------	--	--	--	--	--

<input type="checkbox"/>	2714010	Pneumatic Tires	May result in a finished Seat Height 1/2" greater than stated.			E2211 & E2212	\$ 60
--------------------------	---------	------------------------	--	--	--	---------------	-------

<input type="checkbox"/>	2714020	Pneumatic Tires with Flat-Free Inserts	May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels.			E2211 & E2213	\$ 95
--------------------------	---------	---	---	--	--	---------------	-------

Handrims

<input type="checkbox"/>	2715000	No Handrims
--------------------------	---------	--------------------

<input type="checkbox"/>	2715010	Standard Aluminum Handrims
--------------------------	---------	-----------------------------------

<input type="checkbox"/>	2715020	Plastic Coated Handrims	\$ 150
--------------------------	---------	--------------------------------	--------

<input type="checkbox"/>	2715030	Projection Handrims	Composite Wheels only. 10 Projections, 30", 1-1.8" Length.	\$ 495
--------------------------	---------	----------------------------	--	--------

Wheel Locks & Brakes

User Wheel Locks

<input type="checkbox"/>	2716010	Push to Lock	<input type="checkbox"/>	2716020	Pull to Lock
--------------------------	---------	--------------	--------------------------	---------	--------------

User-activated Unilateral Wheel Lock

Pull to Lock.

<input type="checkbox"/>	1016310	Left Activation	<input type="checkbox"/>	1016311	Right Activation	\$ 795
--------------------------	---------	-----------------	--------------------------	---------	------------------	--------

Wheel Lock Options

<input type="checkbox"/>	2716100	Wheel Lock Extensions	E0961	\$ 80
--------------------------	---------	------------------------------	-------	-------

<input type="checkbox"/>	2716200	Hill Holders	Requires Pneumatic Large Tires with Flat-Free Inserts.	E0974	\$ 450
--------------------------	---------	---------------------	--	-------	--------



BENTLEY LT-R

Order Form - USA
Prices in U.S. Dollars. January 1, 2024

Anti-Tippers

Front Anti-Tippers Required Safety Option. Must choose one type.

- 2717010 Standard E0971 \$ 125
- 2717015 Spring-loaded Dampens the forward-tip movement of the chair during transfers and allows for improved ground clearance. Not removable. E0971 \$ 200

Back Posts

Back Post Style

- 2718010 Stroller Handle Back

Back Post Height Measured from the Seat Pan to the Handle.

- 2718425 Short (22")
- 2718431 Tall (25")

Back Post Folding Fold-down Back is not compatible with Reclining Backrest.

- 2719220 Fold-down

Seat to Back Angle Seat Angles above 98° may cause interference at Seat Depths above 19".

- 2719090 Set at 90°
- 2719098 Set at 98°
- 2719106 Set at 106°
- 2719114 Set at 114°
- 2719122 Set at 122°

Back Options

Angle/Height Adjustable Stroller Handle Not compatible with Headrest Mounting Plate.

- 2719010 Button on Outside \$ 150
- 2719011 Button on Inside

Endomorph Back Backrest will be narrower than Seat by total of 2x the Offset selected below. Minimum width of Backrest is 12".

- 2718102 1" Offset per side \$ 395
- 2718103 1.5" Offset per side
- 2718104 2" Offset per side

10" Backrest Extension Not compatible with Angle/Height Adjustable Stroller Handle.

- 2719210 E0966 \$ 150

Headrest Mounting Plate

Not compatible with Backrest Extension, Angle/Height Adjustable Stroller Handle. Not required if mounting to a solid Back.

- 2719212 \$ 95

Back Upholstery

- 2720000 **No Back Upholstery**

Standard Upholstery

- 2720418 Short (18" Height) E2611 / E2612 \$ 195
- 2720421 Regular (21" Height)
- 2720424 Tall (24" Height)

Tension Adjustable Upholstery Height of Upholstery is 16".

- 2720450 E2611 / E2612 \$ 295

Position of Upholstery Select the desired position of the top of the Upholstery relative to the Seat Pan.

		22" Back Post									
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"		
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>					
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>					
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>					

		25" Back Post									
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"		
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Tall							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Tilt Activation

Attendant Tilt only Select trigger mounting location below.

- 2723011 Right Back Post
- 2723012 Left Back Post

User Tilt only Full Length Arm Pads only. Select trigger mounting location below.

- 2723021 Right Armrest \$ 100
- 2723022 Left Armrest

Attendant & User Tilt Full Length Arm Pads only. Select trigger mounting locations below.

- 2723031 Right Back Post \$ 195
- 2723032 Left Back Post
- 2723041 Right Armrest
- 2723042 Left Armrest

- 2723110 **Tilt Trigger Extension** \$ 50

See next page for more options >

Power Tilt 250 lbs Weight Capacity. 20° tilt range. Includes 2 batteries for off-chair charging. pdgmobility.com/power-tilt		E1002
Minimum Seat Width is 15". Power Tilt is compatible with either O2 Holder or IV Pole on the same chair, not both. Not available for 12" STF Height.		
<input type="checkbox"/> 2723400	Standard Power Tilt	\$ 2,295
<input type="checkbox"/> 2723410	Power Tilt with Manual Tilt Override	Allows for instant manual release of tilt in case of power tilt failure. \$ 2,795
Select Control Switch	Toggle Switch	<input type="checkbox"/> 2723421 Installed on Right side <input type="checkbox"/> 2723422 Installed on Left side
	Button Switch	<input type="checkbox"/> 2723431 Installed on Right side <input type="checkbox"/> 2723432 Installed on Left side
	Stealth Egg Switch	<input type="checkbox"/> 2723440 Not installed
	Optional for Egg Switch:	<input type="checkbox"/> 2723445 Swing Away Mounting Hardware Not compatible with all Headrest makes and models. E1028 \$ 625
Specify Headrest to be used: _____		

Recline		
<input type="checkbox"/> 2723310	Reclining Backrest (Manual)	E1225 \$ 895
30° maximum recline. Not compatible with Backrest Extension, Fold-down Back, Tension Adjustable Upholstery. Requires longer Frame. Maximum Seat Depth is 19". Finished Back Post Height is 3" taller than selected Back Post Height.		

Front Rigging		
<input type="checkbox"/> 2725000	No Front Rigging	
Easy-Out™ Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame. pdgmobility.com/easy-out		
<input type="checkbox"/> 2725410	Fixed 70°	Not available for 12" STF Height. \$ 95
<input type="checkbox"/> 2725420	Universal Angle Adjustable (0°-110°)	Set at angle: _____ (use increments of 10°) 12" STF Height is only compatible with 0°-60°. K0108 \$ 295
Adjust angle by loosening one bolt. Tool-free leg-length adjustment.		
<input type="checkbox"/> 2725110	Standard Plug In 70° Front Rigging	Legacy Design. Not compatible with Easy-Out.

Footplates		
<input type="checkbox"/> 2726010	Composite Footplates	
<input type="checkbox"/> 2726020	Aluminum Angle Adjustable Footplates	K0040 \$ 200
<input type="checkbox"/> 2726030	One Piece Fixed Footplate	Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable. \$ 350
<input type="checkbox"/> 2726040	One Piece Flip-up Footplate	Requires selection of Composite or Aluminum Angle Adjustable. \$ 395
<input type="checkbox"/> 2726050	Interlocking Angle Adjustable Footplates	Available for Seat Widths up to 21" only. K0040 \$ 200

Footrest/Legrest Accessories		
<input type="checkbox"/> 2727012	Calf Pads	Includes reinforced premium Calf Pad mount. \$ 200
<input type="checkbox"/> 2727032	Gel Knee Pads	Compatible with Universal Angle Adjustable Front Rigging only. \$ 200
<input type="checkbox"/> 2727040	Length Adjustable Heel Loops	E0951 \$ 50
<input type="checkbox"/> 2727050	Bodypoint Padded Calf Strap	K0038 \$ 195

Front Rigging Extension Tubes		
Regular Extension Tubes		
<input type="checkbox"/> 2728010	Short	14"-18" for Easy-Out Front Rigging 13"-17" for Standard Plug In Front Rigging
<input type="checkbox"/> 2728020	Long	18"-24" for Easy-Out Front Rigging 15"-19" for Standard Plug In Front Rigging \$ 75
	Lower Leg Length	If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.
Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.		

Miscellaneous Options		
Spoke Guards Made from durable and lightweight aluminum. K0065 \$ 195		
<input type="checkbox"/> 2730020	Black	<input type="checkbox"/> 2730021 Frame Color

See next page for more options >



BENTLEY LT-R

Order Form - USA

Prices in U.S. Dollars. January 1, 2024

<input type="checkbox"/>	2730050	Protective Lower Padding					\$	795
		Stealth Flip Down Abductor	Requires Solid Seat Pan.					
						E0957 & E1028	\$	650
<input type="checkbox"/>	2730060	Small Pad	3"x4"x3"x2" size	<input type="checkbox"/>	2730062	Large Pad	5"x5"x4"x3" size	
<input type="checkbox"/>	2730080	Gel Hip Guides	Reduces available Seat Width by 2".					
						E0956 & E1028	\$	495
<input type="checkbox"/>	2730110	IV Pole						
						E0776	\$	295
<input type="checkbox"/>	2730120	O2 Holder						
						E2208	\$	395
<input type="checkbox"/>	2730200	Q'Strain Lap Belt	Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" Belt Length.					
							\$	245
<input type="checkbox"/>	2730222	Wheelchair Backpack	Black with PDG Logo. 17"x14" size.					
							\$	145
<input type="checkbox"/>	2730170	Adhesive Back Hook Strip, 2" x 30"	For Cushion mounting.					
							\$	10

On-Chair Seating & Positioning

Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. pdgmobility.com/seating

Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions. Use of HCPCS codes does not guarantee coverage or payment.

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.