

### Customer and Order Information

Purchase Order *	Tag	v.0401
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in

### Base Chair

#### Fuze T50n 50° Manual Tilt Wheelchair

- 2601000 Fuze T50 50° Manual Tilt Wheelchair E1161 \$ 3,150
- 2601008 Narrow Frame Upgrade Reduces overall chair width to seat width plus only 4". K0108 \$ 200

### Frame Type & Weight Capacity

#### Fixed Width Frame Select a Fixed Seat Width below.

- 2602010 300 lbs Weight Capacity
- 2602011 400 lbs Weight Capacity Not compatible with Transit Tie Down System.
- 2602025 **Adjustable Width Frame** 250 lbs Weight Capacity. Select an Adjustable Seat Width below. \$ 100

### Frame Options

- 2602222 **Take-Apart Frame** The Take-Apart frame allows attendants to take the Fuze apart in seconds without tools. [pdgmobility.com/take-apart](http://pdgmobility.com/take-apart) \$ 600  
 This feature reduces the weight of the heaviest piece of the wheelchair to 29 lbs or less. Perfect for transporting it in the trunk of a car or SUV.  
 Not compatible with Power Tilt, Drum Brakes, Ventilator Tray, Ventilator Battery Tray, DRB, User Wheel Locks.  
 Take-Apart rear crossbar is not width adjustable. PDG will provide one no-charge crossbar in a different width over the lifetime of each chair ordered with adjustable width frame (can be ordered with the chair or later). Note: Crossbar for the T50n is 4" shorter than the Seat Width.



- 2602300 **Transit Tie Down System (TTDS)**  
 TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details.  
 Tested to max. Weight Capacity of 300 lbs.

### Seat Width Overall chair width is +4.0" for 12" RW

#### Fixed Seat Widths Order with Fixed Frame.

- |   |   |              |
|---|---|--------------|
| <input type="checkbox"/> 2603016 16" Seat Width | <input type="checkbox"/> 2603017 17" Seat Width |              |
| <input type="checkbox"/> 2603018 18" Seat Width | <input type="checkbox"/> 2603019 19" Seat Width |              |
| <input type="checkbox"/> 2603020 20" Seat Width | <input type="checkbox"/> 2603021 21" Seat Width | E2201 \$ 375 |
| <input type="checkbox"/> 2603022 22" Seat Width | <input type="checkbox"/> 2603023 23" Seat Width | E2201 \$ 375 |
| <input type="checkbox"/> 2603024 24" Seat Width |   | E2202 \$ 575 |

#### Adjustable Seat Widths Order with Adjustable Frame.

- |   |   |              |
|---|---|--------------|
| <input type="checkbox"/> 2603216 16"-19" Range - Set at 16" | <input type="checkbox"/> 2603217 16"-19" Range - Set at 17" |              |
| <input type="checkbox"/> 2603218 16"-19" Range - Set at 18" | <input type="checkbox"/> 2603219 16"-19" Range - Set at 19" |              |
| <input type="checkbox"/> 2603319 19"-20" Range - Set at 19" |   |              |
| <input type="checkbox"/> 2603320 19"-20" Range - Set at 20" |   | E2201 \$ 375 |



# FUZE T50n

Order Form - USA

Prices in U.S. Dollars. January 1, 2024

## Seat Depth

### Standard Seat Depths (Adjustable)

- |   |   |   |              |
|---|---|---|--------------|
| <input type="checkbox"/> 2604015 15" Seat Depth | <input type="checkbox"/> 2604016 16" Seat Depth | <input type="checkbox"/> 2604017 17" Seat Depth |              |
| <input type="checkbox"/> 2604018 18" Seat Depth | <input type="checkbox"/> 2604019 19" Seat Depth |   |              |
| <input type="checkbox"/> 2604020 20" Seat Depth |   |   | E2203 \$ 500 |

### Custom Seat Depths Requires Stroller Handle Back - Straight. Not compatible with Reclining or Fold-down Back.

- |   |              |
|---|--------------|
| <input type="checkbox"/> 2604021 21" Seat Depth | E2203 \$ 525 |
| <input type="checkbox"/> 2604022 22" Seat Depth | E2204 \$ 800 |

## Frame Color

### PDG Premium Colors View color swatches on PDG's website. [pdgmobility.com/colors](http://pdgmobility.com/colors)

- |   |   |   |
|---|---|---|
| <input type="checkbox"/> 2605001 Black Shimmer  | <input type="checkbox"/> 2605015 Sky Blue           | <input type="checkbox"/> 2605016 Poppy Red        |
| <input type="checkbox"/> 2605003 Fine Texture Black   | <input type="checkbox"/> 2605010 Lollypop Blue      | <input type="checkbox"/> 2605011 Lollypop Red     |
| <input type="checkbox"/> 2605017 Arctic Grey  | <input type="checkbox"/> 2605018 Tropical Teal      | <input type="checkbox"/> 2605019 Tangerine Orange |
| <input type="checkbox"/> 2605020 Olive Green  | <input type="checkbox"/> 2605021 Lush Lilac         | <input type="checkbox"/> 2605012 Bubble Gum Pink  |
| <input type="checkbox"/> 2605022 Leaf Green   | <input type="checkbox"/> 2605023 Boysenberry Purple | <input type="checkbox"/> 2605014 Daffodil Yellow  |
| <input type="checkbox"/> 2605200 <b>Custom Color</b> <small>Contact PDG for specifications.</small> |   | \$ 300  |

## Armrests

- |   |  |
|---|--|
| <input type="checkbox"/> 2606000 <b>No Armrests</b> |  |
|---|--|

### Height Adjustable Armrests (Single Post)

E0973 \$ 345

- |  |  |
|--|--|
| <input type="checkbox"/> 2606114 10"-14" Range | <input type="checkbox"/> 2606116 12"-16" Range |
|--|--|

- |   |
|---|
| <input type="checkbox"/> 2606100 <b>Fixed Height Armrests (12")</b> <small>Full Length Armpad only.</small> |
|---|



### Armrest Height If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.

Seat Pan to top of Armrest: \_\_\_\_\_ Does not include height of Armpad



### Armrest Position Indicate preference for Armrest position. If left blank, the default setting will be used.

In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

- |                                      |  |
|--------------------------------------|--|
| <input type="checkbox"/> 2606402 Mid | <input type="checkbox"/> 2606403 Front |
|--------------------------------------|--|

### Locking Flip-up Armrests Not recommended for users with excessive strength or tone. Not compatible with Fold-down Back.

E0973 \$ 345

- |   |   |   |
|---|---|---|
| <input type="checkbox"/> 2606310 Set at 10" | <input type="checkbox"/> 2606311 Set at 11" | <input type="checkbox"/> 2606312 Set at 12" |
| <input type="checkbox"/> 2606313 Set at 13" | <input type="checkbox"/> 2606314 Set at 14" |   |

## Armpads

### Regular Armpads

- |   |   |
|---|---|
| <input type="checkbox"/> 2607100 Desk Length - 2" x 10" | <input type="checkbox"/> 2607101 Full Length - 2" x 14" |
|---|---|

### Waterfall Armpads

\$ 115

- |   |   |
|---|---|
| <input type="checkbox"/> 2607105 Desk Length - 2" x 10" | <input type="checkbox"/> 2607106 Full Length - 2" x 14" |
|---|---|

### Wide Armpads May interfere with wheels in some tilt configurations.

\$ 200

- |   |   |
|---|---|
| <input type="checkbox"/> 2607110 Desk Length - 3" x 10" | <input type="checkbox"/> 2607111 Full Length - 3" x 12" |
|---|---|

### Gel Armpads May interfere with wheels in some tilt configurations.

\$ 135

- |   |   |
|---|---|
| <input type="checkbox"/> 2607120 Full Length - 2" x 12" | <input type="checkbox"/> 2607121 Full Length - 3.5" x 12" |
|---|---|

### Contour Arm Trough with Hand Pad

E2209 \$ 125

Not compatible with Lap Tray, Locking Flip-up Armrests. May interfere with wheels in some tilt configurations.

- |   |  |
|---|--|
| <input type="checkbox"/> 2607130 Right (4.5" x 18") | <input type="checkbox"/> 2607131 Left (4.5" x 18") |
|---|--|

## Armrest Options

### Armrest Wideners Move Armrests outward by up to 1" on each side. Seat size does not change.

- |  |        |
|--|--------|
| <input type="checkbox"/> 2607210 1" Total Additional Width | \$ 95  |
| <input type="checkbox"/> 2607211 2" Total Additional Width | \$ 195 |

### Seat Pan

<input type="checkbox"/> 2608000	<b>No Seat Pan</b>		
<input type="checkbox"/> 2608001	<b>Solid Aluminum Seat Pan</b>	E2231	\$ 165

### Seat Height

**Don't forget** Select desired **Seat-to-Floor Height** STF Height is adjustable.

Caster size	12" Rear Wheels				
	16"	17"	18"	19"	20"
5"	<input type="checkbox"/>	<input type="checkbox"/>			
6"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
8"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
6" Soft Roll	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

### Casters

**Urethane Casters**

<input type="checkbox"/> 2611051	5" x 1.25" Caster size	<input type="checkbox"/> 2611061	6" x 1.25" Caster size	<input type="checkbox"/> 2611081	8" x 1.25" Caster size	
<input type="checkbox"/> 2611063	6" x 1.50" Caster size					\$ 85
<input type="checkbox"/> 2611084	8" x 2.00" Caster size					\$ 135

**Pneumatic Casters**

<input type="checkbox"/> 2611161	6" x 1.25" Caster size	E2214 & E2215	\$ 110
<input type="checkbox"/> 2611181	8" x 1.25" Caster size	E2214 & E2215	\$ 120
<input type="checkbox"/> 2611184	8" x 2.00" Caster size	E2214 & E2215	\$ 150

2611260 **Frog Legs 6" x 1.50" Aluminum Soft-Roll Casters** \$ 275

### Rear Wheels

2612012 **12" Composite Wheels**

**Rear Wheel Position** Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.  
In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

<input type="checkbox"/> 2612401	Rear More stable	<input type="checkbox"/> 2612402	Mid	<input type="checkbox"/> 2612403	Front Less stable
----------------------------------	------------------	----------------------------------	-----	----------------------------------	-------------------

### Axle Type

<input type="checkbox"/> 2613000	<b>Quick Release Rear Axles</b>
<input type="checkbox"/> 2613010	<b>Fixed Rear Axles</b>

### Rear Tires

<input type="checkbox"/> 2614000	<b>Urethane Tires</b>		
<input type="checkbox"/> 2614010	<b>Pneumatic Tires</b> May result in a finished Seat Height 1/2" greater than stated.	E2211 & E2212	\$ 60
<input type="checkbox"/> 2614020	<b>Pneumatic Tires with Flat-Free Inserts</b> May result in a finished Seat Height 1/2" greater than stated.	E2211 & E2213	\$ 95

### Handrims

2615000 **No Handrims**

### Wheel Locks & Brakes

**User Wheel Locks** Not compatible with Take-Apart Frame. \$ 150

<input type="checkbox"/> 2616010	Push to Lock	<input type="checkbox"/> 2616020	Pull to Lock
----------------------------------	--------------	----------------------------------	--------------

2616030 **Attendant Wheel Locks (Push to Lock)**

2616005 **Foot Operated Wheel Lock** Rear-mounted unilateral single-activation lock. Not compatible with Ventilator Tray. \$ 595

2616300 **Drum Brakes** Available with Urethane Tires only. Compatible with 17.5", 18.5", 19.5" STF Height only. E2228 \$ 795

**User-activated Unilateral Wheel Lock** Pull to Lock. Available up to 22" Seat Width only. \$ 795

<input type="checkbox"/> 2616310	Left Activation	<input type="checkbox"/> 2616311	Right Activation
----------------------------------	-----------------	----------------------------------	------------------



# FUZE T50n

**Order Form - USA**  
Prices in U.S. Dollars. January 1, 2024

## Wheel Lock Options

<input type="checkbox"/>	2616100	<b>Wheel Lock Extensions</b>	E0961	\$	80
--------------------------	---------	------------------------------	-------	----	----

## Anti-Tippers

<input checked="" type="checkbox"/>	2617010	<b>Rear Anti-Tippers</b> Required Safety Option.	E0971	\$	125
-------------------------------------	---------	--	-------	----	-----

## Back Posts

<b>Back Post Style</b>					
<input type="checkbox"/>	2618100	Stroller Handle Back	<input type="checkbox"/>	2618105	Cane Handle Back Not compatible with Headrest Mounting Plate.
<b>Back Post Type</b>					
<input type="checkbox"/>	2618200	Straight	<input type="checkbox"/>	2618205	8° Bend
<b>Back Post Height</b> Measured from the Seat Pan to the Handle.					
<input type="checkbox"/>	2618420	Short (17")	<input type="checkbox"/>	2618425	Tall (23")
<b>Back Post Folding</b> Fold-down Back is not compatible with Dynamic Rocker Back Interface, Reclining Backrest.					
<input type="checkbox"/>	2619220	Fold-down	<input type="checkbox"/>	2619221	Non-folding
<b>Seat to Back Angle</b>					
<input type="checkbox"/>	2619090	Set at 90°	<input type="checkbox"/>	2619096	Set at 98°
<input type="checkbox"/>	2619108	Set at 114°	<input type="checkbox"/>	2619114	Set at 122°
<input type="checkbox"/>			<input type="checkbox"/>	2619102	Set at 106°

## Back Options

<input type="checkbox"/>	2619010	<b>Angle/Height Adjustable Stroller Handle</b> Not compatible with Headrest Mounting Plate or Cane Handle Back.		\$	150
<input type="checkbox"/>	2619010	Button on Outside	<input type="checkbox"/>	2619011	Button on Inside
<input type="checkbox"/>	2619050	<b>Removable Cane Handle Push Bar</b> Not compatible with Stroller Handle Back.		\$	140
<b>Endomorph Back</b> Backrest will be narrower than Seat by total of 2x the Offset selected below. Minimum width of Backrest is 12".					
<input type="checkbox"/>	2618102	1" Offset per side	<input type="checkbox"/>	2618103	1.5" Offset per side
<input type="checkbox"/>			<input type="checkbox"/>	2618104	2" Offset per side
<input type="checkbox"/>	2619212	<b>Headrest Mounting Plate</b> Not compatible with Angle/Height Adjustable Stroller Handle, Cane Handle Back. Not required if mounting to a solid Back.		\$	95
<input type="checkbox"/>	2619240	<b>Seating Dynamics Dynamic Rocker Back Interface</b>		E2398	\$ 2,995
Not compatible with Take-Apart Frame, Fold-down Back, 30° Recline, Standard Upholstery. Only available with Stroller Handle - 8° Bend. Finished Back Post Height is 7" taller than selected Back Post Height. Maximum compatible Seat Depth is 20". 90° and 98° Seat Angle only.					

## Back Upholstery

<input type="checkbox"/>	2620000	<b>No Back Upholstery</b>																																																																																			
<b>Standard Upholstery</b> Not available with Cane Handle Back. E2611 / E2612 \$ 195																																																																																					
<input type="checkbox"/>	2620418	Short (18" Height)	<input type="checkbox"/>	2620421	Regular (21" Height)																																																																																
<input type="checkbox"/>	2620450	<b>Tension Adjustable Upholstery</b> Height of Upholstery is 16". Not available with Cane Handle Back. E2611 / E2612 \$ 295																																																																																			
<b>Position of Upholstery</b> Select the desired position of the top of the Upholstery relative to the Seat Pan.																																																																																					
<table border="1" style="display: inline-table; margin-right: 20px;"> <thead> <tr> <th colspan="10">17" Back Post</th> </tr> <tr> <th>Style</th> <th>17"</th> <th>18"</th> <th>19"</th> <th>20"</th> <th>21"</th> <th>22"</th> <th>23"</th> <th>24"</th> <th>25"</th> </tr> </thead> <tbody> <tr> <td>Short</td> <td><input type="checkbox"/></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> </tbody> </table> <table border="1" style="display: inline-table;"> <thead> <tr> <th colspan="10">23" Back Post</th> </tr> <tr> <th>Style</th> <th>17"</th> <th>18"</th> <th>19"</th> <th>20"</th> <th>21"</th> <th>22"</th> <th>23"</th> <th>24"</th> <th>25"</th> </tr> </thead> <tbody> <tr> <td>Short</td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td>Regular</td> <td></td> <td></td> <td></td> <td></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td>Tension Adj</td> <td></td> <td></td> <td></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </tbody> </table>						17" Back Post										Style	17"	18"	19"	20"	21"	22"	23"	24"	25"	Short	<input type="checkbox"/>									23" Back Post										Style	17"	18"	19"	20"	21"	22"	23"	24"	25"	Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Regular					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Tension Adj				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
17" Back Post																																																																																					
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"																																																																												
Short	<input type="checkbox"/>																																																																																				
23" Back Post																																																																																					
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"																																																																												
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																																																												
Regular					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																																																												
Tension Adj				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																																																												

## Tilt Activation

<b>Attendant Tilt only</b> Select trigger mounting location below.					
<input type="checkbox"/>	2623011	Right Back Post	<input type="checkbox"/>	2623012	Left Back Post

See next page for more options >



# FUZE T50n

Order Form - USA

Prices in U.S. Dollars. January 1, 2024

<b>Power Tilt</b> 40° tilt range. Includes 2 batteries for off-chair charging. <a href="http://pdgmobility.com/power-tilt">pdgmobility.com/power-tilt</a>		E1002
<input type="checkbox"/>	2623400 Standard Power Tilt (250 lbs Weight Capacity)	\$ 2,295
<input type="checkbox"/>	2623401 Standard Power Tilt (350 lbs Weight Capacity) Requires Fixed Frame.	\$ 2,595
<input type="checkbox"/>	2623410 Power Tilt with Manual Tilt Override (250 lbs Weight Capacity) Allows for instant manual release of tilt in case of power tilt failure.	\$ 2,795
<input type="checkbox"/>	2623411 Power Tilt with Manual Tilt Override (350 lbs Weight Capacity) Allows for instant manual release of tilt in case of power tilt failure. Requires Fixed Frame.	\$ 3,095
Select Control Switch	Toggle Switch <input type="checkbox"/> 2623421 Installed on Right side <input type="checkbox"/> 2623422 Installed on Left side	
	Button Switch <input type="checkbox"/> 2623431 Installed on Right side <input type="checkbox"/> 2623432 Installed on Left side	
	Stealth Egg Switch <input type="checkbox"/> 2623440 Not installed	
	Optional for Egg Switch: <input type="checkbox"/> 2623445 Swing Away Mounting Hardware Not compatible with all Headrest makes and models.	E1028 \$ 625
	Specify Headrest to be used: _____	
<b>Tilt Limit</b> Optional. If no tilt limit is selected, the full tilt range will be available.		
<input type="checkbox"/>	2623220 Set at 20°	
<input type="checkbox"/>	2623230 Set at 30°	
<input type="checkbox"/>	2623240 Set at 40°	

<b>Recline</b>		
<input type="checkbox"/>	2623310 <b>Reclining Backrest (Manual)</b>	E1225 \$ 895
30° maximum recline. Not compatible with Cane Handle Back, Fold-down Back. Not available with 22" Seat Depth. Finished Back Post Height is 3" taller than selected Back Post Height.		

<b>Front Rigging</b>		
<input type="checkbox"/>	2625000 <b>No Front Rigging</b>	
<input type="checkbox"/>	2625001 <b>Sector Blocks only</b>	
<b>Easy-Out™ Front Rigging</b> Effortless insertion and removal. No Hanger Bracket remains on frame. <a href="http://pdgmobility.com/easy-out">pdgmobility.com/easy-out</a>		
<input type="checkbox"/>	2625410 Fixed 70°	\$ 95
<input type="checkbox"/>	2625420 Universal Angle Adjustable (0°-110°) Set at angle: _____ (use increments of 10°)	\$ 295
Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings.		
<b>Swing Away Front Rigging</b> Legacy Design. Not compatible with Easy-Out. Requires pre-tilt of 5°.		
<input type="checkbox"/>	2625130 Swing Away 60° <input type="checkbox"/> 2625140 Swing Away 70°	
<input type="checkbox"/>	2625160 Swing Away 90° Includes 11"-16" Extension Tubes. Not compatible with 8" Casters.	\$ 75
<input type="checkbox"/>	2625200 Swing Away Elevating Legrests Includes Calf Pads.	E0990 \$ 415
<input type="checkbox"/>	2625300 1" Offset Hangers for Swing Away	\$ 75
<input type="checkbox"/>	2625110 <b>Standard Plug In 70° Front Rigging</b> Legacy Design. Not compatible with Easy-Out.	

<b>Footplates</b>		
<input type="checkbox"/>	2626010 <b>Composite Footplates</b>	
<input type="checkbox"/>	2626020 <b>Aluminum Angle Adjustable Footplates</b>	K0040 \$ 200
<input type="checkbox"/>	2626030 <b>One Piece Fixed Footplate</b> Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.	\$ 350
<input type="checkbox"/>	2626040 <b>One Piece Flip-up Footplate</b> Requires selection of Composite or Aluminum Angle Adjustable.	\$ 395
<input type="checkbox"/>	2626050 <b>Interlocking Angle Adjustable Footplates</b> Available for Seat Widths up to 21" only.	K0040 \$ 200

<b>Footrest/Legrest Accessories</b>		
<b>Calf Pad</b> Not compatible with Gel Padding.		
<input type="checkbox"/>	2627012 For Easy-Out™ Front Rigging Includes reinforced premium Calf Pad mount.	\$ 200
<input type="checkbox"/>	2627010 For Swing Away Front Rigging	\$ 135
<input type="checkbox"/>	2627020 <b>Gel Padding</b> Not compatible with Calf Pad or Easy-Out Front Rigging.	\$ 125
<b>Gel Knee Pad</b>		
<input type="checkbox"/>	2627032 For Easy-Out™ Angle Adjustable Front Rigging	\$ 200
<input type="checkbox"/>	2627030 For Swing Away ELR Front Rigging	\$ 100
<input type="checkbox"/>	2627040 <b>Length Adjustable Heel Loops</b>	E0951 \$ 50
<input type="checkbox"/>	2627050 <b>Bodypoint Padded Calf Strap</b> Sizes greater than 22" are not padded.	K0038 \$ 195



# FUZE T50n

**Order Form - USA**

Prices in U.S. Dollars. January 1, 2024

## Front Rigging Extension Tubes

**Regular Extension Tubes** No selection required for 90° Front Rigging

- 2628010 Short 14"-18" for Easy-Out Front Rigging 13"-17" for Swing Away and Standard Plug In Front Rigging
- 2628020 Long 18"-24" for Easy-Out Front Rigging 15"-19" for Swing Away and Standard Plug In Front Rigging

\$ 75



**Lower Leg Length** If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.

Seat Pan to Foot Plates: \_\_\_\_\_ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.

## Miscellaneous Options

**Stealth Flip Down Abductor** Requires Solid Seat Pan.

E0957 & E1028 \$ 650

- 2630060 Small Pad 3"x4"x3"x2" size  2630062 Large Pad 5"x5"x4"x3" size

2630070 **QR Bracket - Accessory Mount**

E1028 \$ 200

2630080 **Gel Hip Guides** Reduces available Seat Width by 2".

E0956 & E1028 \$ 495

2630090 **Hip Belt Mount** Provides improved mounting options for Bodypoint Hip Belt or other third-party belts.

E1028 \$ 55

2630150 **Ventilator Tray** Reduces ground clearance. Chair may be unstable when unoccupied. Minimum Seat Width is 20".  
Not compatible with Power Tilt, Foot Operated Wheel Locks.

E1029 \$ 295

2630152 **Ventilator Bracket** Compatible with LTV1 (laptop) ventilator.

E1029 \$ 295

2630160 **Clear Lap Tray (Slide On)** Available for Seat Widths of 17"-22" only. Not compatible with Contour Arm Trough.

E0950 \$ 490

2630200 **Q'Straint Lap Belt** Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" Belt Length.

\$ 245

2630222 **Wheelchair Backpack** Black with PDG Logo. 17"x14" size.

\$ 145

2630170 **Adhesive Back Hook Strip, 2" x 30"** For Cushion mounting.

\$ 10

## On-Chair Seating & Positioning

Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. [pdgmobility.com/seating](http://pdgmobility.com/seating)

## Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions. Use of HCPCS codes does not guarantee coverage or payment.

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.