



# STELLAR GL

Order Form - USA

Prices in U.S. Dollars. January 1, 2024

Customer and Order Information		
Purchase Order *	Tag	v.0401
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in <span style="float: right;">*required</span>

Base Chair		
<input checked="" type="checkbox"/> 1901000	<b>Stellar GL 45° Manual Tilt Wheelchair</b>	E1161 \$ 3,150
<input type="checkbox"/> Also available: the lightweight <b>Stellar GLT</b> with aluminum frame. See separate order form. <a href="http://pdgmobility.com/glt">pdgmobility.com/glt</a>		

Weight Capacity		
<b>Chair Weight Capacity</b>		
<input type="checkbox"/> 1902000	250 lbs Capacity Available for 14" to 22" Seat Width.	
<input type="checkbox"/> 1902050	300 lbs Capacity Available for 14" to 22" Seat Width.	\$ 200
<input type="checkbox"/> 1905006	450 lbs Capacity Available for 20" to 32" Seat Width. The <b>Stellar HD</b> offers superior strength and performance. See separate order form.	\$ 500
<b>User Weight Range</b> Required for setup of tilt uplift force.		
<input type="checkbox"/> 1902430	Up to 75 lbs	
<input type="checkbox"/> 1902431	75 lbs to 100 lbs	
<input type="checkbox"/> 1902432	100 lbs to 115 lbs	
<input type="checkbox"/> 1902433	115 lbs to 150 lbs	
<input type="checkbox"/> 1902434	150 lbs to 200 lbs	
<input type="checkbox"/> 1902435	200 lbs to 250 lbs	
<input type="checkbox"/> 1902436	250 lbs to 300 lbs	
<input type="checkbox"/> 1902437	300 lbs to 375 lbs	
<input type="checkbox"/> 1902438	375 lbs to 450 lbs	

Frame Options		
<input type="checkbox"/> 1902200	<b>Extended Base with Rear Floor Jacks</b>	K0108 \$ 900
Lifts Rear Wheels off ground for increased stability. Limited ground clearance. Not compatible with Attendant Wheel Lock, HD Anti-Tipper, IV Pole, O2 Holder.		
<input type="checkbox"/> 1902300	<b>Transit Tie Down System (TTDS)</b>	
TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Not compatible with Seat Widths greater than 22". Tested to max. Weight Capacity of 300 lbs.		

Seat Width		
Overall chair width is +9.5" with 12/16" Rear Wheels and +10.5" with 20/22/24" Rear Wheels.		
<b>Standard Seat Widths</b>		
<input type="checkbox"/> 1903014	14" Seat Width	
<input type="checkbox"/> 1903015	15" Seat Width	
<input type="checkbox"/> 1903016	16" Seat Width	
<input type="checkbox"/> 1903017	17" Seat Width	
<input type="checkbox"/> 1903018	18" Seat Width	
<input type="checkbox"/> 1903019	19" Seat Width	
<input type="checkbox"/> 1903020	20" Seat Width	
<input type="checkbox"/> 1903021	21" Seat Width	
<input type="checkbox"/> 1903022	22" Seat Width	E2201 \$ 375
<b>Custom Seat Widths</b> Requires 450 lbs Weight Capacity Frame. Extended delivery lead times may apply.		
<input type="checkbox"/> 1903023	23" Seat Width	E2201 \$ 375
<input type="checkbox"/> 1903024	24" Seat Width	E2202 \$ 575
<input type="checkbox"/> 1903025	25" Seat Width	E2202 \$ 575
<input type="checkbox"/> 1903026	26" Seat Width	E2202 \$ 575
<input type="checkbox"/> 1903027	27" Seat Width	E2202 \$ 575
<input type="checkbox"/> 1903028	28" Seat Width	
<input type="checkbox"/> 1903029	29" Seat Width	
<input type="checkbox"/> 1903030	30" Seat Width	K0108 \$ 1,100
<input type="checkbox"/> 1903031	31" Seat Width	
<input type="checkbox"/> 1903032	32" Seat Width	K0108 \$ 1,100

Seat Depth		
<b>Standard Seat Depths (Adjustable)</b>		
<input type="checkbox"/> 1904016	16" Seat Depth	
<input type="checkbox"/> 1904017	17" Seat Depth	
<input type="checkbox"/> 1904018	18" Seat Depth	
<input type="checkbox"/> 1904019	19" Seat Depth	
<input type="checkbox"/> 1904020	20" Seat Depth	E2203 \$ 500
<b>Custom Seat Depths</b>		
<input type="checkbox"/> 1904021	21" Seat Depth	E2203 \$ 525
<input type="checkbox"/> 1904022	22" Seat Depth	E2204 \$ 800

Frame Color			
<b>PDG Premium Colors</b> View color swatches on PDG's website. <a href="http://pdgmobility.com/colors">pdgmobility.com/colors</a>			
<input type="checkbox"/> 1905001	Black Shimmer	<input type="checkbox"/> 1905015	Sky Blue
<input type="checkbox"/> 1905003	Fine Texture Black	<input type="checkbox"/> 1905010	Lollypop Blue
<input type="checkbox"/> 1905017	Arctic Grey	<input type="checkbox"/> 1905018	Tropical Teal
<input type="checkbox"/> 1905020	Olive Green	<input type="checkbox"/> 1905021	Lush Lilac
<input type="checkbox"/> 1905022	Leaf Green	<input type="checkbox"/> 1905023	Boysenberry Purple
<input type="checkbox"/> 1905016	Poppy Red	<input type="checkbox"/> 1905011	Lollypop Red
		<input type="checkbox"/> 1905019	Tangerine Orange
		<input type="checkbox"/> 1905012	Bubble Gum Pink
		<input type="checkbox"/> 1905014	Daffodil Yellow
<input type="checkbox"/> 1905200	<b>Custom Color</b>	Contact PDG for specifications.	\$ 300

Armrests			
<input type="checkbox"/> 1906000	<b>No Armrests</b>		
<b>Single Post Height Adjustable Armrests</b>			E0973 \$ 345
<input type="checkbox"/> 1906111	8"-11" Range	<input type="checkbox"/> 1906114	10"-14" Range
		<input type="checkbox"/> 1906116	12"-16" Range
<b>Dual Post Height Adjustable Armrests</b>			E0973 \$ 345
<input type="checkbox"/> 1906211	8"-11" Range	<input type="checkbox"/> 1906214	10"-14" Range
		<input type="checkbox"/> 1906216	12"-16" Range
<b>Fixed Height Armrests (12")</b> Full Length Armpad only.			
<input type="checkbox"/> 1906100	Single Post	<input type="checkbox"/> 1906200	Dual Post
	<b>Armrest Height</b> If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used. Seat Pan to top of Armrest: _____ Does not include height of Armpad		
	<b>Armrest Position</b> Indicate preference for Armrest position. If left blank, the default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.		
<input type="checkbox"/> 1906401	Rear	<input type="checkbox"/> 1906402	Mid
		<input type="checkbox"/> 1906403	Front
<b>Locking Flip-up</b> Not recommended for users with excessive strength or tone. Not compatible with Fold-down Back, Clothing Guards, Back Upholstery. Available with Full Length and Desk Length Armpads.			E0973 \$ 345
<input type="checkbox"/> 1906310	Set at 10"	<input type="checkbox"/> 1906311	Set at 11"
		<input type="checkbox"/> 1906312	Set at 12"
<input type="checkbox"/> 1906313	Set at 13"	<input type="checkbox"/> 1906314	Set at 14"

Armpads			
<b>Regular Armpads</b>			
<input type="checkbox"/> 1907100	Desk Length - 2" x 10"	<input type="checkbox"/> 1907101	Full Length - 2" x 14"
<b>Waterfall Armpads</b> Not compatible with Dual Post Armrests			\$ 115
<input type="checkbox"/> 1907105	Desk Length - 2" x 10"	<input type="checkbox"/> 1907106	Full Length - 2" x 14"
<b>Wide Armpads</b> May interfere with wheels in some tilt configurations.			\$ 200
<input type="checkbox"/> 1907110	Desk Length - 3" x 10"	<input type="checkbox"/> 1907111	Full Length - 3" x 12"
<b>Gel Armpads</b> May interfere with wheels in some tilt configurations.			\$ 135
<input type="checkbox"/> 1907120	Full Length - 2" x 12"	<input type="checkbox"/> 1907121	Full Length - 3.5" x 12"
<b>Contour Arm Trough with Hand Pad</b> May interfere with wheels in some tilt configurations.			E2209 \$ 125
<input type="checkbox"/> 1907130	Right (4.5" x 18")	<input type="checkbox"/> 1907131	Left (4.5" x 18")

Armrest Options			
<b>Armrest Wideners</b> Move Armrests outward by up to 1" on each side. Seat size does not change. May cause tilt interference.			
<input type="checkbox"/> 1907210	1" Total Additional Width		\$ 95
<input type="checkbox"/> 1907211	2" Total Additional Width		\$ 195
<b>Clothing Guards</b> Not compatible with Flip-up Armrests.			
<input type="checkbox"/> 1907310	For Single Post Armrests	Not compatible with 8"-11" Height range.	
<input type="checkbox"/> 1907320	For Dual Post Armrests		
<b>Armrest Extensions</b> Not compatible with Flip-up Armrests, Desk Length Armpads.			\$ 75
<input type="checkbox"/> 1907400	3" Extension - Forward		<input type="checkbox"/> 1907401
			3" Extension - Rearward

### Seat Pan

1908000 **No Seat Pan**

1908001 **Solid Aluminum Seat Pan** E2231 \$ 165

### Seat Height

**Don't forget** **Select desired Seat-to-Floor Height** STF Height is adjustable within the range of the chosen Subframe Plate (Low/High).

6" Casters	Low							High						
	RW	13"	14"	15"	16"	17"	18"	RW	15"	16"	17"	18"	19"	20"
	12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
16"				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
24"				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

8" Casters	Low							High						
	RW	13"	14"	15"	16"	17"	18"	RW	15"	16"	17"	18"	19"	20"
	16"				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
24"				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

### Casters

**Urethane Casters** 8" Casters are compatible only with 16" and 24" Rear Wheels.

1911061 6" x 1.25" Caster size  1911081 8" x 1.25" Caster size

1911063 6" x 1.50" Caster size \$ 85

1911084 8" x 2.00" Caster size \$ 135

**Pneumatic Casters** 8" Casters are compatible only with 16" and 24" Rear Wheels.

1911161 6" x 1.25" Caster size E2214 & E2215 \$ 110

1911181 8" x 1.25" Caster size E2214 & E2215 \$ 120

1911184 8" x 2.00" Caster size E2214 & E2215 \$ 150

1911260 **Frog Legs 6" x 1.50" Aluminum Soft-Roll Casters** \$ 275

1911400 **Frog Legs Suspension Forks** Only compatible with 16 and 24" Rear Wheels. Requires 6" Soft-Roll Caster. E1015 \$ 595

### Rear Wheels

**Composite Wheels**

1912012 12" Wheel size  1912016 16" Wheel size  1912020 20" Wheel size

1912022 22" Wheel size  1912024 24" Wheel size

**Spoke Wheels** \$ 75

1912120 20" Wheel size  1912122 22" Wheel size  1912124 24" Wheel size

1912224 **24" Heavy Duty Spoke Wheels** \$ 350

**Rear Wheel Position** Indicate Rear Wheel positioning preference. If left blank, the default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

1912401 Rear  1912402 Mid  1912403 Front

### Axle Type

1913000 **Quick Release Rear Axles**

1913010 **Fixed Rear Axles**

### Rear Tires

1914000 **Urethane Tires**

1914010 **Pneumatic Tires** May result in a finished Seat Height 1/2" greater than stated. E2211 & E2212 \$ 60

1914020 **Pneumatic Tires with Flat-Free Inserts** May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels. E2211 & E2213 \$ 95



Handrims		Handrims are only available with 20", 22", and 24" Rear Wheels.
<input type="checkbox"/>	1915000 <b>No Handrims</b>	
<input type="checkbox"/>	1915010 <b>Standard Aluminum Handrims</b>	
<input type="checkbox"/>	1915020 <b>Plastic Coated Handrims</b> Not compatible with Drum Brakes.	\$ 150
<input type="checkbox"/>	1915030 <b>Projection Handrims</b> Composite Wheels only. 10 Projections, 30°, 1-1.8" Length.	\$ 495

Wheel Locks & Brakes			
<b>User Wheel Locks</b>			
<input type="checkbox"/>	1916010 Push to Lock	<input type="checkbox"/>	1916020 Pull to Lock
<input type="checkbox"/>	1916030 <b>Attendant Wheel Locks (Push to Lock)</b>		\$ 85
<input type="checkbox"/>	1916040 <b>User &amp; Attendant Wheel Locks</b> Not compatible with Drum Brakes.		\$ 125
<input type="checkbox"/>	1916300 <b>Drum Brakes</b> Available with 12" Composite or 20" Spoke Rear Wheels. Urethane Tires only.	E2228	\$ 795
<b>User-activated Unilateral Wheel Lock</b> Pull to Lock. Available up to 22" Seat Width only.			\$ 795
<input type="checkbox"/>	1916310 Left Activation	<input type="checkbox"/>	1916311 Right Activation

Wheel Lock Options			
<input type="checkbox"/>	1916100 <b>Wheel Lock Extensions</b> Not compatible with Attendant Wheel Lock.	E0961	\$ 80
<input type="checkbox"/>	1916200 <b>Hill Holders</b> Requires Pneumatic Rear Tires with Flat-Free Inserts.	E0974	\$ 450

Anti-Tippers			
<input checked="" type="checkbox"/>	1917010 <b>Rear Anti-Tippers</b> Required Safety Option.	E0971	\$ 125
<input type="checkbox"/>	1917020 <b>Front Anti-Tippers</b>	E0971	\$ 795
<input type="checkbox"/>	1917030 <b>Heavy Duty Rear Anti-Tippers</b> Provide additional rearward stability for users with excessive movement. Not compatible with Attendant Wheel Locks, Rear Floor Jacks, IV Pole, O2 Holder.	E0971	\$ 895

Back Posts					
<b>Back Post Style</b>					
<input type="checkbox"/>	1918010 Stroller Handle Back	<input type="checkbox"/>	1918020 Cane Handle Back	Seat Widths of up to 22" only. May reduce tilt range.	
<b>Back Post Height</b> Measured from the Seat Pan to the Handle.					
<input type="checkbox"/>	1918425 Short (22")	<input type="checkbox"/>	1918431 Tall (25")		
<b>Back Post Folding</b> Folding Back is not compatible with Dynamic Rocker Back Interface, Reclining Backrest.					
<input type="checkbox"/>	1919220 Fold-down	<input type="checkbox"/>	1919221 Non-folding		
<b>Seat to Back Angle</b>					
<input type="checkbox"/>	1919090 Set at 90°	<input type="checkbox"/>	1919098 Set at 98°	<input type="checkbox"/>	1919106 Set at 106°
<input type="checkbox"/>	1919114 Set at 114°	<input type="checkbox"/>	1919122 Set at 122°		

Back Options					
<b>Angle/Height Adjustable Stroller Handle</b> Not compatible with Headrest Mounting Plate or Cane Handle Back.		\$ 150			
<input type="checkbox"/>	1919010 Button on Outside	<input type="checkbox"/>	1919011 Button on Inside		
<b>Endomorph Back</b> Backrest will be narrower than Seat by total of 2x the Offset selected below. Minimum width of Backrest is 12".		\$ 395			
<input type="checkbox"/>	1918102 1" Offset per side	<input type="checkbox"/>	1918103 1.5" Offset per side	<input type="checkbox"/>	1918104 2" Offset per side
<input type="checkbox"/>	1919210 <b>10" Backrest Extension</b> Not compatible with Angle/Height Adjustable Stroller Handle.	E0966	\$ 150		
<input type="checkbox"/>	1919212 <b>Headrest Mounting Plate</b> Not compatible with Backrest Extension, Angle/Height Adjustable Stroller Handle, Cane Handle Back. Not required if mounting to a solid Back.		\$ 95		
<input type="checkbox"/>	1919230 <b>Spring Back Suspension</b> Not compatible with Reclining Back.		\$ 105		
<input type="checkbox"/>	1919240 <b>Seating Dynamics Dynamic Rocker Back Interface</b> Not compatible with Fold-down Back, Spring Back Suspension, 30° Recline, Cane Handle, Standard Upholstery. Finished Back Post Height is 4" taller than selected Back Post Height.	E2398	\$ 2,995		



### Back Upholstery

1920000 **No Back Upholstery**

**Standard Upholstery** Not available with Cane Handle Back. E2611 / E2612 \$ 195

1920418 Short (18" Height)       1920421 Regular (21" Height)       1920424 Tall (24" Height)

**Tension Adjustable Upholstery** Height of Upholstery is 16". Not available with Cane Handle Back.

1920450 For 14" to 22" Seat Width E2611 / E2612 \$ 295

1920451 For 24" to 32" Seath Width E2611 / E2612 \$ 375

**Position of Upholstery** Select the desired position of the top of the Upholstery relative to the Seat Pan.

22" Back Post									
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

25" Back Post									
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tall								<input type="checkbox"/>	<input type="checkbox"/>
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

### Tilt Activation

**Attendant Tilt only** Select trigger mounting location below. Chairs with 450 lbs Weight Capacity require two tilt triggers (L and R) for tilt activation.

1923011 Right Back Post       1923012 Left Back Post

**User Tilt only** Not available for 450 lbs Weight Capacity. Full Length Arm Pads only. Select trigger mounting location below. \$ 100

1923021 Right Armrest       1923022 Left Armrest

**Attendant & User Tilt** Not available for 450 lbs Weight Capacity. Full Length Arm Pads only. Select trigger mounting locations below. \$ 195

1923031 Right Back Post       1923032 Left Back Post

1923041 Right Armrest       1923042 Left Armrest

1923110 **Tilt Trigger Extension** \$ 50

**Tilt Limit** Optional. If no tilt limit is selected, the full tilt range will be available.

1923220 Set at 20°       1923230 Set at 30°       1923240 Set at 40°

**Power Tilt** 250 lbs Weight Capacity. 45° tilt range. Includes 2 batteries for off-chair charging. [pdgmobility.com/power-tilt](http://pdgmobility.com/power-tilt) E1002

Minimum Seat Width is 16". Power Tilt is compatible with either O2 Holder or IV Pole on the same chair, not both.

1923400 Standard Power Tilt \$ 2,295

1923410 Power Tilt with Manual Tilt Override Allows for instant manual release of tilt in case of power tilt failure. \$ 2,795

Toggle Switch       1923421 Installed on Right side       1923422 Installed on Left side

Button Switch       1923431 Installed on Right side       1923432 Installed on Left side

Stealth Egg Switch       1923440 Not installed

Optional for Egg Switch:       1923445 Swing Away Mounting Hardware Not compatible with all Headrest makes and models. E1028 \$ 625

*Specify Headrest to be used: \_\_\_\_\_*

### Recline

1923310 **Reclining Backrest (Manual)** \$ 895

30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold-down Back, Tension Adjustable Upholstery. Finished Back Post Height is 3" taller than selected Back Post Height. Limits Tilt by up to 20°.

### Front Rigging

1925000 **No Front Rigging**       1925001 **Sector Blocks only**

**Easy-Out™ Front Rigging** Effortless insertion and removal. No Hanger Bracket remains on frame. [pdgmobility.com/easy-out](http://pdgmobility.com/easy-out)

1925410 Fixed 70° \$ 95

1925420 Universal Angle Adjustable (0°-110°) Set at angle: \_\_\_\_\_ (use increments of 10°) \$ 295


Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings.

See next page for more options >

<b>Swing Away Front Rigging</b> Legacy Design. Not compatible with Easy-Out.			
<input type="checkbox"/>	1925130 Swing Away 60°	<input type="checkbox"/>	1925140 Swing Away 70°
<input type="checkbox"/>	1925200 Swing Away Elevating Legrests Includes Calf Pads.		E0990 \$ 415
<input type="checkbox"/>	1925300 1" Offset Hangers for Swing Away		\$ 75
<input type="checkbox"/>	1925110 <b>Standard Plug In 70° Front Rigging</b> Legacy Design. Not compatible with Easy-Out.		

Footplates			
<input type="checkbox"/>	1926010 <b>Composite Footplates</b>		
<input type="checkbox"/>	1926020 <b>Aluminum Angle Adjustable Footplates</b>	K0040	\$ 200
<input type="checkbox"/>	1926030 <b>One Piece Fixed Footplate</b> Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.		\$ 350
<input type="checkbox"/>	1926040 <b>One Piece Flip-up Footplate</b> Requires selection of Composite or Aluminum Angle Adjustable.		\$ 395
<input type="checkbox"/>	1926050 <b>Interlocking Angle Adjustable Footplates</b> Available for Seat Widths up to 21" only.	K0040	\$ 200

Footrest/Legrest Accessories			
<b>Calf Pads</b> Not compatible with Gel Padding.			
<input type="checkbox"/>	1927012 For Easy-Out™ Front Rigging Includes reinforced premium Calf Pad mount.		\$ 200
<input type="checkbox"/>	1927010 For Swing Away Front Rigging		\$ 135
<input type="checkbox"/>	1927020 <b>Gel Padding</b> Not compatible with Calf Pad or Easy-Out Front Rigging.		\$ 125
<b>Gel Knee Pads</b>			
<input type="checkbox"/>	1927032 For Easy-Out™ Angle Adjustable Front Rigging		\$ 200
<input type="checkbox"/>	1927030 For Swing Away ELR Front Rigging		\$ 100
<input type="checkbox"/>	1927040 <b>Length Adjustable Heel Loops</b>	E0951	\$ 50
<input type="checkbox"/>	1927050 <b>Bodypoint Padded Calf Strap</b> Sizes greater than 22" are not padded.	K0038	\$ 195
<b>Miller's Dynamic Footrest Mounts</b> Not compatible with Easy-Out Front Rigging.			
<input type="checkbox"/>	1927060 For Standard Plug In Front Rigging		\$ 445
<input type="checkbox"/>	1927065 For Swing Away 70° Front Rigging		\$ 575

Front Rigging Extension Tubes			
<b>Regular Extension Tubes</b>			
<input type="checkbox"/>	1928010 Short 14"-18" for Easy-Out Front Rigging 13"-17" for Swing Away and Standard Plug In Front Rigging		
<input type="checkbox"/>	1928020 Long 18"-24" for Easy-Out Front Rigging 15"-19" for Swing Away and Standard Plug In Front Rigging		\$ 75
<b>Extra Strength Extension Tubes (Solid)</b> Available on Plug In Front Rigging only. Not compatible with Swing Away or Easy-Out Front Rigging.			\$ 100
<input type="checkbox"/>	1928110 Standard (13"-17")	<input type="checkbox"/>	1928120 Long (15"-19")
	<b>Lower Leg Length</b> If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.		
	Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.		

Miscellaneous Options			
<b>Spoke Guards</b> Available for 20"-24" Composite Rear Wheels. Made from durable and lightweight aluminum.		K0065	\$ 195
<input type="checkbox"/>	1930020 Black	<input type="checkbox"/>	1930021 Frame Color
<input type="checkbox"/>	1930050 <b>Protective Lower Padding</b>		\$ 795
<b>Stealth Flip Down Abductor</b> Requires Solid Seat Pan.		E0957 & E1028	\$ 650
<input type="checkbox"/>	1930060 Small Pad 3"x4"x3"x2" size	<input type="checkbox"/>	1930062 Large Pad 5"x5"x4"x3" size
<input type="checkbox"/>	1930080 <b>Gel Hip Guides</b> Reduces available Seat Width by 2".	E0956 & E1028	\$ 495
<input type="checkbox"/>	1930110 <b>IV Pole</b> Not compatible with HD Anti Tipper, Attendant Wheel Locks.	E0776	\$ 295
<input type="checkbox"/>	1930120 <b>O2 Holder</b> Not compatible with HD Anti Tipper, Attendant Wheel Locks.	E2208	\$ 395
<input type="checkbox"/>	1930200 <b>Q'Strain Lap Belt</b> Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" Belt Length.		\$ 245
<input type="checkbox"/>	1930222 <b>Wheelchair Backpack</b> Black with PDG Logo. 17"x14" size.		\$ 145
<input type="checkbox"/>	1930170 <b>Adhesive Back Hook Strip, 2" x 30"</b> For Cushion mounting.		\$ 10



# STELLAR GL

**Order Form - USA**  
Prices in U.S. Dollars. January 1, 2024

## On-Chair Seating & Positioning



Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. [pdgmobility.com/seating](https://pdgmobility.com/seating)

## Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions. Use of HCPCS codes does not guarantee coverage or payment.

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.