

# BENTLEY LT-R



The *lightweight Bentley LT-R* promotes user independence by offering *improved self-propulsion* and *superior maneuverability* through a reverse wheel layout. It tilts up to 20° from the front of the seat (no knee-rise during tilt).



Hand Propulsion



Foot Propulsion



Self-Tilt



Desk/Table Accessible



Transit Tie Downs

**Bentley LT-R** 20° lightweight tilt chair offering unmatched self-propulsion options

visit [pdgmobility.com/blt-r](https://pdgmobility.com/blt-r) for more information



Drive wheels remain in an ergonomic position throughout the tilt range, providing superior access for hand-propelling

Super-low 12" STF Height available for easier foot-propelling



DYNAMIC TILT  
Up to 20°



WEIGHT CAPACITY  
250 lbs



SEAT WIDTH  
14" to 22"



SEAT DEPTH  
16" to 22"



SEAT HEIGHT  
12" to 19"



BACK SEAT ANGLE  
90° to 122°



OVERALL WIDTH  
10.5"



CHAIR WEIGHT  
34 lbs



# BENTLEY LT-R



visit [pdgmobility.com/blt-r](https://pdgmobility.com/blt-r)  
for more information

## *Bentley exclusive features*

- ✓ Foot and hand propel in tilt position with no knee-rise
- ✓ Short wheelbase for maneuverability
- ✓ Adjustable axle position to move center of gravity



## *Available options*

- ✓ Transit tie downs WC19 & ISO7176-19
- ✓ 30° manual dynamic recline
- ✓ Power tilt
- ✓ Easy-Out front rigging

*“ the Bentley LT-R  
complements PDG’s 20° tilt  
chair lineup by supporting  
additional ergonomic hand  
and foot positions for  
self-propulsion ”*

