

# STELLAR GL



The *Stellar GL* is a versatile 45° low knee-rise tilt chair. It promotes self propulsion, table access, and *participation in mobility related activities of daily living*. The rigid frame provides superior stability and maneuverability.



Hand Propulsion



Foot Propulsion



Self-Tilt



Desk/Table Accessible



Transit Tie Downs

## Exclusive features

- ✓ Foot and hand propel in tilt position with low knee-rise
- ✓ Short wheelbase for maneuverability
- ✓ Simple center of gravity adjustments

## Available options *See order form*

- ✓ Transit tie downs WC-19 & ISO-7176-19
- ✓ 30° manual dynamic recline
- ✓ Power tilt



**Stellar GLT** 45° tilt chair with a lightweight aluminum frame

**Stellar GL** 45° tilt chair with a rugged steel frame

The lightweight design of the *Stellar GLT* combines the benefits of 45° tilt with a high degree of mobility independence.

**YOUR PICK**

The rugged steel frame of the standard *Stellar GL* makes it a great choice for heavier users or those with high tone/agitation.

# STELLAR GL



## Stellar GLT



DYNAMIC TILT  
Up to 45°



WEIGHT CAPACITY  
250 lbs



SEAT WIDTH  
14" to 22"



SEAT DEPTH  
16" to 22"



SEAT HEIGHT  
13" to 20"



BACK SEAT ANGLE  
90° to 122°



OVERALL WIDTH  
+9.5" to 10.5"



CHAIR WEIGHT  
34 lbs

visit [pdgmobility.com/glt](https://pdgmobility.com/glt)  
for more information



## Stellar GL

visit [pdgmobility.com/gl](https://pdgmobility.com/gl)  
for more information



DYNAMIC TILT  
Up to 45°



WEIGHT CAPACITY  
Up to 450 lbs



SEAT WIDTH  
14" to 24"



SEAT DEPTH  
16" to 22"



SEAT HEIGHT  
13" to 20"



BACK SEAT ANGLE  
90° to 122°



OVERALL WIDTH  
+9.5" to 10.5"



CHAIR WEIGHT  
51 lbs

