

# STELLAR HD



The *Stellar HD* is a bariatric wheelchair offering up to 45° of tilt and 600 lbs weight capacity. It features heavy-gauge *steel reinforcement of the frame and other components*. Multiple gas struts ensure *effortless tilt operation*.



Hand Propulsion



Foot Propulsion



Self-Tilt



Desk/Table Accessible



Transit Tie Downs

## Stellar HD Bariatric tilt chair up to 600 lbs

visit [pdgmobility.com/hd](https://pdgmobility.com/hd) for more information

Available options See order form

- ✓ 30° manual dynamic recline
- ✓ Dynamic Rocker Back



DYNAMIC TILT  
Up to 45°



WEIGHT CAPACITY  
Up to 600 lbs



SEAT WIDTH  
22" to 32"



SEAT DEPTH  
18" to 22"



SEAT HEIGHT  
16" to 18"



BACK SEAT ANGLE  
90°



OVERALL WIDTH  
+11.0"



CHAIR WEIGHT  
62 lbs

