



BENTLEY LT

Order Form - Distributor

Prices in Canadian Dollars. January 1, 2024

Customer and Order Information

| | | |
|--|------------------|--------------------|
| Purchase Order * | Tag | v.0401 |
| Dealer/Provider Name * | | |
| Purchaser Name and Contact Information * | | |
| Shipping Address * | | |
| User Information (if available) | Weight _____ lbs | Height _____ ft/in |
| | | *required |

Base Chair

Bentley LT 20° Manual Tilt Wheelchair

- 2401000 Bentley 20° Manual Tilt Wheelchair \$ 4,400
 - 2401009 Lightweight Aluminum Frame Upgrade Aluminum Frame is not recommended for agitated users. \$ 300
- i** Choose the standard **Bentley** with a steel frame for agitated or heavier users. See separate order form.

Weight Capacity

- 2402000 **250 lbs Chair Weight Capacity**
- User Weight Range** Required for setup of tilt uplift force.
- 2402430 Up to 75 lbs
 - 2402431 75 lbs to 100 lbs
 - 2402432 100 lbs to 115 lbs
 - 2402433 115 lbs to 150 lbs
 - 2402434 150 lbs to 200 lbs
 - 2402435 200 lbs to 250 lbs

Frame Options

- 2402300 **Transit Tie Down System (TTDS)**
TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details.

Seat Width

Overall chair width is +9.5" with 12/16" Rear Wheels and +10.5" with 20/22/24" Rear Wheels.

Standard Seat Widths

- 2403014 14" Seat Width
- 2403015 15" Seat Width
- 2403016 16" Seat Width
- 2403017 17" Seat Width
- 2403018 18" Seat Width
- 2403019 19" Seat Width
- 2403020 20" Seat Width
- 2403021 21" Seat Width
- 2403022 22" Seat Width

Seat Depth

Standard Seat Depths (Adjustable)

- 2404016 16" Seat Depth
- 2404017 17" Seat Depth
- 2404018 18" Seat Depth
- 2404019 19" Seat Depth
- 2404020 20" Seat Depth

Custom Seat Depths

- 2404021 21" Seat Depth
 - 2404022 22" Seat Depth
- \$ 400

Frame Color

PDG Premium Colors

 View color swatches on PDG's website. pdgmobility.com/colors

- 2405001 Black Shimmer
 - 2405003 Fine Texture Black
 - 2405017 Arctic Grey
 - 2405020 Olive Green
 - 2405022 Leaf Green
 - 2405001 Sky Blue
 - 2405010 Lollypop Blue
 - 2405018 Tropical Teal
 - 2405021 Lush Lilac
 - 2405023 Boysenberry Purple
 - 2405016 Poppy Red
 - 2405011 Lollypop Red
 - 2405019 Tangerine Orange
 - 2405012 Bubble Gum Pink
 - 2405014 Daffodil Yellow
- 2405200 **Custom Color** Contact PDG for specifications. \$ 300

Armrests

2406000 **No Armrests**

Height Adjustable Armrests (Single Post) \$ 295

2406111 8"-11" Range 2406114 10"-14" Range 2406116 12"-16" Range

2406100 **Fixed Height Armrests (12")** Full Length Armpad only.

Armrest Height If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.
 Seat Pan to top of Armrest: _____ Does not include height of Armpad

Armrest Position Indicate preference for Armrest position. If left blank, the default setting will be used.
 In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

2406401 Rear 2406402 Mid 2406403 Front

Armpads

Regular Armpads

2407100 Desk Length - 2" x 10" 2407101 Full Length - 2" x 14"

Waterfall Armpads \$ 115

2407105 Desk Length - 2" x 10" 2407106 Full Length - 2" x 14"

Wide Armpads May interfere with wheels in some tilt configurations. \$ 200

2407110 Desk Length - 3" x 10" 2407111 Full Length - 3" x 12"

Gel Armpads May interfere with wheels in some tilt configurations. \$ 135

2407120 Full Length - 2" x 12" 2407121 Full Length - 3.5" x 12"

Contour Arm Trough with Hand Pad May interfere with wheels in some tilt configurations. \$ 125

2407130 Right (4.5" x 18") 2407131 Left (4.5" x 18")

Armrest Options

Armrest Wideners Move Armrests outward by up to 1" on each side. Seat size does not change. May cause tilt interference.

2407210 1" Total Additional Width \$ 95

2407211 2" Total Additional Width \$ 195

2407310 **Clothing Guards** Not compatible with 8"-11" Height range. \$ 110

Armrest Extensions Not compatible with Desk Length Armpads. \$ 75

2407400 3" Extension - Forward 2407401 3" Extension - Rearward

Seat Pan

2408000 **No Seat Pan**

2408001 **Solid Aluminum Seat Pan** \$ 135

Seat Height

Don't forget **Select desired Seat-to-Floor Height** STF Height is adjustable within the range of the chosen Subframe Plate (Low/High).

| Casters | Low | | | | | | High | | | | | | |
|------------|-----|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|
| | RW | 13" | 14" | 15" | 16" | 17" | RW | 15" | 16" | 17" | 18" | 19" | 20" |
| 6" Casters | 12" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | 12" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| | 16" | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | 16" | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| | 20" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | 20" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| | 22" | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | 22" | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| | 24" | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | 24" | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 8" Casters | RW | 13" | 14" | 15" | 16" | 17" | RW | 15" | 16" | 17" | 18" | 19" | 20" |
| | 16" | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | 16" | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| | 24" | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | 24" | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |

Casters

Urethane Casters 8" Casters are compatible only with 16" and 24" Rear Wheels.

| | | | | |
|----------------------------------|------------------------|----------------------------------|------------------------|--------|
| <input type="checkbox"/> 2411061 | 6" x 1.25" Caster size | <input type="checkbox"/> 2411081 | 8" x 1.25" Caster size | |
| <input type="checkbox"/> 2411063 | 6" x 1.50" Caster size | | | \$ 85 |
| <input type="checkbox"/> 2411084 | 8" x 2.00" Caster size | | | \$ 135 |

Pneumatic Casters 8" Casters are compatible only with 16" and 24" Rear Wheels.

| | | | | |
|----------------------------------|------------------------|--|--|--------|
| <input type="checkbox"/> 2411161 | 6" x 1.25" Caster size | | | \$ 110 |
| <input type="checkbox"/> 2411181 | 8" x 1.25" Caster size | | | \$ 120 |
| <input type="checkbox"/> 2411184 | 8" x 2.00" Caster size | | | \$ 150 |

| | | | | |
|----------------------------------|--|--|--|--------|
| <input type="checkbox"/> 2411260 | Frog Legs 6" x 1.50" Aluminum Soft-Roll Casters | | | \$ 275 |
|----------------------------------|--|--|--|--------|

| | | | | |
|----------------------------------|-----------------------------------|--|--|--------|
| <input type="checkbox"/> 2411400 | Frog Legs Suspension Forks | Only compatible with 16 and 24" Rear Wheels. Requires 6" Soft-Roll Caster. | | \$ 695 |
|----------------------------------|-----------------------------------|--|--|--------|

Rear Wheels

Composite Wheels

| | | | | | |
|----------------------------------|----------------|----------------------------------|----------------|----------------------------------|----------------|
| <input type="checkbox"/> 2412012 | 12" Wheel size | <input type="checkbox"/> 2412016 | 16" Wheel size | <input type="checkbox"/> 2412020 | 20" Wheel size |
| <input type="checkbox"/> 2412022 | 22" Wheel size | <input type="checkbox"/> 2412024 | 24" Wheel size | | |

Spoke Wheels

| | | | | | |
|----------------------------------|----------------|----------------------------------|----------------|----------------------------------|----------------|
| <input type="checkbox"/> 2412120 | 20" Wheel size | <input type="checkbox"/> 2412122 | 22" Wheel size | <input type="checkbox"/> 2412124 | 24" Wheel size |
|----------------------------------|----------------|----------------------------------|----------------|----------------------------------|----------------|



Rear Wheel Position Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.

In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

| | | | | | |
|----------------------------------|------|----------------------------------|-----|----------------------------------|-------|
| <input type="checkbox"/> 2412401 | Rear | <input type="checkbox"/> 2412402 | Mid | <input type="checkbox"/> 2412403 | Front |
|----------------------------------|------|----------------------------------|-----|----------------------------------|-------|

Axle Type

| | |
|----------------------------------|---------------------------------|
| <input type="checkbox"/> 2413000 | Quick Release Rear Axles |
| <input type="checkbox"/> 2413010 | Fixed Rear Axles |

Rear Tires

| | | | |
|----------------------------------|---|---|-------|
| <input type="checkbox"/> 2414000 | Urethane Tires | | |
| <input type="checkbox"/> 2414010 | Pneumatic Tires | May result in a finished Seat Height 1/2" greater than stated. | \$ 60 |
| <input type="checkbox"/> 2414020 | Pneumatic Tires with Flat-Free Inserts | May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels. | \$ 95 |

Handrims

Handrims are only available with 20", 22", and 24" Rear Wheels.

| | | | |
|----------------------------------|-----------------------------------|--|--------|
| <input type="checkbox"/> 2415000 | No Handrims | | |
| <input type="checkbox"/> 2415010 | Standard Aluminum Handrims | | |
| <input type="checkbox"/> 2415020 | Plastic Coated Handrims | Not compatible with Drum Brakes. | \$ 150 |
| <input type="checkbox"/> 2415030 | Projection Handrims | Composite Wheels only. 10 Projections, 30°, 1-1.8" Length. | \$ 595 |

Wheel Locks & Brakes

User Wheel Locks

| | | | | |
|----------------------------------|---|---|------------------|--------|
| <input type="checkbox"/> 2416010 | Push to Lock | <input type="checkbox"/> 2416020 | Pull to Lock | |
| <input type="checkbox"/> 2416030 | Attendant Wheel Locks (Push to Lock) | | | \$ 85 |
| <input type="checkbox"/> 2416040 | User & Attendant Wheel Locks | Not compatible with Drum Brakes. | | \$ 125 |
| <input type="checkbox"/> 2416300 | Drum Brakes | Available with 12" Composite or 20" Spoke Rear Wheels. Urethane Tires only. | | \$ 995 |
| | User-activated Unilateral Wheel Lock | Pull to Lock. | | \$ 995 |
| <input type="checkbox"/> 2416310 | Left Activation | <input type="checkbox"/> 2416311 | Right Activation | |



| Wheel Lock Options | | |
|--------------------------|--------------------------------------|--|
| <input type="checkbox"/> | 2416100 Wheel Lock Extensions | Not compatible with Attendant Wheel Lock. \$ 60 |
| <input type="checkbox"/> | 2416200 Hill Holders | Requires Pneumatic Rear Tires with Flat-Free Inserts. \$ 450 |

| Anti-Tippers | | |
|-------------------------------------|----------------------------------|-------------------------------|
| <input checked="" type="checkbox"/> | 2417010 Rear Anti-Tippers | Required Safety Option. \$ 75 |

| Back Posts | | |
|--|---------|----------------------|
| Back Post Style | | |
| <input checked="" type="checkbox"/> | 2418010 | Stroller Handle Back |
| Back Post Height Measured from the Seat Pan to the Handle. | | |
| <input type="checkbox"/> | 2418425 | Short (22") |
| <input type="checkbox"/> | 2418431 | Tall (25") |
| Back Post Folding Fold-down Back is not compatible with Reclining Backrest. | | |
| <input checked="" type="checkbox"/> | 2419220 | Fold-down |
| Seat to Back Angle | | |
| <input type="checkbox"/> | 2419090 | Set at 90° |
| <input type="checkbox"/> | 2419098 | Set at 98° |
| <input type="checkbox"/> | 2419106 | Set at 106° |
| <input type="checkbox"/> | 2419114 | Set at 114° |
| <input type="checkbox"/> | 2419122 | Set at 122° |

| Back Options | | |
|--|---------|---|
| Angle/Height Adjustable Stroller Handle Not compatible with Headrest Mounting Plate. \$ 150 | | |
| <input type="checkbox"/> | 2419010 | Button on Outside |
| <input type="checkbox"/> | 2419011 | Button on Inside |
| Endomorph Back Backrest will be narrower than Seat by total of 2x the Offset selected below. Minimum width of Backrest is 12". \$ 395 | | |
| <input type="checkbox"/> | 2418102 | 1" Offset per side |
| <input type="checkbox"/> | 2418103 | 1.5" Offset per side |
| <input type="checkbox"/> | 2418104 | 2" Offset per side |
| <input type="checkbox"/> | 2419210 | 10" Backrest Extension Not compatible with Angle/Height Adjustable Stroller Handle. \$ 195 |
| <input type="checkbox"/> | 2419212 | Headrest Mounting Plate \$ 95 |
| Not compatible with Backrest Extension, Angle/Height Adjustable Stroller Handle. Not required if mounting to a solid Back. | | |

| Back Upholstery | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
|---|--------------------------|--|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|-------|---------------|--|--|--|--|--|--|--|--|-----|-----|-----|-----|-----|-----|-----|-----|-----|-------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|---------|--|--|--|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|-------------|--|--|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|-------------|--|--|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|
| <input type="checkbox"/> | 2420000 | No Back Upholstery | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Standard Upholstery \$ 265 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | 2420418 | Short (18" Height) | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | 2420421 | Regular (21" Height) | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | 2420424 | Tall (24" Height) | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | 2420450 | Tension Adjustable Upholstery Height of Upholstery is 16". \$ 365 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Position of Upholstery Select the desired position of the top of the Upholstery relative to the Seat Pan. | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| <table border="1"> <thead> <tr> <th rowspan="2">Style</th> <th colspan="9">22" Back Post</th> </tr> <tr> <th>17"</th> <th>18"</th> <th>19"</th> <th>20"</th> <th>21"</th> <th>22"</th> <th>23"</th> <th>24"</th> <th>25"</th> </tr> </thead> <tbody> <tr> <td>Short</td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Regular</td> <td></td> <td></td> <td></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Tension Adj</td> <td></td> <td></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td></td> <td></td> <td></td> </tr> </tbody> </table> | | | | | | | | | | Style | 22" Back Post | | | | | | | | | 17" | 18" | 19" | 20" | 21" | 22" | 23" | 24" | 25" | Short | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | | Regular | | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | | Tension Adj | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | | | | | | | | | | | |
| Style | 22" Back Post | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 17" | 18" | 19" | 20" | 21" | 22" | 23" | 24" | 25" | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Short | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Regular | | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Tension Adj | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| <table border="1"> <thead> <tr> <th rowspan="2">Style</th> <th colspan="9">25" Back Post</th> </tr> <tr> <th>17"</th> <th>18"</th> <th>19"</th> <th>20"</th> <th>21"</th> <th>22"</th> <th>23"</th> <th>24"</th> <th>25"</th> </tr> </thead> <tbody> <tr> <td>Short</td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td>Regular</td> <td></td> <td></td> <td></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td>Tall</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td>Tension Adj</td> <td></td> <td></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </tbody> </table> | | | | | | | | | | Style | 25" Back Post | | | | | | | | | 17" | 18" | 19" | 20" | 21" | 22" | 23" | 24" | 25" | Short | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | Regular | | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | Tall | | | | | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | Tension Adj | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| Style | 25" Back Post | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 17" | 18" | 19" | 20" | 21" | 22" | 23" | 24" | 25" | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Short | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Regular | | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Tall | | | | | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Tension Adj | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |

| Tilt Activation | | |
|---|---------|-----------------|
| Attendant Tilt only Select trigger mounting location below. | | |
| <input type="checkbox"/> | 2423011 | Right Back Post |
| <input type="checkbox"/> | 2423012 | Left Back Post |
| User Tilt only Full Length Arm Pads only. Select trigger mounting location below. \$ 100 | | |
| <input type="checkbox"/> | 2423021 | Right Armrest |
| <input type="checkbox"/> | 2423022 | Left Armrest |
| Attendant & User Tilt Full Length Arm Pads only. Select trigger mounting locations below. \$ 195 | | |
| <input type="checkbox"/> | 2423031 | Right Back Post |
| <input type="checkbox"/> | 2423032 | Left Back Post |
| <input type="checkbox"/> | 2423041 | Right Armrest |
| <input type="checkbox"/> | 2423042 | Left Armrest |

See next page for more options >



| | | | | | |
|--|---------|--------------------------------------|--|---|----------|
| <input type="checkbox"/> | 2423110 | Tilt Trigger Extension | | | \$ 50 |
| <p>Power Tilt 250 lbs Weight Capacity. 20° tilt range. Includes 2 batteries for off-chair charging. pdgmobility.com/power-tilt Minimum Seat Width is 15". Power Tilt is compatible with either O2 Holder or IV Pole on the same chair, not both.</p> | | | | | |
| <input type="checkbox"/> | 2423400 | Standard Power Tilt | | | \$ 2,095 |
| <input type="checkbox"/> | 2423410 | Power Tilt with Manual Tilt Override | Allows for instant manual release of tilt in case of power tilt failure. | | \$ 2,695 |
| | | Toggle Switch | <input type="checkbox"/> 2423421 Installed on Right side | <input type="checkbox"/> 2423422 Installed on Left side | |
| | | Button Switch | <input type="checkbox"/> 2423431 Installed on Right side | <input type="checkbox"/> 2423432 Installed on Left side | |
| | | Stealth Egg Switch | <input type="checkbox"/> 2423440 Not installed | | |
| | | Optional for Egg Switch: | <input type="checkbox"/> 2423445 Swing Away Mounting Hardware | Not compatible with all Headrest makes and models. | \$ 625 |
| Specify Headrest to be used: _____ | | | | | |

| Recline | | | | | |
|---|---------|------------------------------------|--|--|--------|
| <input type="checkbox"/> | 2423310 | Reclining Backrest (Manual) | | | \$ 895 |
| 30° maximum recline. Not compatible with Backrest Extension, Fold-down Back, Tension Adjustable Upholstery, Dynamic Rocker Back. Requires longer Frame. Maximum Seat Depth is 19". Finished Back Post Height is 3" taller than selected Back Post Height. | | | | | |

| Front Rigging | | | | | |
|--|---------|---|--|--|--------|
| <input type="checkbox"/> | 2425000 | No Front Rigging | <input type="checkbox"/> 2425001 Sector Blocks only | | |
| <p>Easy-Out™ Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame. pdgmobility.com/easy-out</p> | | | | | |
| <input type="checkbox"/> | 2425410 | Fixed 70° | | | \$ 95 |
| <input type="checkbox"/> | 2425420 | Universal Angle Adjustable (0°-110°) | Set at angle: _____ (use increments of 10°) | | \$ 295 |
| Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings. | | | | | |
| <p>Swing Away Front Rigging Legacy Design. Not compatible with Easy-Out.</p> | | | | | |
| <input type="checkbox"/> | 2425130 | Swing Away 60° | <input type="checkbox"/> 2425140 Swing Away 70° | | |
| <input type="checkbox"/> | 2425200 | Swing Away Elevating Legrests | Includes Calf Pads. | | \$ 375 |
| <input type="checkbox"/> | 2425300 | 1" Offset Hangers for Swing Away | | | \$ 75 |
| <input type="checkbox"/> | 2425110 | Standard Plug In 70° Front Rigging | Legacy Design. Not compatible with Easy-Out. | | |

| Footplates | | | | | |
|--------------------------|---------|---|--|--|--------|
| <input type="checkbox"/> | 2426010 | Composite Footplates | | | |
| <input type="checkbox"/> | 2426020 | Aluminum Angle Adjustable Footplates | | | \$ 165 |
| <input type="checkbox"/> | 2426030 | One Piece Fixed Footplate | Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable. | | \$ 350 |
| <input type="checkbox"/> | 2426040 | One Piece Flip-up Footplate | Requires selection of Composite or Aluminum Angle Adjustable. | | \$ 395 |
| <input type="checkbox"/> | 2426050 | Interlocking Angle Adjustable Footplates | Available for Seat Widths up to 21" only. | | \$ 195 |

| Footrest/Legrest Accessories | | | | | |
|--|---------|--|---|--|--------|
| <p>Calf Pads Not compatible with Gel Padding.</p> | | | | | |
| <input type="checkbox"/> | 2427012 | For Easy-Out™ Front Rigging | Includes reinforced premium Calf Pad mount. | | \$ 200 |
| <input type="checkbox"/> | 2427010 | For Swing Away Front Rigging | | | \$ 135 |
| <input type="checkbox"/> | 2427020 | Gel Padding | Not compatible with Calf Pad or Easy-Out Front Rigging. | | \$ 125 |
| <p>Gel Knee Pads</p> | | | | | |
| <input type="checkbox"/> | 2427032 | For Easy-Out™ Angle Adjustable Front Rigging | | | \$ 200 |
| <input type="checkbox"/> | 2427030 | For Swing Away ELR Front Rigging | | | \$ 100 |
| <input type="checkbox"/> | 2427040 | Length Adjustable Heel Loops | | | \$ 50 |
| <input type="checkbox"/> | 2427050 | Bodypoint Padded Calf Strap | | | \$ 195 |

| Front Rigging Extension Tubes | | | | | |
|---------------------------------------|---------|-------|------------------------------------|---|-------|
| <p>Regular Extension Tubes</p> | | | | | |
| <input type="checkbox"/> | 2428010 | Short | 14"-18" for Easy-Out Front Rigging | 13"-17" for Swing Away and Standard Plug In Front Rigging | |
| <input type="checkbox"/> | 2428020 | Long | 18"-24" for Easy-Out Front Rigging | 15"-19" for Swing Away and Standard Plug In Front Rigging | \$ 75 |

See next page for more options >



BENTLEY LT

Order Form - Distributor

Prices in Canadian Dollars. January 1, 2024



Lower Leg Length If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.

Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.

Miscellaneous Options

One-Arm Drive Available with 22" and 24" Composite Rear Wheels. Minimum 17" STF and 18" STF, respectively. Limits Tilt. Not compatible with User & Attendant Wheel Locks, Drum Brakes, Quick Release Rear Axles, Pneumatic Tires. \$ 3,500

2430010 Installed on Right side 2430011 Installed on Left side

2430012 Large Hand Setup 24" Rear Wheels with 22" Driving Handrim.

Spoke Guards Available for 20"-24" Composite Rear Wheels. Made from durable and lightweight aluminum. \$ 235

2430020 Black 2430021 Frame Color

2430050 **Protective Lower Padding** \$ 995

Stealth Flip Down Abductor Requires Solid Seat Pan. \$ 825

2430060 Small Pad 3"x4"x3"x2" size 2430062 Large Pad 5"x5"x4"x3" size

2430080 **Gel Hip Guides** Reduces available Seat Width by 2". \$ 595

2430110 **IV Pole** Not compatible with Attendant Wheel Locks. Requires frame modification and cannot be added to chair later. \$ 345

2430120 **O2 Holder** Not compatible with Attendant Wheel Locks. Requires frame modification and cannot be added to chair later. \$ 495

2430200 **Q'Straint Lap Belt** Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" Belt Length. \$ 295

2430222 **Wheelchair Backpack** Black with PDG Logo. 17"x14" size. \$ 165

2430170 **Adhesive Back Hook Strip, 2" x 30"** For Cushion mounting. \$ 10

On-Chair Seating & Positioning

Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. pdgmobility.com/seating

Order Notes

Please send your completed order form by email to info@pdgmobility.com.