



### Customer and Order Information

|  |                  |                    |
|--|------------------|--------------------|
| Purchase Order *                         | Tag              | v.0101             |
| Dealer/Provider Name *                   |                  |                    |
| Purchaser Name and Contact Information * |                  |                    |
| Shipping Address *                       |                  |                    |
| User Information (if available)          | Weight _____ lbs | Height _____ ft/in |
|  |                  | *required          |

### Base Chair

|   |  |          |
|---|--|----------|
| <input checked="" type="checkbox"/> 1101000 | <b>Eclipse Bariatric Manual Wheelchair</b> | \$ 3,600 |
|---|--|----------|

### Weight Capacity

|                                  |                   |  |
|----------------------------------|-------------------|--|
| <b>Chair Weight Capacity</b>     |                   |  |
| <input type="checkbox"/> 1102000 | 600 lbs Capacity  |  |
| <input type="checkbox"/> 1102100 | 1000 lbs Capacity | Not compatible with Transit Tie Down System. |
|                                  |                   | \$ 300                                       |

### Frame Options

|  |                                       |        |
|--|---------------------------------------|--------|
| <input type="checkbox"/> 1102300   | <b>Transit Tie Down System (TTDS)</b> | \$ 395 |
| TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Tested to max. Weight Capacity of 600 lbs. Not compatible with O2 holder or IV pole. |                                       |        |

### Seat Width

Overall chair width is +8"

**Seat Widths** Seat Width is adjustable with the optional Adjustment Kit (sold separately).

|                                  |                                     |  |   |          |
|----------------------------------|-------------------------------------|--|---|----------|
| <input type="checkbox"/> 1103520 | 20" Seat Width                      | <input type="checkbox"/> 1103522         | 22" Seat Width                                    | \$ 675   |
| <input type="checkbox"/> 1103524 | 24" Seat Width                      | <input type="checkbox"/> 1103526         | 26" Seat Width                                    | \$ 675   |
| <input type="checkbox"/> 1103528 | 28" Seat Width                      | <input type="checkbox"/> 1103530         | 30" Seat Width                                    | \$ 875   |
| <input type="checkbox"/> 1103532 | 32" Seat Width                      | <input type="checkbox"/> 1103536         | 36" Seat Width                                    | \$ 1,100 |
| <input type="checkbox"/> 1103540 | 40" Seat Width                      | <input type="checkbox"/> 1103544         | 44" Seat Width                                    | \$ 1,300 |
| <input type="checkbox"/> 1103548 | 48" Seat Width                      |  |   | \$ 1,300 |
| <input type="checkbox"/> 1103599 | Unlisted Seat Width (between sizes) | Seat Width: _____ (use increments of 1") | Add this fee to the price of the next-bigger size | \$ 100   |

**i** Seat Widths greater than 48" are available for Custom Chairs

### Seat Depth

**Seat Depths (Adjustable)**

|                                  |                |                                  |                |                                  |                |
|----------------------------------|----------------|----------------------------------|----------------|----------------------------------|----------------|
| <input type="checkbox"/> 1104016 | 16" Seat Depth | <input type="checkbox"/> 1104017 | 17" Seat Depth | <input type="checkbox"/> 1104018 | 18" Seat Depth |
| <input type="checkbox"/> 1104019 | 19" Seat Depth | <input type="checkbox"/> 1104020 | 20" Seat Depth | <input type="checkbox"/> 1104021 | 21" Seat Depth |
| <input type="checkbox"/> 1104022 | 22" Seat Depth |                                  |                |                                  |                |

**Custom Seat Depths**


|                                  |                |                                  |                |                          |        |
|----------------------------------|----------------|----------------------------------|----------------|--------------------------|--------|
| <input type="checkbox"/> 1104023 | 23" Seat Depth | <input type="checkbox"/> 1104024 | 24" Seat Depth | Not available with TTDS. | \$ 400 |
|----------------------------------|----------------|----------------------------------|----------------|--------------------------|--------|

### Frame Color

**PDG Premium Colors** View color swatches on PDG's website. [pdgmobility.com/colors](http://pdgmobility.com/colors)

|                                  |                    |                                  |                    |                                  |                  |
|----------------------------------|--------------------|----------------------------------|--------------------|----------------------------------|------------------|
| <input type="checkbox"/> 1105001 | Black Shimmer      | <input type="checkbox"/> 1105015 | Sky Blue           | <input type="checkbox"/> 1105016 | Poppy Red        |
| <input type="checkbox"/> 1105003 | Fine Texture Black | <input type="checkbox"/> 1105010 | Lollypop Blue      | <input type="checkbox"/> 1105011 | Lollypop Red     |
| <input type="checkbox"/> 1105017 | Arctic Grey        | <input type="checkbox"/> 1105018 | Tropical Teal      | <input type="checkbox"/> 1105019 | Tangerine Orange |
| <input type="checkbox"/> 1105020 | Olive Green        | <input type="checkbox"/> 1105021 | Lush Lilac         | <input type="checkbox"/> 1105012 | Bubble Gum Pink  |
| <input type="checkbox"/> 1105022 | Leaf Green         | <input type="checkbox"/> 1105023 | Boysenberry Purple | <input type="checkbox"/> 1105014 | Daffodil Yellow  |

1105200 **Custom Color** Contact PDG for specifications. \$ 300

| Armrests  |  |
|---|--|
| <input type="checkbox"/> 1106000  | <b>No Armrests</b>   |
| <b>Dual Post Height Adjustable Armrests</b> <span style="float: right;">\$ 295</span> |  |
| <input type="checkbox"/> 1106214 10"-14"  | <input type="checkbox"/> 1106216 12"-16" Not compatible with Desk Length Armpads.  |
| <input type="checkbox"/> 1106200  | <b>Dual Post Fixed Height Armrests (12")</b>   |
|       | <b>Armrest Height</b> If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.<br>Seat Pan to top of Armrest: _____ Does not include height of Armpad |


| Armpads  |   |
|--|---|
| <b>Regular Armpads</b>   |   |
| <input type="checkbox"/> 1107100 Desk Length - 2" x 10"  | <input type="checkbox"/> 1107101 Full Length - 2" x 14"   |
| <b>Wide Armpads</b> May interfere with wheels in some tilt configurations. <span style="float: right;">\$ 200</span> |   |
| <input type="checkbox"/> 1107110 Desk Length - 3" x 10"  | <input type="checkbox"/> 1107111 Full Length - 3" x 12"   |
| <b>Gel Armpads</b> May interfere with wheels in some tilt configurations. <span style="float: right;">\$ 135</span>  |   |
| <input type="checkbox"/> 1107120 Full Length - 2" x 12"  | <input type="checkbox"/> 1107121 Full Length - 3.5" x 12" |

| Armrest Options                  |   |
|----------------------------------|---|
| <input type="checkbox"/> 1107250 | <b>Forward Guards</b> <span style="float: right;">\$ 200</span>   |
| <input type="checkbox"/> 1107400 | <b>Armrest Extensions - Forward</b> Not compatible with Desk Length Armpads. <span style="float: right;">\$ 75</span> |

| Seat Pan                                    |  |
|---|--|
| <input checked="" type="checkbox"/> 1108001 | <b>Solid Aluminum Seat Pan</b> <span style="float: right;">\$ 165</span> |

| Seat Height  |  |
|--|--|
| <b>Front Seat Height</b> Front and Rear Height can be max. 1" apart.                     |  |
| <input type="checkbox"/> 1110014 14" STF Front   | <input type="checkbox"/> 1110015 15" STF Front |
| <input type="checkbox"/> 1110017 17" STF Front   | <input type="checkbox"/> 1110018 18" STF Front |
| <input type="checkbox"/> 1110019 19" STF Front   | <input type="checkbox"/> 1110020 20" STF Front |
| <input type="checkbox"/> 1110016 16" STF Front <span style="float: right;">\$ 200</span> |  |
| <b>Rear Seat Height</b> Front and Rear Height can be max. 1" apart.                      |  |
| <input type="checkbox"/> 1110114 14" STF Rear  | <input type="checkbox"/> 1110115 15" STF Rear  |
| <input type="checkbox"/> 1110117 17" STF Rear  | <input type="checkbox"/> 1110118 18" STF Rear  |
| <input type="checkbox"/> 1110119 19" STF Rear  | <input type="checkbox"/> 1110120 20" STF Rear  |
| <input type="checkbox"/> 1110116 16" STF Rear <span style="float: right;">\$ 200</span>  |  |

| Casters                                     |  |
|---|--|
| <input checked="" type="checkbox"/> 1111050 | <b>5" Casters With Heavy Duty Steel Fork</b> |

| Rear Wheels  |  |
|--|--|
| <input checked="" type="checkbox"/> 1112224  | <b>24" Heavy Duty Spoke Wheels</b>   |
|  | <b>Rear Wheel Position</b> Positions 3 and 4 not compatible with TTDS.<br>In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference. |
| <input type="checkbox"/> 1113101 Set at 1 (Rear)                                   | <input type="checkbox"/> 1113102 Set at 2  |
| <input type="checkbox"/> 1113103 Set at 3  | <input type="checkbox"/> 1113104 Set at 4 (Front)  |

| Axle Type                                   |                                 |
|---|---------------------------------|
| <input checked="" type="checkbox"/> 1113000 | <b>Quick Release Rear Axles</b> |

| Rear Tires                                  |                       |
|---|-----------------------|
| <input checked="" type="checkbox"/> 1114000 | <b>Urethane Tires</b> |



**Handrims**

1115000 **No Handrims**

1115010 **Standard Aluminum Handrims**

1115020 **Plastic Coated Handrims** \$ 250

**Wheel Locks**

**User Wheel Locks**

1116010 Push to Lock  1116020 Pull to Lock

**Wheel Lock Options**

1116100 **Wheel Lock Extensions** \$ 60

**Anti-Tippers**

1117010 **Rear Anti-Tippers** Required Safety Option. \$ 75

**Back Posts**

**Back Post Style**

1118030 Stroller Handle Back  1118040 Cane Handle Back Not compatible with Headrest Mounting Plate.

**Back Post Height** Measured from the Seat Pan to the Handle.

1118425 Short (22")  1118431 Tall (25")

**Seat to Back Angle**

1119085 Set at 85°  1119090 Set at 90°

1119095 Set at 95°  1119100 Set at 100°

**Back Options**

**Endomorph Back** Backrest will be narrower than Seat by total of 2x the Offset selected below. \$ 395

1118102 1" Offset per side  1118103 1.5" Offset per side  1118104 2" Offset per side

1119210 **10" Backrest Extension** \$ 150

1119212 **Headrest Mounting Plate** Not required if mounting to a solid Back. \$ 95

**Back Upholstery**

1120000 **No Back Upholstery**

**Tension Adjustable** \$ 525

1120421 Regular (18" Height)  1120424 Tall (22" Height)

**Position of Upholstery** Select the desired position of the top of the Upholstery relative to the Seat Pan.

| 22" Back Post |                          |                          |                          |                          |                          |                          |     |     |     |
|---------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|-----|-----|-----|
| Style         | 17"                      | 18"                      | 19"                      | 20"                      | 21"                      | 22"                      | 23" | 24" | 25" |
| Regular       | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |     |     |     |
| Tall          |                          |                          |                          | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |     |     |     |

| 25" Back Post |                          |                          |                          |                          |                          |                          |                          |                          |                          |
|---------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|
| Style         | 17"                      | 18"                      | 19"                      | 20"                      | 21"                      | 22"                      | 23"                      | 24"                      | 25"                      |
| Regular       | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| Tall          |                          |                          |                          |                          | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |

**Recline**

1123310 **Reclining Backrest (Manual)** \$ 895

30° maximum recline. Not compatible with Backrest Extension. Reclining Backrest is 1" narrower than Seat Width. Back Frame is not removable. Finished Back Post Height is 3" taller than selected Back Post Height. Not recommended for user weight above 450 lbs. Maximum Seat Depth is 23".

### Front Rigging

|  |         |   |                          |                         |                           |        |
|--|---------|---|--------------------------|-------------------------|---------------------------|--------|
| <input type="checkbox"/>   | 1125000 | <b>No Front Rigging</b>                   | <input type="checkbox"/> | 1125001                 | <b>Sector Blocks only</b> |        |
| <input type="checkbox"/>   | 1125002 | <b>Center Mount Front Rigging Bracket</b> |                          |                         |                           | \$ 300 |
| <b>PDG Easy-Out™ Heavy Duty Front Rigging</b> Effortless insertion and removal. No Hanger Bracket remains on frame. <a href="http://pdgmobility.com/easy-out">pdgmobility.com/easy-out</a> |         |   |                          |                         |                           |        |
| <input type="checkbox"/>   | 1125410 | Fixed 60°                                 |                          |                         |                           | \$ 95  |
| <input type="checkbox"/>   | 1125420 | Universal Angle Adjustable (0°-60°)       | Set at angle: _____      | (use increments of 10°) |                           | \$ 295 |
| Adjust angle by loosening one bolt. Tool-free leg-length adjustment.   |         |   |                          |                         |                           |        |
| <b>Swing Away Front Rigging</b> Legacy Design. Not compatible with Easy-Out.   |         |   |                          |                         |                           |        |
| <input type="checkbox"/>   | 1125140 | Swing Away 60°                            |                          |                         |                           |        |
| <input type="checkbox"/>   | 1125200 | Swing Away Elevating Legrests             | Includes Calf Pads.      |                         |                           | \$ 375 |
| <input type="checkbox"/>   | 1125300 | 1" Offset Hangers for Swing Away          |                          |                         |                           | \$ 75  |


### Footplates

|                                     |         |   |  |  |  |        |
|-------------------------------------|---------|---|--|--|--|--------|
| <input checked="" type="checkbox"/> | 1126020 | <b>Aluminum Angle Adjustable Footplates</b> |  |  |  | \$ 185 |
|-------------------------------------|---------|---|--|--|--|--------|

### Footrest/Legrest Accessories

|                          |         |  |  |  |  |        |
|--------------------------|---------|--|--|--|--|--------|
| <input type="checkbox"/> | 1127015 | <b>Oversized 9" x 10" Calf Pad</b>           | Not compatible with Seat Width of less than 25". |  |  | \$ 695 |
| <b>Calf Pads</b>         |         |  |  |  |  |        |
| <input type="checkbox"/> | 1127012 | For Easy-Out™ Front Rigging                  | Includes reinforced premium Calf Pad mount.      | <b>Currently not available for Easy-Out.</b> |  | \$ 200 |
| <input type="checkbox"/> | 1127010 | For Swing Away Front Rigging                 |  |  |  | \$ 135 |
| <b>Gel Knee Pads</b>     |         |  |  |  |  |        |
| <input type="checkbox"/> | 1127032 | For Easy-Out™ Angle Adjustable Front Rigging |  |  |  | \$ 200 |
| <input type="checkbox"/> | 1127030 | For Swing Away ELR Front Rigging             |  |  |  | \$ 100 |
| <input type="checkbox"/> | 1127040 | <b>Length Adjustable Heel Loops</b>          |  |  |  | \$ 50  |
| <input type="checkbox"/> | 1127051 | <b>Calf Strap</b>                            | Available up to maximum size of 32".             |  |  | \$ 195 |

### Front Rigging Extension Tubes

|  |  |       |                                    |                                      |  |       |
|--|--|-------|------------------------------------|--------------------------------------|--|-------|
| <b>Regular Extension Tubes</b>   |  |       |                                    |                                      |  |       |
| <input type="checkbox"/>   | 1128010  | Short | 14"-18" for Easy-Out Front Rigging | 13"-17" for Swing Away Front Rigging |  |       |
| <input type="checkbox"/>   | 1128020  | Long  | 18"-24" for Easy-Out Front Rigging | 15"-19" for Swing Away Front Rigging |  | \$ 75 |
|  | <b>Lower Leg Length</b> If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.<br>Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion. |       |                                    |                                      |  |       |

### Miscellaneous Options

|                          |         |   |  |  |  |        |
|--------------------------|---------|---|--|--|--|--------|
| <input type="checkbox"/> | 1130110 | <b>IV Pole</b>                            | Not compatible with TTDS.  |  |  | \$ 345 |
| <input type="checkbox"/> | 1130120 | <b>O2 Holder</b>                          | Not compatible with TTDS.  |  |  | \$ 495 |
| <input type="checkbox"/> | 1130200 | <b>Q'Straint Lap Belt</b>                 | Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" Belt Length. |  |  | \$ 295 |
| <input type="checkbox"/> | 1130222 | <b>Wheelchair Backpack</b>                | Black with PDG Logo. 17"x14" size.   |  |  | \$ 165 |
| <input type="checkbox"/> | 1130170 | <b>Adhesive Back Hook Strip, 2" x 30"</b> | For Cushion mounting.  |  |  | \$ 10  |

### On-Chair Seating & Positioning

|   |   |  |  |  |  |  |
|---|---|--|--|--|--|--|
|  | Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. <a href="http://pdgmobility.com/seating">pdgmobility.com/seating</a> |  |  |  |  |  |
|---|---|--|--|--|--|--|

### Order Notes

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com).