



FUZE T20

Order Form - Distributor

Prices in Canadian Dollars. January 1, 2024

Customer and Order Information

Purchase Order *	Tag	v.0401
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in

Base Chair

<input checked="" type="checkbox"/> 1601000	Fuze T20 20° Manual Tilt Wheelchair	\$ 4,400
i	Also available: the lightweight Bentley LT with aluminum frame. See separate order form. pdgmobility.com/blt	

Frame Type & Weight Capacity

Frame Type		
<input type="checkbox"/> 1602010	Fixed Width Frame 250 lbs Weight Capacity. Select a Fixed Seat Width below.	
<input type="checkbox"/> 1602025	Adjustable Width Frame 250 lbs Weight Capacity. Select an Adjustable Seat Width below.	\$ 100
User Weight Range Required for setup of tilt uplift force.		
<input type="checkbox"/> 1602406	Up to 150 lbs	<input type="checkbox"/> 1602410 150 lbs to 200 lbs
		<input type="checkbox"/> 1602412 200 lbs to 250 lbs

Frame Options

i	The Transit Tie Down System (TTDS) is no longer available for the Fuze T20. We apologize for the inconvenience. Consider the lightweight Bentley LT with aluminum frame as an alternative. See separate order form.	
----------	---	--

Seat Width

Overall chair width is +9.5" for 12/16" RW and +10.5" for 20/22/24" RW

Fixed Seat Widths Order with Fixed Frame.		
<input type="checkbox"/> 1603015	15" Seat Width	<input type="checkbox"/> 1603016 16" Seat Width
<input type="checkbox"/> 1603017	17" Seat Width	<input type="checkbox"/> 1603018 18" Seat Width
<input type="checkbox"/> 1603019	19" Seat Width	<input type="checkbox"/> 1603020 20" Seat Width
<input type="checkbox"/> 1603021	21" Seat Width	<input type="checkbox"/> 1603022 22" Seat Width
Adjustable Seat Widths Order with Adjustable Frame.		
<input type="checkbox"/> 1603113	13"-15" Range - Set at 13"	<input type="checkbox"/> 1603114 13"-15" Range - Set at 14"
<input type="checkbox"/> 1603115	13"-15" Range - Set at 15"	<input type="checkbox"/> 1603215 15"-20" Range - Set at 15"
<input type="checkbox"/> 1603216	15"-20" Range - Set at 16"	<input type="checkbox"/> 1603217 15"-20" Range - Set at 17"
<input type="checkbox"/> 1603218	15"-20" Range - Set at 18"	<input type="checkbox"/> 1603219 15"-20" Range - Set at 19"
<input type="checkbox"/> 1603220	15"-20" Range - Set at 20"	<input type="checkbox"/> 1603321 21"-22" Range - Set at 21"
<input type="checkbox"/> 1603322	21"-22" Range - Set at 22"	

Seat Depth

Standard Seat Depths (Adjustable)		
<input type="checkbox"/> 1604015	15" Seat Depth	<input type="checkbox"/> 1604016 16" Seat Depth
<input type="checkbox"/> 1604017	17" Seat Depth	<input type="checkbox"/> 1604018 18" Seat Depth
<input type="checkbox"/> 1604019	19" Seat Depth	<input type="checkbox"/> 1604020 20" Seat Depth

Frame Color

PDG Premium Colors View color swatches on PDG's website. pdgmobility.com/colors		
<input type="checkbox"/> 1605001	Black Shimmer	<input type="checkbox"/> 1605015 Sky Blue
<input type="checkbox"/> 1605003	Fine Texture Black	<input type="checkbox"/> 1605016 Poppy Red
<input type="checkbox"/> 1605017	Arctic Grey	<input type="checkbox"/> 1605011 Lollypop Red
<input type="checkbox"/> 1605020	Olive Green	<input type="checkbox"/> 1605018 Tropical Teal
<input type="checkbox"/> 1605019	Tangerine Orange	<input type="checkbox"/> 1605021 Lush Lilac
<input type="checkbox"/> 1605022	Leaf Green	<input type="checkbox"/> 1605012 Bubble Gum Pink
<input type="checkbox"/> 1605023	Boysenberry Purple	<input type="checkbox"/> 1605014 Daffodil Yellow
<input type="checkbox"/> 1605200	Custom Color Contact PDG for specifications.	\$ 300

Armrests

1606000 **No Armrests**

Height Adjustable Armrests (Single Post) 8"-11" not compatible with 22" and 24" Rear Wheels. \$ 295

1606111 8"-11" Range 1606114 10"-14" 1606116 12"-16" Range

1606100 **Fixed Height Armrests (12")** Full Length Armpad only.

Armrest Height If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.
 Seat Pan to top of Armrest: _____ Does not include height of Armpad

Locking Flip-up Armrests Not recommended for users with excessive strength or tone. Not compatible with Fold-down Back. \$ 295

1606310 Set at 10" 1606311 Set at 11" 1606312 Set at 12"
 1606313 Set at 13" 1606314 Set at 14"

Armpads

Regular Armpads

1607100 Desk Length - 2" x 10" 1607101 Full Length - 2" x 14"

Waterfall Armpads \$ 115

1607105 Desk Length - 2" x 10" 1607106 Full Length - 2" x 14"

Wide Armpads May interfere with wheels in some tilt configurations. \$ 200

1607110 Desk Length - 3" x 10" 1607111 Full Length - 3" x 12"

Gel Armpads May interfere with wheels in some tilt configurations. \$ 135

1607120 Full Length - 2" x 12" 1607121 Full Length - 3.5" x 12"

Contour Arm Trough with Hand Pad May interfere with wheels in some tilt configurations. \$ 125

1607130 Right (4.5" x 18") 1607131 Left (4.5" x 18")

Armrest Options

Armrest Wideners Move Armrests outward by up to 1" on each side. Seat size does not change.

1607210 1" Total Additional Width \$ 95
 1607211 2" Total Additional Width \$ 195

Seat Pan

1608000 **No Seat Pan**

1608001 **Solid Aluminum Seat Pan** \$ 135

Seat Height

Don't forget **Select desired Seat-to-Floor Height** STF Height is adjustable.

RW	13"	14"	15"	16"	17"	18"	19"
12"	<input type="checkbox"/>						
16"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				
20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

RW	13"	14"	15"	16"	17"	18"	19"
12"							
16"		<input type="checkbox"/>	<input type="checkbox"/>				
20"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

RW	13"	14"	15"	16"	17"	18"	19"
12"							
16"			<input type="checkbox"/>				
20"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
22"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

RW	13"	14"	15"	16"	17"	18"	19"
12"							
16"							
20"				<input type="checkbox"/>	<input type="checkbox"/>		
22"				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
24"				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Casters

Urethane Casters

- | | | | | | | |
|----------------------------------|------------------------|----------------------------------|------------------------|----------------------------------|------------------------|--------|
| <input type="checkbox"/> 1611051 | 5" x 1.25" Caster size | <input type="checkbox"/> 1611061 | 6" x 1.25" Caster size | <input type="checkbox"/> 1611081 | 8" x 1.25" Caster size | |
| <input type="checkbox"/> 1611063 | 6" x 1.50" Caster size | | | | | \$ 85 |
| <input type="checkbox"/> 1611084 | 8" x 2.00" Caster size | | | | | \$ 135 |

Pneumatic Casters

- | | | | | | | |
|----------------------------------|------------------------|--|--|--|--|--------|
| <input type="checkbox"/> 1611161 | 6" x 1.25" Caster size | | | | | \$ 110 |
| <input type="checkbox"/> 1611181 | 8" x 1.25" Caster size | | | | | \$ 120 |
| <input type="checkbox"/> 1611184 | 8" x 2.00" Caster size | | | | | \$ 150 |

Rear Wheels

Composite Wheels

- | | | | | | | |
|----------------------------------|----------------|----------------------------------|----------------|----------------------------------|----------------|--|
| <input type="checkbox"/> 1612012 | 12" Wheel size | <input type="checkbox"/> 1612016 | 16" Wheel size | <input type="checkbox"/> 1612020 | 20" Wheel size | |
| <input type="checkbox"/> 1612022 | 22" Wheel size | <input type="checkbox"/> 1612024 | 24" Wheel size | | | |

Spoke Wheels

- | | | | | | | |
|----------------------------------|----------------|----------------------------------|----------------|----------------------------------|----------------|-------|
| <input type="checkbox"/> 1612120 | 20" Wheel size | <input type="checkbox"/> 1612122 | 22" Wheel size | <input type="checkbox"/> 1612124 | 24" Wheel size | \$ 75 |
|----------------------------------|----------------|----------------------------------|----------------|----------------------------------|----------------|-------|

- | | | | | | | |
|----------------------------------|------------------------------------|--|--|--|--|--------|
| <input type="checkbox"/> 1612224 | 24" Heavy Duty Spoke Wheels | | | | | \$ 350 |
|----------------------------------|------------------------------------|--|--|--|--|--------|



Rear Wheel Position

Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.

In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

- | | | | | | | | |
|----------------------------------|------|-------------|----------------------------------|-----|----------------------------------|-------|-------------|
| <input type="checkbox"/> 1612401 | Rear | More stable | <input type="checkbox"/> 1612402 | Mid | <input type="checkbox"/> 1612403 | Front | Less stable |
|----------------------------------|------|-------------|----------------------------------|-----|----------------------------------|-------|-------------|

Axle Type

- | | |
|----------------------------------|---------------------------------|
| <input type="checkbox"/> 1613000 | Quick Release Rear Axles |
| <input type="checkbox"/> 1613010 | Fixed Rear Axles |

Rear Tires

- | | | | |
|----------------------------------|---|---|-------|
| <input type="checkbox"/> 1614000 | Urethane Tires | | |
| <input type="checkbox"/> 1614010 | Pneumatic Tires | May result in a finished Seat Height 1/2" greater than stated. | \$ 60 |
| <input type="checkbox"/> 1614020 | Pneumatic Tires with Flat-Free Inserts | May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels. | \$ 95 |

Handrims

Handrims are available only with 20", 22", and 24" Rear Wheels.

- | | | | |
|----------------------------------|--------------------------------|--|--------|
| <input type="checkbox"/> 1615000 | No Handrims | | |
| <input type="checkbox"/> 1615010 | Standard Handrims | | |
| <input type="checkbox"/> 1615020 | Plastic Coated Handrims | Not compatible with Drum Brakes. | \$ 150 |
| <input type="checkbox"/> 1615030 | Projection Handrims | Composite Wheels only. 10 Projections, 30°, 1-1.8" Length. | \$ 595 |

Wheel Locks & Brakes

User Wheel Locks

- | | | | | |
|---|---|---|------------------|---------------|
| <input type="checkbox"/> 1616010 | Push to Lock | <input type="checkbox"/> 1616020 | Pull to Lock | |
| <input type="checkbox"/> 1616030 | Attendant Wheel Locks (Push to Lock) | Not compatible with Drum Brakes, Wheel Lock Extensions. | | \$ 85 |
| <input type="checkbox"/> 1616040 | User & Attendant Wheel Locks | Not compatible with Drum Brakes. | | \$ 125 |
| <input type="checkbox"/> 1616300 | Drum Brakes | Available with 12" Composite or 20" Spoke Rear Wheels. Urethane Tires only. | | \$ 995 |
| User-activated Unilateral Wheel Lock | | | | Pull to Lock. |
| <input type="checkbox"/> 1616310 | Left Activation | <input type="checkbox"/> 1616311 | Right Activation | \$ 995 |

Wheel Lock Options

- | | | | |
|----------------------------------|------------------------------|---|-------|
| <input type="checkbox"/> 1616100 | Wheel Lock Extensions | Not compatible with Attendant Wheel Lock. | \$ 60 |
|----------------------------------|------------------------------|---|-------|

Anti-Tippers

- | | | | |
|---|--------------------------|-------------------------|-------|
| <input checked="" type="checkbox"/> 1617010 | Rear Anti-Tippers | Required Safety Option. | \$ 75 |
|---|--------------------------|-------------------------|-------|



FUZE T20

Order Form - Distributor

Prices in Canadian Dollars. January 1, 2024

Back Posts

Back Post Style

- 1618100 Stroller Handle Back 1618105 Cane Handle Back Not compatible with Headrest Mounting Plate.

Back Post Type

- 1618200 Straight

Back Post Height Measured from the Seat Pan to the Handle.

- 1618420 Short (17") 1618425 Tall (23")

Back Post Folding

- 1619220 Fold-down 1619221 Non-folding

Seat to Back Angle

- 1619090 Set at 90° 1619096 Set at 96° 1619102 Set at 102°
 1619108 Set at 108° 1619114 Set at 114° 1619120 Set at 120°

Back Options

Angle/Height Adjustable Stroller Handle Not compatible with Headrest Mounting Plate or Cane Handle Back.

\$ 150

- 1619010 Button on Outside 1619011 Button on Inside

Removable Cane Handle Push Bar Not compatible with Stroller Handle Back.

\$ 140

Endomorph Back Backrest will be narrower than Seat by total of 2x the Offset selected below. Minimum width of Backrest is 12".

\$ 395

- 1618102 1" Offset per side 1618103 1.5" Offset per side 1618104 2" Offset per side

Headrest Mounting Plate Not compatible with Angle/Height Adjustable Stroller Handle, Cane Handle Back. Not required if mounting to a solid Back.

\$ 95

Back Upholstery

No Back Upholstery

- 1620000 **No Back Upholstery**

Standard Upholstery Not available with Cane Handle Back.

\$ 265

- 1620418 Short (18" Height) 1620421 Regular (21" Height)

Tension Adjustable Upholstery Height of Upholstery is 16". Not available with Cane Handle Back.

\$ 365

Position of Upholstery Select the desired position of the top of the Upholstery relative to the Seat Pan.

		17" Back Post								
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"	
Short	<input type="checkbox"/>									

		23" Back Post								
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"	
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Regular					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Tension Adj				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

Tilt Activation

Attendant Tilt only Select trigger mounting location below.

- 1623011 Right Back Post 1623012 Left Back Post

Front Rigging

No Front Rigging 1625001 Sector Blocks only

- 1625000 **No Front Rigging** 1625001 **Sector Blocks only**

Easy-Out™ Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame. pdgmobility.com/easy-out

- 1625410 Fixed 70° \$ 95

Universal Angle Adjustable (0°-110°) Set at angle: _____ (use increments of 10°) Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings.

- 1625420 **Universal Angle Adjustable (0°-110°)** \$ 295

Swing Away Front Rigging Legacy Design. Not compatible with Easy-Out.

- 1625130 Swing Away 60° 1625140 Swing Away 70°

Swing Away 90° Includes 11"-16" Extension Tubes. Not compatible with 8" Casters.

- 1625160 **Swing Away 90°** \$ 75

Swing Away Elevating Legrests Includes Calf Pads.

- 1625200 **Swing Away Elevating Legrests** \$ 375

1" Offset Hangers for Swing Away

- 1625300 **1" Offset Hangers for Swing Away** \$ 75

Standard Plug In 70° Front Rigging Legacy Design. Not compatible with Easy-Out.

- 1625110 **Standard Plug In 70° Front Rigging**


Footplates

<input type="checkbox"/>	1626010	Composite Footplates		
<input type="checkbox"/>	1626020	Aluminum Angle Adjustable Footplates		\$ 165
<input type="checkbox"/>	1626030	One Piece Fixed Footplate	Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.	\$ 350
<input type="checkbox"/>	1626040	One Piece Flip-up Footplate	Requires selection of Composite or Aluminum Angle Adjustable.	\$ 395
<input type="checkbox"/>	1626050	Interlocking Angle Adjustable Foot Plates	Available for Seat Widths up to 21" only.	\$ 195

Footrest/Legrest Accessories

Calf Pad Not compatible with Gel Padding.				
<input type="checkbox"/>	1627012	For Easy-Out™ Front Rigging	Includes reinforced premium Calf Pad mount.	\$ 200
<input type="checkbox"/>	1627010	For Swing Away Front Rigging		\$ 135
<input type="checkbox"/>	1627020	Gel Padding	Not compatible with Calf Pad or Easy-Out Front Rigging.	\$ 125
Gel Knee Pad				
<input type="checkbox"/>	1627032	For Easy-Out™ Angle Adjustable Front Rigging		\$ 200
<input type="checkbox"/>	1627030	For Swing Away ELR Front Rigging		\$ 100
<input type="checkbox"/>	1627040	Length Adjustable Heel Loops		\$ 50
<input type="checkbox"/>	1627050	Bodypoint Padded Calf Strap		\$ 195

Front Rigging Extension Tubes

Regular Extension Tubes No selection required for 90° Front Rigging				
<input type="checkbox"/>	1628010	Short	14"-18" for Easy-Out Front Rigging 13"-17" for Swing Away and Standard Plug In Front Rigging	
<input type="checkbox"/>	1628020	Long	18"-24" for Easy-Out Front Rigging 15"-19" for Swing Away and Standard Plug In Front Rigging	\$ 75
	Lower Leg Length If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.			
	Seat Pan to Foot Plates: _____		Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.	

Miscellaneous Options

Spoke Guards Available for 20"-24" Composite Rear Wheels. Made from durable and lightweight aluminum.				
<input type="checkbox"/>	1630020	Black	<input type="checkbox"/>	1630021 Frame Color
Stealth Flip Down Abductor Requires Solid Seat Pan.				
<input type="checkbox"/>	1630060	Small Pad	3"x4"x3"x2" size	<input type="checkbox"/>
<input type="checkbox"/>	1630062	Large Pad	5"x5"x4"x3" size	
<input type="checkbox"/>	1630080	Gel Hip Guides	Reduces available Seat Width by 2".	\$ 595
<input type="checkbox"/>	1630090	Hip Belt Mount	Provides improved mounting options for Bodypoint Hip Belt or other third-party belts.	\$ 55
IV Pole Requires minimum 15" STF Height.				
<input type="checkbox"/>	1630110	Installed on Left Side	<input type="checkbox"/>	1630111 Installed on Right Side
<input type="checkbox"/>	1630120	O2 Holder	15" STF and higher.	\$ 495
<input type="checkbox"/>	1630160	Clear Lap Tray (Slide On)	Available for Seat Widths of 17"-22" only. Not compatible with Contour Arm Trough.	\$ 490
<input type="checkbox"/>	1630200	Q'Strain Lap Belt	Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" Belt Length.	\$ 295
<input type="checkbox"/>	1630222	Wheelchair Backpack	Black with PDG Logo. 17"x14" size.	\$ 165
<input type="checkbox"/>	1630170	Adhesive Back Hook Strip, 2" x 30"	For Cushion mounting.	\$ 10

On-Chair Seating & Positioning

	Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. pdgmobility.com/seating
---	---

Order Notes

Please send your completed order form by email to info@pdgmobility.com.