



BENTLEY LT

Order Form - USA

Prices in U.S. Dollars. January 1, 2024

Customer and Order Information

Purchase Order *	Tag	v.0401
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in
		*required

Base Chair

Bentley LT 20° Manual Tilt Wheelchair

- | | | | | |
|-------------------------------------|--|--|-------|----------|
| <input checked="" type="checkbox"/> | 2401000 | Bentley 20° Manual Tilt Wheelchair | E1161 | \$ 3,150 |
| <input checked="" type="checkbox"/> | 2401009 | Lightweight Aluminum Frame Upgrade Aluminum Frame is not recommended for agitated users. | K0108 | \$ 300 |
| <input type="checkbox"/> | i Choose the standard Bentley with a steel frame for agitated or heavier users. See separate order form. | | | |

Weight Capacity

- | | | | | |
|---|---------|--------------------------------------|--------------------------|----------------------------|
| <input checked="" type="checkbox"/> | 2402000 | 250 lbs Chair Weight Capacity | | |
| User Weight Range Required for setup of tilt uplift force. | | | | |
| <input type="checkbox"/> | 2402430 | Up to 75 lbs | <input type="checkbox"/> | 2402431 75 lbs to 100 lbs |
| <input type="checkbox"/> | 2402432 | 100 lbs to 115 lbs | <input type="checkbox"/> | 2402433 115 lbs to 150 lbs |
| <input type="checkbox"/> | 2402434 | 150 lbs to 200 lbs | <input type="checkbox"/> | 2402435 200 lbs to 250 lbs |

Frame Options

- | | | |
|---|---------|---------------------------------------|
| <input type="checkbox"/> | 2402300 | Transit Tie Down System (TTDS) |
| TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. | | |

Seat Width

Overall chair width is +9.5" with 12/16" Rear Wheels and +10.5" with 20/22/24" Rear Wheels.

Standard Seat Widths

- | | | | | | | | | | |
|--------------------------|---------|----------------|--------------------------|---------|----------------|--------------------------|---------|----------------|--------|
| <input type="checkbox"/> | 2403014 | 14" Seat Width | <input type="checkbox"/> | 2403015 | 15" Seat Width | <input type="checkbox"/> | 2403016 | 16" Seat Width | |
| <input type="checkbox"/> | 2403017 | 17" Seat Width | <input type="checkbox"/> | 2403018 | 18" Seat Width | <input type="checkbox"/> | 2403019 | 19" Seat Width | |
| <input type="checkbox"/> | 2403020 | 20" Seat Width | <input type="checkbox"/> | 2403021 | 21" Seat Width | <input type="checkbox"/> | 2403022 | 22" Seat Width | |
| | | | | | | | | E2201 | \$ 375 |

Seat Depth

Standard Seat Depths (Adjustable)

- | | | | | | | | |
|--------------------------|---------|----------------|--------------------------|---------|----------------|-------|--------|
| <input type="checkbox"/> | 2404016 | 16" Seat Depth | <input type="checkbox"/> | 2404017 | 17" Seat Depth | | |
| <input type="checkbox"/> | 2404018 | 18" Seat Depth | <input type="checkbox"/> | 2404019 | 19" Seat Depth | | |
| <input type="checkbox"/> | 2404020 | 20" Seat Depth | | | | | |
| | | | | | | E2203 | \$ 500 |

Custom Seat Depths

- | | | | | |
|--------------------------|---------|----------------|-------|--------|
| <input type="checkbox"/> | 2404021 | 21" Seat Depth | E2203 | \$ 525 |
| <input type="checkbox"/> | 2404022 | 22" Seat Depth | E2204 | \$ 800 |

Frame Color

PDG Premium Colors

View color swatches on PDG's website. pdgmobility.com/colors

- | | | | | | | | | |
|--------------------------|---------|---|--------------------------|---------|--------------------|--------------------------|---------|------------------|
| <input type="checkbox"/> | 2405001 | Black Shimmer | <input type="checkbox"/> | 2405015 | Sky Blue | <input type="checkbox"/> | 2405016 | Poppy Red |
| <input type="checkbox"/> | 2405003 | Fine Texture Black | <input type="checkbox"/> | 2405010 | Lollypop Blue | <input type="checkbox"/> | 2405011 | Lollypop Red |
| <input type="checkbox"/> | 2405017 | Arctic Grey | <input type="checkbox"/> | 2405018 | Tropical Teal | <input type="checkbox"/> | 2405019 | Tangerine Orange |
| <input type="checkbox"/> | 2405020 | Olive Green | <input type="checkbox"/> | 2405021 | Lush Lilac | <input type="checkbox"/> | 2405012 | Bubble Gum Pink |
| <input type="checkbox"/> | 2405022 | Leaf Green | <input type="checkbox"/> | 2405023 | Boysenberry Purple | <input type="checkbox"/> | 2405014 | Daffodil Yellow |
| <input type="checkbox"/> | 2405200 | Custom Color Contact PDG for specifications. | | | | | | \$ 300 |

Armrests

2406000 **No Armrests**

Height Adjustable Armrests (Single Post) E0973 \$ 345

2406111 8"-11" Range 2406114 10"-14" Range 2406116 12"-16" Range

2406100 **Fixed Height Armrests (12")** Full Length Armpad only.

Armrest Height If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.
 Seat Pan to top of Armrest: _____ Does not include height of Armpad

Armrest Position Indicate preference for Armrest position. If left blank, the default setting will be used.
 In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

2406401 Rear 2406402 Mid 2406403 Front

Armpads

Regular Armpads

2407100 Desk Length - 2" x 10" 2407101 Full Length - 2" x 14"

Waterfall Armpads \$ 115

2407105 Desk Length - 2" x 10" 2407106 Full Length - 2" x 14"

Wide Armpads May interfere with wheels in some tilt configurations. \$ 200

2407110 Desk Length - 3" x 10" 2407111 Full Length - 3" x 12"

Gel Armpads May interfere with wheels in some tilt configurations. \$ 135

2407120 Full Length - 2" x 12" 2407121 Full Length - 3.5" x 12"

Contour Arm Trough with Hand Pad May interfere with wheels in some tilt configurations. E2209 \$ 125

2407130 Right (4.5" x 18") 2407131 Left (4.5" x 18")

Armrest Options

Armrest Wideners Move Armrests outward by up to 1" on each side. Seat size does not change. May cause tilt interference.

2407210 1" Total Additional Width \$ 95

2407211 2" Total Additional Width \$ 195

2407310 **Clothing Guards** Not compatible with 8"-11" Height range. \$ 110

Armrest Extensions Not compatible with Desk Length Armpads. \$ 75

2407400 3" Extension - Forward 2407401 3" Extension - Rearward

Seat Pan

2408000 **No Seat Pan**

2408001 **Solid Aluminum Seat Pan** E2231 \$ 165

Seat Height

Don't forget **Select desired Seat-to-Floor Height** STF Height is adjustable within the range of the chosen Subframe Plate (Low/High).

Casters	Low						High						
	RW	13"	14"	15"	16"	17"	RW	15"	16"	17"	18"	19"	20"
6" Casters	12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
8" Casters	RW	13"	14"	15"	16"	17"	RW	15"	16"	17"	18"	19"	20"
	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Casters

Urethane Casters 8" Casters are compatible only with 16" and 24" Rear Wheels.

<input type="checkbox"/> 2411061	6" x 1.25" Caster size	<input type="checkbox"/> 2411081	8" x 1.25" Caster size		
<input type="checkbox"/> 2411063	6" x 1.50" Caster size				\$ 85
<input type="checkbox"/> 2411084	8" x 2.00" Caster size				\$ 135

Pneumatic Casters 8" Casters are compatible only with 16" and 24" Rear Wheels.

<input type="checkbox"/> 2411161	6" x 1.25" Caster size			E2214 & E2215	\$ 110
<input type="checkbox"/> 2411181	8" x 1.25" Caster size			E2214 & E2215	\$ 120
<input type="checkbox"/> 2411184	8" x 2.00" Caster size			E2214 & E2215	\$ 150

<input type="checkbox"/> 2411260	Frog Legs 6" x 1.50" Aluminum Soft-Roll Casters				\$ 275
----------------------------------	--	--	--	--	--------

<input type="checkbox"/> 2411400	Frog Legs Suspension Forks	Only compatible with 16 and 24" Rear Wheels. Requires 6" Soft-Roll Caster.		E1015	\$ 595
----------------------------------	-----------------------------------	--	--	-------	--------

Rear Wheels

Composite Wheels

<input type="checkbox"/> 2412012	12" Wheel size	<input type="checkbox"/> 2412016	16" Wheel size	<input type="checkbox"/> 2412020	20" Wheel size
<input type="checkbox"/> 2412022	22" Wheel size	<input type="checkbox"/> 2412024	24" Wheel size		

Spoke Wheels

<input type="checkbox"/> 2412120	20" Wheel size	<input type="checkbox"/> 2412122	22" Wheel size	<input type="checkbox"/> 2412124	24" Wheel size	\$ 75
----------------------------------	----------------	----------------------------------	----------------	----------------------------------	----------------	-------



Rear Wheel Position Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.

In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

<input type="checkbox"/> 2412401	Rear	<input type="checkbox"/> 2412402	Mid	<input type="checkbox"/> 2412403	Front
----------------------------------	------	----------------------------------	-----	----------------------------------	-------

Axle Type

<input type="checkbox"/> 2413000	Quick Release Rear Axles
<input type="checkbox"/> 2413010	Fixed Rear Axles

Rear Tires

<input type="checkbox"/> 2414000	Urethane Tires		
<input type="checkbox"/> 2414010	Pneumatic Tires	May result in a finished Seat Height 1/2" greater than stated.	E2211 & E2212 \$ 60
<input type="checkbox"/> 2414020	Pneumatic Tires with Flat-Free Inserts	May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels.	E2211 & E2213 \$ 95

Handrims

Handrims are only available with 20", 22", and 24" Rear Wheels.

<input type="checkbox"/> 2415000	No Handrims		
<input type="checkbox"/> 2415010	Standard Aluminum Handrims		
<input type="checkbox"/> 2415020	Plastic Coated Handrims	Not compatible with Drum Brakes.	\$ 150
<input type="checkbox"/> 2415030	Projection Handrims	Composite Wheels only. 10 Projections, 30°, 1-1.8" Length.	\$ 495

Wheel Locks & Brakes

User Wheel Locks

<input type="checkbox"/> 2416010	Push to Lock	<input type="checkbox"/> 2416020	Pull to Lock		
<input type="checkbox"/> 2416030	Attendant Wheel Locks (Push to Lock)				\$ 85
<input type="checkbox"/> 2416040	User & Attendant Wheel Locks	Not compatible with Drum Brakes.			\$ 125
<input type="checkbox"/> 2416300	Drum Brakes	Available with 12" Composite or 20" Spoke Rear Wheels. Urethane Tires only.		E2228	\$ 795
	User-activated Unilateral Wheel Lock	Pull to Lock.			\$ 795
<input type="checkbox"/> 2416310	Left Activation	<input type="checkbox"/> 2416311	Right Activation		



Wheel Lock Options			
<input type="checkbox"/>	2416100	Wheel Lock Extensions Not compatible with Attendant Wheel Lock.	E0961 \$ 80
<input type="checkbox"/>	2416200	Hill Holders Requires Pneumatic Rear Tires with Flat-Free Inserts.	E0974 \$ 450

Anti-Tippers			
<input checked="" type="checkbox"/>	2417010	Rear Anti-Tippers Required Safety Option.	E0971 \$ 125

Back Posts			
Back Post Style			
<input checked="" type="checkbox"/>	2418010	Stroller Handle Back	
Back Post Height Measured from the Seat Pan to the Handle.			
<input type="checkbox"/>	2418425	Short (22")	<input type="checkbox"/> 2418431 Tall (25")
Back Post Folding Fold-down Back is not compatible with Reclining Backrest.			
<input checked="" type="checkbox"/>	2419220	Fold-down	
Seat to Back Angle			
<input type="checkbox"/>	2419090	Set at 90°	<input type="checkbox"/> 2419098 Set at 98° <input type="checkbox"/> 2419106 Set at 106°
<input type="checkbox"/>	2419114	Set at 114°	<input type="checkbox"/> 2419122 Set at 122°

Back Options			
Angle/Height Adjustable Stroller Handle Not compatible with Headrest Mounting Plate. \$ 150			
<input type="checkbox"/>	2419010	Button on Outside	<input type="checkbox"/> 2419011 Button on Inside
Endomorph Back Backrest will be narrower than Seat by total of 2x the Offset selected below. Minimum width of Backrest is 12". \$ 395			
<input type="checkbox"/>	2418102	1" Offset per side	<input type="checkbox"/> 2418103 1.5" Offset per side <input type="checkbox"/> 2418104 2" Offset per side
<input type="checkbox"/>	2419210	10" Backrest Extension Not compatible with Angle/Height Adjustable Stroller Handle.	E0966 \$ 150
<input type="checkbox"/>	2419212	Headrest Mounting Plate	\$ 95
Not compatible with Backrest Extension, Angle/Height Adjustable Stroller Handle. Not required if mounting to a solid Back.			

Back Upholstery									
<input type="checkbox"/>	2420000	No Back Upholstery							
Standard Upholstery E2611 / E2612 \$ 195									
<input type="checkbox"/>	2420418	Short (18" Height)	<input type="checkbox"/> 2420421 Regular (21" Height) <input type="checkbox"/> 2420424 Tall (24" Height)						
<input type="checkbox"/>	2420450	Tension Adjustable Upholstery Height of Upholstery is 16".	E2611 / E2612 \$ 295						
Position of Upholstery Select the desired position of the top of the Upholstery relative to the Seat Pan.									
		22" Back Post							
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
		25" Back Post							
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tall							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Tilt Activation			
Attendant Tilt only Select trigger mounting location below.			
<input type="checkbox"/>	2423011	Right Back Post	<input type="checkbox"/> 2423012 Left Back Post
User Tilt only Full Length Arm Pads only. Select trigger mounting location below. \$ 100			
<input type="checkbox"/>	2423021	Right Armrest	<input type="checkbox"/> 2423022 Left Armrest
Attendant & User Tilt Full Length Arm Pads only. Select trigger mounting locations below. \$ 195			
<input type="checkbox"/>	2423031	Right Back Post	<input type="checkbox"/> 2423032 Left Back Post
<input type="checkbox"/>	2423041	Right Armrest	<input type="checkbox"/> 2423042 Left Armrest

See next page for more options >



<input type="checkbox"/>	2423110	Tilt Trigger Extension				\$	50
		Power Tilt	250 lbs Weight Capacity. 20° tilt range. Includes 2 batteries for off-chair charging. pdgmobility.com/power-tilt				E1002
			Minimum Seat Width is 15". Power Tilt is compatible with either O2 Holder or IV Pole on the same chair, not both.				
<input type="checkbox"/>	2423400	Standard Power Tilt				\$	2,295
<input type="checkbox"/>	2423410	Power Tilt with Manual Tilt Override	Allows for instant manual release of tilt in case of power tilt failure.			\$	2,795
		Toggle Switch	<input type="checkbox"/> 2423421 Installed on Right side	<input type="checkbox"/> 2423422 Installed on Left side			
		Button Switch	<input type="checkbox"/> 2423431 Installed on Right side	<input type="checkbox"/> 2423432 Installed on Left side			
		Stealth Egg Switch	<input type="checkbox"/> 2423440 Not installed				
		Optional for Egg Switch:	<input type="checkbox"/> 2423445 Swing Away Mounting Hardware	Not compatible with all Headrest makes and models.		E1028	\$ 625
			Specify Headrest to be used: _____				

Recline							
<input type="checkbox"/>	2423310	Reclining Backrest (Manual)				\$	895
			30° maximum recline. Not compatible with Backrest Extension, Fold-down Back, Tension Adjustable Upholstery, Dynamic Rocker Back. Requires longer Frame. Maximum Seat Depth is 19". Finished Back Post Height is 3" taller than selected Back Post Height.				E1225

Front Rigging							
<input type="checkbox"/>	2425000	No Front Rigging		<input type="checkbox"/> 2425001 Sector Blocks only			
		Easy-Out™ Front Rigging	Effortless insertion and removal. No Hanger Bracket remains on frame. pdgmobility.com/easy-out				
<input type="checkbox"/>	2425410	Fixed 70°				\$	95
<input type="checkbox"/>	2425420	Universal Angle Adjustable (0°-110°)	Set at angle: _____ (use increments of 10°)			\$	295
			Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings.				K0108
		Swing Away Front Rigging	Legacy Design. Not compatible with Easy-Out.				
<input type="checkbox"/>	2425130	Swing Away 60°	<input type="checkbox"/> 2425140 Swing Away 70°				
<input type="checkbox"/>	2425200	Swing Away Elevating Legrests	Includes Calf Pads.			\$	415
<input type="checkbox"/>	2425300	1" Offset Hangers for Swing Away				\$	75
<input type="checkbox"/>	2425110	Standard Plug In 70° Front Rigging	Legacy Design. Not compatible with Easy-Out.				

Footplates							
<input type="checkbox"/>	2426010	Composite Footplates					
<input type="checkbox"/>	2426020	Aluminum Angle Adjustable Footplates				\$	200
							K0040
<input type="checkbox"/>	2426030	One Piece Fixed Footplate	Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.			\$	350
<input type="checkbox"/>	2426040	One Piece Flip-up Footplate	Requires selection of Composite or Aluminum Angle Adjustable.			\$	395
<input type="checkbox"/>	2426050	Interlocking Angle Adjustable Footplates	Available for Seat Widths up to 21" only.			\$	200
							K0040

Footrest/Legrest Accessories							
		Calf Pads	Not compatible with Gel Padding.				
<input type="checkbox"/>	2427012	For Easy-Out™ Front Rigging	Includes reinforced premium Calf Pad mount.			\$	200
<input type="checkbox"/>	2427010	For Swing Away Front Rigging				\$	135
<input type="checkbox"/>	2427020	Gel Padding	Not compatible with Calf Pad or Easy-Out Front Rigging.			\$	125
		Gel Knee Pads					
<input type="checkbox"/>	2427032	For Easy-Out™ Angle Adjustable Front Rigging				\$	200
<input type="checkbox"/>	2427030	For Swing Away ELR Front Rigging				\$	100
<input type="checkbox"/>	2427040	Length Adjustable Heel Loops				\$	50
							E0951
<input type="checkbox"/>	2427050	Bodypoint Padded Calf Strap				\$	195
							K0038

Front Rigging Extension Tubes							
		Regular Extension Tubes					
<input type="checkbox"/>	2428010	Short	14"-18" for Easy-Out Front Rigging	13"-17" for Swing Away and Standard Plug In Front Rigging			
<input type="checkbox"/>	2428020	Long	18"-24" for Easy-Out Front Rigging	15"-19" for Swing Away and Standard Plug In Front Rigging		\$	75

See next page for more options >



BENTLEY LT

Order Form - USA

Prices in U.S. Dollars. January 1, 2024



Lower Leg Length If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.

Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.

Miscellaneous Options

One-Arm Drive Available with 22" and 24" Composite Rear Wheels. Minimum 17" STF and 18" STF, respectively. Limits Tilt. Not compatible with User & Attendant Wheel Locks, Drum Brakes, Quick Release Rear Axles, Pneumatic Tires.		E0958	\$ 3,500
<input type="checkbox"/> 2430010	Installed on Right side	<input type="checkbox"/> 2430011	Installed on Left side
<input type="checkbox"/> 2430012	Large Hand Setup 24" Rear Wheels with 22" Driving Handrim.		
Spoke Guards Available for 20"-24" Composite Rear Wheels. Made from durable and lightweight aluminum.		K0065	\$ 195
<input type="checkbox"/> 2430020	Black	<input type="checkbox"/> 2430021	Frame Color
<input type="checkbox"/> 2430050	Protective Lower Padding		\$ 795
Stealth Flip Down Abductor Requires Solid Seat Pan.		E0957 & E1028	\$ 650
<input type="checkbox"/> 2430060	Small Pad 3"x4"x3"x2" size	<input type="checkbox"/> 2430062	Large Pad 5"x5"x4"x3" size
<input type="checkbox"/> 2430080	Gel Hip Guides Reduces available Seat Width by 2".	E0956 & E1028	\$ 495
<input type="checkbox"/> 2430110	IV Pole Not compatible with Attendant Wheel Locks. Requires frame modification and cannot be added to chair later.	E0776	\$ 295
<input type="checkbox"/> 2430120	O2 Holder Not compatible with Attendant Wheel Locks. Requires frame modification and cannot be added to chair later.	E2208	\$ 395
<input type="checkbox"/> 2430200	Q'Straint Lap Belt Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" Belt Length.		\$ 245
<input type="checkbox"/> 2430222	Wheelchair Backpack Black with PDG Logo. 17"x14" size.		\$ 145
<input type="checkbox"/> 2430170	Adhesive Back Hook Strip, 2" x 30" For Cushion mounting.		\$ 10

On-Chair Seating & Positioning

Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. pdgmobility.com/seating

Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions. Use of HCPCS codes does not guarantee coverage or payment.

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.