



# STELLAR GLT

Order Form - USA

Prices in U.S. Dollars. January 1, 2024

## Customer and Order Information

Purchase Order *	Tag	v.0401
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in
		*required

## Base Chair

**Stellar GLT 45° Manual Tilt Wheelchair**

<input checked="" type="checkbox"/> 2501000	Stellar GL 45° Manual Tilt Wheelchair	E1161	\$ 3,150
<input checked="" type="checkbox"/> 2501009	Lightweight Aluminum Frame Upgrade Aluminum Frame is not recommended for agitated users.	K0108	\$ 300
<b>i</b> Choose the standard <b>Stellar GL</b> with a steel frame for agitated or heavier users. See separate order form.			

## Weight Capacity

2502000 **250 lbs Chair Weight Capacity**

**User Weight Range** Required for setup of tilt uplift force.

<input type="checkbox"/> 1902430	Up to 75 lbs	<input type="checkbox"/> 2502431	75 lbs to 100 lbs	<input type="checkbox"/> 2502432	100 lbs to 115 lbs
<input type="checkbox"/> 1902433	115 lbs to 150 lbs	<input type="checkbox"/> 2502434	150 lbs to 200 lbs	<input type="checkbox"/> 2502435	200 lbs to 250 lbs

## Frame Options

2502300 **Transit Tie Down System (TTDS)**

TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details.

## Seat Width

Overall chair width is +9.5" with 12/16" Rear Wheels and +10.5" with 20/22/24" Rear Wheels.

**Standard Seat Widths**

<input type="checkbox"/> 2503014	14" Seat Width	<input type="checkbox"/> 2503015	15" Seat Width	<input type="checkbox"/> 2503016	16" Seat Width
<input type="checkbox"/> 2503017	17" Seat Width	<input type="checkbox"/> 2503018	18" Seat Width	<input type="checkbox"/> 2503019	19" Seat Width
<input type="checkbox"/> 2503020	20" Seat Width	<input type="checkbox"/> 2503021	21" Seat Width	<input type="checkbox"/> 2503022	22" Seat Width
					E2201 \$ 375

## Seat Depth

**Standard Seat Depths (Adjustable)**

<input type="checkbox"/> 2504016	16" Seat Depth	<input type="checkbox"/> 1904017	17" Seat Depth
<input type="checkbox"/> 1904018	18" Seat Depth	<input type="checkbox"/> 2504019	19" Seat Depth
<input type="checkbox"/> 2504020	20" Seat Depth	E2203 \$ 500	
<b>Custom Seat Depths</b>			
<input type="checkbox"/> 2504021	21" Seat Depth	E2203 \$ 525	
<input type="checkbox"/> 2504022	22" Seat Depth	E2204 \$ 800	

## Frame Color

**PDG Premium Colors** View color swatches on PDG's website. [pdgmobility.com/colors](http://pdgmobility.com/colors)

<input type="checkbox"/> 2505001	Black Shimmer	<input type="checkbox"/> 2505015	Sky Blue	<input type="checkbox"/> 2505016	Poppy Red
<input type="checkbox"/> 2505003	Fine Texture Black	<input type="checkbox"/> 2505010	Lollypop Blue	<input type="checkbox"/> 2505011	Lollypop Red
<input type="checkbox"/> 2505017	Arctic Grey	<input type="checkbox"/> 2505018	Tropical Teal	<input type="checkbox"/> 2505019	Tangerine Orange
<input type="checkbox"/> 2505020	Olive Green	<input type="checkbox"/> 2505021	Lush Lilac	<input type="checkbox"/> 2505012	Bubble Gum Pink
<input type="checkbox"/> 2505022	Leaf Green	<input type="checkbox"/> 2505023	Boysenberry Purple	<input type="checkbox"/> 2505014	Daffodil Yellow
<input type="checkbox"/> 2505200	<b>Custom Color</b> Contact PDG for specifications.				\$ 300

### Armrests

2506000 **No Armrests**

**Height Adjustable Armrests (Single Post)** E0973 \$ 345

2506111 8"-11" Range       2506114 10"-14" Range       2506116 12"-16" Range

2506100 **Fixed Height Armrests (12")** Full Length Armpad only.

**Armrest Height** If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.  
 Seat Pan to top of Armrest: \_\_\_\_\_ Does not include height of Armpad

**Armrest Position** Indicate preference for Armrest position. If left blank, the default setting will be used.  
 In the interest of maintaining chair stably and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

2506401 Rear       2506402 Mid       2506403 Front

### Armpads

**Regular Armpads**

2507100 Desk Length - 2" x 10"       2507101 Full Length - 2" x 14"

**Waterfall Armpads** \$ 115

2507105 Desk Length - 2" x 10"       2507106 Full Length - 2" x 14"

**Wide Armpads** May interfere with wheels in some tilt configurations. \$ 200

2507110 Desk Length - 3" x 10"       2507111 Full Length - 3" x 12"

**Gel Armpads** May interfere with wheels in some tilt configurations. \$ 135

2507120 Full Length - 2" x 12"       2507121 Full Length - 3.5" x 12"

**Contour Arm Trough with Hand Pad** May interfere with wheels in some tilt configurations. E2209 \$ 125

2507130 Right (4.5" x 18")       2507131 Left (4.5" x 18")

### Armrest Options

**Armrest Wideners** Move Armrests outward by up to 1" on each side. Seat size does not change. May cause tilt interference.

2507210 1" Total Additional Width \$ 95

2507211 2" Total Additional Width \$ 195

2507310 **Clothing Guards** Not compatible with 8"-11" Height range. \$ 110

**Armrest Extensions** Not compatible with Desk Length Armpads. \$ 75

2507400 3" Extension - Forward       2507401 3" Extension - Rearward

### Seat Pan

2508000 **No Seat Pan**

2508001 **Solid Aluminum Seat Pan** E2231 \$ 165

### Seat Height

**Don't forget** **Select desired Seat-to-Floor Height** STF Height is adjustable within the range of the chosen Subframe Plate (Low/High).

6" Casters	Low							High						
	RW	13"	14"	15"	16"	17"	18"	RW	15"	16"	17"	18"	19"	20"
12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>


  

8" Casters	Low							High						
	RW	13"	14"	15"	16"	17"	18"	RW	15"	16"	17"	18"	19"	20"
16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

### Casters

<b>Urethane Casters</b> 8" Casters are compatible only with 16" and 24" Rear Wheels.		
<input type="checkbox"/> 2511061 6" x 1.25" Caster size	<input type="checkbox"/> 2511081 8" x 1.25" Caster size	
<input type="checkbox"/> 2511063 6" x 1.50" Caster size		\$ 85
<input type="checkbox"/> 2511084 8" x 2.00" Caster size		\$ 135
<b>Pneumatic Casters</b> 8" Casters are compatible only with 16" and 24" Rear Wheels.		
<input type="checkbox"/> 2511161 6" x 1.25" Caster size		E2214 & E2215 \$ 110
<input type="checkbox"/> 2511181 8" x 1.25" Caster size		E2214 & E2215 \$ 120
<input type="checkbox"/> 2511184 8" x 2.00" Caster size		E2214 & E2215 \$ 150
<input type="checkbox"/> 2511260 <b>Frog Legs 6" x 1.50" Aluminum Soft-Roll Casters</b>		\$ 275
<input type="checkbox"/> 2511400 <b>Frog Legs Suspension Forks</b> Only compatible with 16 and 24" Rear Wheels. Requires 6" Soft-Roll Caster.		E1015 \$ 595

### Rear Wheels

<b>Composite Wheels</b>		
<input type="checkbox"/> 2512012 12" Wheel size	<input type="checkbox"/> 2512016 16" Wheel size	<input type="checkbox"/> 2512020 20" Wheel size
<input type="checkbox"/> 2512022 22" Wheel size	<input type="checkbox"/> 2512024 24" Wheel size	
<b>Spoke Wheels</b> \$ 75		
<input type="checkbox"/> 2512120 20" Wheel size	<input type="checkbox"/> 2512122 22" Wheel size	<input type="checkbox"/> 2512124 24" Wheel size
 <b>Rear Wheel Position</b> Indicate Rear Wheel positioning preference. If left blank, the default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.		
<input type="checkbox"/> 2512401 Rear	<input type="checkbox"/> 2512402 Mid	<input type="checkbox"/> 2512403 Front

### Axle Type

<input type="checkbox"/> 2513000 <b>Quick Release Rear Axles</b>
<input type="checkbox"/> 2513010 <b>Fixed Rear Axles</b>

### Rear Tires

<input type="checkbox"/> 2514000 <b>Urethane Tires</b>		
<input type="checkbox"/> 2514010 <b>Pneumatic Tires</b> May result in a finished Seat Height 1/2" greater than stated.	E2211 & E2212	\$ 60
<input type="checkbox"/> 2514020 <b>Pneumatic Tires with Flat-Free Inserts</b> May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels.	E2211 & E2213	\$ 95

### Handrims

Handrims are only available with 20", 22", and 24" Rear Wheels.

<input type="checkbox"/> 2515000 <b>No Handrims</b>		
<input type="checkbox"/> 2515010 <b>Standard Aluminum Handrims</b>		
<input type="checkbox"/> 2515020 <b>Plastic Coated Handrims</b> Not compatible with Drum Brakes.		\$ 135
<input type="checkbox"/> 2515030 <b>Projection Handrims</b> Composite Wheels only. 10 Projections, 30°, 1-1.8" Length.		\$ 495

### Wheel Locks & Brakes

<b>User Wheel Locks</b>		
<input type="checkbox"/> 2516010 Push to Lock	<input type="checkbox"/> 2516020 Pull to Lock	
<input type="checkbox"/> 2516030 <b>Attendant Wheel Locks (Push to Lock)</b>		\$ 85
<input type="checkbox"/> 2516040 <b>User &amp; Attendant Wheel Locks</b> Not compatible with Drum Brakes.		\$ 125
<input type="checkbox"/> 2516300 <b>Drum Brakes</b> Available with 12" Composite or 20" Spoke Rear Wheels. Urethane Tires only.	E2228	\$ 795
<b>User-activated Unilateral Wheel Lock</b> Pull to Lock. \$ 795		
<input type="checkbox"/> 2516310 Left Activation	<input type="checkbox"/> 2516311 Right Activation	



Wheel Lock Options			
<input type="checkbox"/>	2516100	<b>Wheel Lock Extensions</b> Not compatible with Attendant Wheel Lock.	E0961 \$ 80
<input type="checkbox"/>	2516200	<b>Hill Holders</b> Requires Pneumatic Rear Tires with Flat-Free Inserts.	E0974 \$ 450

Anti-Tippers			
<input checked="" type="checkbox"/>	2517010	<b>Rear Anti-Tippers</b> Required Safety Option.	E0971 \$ 125

Back Posts			
<b>Back Post Style</b>			
<input checked="" type="checkbox"/>	2518010	Stroller Handle Back	
<b>Back Post Height</b> Measured from the Seat Pan to the Handle.			
<input type="checkbox"/>	2518425	Short (22")	<input type="checkbox"/> 2518431 Tall (25")
<b>Back Post Folding</b> Fold-down Back is not compatible with Reclining Backrest.			
<input checked="" type="checkbox"/>	2519220	Fold-down	
<b>Seat to Back Angle</b>			
<input type="checkbox"/>	2519090	Set at 90°	<input type="checkbox"/> 2519098 Set at 98° <input type="checkbox"/> 2519106 Set at 106°
<input type="checkbox"/>	2519114	Set at 114°	<input type="checkbox"/> 2519122 Set at 122°

Back Options			
<b>Angle/Height Adjustable Stroller Handle</b> Not compatible with Headrest Mounting Plate. \$ 150			
<input type="checkbox"/>	2519010	Button on Outside	<input type="checkbox"/> 2519011 Button on Inside
<b>Endomorph Back</b> Backrest will be narrower than Seat by total of 2x the Offset selected below. Minimum width of Backrest is 12". \$ 395			
<input type="checkbox"/>	2518102	1" Offset per side	<input type="checkbox"/> 2518103 1.5" Offset per side <input type="checkbox"/> 2518104 2" Offset per side
<input type="checkbox"/>	2519210	<b>10" Backrest Extension</b> Not compatible with Angle/Height Adjustable Stroller Handle.	E0966 \$ 150
<input type="checkbox"/>	2519212	<b>Headrest Mounting Plate</b>	\$ 95
Not compatible with Backrest Extension, Angle/Height Adjustable Stroller Handle. Not required if mounting to a solid Back.			

Back Upholstery																						
<input type="checkbox"/>	2520000	<b>No Back Upholstery</b>																				
<b>Standard Upholstery</b> E2611 / E2612 \$ 195																						
<input type="checkbox"/>	2520418	Short (18" Height)	<input type="checkbox"/> 2520421 Regular (21" Height) <input type="checkbox"/> 2520424 Tall (24" Height)																			
<input type="checkbox"/>	2520450	<b>Tension Adjustable Upholstery</b> Height of Upholstery is 16".	E2611 / E2612 \$ 295																			
<b>Position of Upholstery</b> Select the desired position of the top of the Upholstery relative to the Seat Pan.																						
		22" Back Post									25" Back Post											
	<b>Style</b>	17"	18"	19"	20"	21"	22"	23"	24"	25"		<b>Style</b>	17"	18"	19"	20"	21"	22"	23"	24"	25"	
	Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>					Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>					Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>					Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Tilt Activation			
<b>Attendant Tilt only</b> Select trigger mounting location below.			
<input type="checkbox"/>	2523011	Right Back Post	<input type="checkbox"/> 2523012 Left Back Post
<b>User Tilt only</b> Not available on widths greater than 22". Full Length Arm Pads only. Select trigger mounting location below. \$ 100			
<input type="checkbox"/>	2523021	Right Armrest	<input type="checkbox"/> 2523022 Left Armrest
<b>Attendant &amp; User Tilt</b> Not available on widths greater than 22". Full Length Arm Pads only. Select trigger mounting locations below. \$ 195			
<input type="checkbox"/>	2523031	Right Back Post	<input type="checkbox"/> 2523032 Left Back Post
<input type="checkbox"/>	2523041	Right Armrest	<input type="checkbox"/> 2523042 Left Armrest

See next page for more options >



<input type="checkbox"/>	2523110	<b>Tilt Trigger Extension</b>				\$	50	
<b>Tilt Limit</b> Optional. If no tilt limit is selected, the full tilt range will be available.								
<input type="checkbox"/>	2523220	Set at 20°	<input type="checkbox"/>	2523230	Set at 30°	<input type="checkbox"/>	2523240	Set at 40°
<b>Power Tilt</b> 250 lbs Weight Capacity. 45° tilt range. Includes 2 batteries for off-chair charging. <a href="http://pdgmobility.com/power-tilt">pdgmobility.com/power-tilt</a> <span style="float: right;">E1002</span> Minimum Seat Width is 16". Power Tilt is compatible with either O2 Holder or IV Pole on the same chair, not both.								
<input type="checkbox"/>	2523400	Standard Power Tilt				\$	2,295	
<input type="checkbox"/>	2523410	Power Tilt with Manual Tilt Override	Allows for instant manual release of tilt in case of power tilt failure.			\$	2,795	
Select Control Switch		Toggle Switch	<input type="checkbox"/>	2523421	Installed on Right side	<input type="checkbox"/>	2523422	Installed on Left side
		Button Switch	<input type="checkbox"/>	2523431	Installed on Right side	<input type="checkbox"/>	2523432	Installed on Left side
		Stealth Egg Switch	<input type="checkbox"/>	2523440	Not installed			
		Optional for Egg Switch:	<input type="checkbox"/>	2523445	Swing Away Mounting Hardware	Not compatible with all Headrest makes and models.	E1028	\$ 625
Specify Headrest to be used: _____								

Recline							
<input type="checkbox"/>	2523310	<b>Reclining Backrest (Manual)</b>				\$	895
30° maximum recline. Not compatible with Backrest Extension, Fold-down Back, Tension Adjustable Upholstery. Finished Back Post Height is 3" taller than selected Back Post Height. Limits Tilt by up to 20°.							

Front Rigging							
<input type="checkbox"/>	2525000	<b>No Front Rigging</b>	<input type="checkbox"/>	2525001	<b>Sector Blocks only</b>		
<b>Easy-Out™ Front Rigging</b> Effortless insertion and removal. No Hanger Bracket remains on frame. <a href="http://pdgmobility.com/easy-out">pdgmobility.com/easy-out</a>							
<input type="checkbox"/>	2525410	Fixed 70°				\$	95
<input type="checkbox"/>	2525420	Universal Angle Adjustable (0°-110°)	Set at angle: _____ (use increments of 10°)			K0108	\$ 295
Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings.							
<b>Swing Away Front Rigging</b> Legacy Design. Not compatible with Easy-Out.							
<input type="checkbox"/>	2525130	Swing Away 60°	<input type="checkbox"/>	2525140	Swing Away 70°		
<input type="checkbox"/>	2525200	Swing Away Elevating Legrests	Includes Calf Pads.			E0990	\$ 415
<input type="checkbox"/>	2525300	1" Offset Hangers for Swing Away					\$ 75
<input type="checkbox"/>	2525110	<b>Standard Plug In 70° Front Rigging</b>	Legacy Design. Not compatible with Easy-Out.				

Footplates							
<input type="checkbox"/>	2526010	<b>Composite Footplates</b>					
<input type="checkbox"/>	2526020	<b>Aluminum Angle Adjustable Footplates</b>				K0040	\$ 200
<input type="checkbox"/>	2526030	<b>One Piece Fixed Footplate</b>	Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.				\$ 350
<input type="checkbox"/>	2526040	<b>One Piece Flip-up Footplate</b>	Requires selection of Composite or Aluminum Angle Adjustable.				\$ 395
<input type="checkbox"/>	2526050	<b>Interlocking Angle Adjustable Footplates</b>	Available for Seat Widths up to 21" only.			K0040	\$ 200

Footrest/Legrest Accessories							
<b>Calf Pads</b> Not compatible with Gel Padding.							
<input type="checkbox"/>	2527012	For Easy-Out™ Front Rigging	Includes reinforced premium Calf Pad mount.			\$	200
<input type="checkbox"/>	2527010	For Swing Away Front Rigging				\$	135
<input type="checkbox"/>	2527020	<b>Gel Padding</b>	Not compatible with Calf Pad or Easy-Out Front Rigging.			\$	125
<b>Gel Knee Pads</b>							
<input type="checkbox"/>	2527032	For Easy-Out™ Angle Adjustable Front Rigging				\$	200
<input type="checkbox"/>	2527030	For Swing Away ELR Front Rigging				\$	100
<input type="checkbox"/>	2527040	<b>Length Adjustable Heel Loops</b>				E0951	\$ 50
<input type="checkbox"/>	2527050	<b>Bodypoint Padded Calf Strap</b>				K0038	\$ 195



### Front Rigging Extension Tubes

#### Regular Extension Tubes

- 2528010 Short 14"-18" for Easy-Out Front Rigging 13"-17" for Swing Away and Standard Plug In Front Rigging
- 2528020 Long 18"-24" for Easy-Out Front Rigging 15"-19" for Swing Away and Standard Plug In Front Rigging \$ 75



**Lower Leg Length** If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.

Seat Pan to Foot Plates: \_\_\_\_\_ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.

### Miscellaneous Options

- 2530020 Black  2530021 Frame Color K0065 \$ 195
- 2530060 Small Pad 3"x4"x3"x2" size  2530062 Large Pad 5"x5"x4"x3" size E0957 & E1028 \$ 650
- 2530080 **Gel Hip Guides** Reduces available Seat Width by 2". E0956 & E1028 \$ 495
- 2530110 **IV Pole** Not compatible with Attendant Wheel Locks. Requires frame modification and cannot be added to chair later. E0776 \$ 295
- 2530120 **O2 Holder** Not compatible with Attendant Wheel Locks. Requires frame modification and cannot be added to chair later. E2208 \$ 395
- 2530200 **Q'Straint Lap Belt** Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" Belt Length. \$ 245
- 2530222 **Wheelchair Backpack** Black with PDG Logo. 17"x14" size. \$ 145
- 2530170 **Adhesive Back Hook Strip, 2" x 30"** For Cushion mounting. \$ 10

### On-Chair Seating & Positioning

Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. [pdgmobility.com/seating](http://pdgmobility.com/seating)

### Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions. Use of HCPCS codes does not guarantee coverage or payment.

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.