

# Bentley Manual Tilt Wheelchair

Ideal for frail individuals with limited stamina who wish to maintain independent mobility.

Optional long-frame version is ideal for individuals with high agitation.



- Tilt's 20 degrees, rotating about the front of the seat – ideal for foot propulsion while in-tilt
- Improves posture and prevents sliding forward out of the wheelchair
- Excellent access to rear wheels for effortless propulsion
- Adjusts to fit individual need
- Lots of options available



## Ideal for foot propulsion

- Seat tilts 20 degrees, rotating about the front of the seat pan so user can foot propel while in tilt
- Ample clear space under the front of the seat ensures excellent heel clearance
- Front casters are located further apart compared to most wheelchairs, allowing extra room for feet
- Seat height adjusts 13" (33cm) – 20" (50.8cm)




## Plug-in Front Rigging Option

- Choose traditional swing-away front rigging or plug-in style to ensure legs have ample clear space for foot propulsion



# Rear Wheel Adjusts to Individual Need



- Rear wheel is infinitely adjustable
- Wheel lock and anti-tipper adjust automatically with rear wheel
- Ensures excellent access and effortless propulsion for frail individuals



# Adjusts to Fit Individual Need

- Choose tilt activator location for Care-giver or user access




- Seat width 14" (35.6cm) – 32" (81.3cm)
- Seat depth 16" (40.6cm) – 22" (55.9cm)
- Seat height 13" (33cm) – 20" (50.8cm)
- Backrest angle 90 ° - 120°
- Weight limit 250 lb (113.4kg), option to 450 lb (203 kg)



# High Agitation Option

- Frame extensions ensure extra stability for high agitation applications
- Choose one of two anti-tipping options for extra stability



# High Agitation Padding



- Optional padding available for high agitation applications ensures user comfort and safety



# Thanks for taking a look at the Bentley

Visit [www.PDGmobility.com](http://www.PDGmobility.com) and  
have a look at Bentley Case  
Histories.

