

<b>PO#</b>	<b>Tag:</b>	<b>Dealer Name:</b>	<b>Dealer Contact:</b>
<b>Mailing Address:</b>		<b>Shipping Address:</b>	
<b>City:</b>	<b>Province:</b>	<b>Postal Code:</b>	<b>City:</b>
<b>Tel:</b>	<b>Fax:</b>		<b>Tel:</b>
<b>BENTLEY FRAME TYPE - 20 DEGREES OF TILT</b>		<b>Confirm Order Via E-mail:</b> <input type="checkbox"/> YES <input type="checkbox"/> NO <b>Send PDG e-News:</b> <input type="checkbox"/> YES <input type="checkbox"/> NO <b>E-Mail:</b>	
<input checked="" type="checkbox"/> 50330	Bentley: 20 Degree Manual Tilt Wheelchair	\$3,995	

<b>CLIENT INFORMATION - (Ensures optimal set-up of chair)</b>		
<input checked="" type="checkbox"/> 50545	Weight (lbs.) _____ Height (ft./in.) _____	NC

FRAME OPTIONS		
<input checked="" type="checkbox"/> 50331	Standard Base & Seat Frame Max Weight Capacity 250 lbs.	NC
<input type="checkbox"/> 50546	Transit Tie-Down System (TTDS) <small>(TTDS includes 4 factory installed transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks or Extended Base.</small>	NC
<input type="checkbox"/> 50547	Heavy Duty Package (Double Gas Strut) <small>Max Weight Capacity 450 lbs. Requires 20" composite or 24" heavy duty spoked wheel. Fixed axles only. Order as "add on" for seat widths up to 22". Heavy duty package included on seat width 23" or wider.</small>	\$300
<input type="checkbox"/> 51304	Extended Base with Rear Floor Jacks <small>Lifts rear wheels off ground for increased stability. Not compatible with 22" rear wheel, Attendant Wheel Locks, Anti tipper, IV Pole, O2 Holder. Limited ground clearance.</small>	\$600

SEAT WIDTH - (Inches)		
<input type="checkbox"/> 50549	14" Seat Width	NC
<input type="checkbox"/> 50550	15" Seat Width	NC
<input type="checkbox"/> 50551	16" Seat Width	NC
<input type="checkbox"/> 50552	17" Seat Width	NC
<input type="checkbox"/> 50553	18" Seat Width	NC
<input type="checkbox"/> 50554	19" Seat Width	NC
<input type="checkbox"/> 50555	20" Seat Width	NC
<input type="checkbox"/> 50556	21" Seat Width	\$195
<input type="checkbox"/> 50557	22" Seat Width	\$195

CUSTOM SEAT WIDTH - (Requires 4 Week Delivery)		
<small>Max weight capacity 450 lbs. Requires 20" composite or 24" heavy duty spoked wheel. Fixed axles only.</small>		
<input type="checkbox"/> 50558	23" Seat Width	\$600
<input type="checkbox"/> 50559	24" Seat Width	\$600
<input type="checkbox"/> 50560	25" Seat Width	\$600
<input type="checkbox"/> 50561	26" Seat Width	\$600
<input type="checkbox"/> 50562	27" Seat Width	\$600
<input type="checkbox"/> 50563	28" Seat Width	\$900
<input type="checkbox"/> 50564	29" Seat Width	\$900
<input type="checkbox"/> 50565	30" Seat Width	\$900
<input type="checkbox"/> 50566	31" Seat Width	\$900
<input type="checkbox"/> 50567	32" Seat Width	\$900

SEAT DEPTH - (Adjustable)		
<input type="checkbox"/> 50568	16" Seat Depth	NC
<input type="checkbox"/> 50569	17" Seat Depth	NC
<input type="checkbox"/> 50570	18" Seat Depth	NC
<input type="checkbox"/> 50571	19" Seat Depth	NC
<input type="checkbox"/> 50572	20" Seat Depth	NC
<input type="checkbox"/> 50573	21" Seat Depth	NC
<input type="checkbox"/> 50574	22" Seat Depth	NC

FRAME COLOR		
<input type="checkbox"/> 50575	Blue Shimmer	NC
<input type="checkbox"/> 50576	Tracer Red	NC
<input type="checkbox"/> 50577	Black Shimmer	NC
<input type="checkbox"/> 50578	Burnt Copper	NC
<input type="checkbox"/> 50579	Misty Purple	NC
<input type="checkbox"/> 50580	Forest Green	NC
<input type="checkbox"/> 51790	Custom Color: _____ (Contact PDG for Specification)	\$300

ARM TYPE		
<input type="checkbox"/> 50582	No Armrests	STD
<input type="checkbox"/> 50583	Dual Post Fixed Height 12" Full length pad only.	NC
<input type="checkbox"/> 50584	Dual Post Adjustable Height Armrest <small>Option: 8" - 11" <input type="checkbox"/> 10" - 14" <input type="checkbox"/> 12" - 16" <input type="checkbox"/></small>	\$295
<input type="checkbox"/> 50587	3" Armrest Extension - Forward <small>Not compatible with flip-up armrests.</small>	\$59
<input type="checkbox"/> 50588	3" Armrest Extension - Rearward <small>Not compatible with flip-up armrests.</small>	\$59
<input type="checkbox"/> 50589	Locking Flip-up Armrests <small>Not recommended for clients with excessive strength or tone. Not compatible with fold down back, clothing guards, back upholstery.</small>	\$295
<small>Set At: 10" <input type="checkbox"/> 11" <input type="checkbox"/> 12" <input type="checkbox"/> 13" <input type="checkbox"/> 14" <input type="checkbox"/></small>		

ARMPAD OPTIONS		
<input type="checkbox"/> 50594	Desk Length Armpad, 2" X 10" (pair)	NC
<input type="checkbox"/> 50595	Full Length Armpad, 2" X 12" (pair)	NC
<input type="checkbox"/> 50596	3" Wide Arm Pads <small>May interfere with wheels in some tilt configuration.</small>	\$49
<input type="checkbox"/> 50597	Full Length 3"x12" (pair)	
<input type="checkbox"/> 50598	Desk Length 3"x12" (pair)	
<input type="checkbox"/> 50599	Gel Armpad <small>May interfere with wheels in some tilt configuration.</small>	\$115
<input type="checkbox"/> 50598	Full Length 2"x12" (pair)	
<input type="checkbox"/> 50599	Full Length 3.5"x12" (pair)	
<input type="checkbox"/> 50600	Contour Arm Trough <small>May interfere with wheels in some tilt configuration.</small>	\$85
<input type="checkbox"/> 50600	Right 4.5"x13"	
<input type="checkbox"/> 50601	Left 4.5"x13"	
<input type="checkbox"/> 50602	Flat Hand Pan <small>Requires Contour Arm Trough.</small>	\$55
<input type="checkbox"/> 50602	Right 4.5"x5"	
<input type="checkbox"/> 50603	Left 4.5"x5"	
<input type="checkbox"/> 50604	Clothing Guards (Black Plastic) <small>Not compatible with flip-up armrests.</small>	\$205
<input type="checkbox"/> 50605	Armrest wideners (pr) <small>Not compatible with Clothing Guards.</small>	\$375

SEAT FRAME OPTIONS		
<input type="checkbox"/> 50606	No Seat Pan	STD
<input type="checkbox"/> 50607	Solid Seat Pan	\$140

<b>PO#</b>	<b>Tag:</b>	<b>Dealer Name:</b>	<b>Dealer Contact:</b>
------------	-------------	---------------------	------------------------

SEAT CUSHION OPTIONS - (The Comfort Company)		
<input type="checkbox"/> 50608	No Seat Cushion	STD
<input type="checkbox"/> 52101	Cushion M2 Zero Elevation Stretch-Air	\$577
Width: 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20" <input type="checkbox"/>		
Depth: 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20" <input type="checkbox"/>		
<input type="checkbox"/> 52116	Spare Cushion Cover	\$109
<input type="checkbox"/> 52110	Cushion M2 Zero Elevation Stretch-Air (Extra Wide)	\$677
Width: 22" <input type="checkbox"/> 24" <input type="checkbox"/>		
Depth: 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20" <input type="checkbox"/>		
<input type="checkbox"/> 52125	Spare Cushion Cover (Extra Wide)	\$109
<input type="checkbox"/>	Incontinent Cushion Liner	\$85

CASTER OPTIONS		
<input type="checkbox"/> 50615	6" x 1-1/4" Urethane	NC
<input type="checkbox"/> 50616	6" x 1-1/4" Pneumatic	\$110
<input type="checkbox"/> 51311	6" x 1-3/4" Urethane	\$75
<input type="checkbox"/> 50617	8" x 1-1/4" Urethane	NC
<input type="checkbox"/> 50618	8" x 1-1/4" Pneumatic	\$110
<input type="checkbox"/> 50619	8" x 2" Urethane Caster	\$135
<input type="checkbox"/> 50620	8" x 2" Pneumatic Caster	\$135

8" Casters - Compatible only with 16" and 24" rear wheels.

REAR WHEEL OPTIONS		
<input type="checkbox"/> 50621	12" Composite Wheel	NC
<input type="checkbox"/> 50622	16" Composite Wheel	NC
<input type="checkbox"/> 50623	20" Composite Wheel	Required for Heavy Duty Package. NC
<input type="checkbox"/> 51312	20" Spoke Wheel	Required for Heavy Duty Package. \$75
<input type="checkbox"/> 51461	22" Composite Wheel	NC
<input type="checkbox"/> 51460	22" Spoke Wheel	\$75
<input type="checkbox"/> 50624	24" Composite Wheel	NC
<input type="checkbox"/> 51313	24" Spoke Wheel	\$75
<input type="checkbox"/> 50625	24" Spoked Wheel - Heavy Duty	\$350

Required for Heavy Duty Package. Not compatible with Pneumatic rear tires with Flat Free Inserts.

SEAT HEIGHT OPTIONS																																				
Rear Wheel	6" Casters					<p><b>(CIRCLE) desired seat-to-floor height in chart.</b></p> <table border="1" style="margin: 10px auto; border-collapse: collapse;"> <thead> <tr> <th style="text-align: center;">Rear Wheel</th> <th colspan="5" style="text-align: center;">8" Casters</th> </tr> </thead> <tbody> <tr> <td style="text-align: center;">16"</td> <td style="text-align: center;">15"</td> <td style="text-align: center;">16"</td> <td style="text-align: center;">17"</td> <td style="text-align: center;">18"</td> <td style="text-align: center;">nc</td> </tr> <tr> <td style="text-align: center;">18"</td> <td style="text-align: center;">18"</td> <td style="text-align: center;">19"</td> <td style="text-align: center;">20"</td> <td style="text-align: center;">\$24</td> <td></td> </tr> <tr> <td style="text-align: center;">24"</td> <td style="text-align: center;">15"</td> <td style="text-align: center;">16"</td> <td style="text-align: center;">17"</td> <td style="text-align: center;">18"</td> <td style="text-align: center;">nc</td> </tr> <tr> <td></td> <td style="text-align: center;">19"</td> <td style="text-align: center;">20"</td> <td style="text-align: center;">\$24</td> <td></td> <td></td> </tr> </tbody> </table>	Rear Wheel	8" Casters					16"	15"	16"	17"	18"	nc	18"	18"	19"	20"	\$24		24"	15"	16"	17"	18"	nc		19"	20"	\$24		
Rear Wheel	8" Casters																																			
16"	15"	16"	17"	18"	nc																															
18"	18"	19"	20"	\$24																																
24"	15"	16"	17"	18"	nc																															
	19"	20"	\$24																																	
12"	13"	14"	15"	16"	nc																															
			17"	18"	\$24																															
16"	15"	16"	17"	18"	nc																															
			19"	20"	\$24																															
20"	13"	14"	15"	16"	nc																															
			17"	18"	\$24																															
22"	14"	15"	16"	17"	nc																															
			18"	19"	\$24																															
24"	15"	16"	17"	18"	nc																															
			19"	20"	\$24																															

REAR TIRE OPTIONS		
<input type="checkbox"/> 50637	Urethane Tires	STD
<input type="checkbox"/> 50638	Pneumatic Tires	NC
Pneumatic tires may result in a finished seat height 1/2" higher than stated.		
<input type="checkbox"/> 50639	Pneumatic Tires With Flat-Free Inserts	\$200
Pneumatic tires may result in a finished seat height 1/2" higher than stated. Flat Free Inserts not compatible with Spoked Wheels.		

AXLE TYPE		
<input type="checkbox"/> 50640	Fixed Rear Axle	Required for Heavy Duty Package. STD
<input type="checkbox"/> 50641	Quick Release Rear Axles	NC

HAND RIM OPTIONS		
<input type="checkbox"/> 51632	Omit Hand Rim	20", 22" or 24" rear wheels only. NC
<input type="checkbox"/> 50642	Standard Hand Rims	20", 22" or 24" rear wheels only. NC
<input type="checkbox"/> 50643	Plastic Coated Hand Rims	Require a 4 to 6 week lead time. \$95
<input type="checkbox"/> 50644	Plastic Coated Projection Hand Rims	20", 22" or 24" rear wheels only. \$225
<input type="checkbox"/> 50645	Number of Projections - 8	
<input type="checkbox"/> 50649	Projection Length - 1"	
<input type="checkbox"/> 50650	Select Projection Angle: 0° <input type="checkbox"/> 45° <input type="checkbox"/>	
<input type="checkbox"/> 50651		

WHEEL LOCKS		
<input type="checkbox"/> 50653	Push To Lock Wheel Locks	STD
<input type="checkbox"/> 50654	Pull To Lock Wheel Locks	NC
<input type="checkbox"/> 50655	Attendant Wheel Lock (Push To Lock)	\$85
Not compatible with Transit Tie Down System.		
<input type="checkbox"/> 50656	User Wheel Locks And Attendant Wheel Locks	\$150
Available only on 20" and 24" rear wheels. Not compatible with Transit Tie Down System.		
<input type="checkbox"/> 50657	Wheel Lock Extensions (Pair)	\$60
Not compatible with Attendant wheel locks.		
<input type="checkbox"/> 50658	Hill Holders	\$350

ANTI-TIP OPTIONS		
<input checked="" type="checkbox"/> 50659	Rear Anti-Tippers (Pair)	Required safety option. \$99
<input type="checkbox"/> 50660	Front Anti-Tippers (Pair)	\$300
Not compatible with 13" Seat to Floor Height.		
<input type="checkbox"/> 50661	Extended Base - with Heavy Duty Rear Stabilizers.	\$600
Provides additional rearward stability for clients with excessive movement. Not compatible with Attendant Wheel Locks, Rear Floor Jacks, IV Pole, O2 Holder.		

BACK POST STYLE		
<input type="checkbox"/> 50662	Push Bar Back - 25" High	STD
<input type="checkbox"/> 50663	Push Bar Back - 31" High	\$35
<input type="checkbox"/> 50664	Cane Handle Back - 25" High	\$35
Seat widths 14" - 22" only. Not compatible with headrest mounting plate.		
<input type="checkbox"/> 50665	Cane Handle Back - 31" High	\$35
Seat widths 14" - 22" only. Not compatible with headrest mounting plate.		

BACK OPTIONS		
<input type="checkbox"/> 50666	Adjustable - Height / Angle Push Bar Handle	\$395
Not compatible with headrest mounting plate or cane Handle Back.		
<input type="checkbox"/> 50667	Seat To Back Angle:	STD
Set At: 90° <input type="checkbox"/> 98° <input type="checkbox"/> 106° <input type="checkbox"/> 114° <input type="checkbox"/> 122° <input type="checkbox"/>		
<input type="checkbox"/> 50673	Backrest Extension 10" High	\$150
Not compatible with adjustable height / angle push bar handle.		
<input type="checkbox"/> 50674	Fold Down Back	NC
Not compatible with Reclining Back.		
<input checked="" type="checkbox"/> 50675	Spring Back Suspension	STD
Not compatible with Reclining Back.		

BACK UPHOLSTERY - (Back height measured from seat pan to top of upholstery)		
<input type="checkbox"/> 50676	No Upholstery	STD
<input type="checkbox"/> 50677	Standard Back Upholstery	\$125
Set Height At: 21" <input type="checkbox"/> 22" <input type="checkbox"/> 23" <input type="checkbox"/>		
<input type="checkbox"/> 50681	Short Upholstery - (25" frame only)	\$125
Set Height At: 18" <input type="checkbox"/> 19" <input type="checkbox"/> 20" <input type="checkbox"/>		
<input type="checkbox"/> 50684	Tall Upholstery - (31" frame only)	\$125
Set Height At: 24" <input type="checkbox"/> 25" <input type="checkbox"/> 26" <input type="checkbox"/>		
<input type="checkbox"/> 50688	Adjustable Tension Upholstery	\$250
Set Height At: 21" <input type="checkbox"/> 22" <input type="checkbox"/> 23" <input type="checkbox"/>		

<b>PO#</b>	<b>Tag:</b>	<b>Dealer Name:</b>	<b>Dealer Contact:</b>
------------	-------------	---------------------	------------------------

BACK OPTIONS A - (The Comfort Company)		
<input type="checkbox"/>	52064 Acta-Back Standard Stretch Air	
	Height: 16" <input type="checkbox"/> 18" <input type="checkbox"/>	
	Width: 15"- 17" <input type="checkbox"/> 17"- 19" <input type="checkbox"/> 19"- 21" <input type="checkbox"/>	\$815
	Width: 21"- 23" <input type="checkbox"/>	\$898
<input type="checkbox"/>	52066 Acta-Back Standard Stretch Air (Extra Tall-20")	
	Width: 15"- 17" <input type="checkbox"/> 17"- 19" <input type="checkbox"/> 19"- 21" <input type="checkbox"/>	\$898
	Width: 21"- 23" <input type="checkbox"/> 23"- 25" <input type="checkbox"/>	\$980

BACK OPTIONS B - (The Comfort Company)		
<input type="checkbox"/>	52081 Acta-Back Standard Stretch Air DEEP	
	Height: 16" <input type="checkbox"/> 18" <input type="checkbox"/>	
	Width: 15"- 17" <input type="checkbox"/> 17"- 19" <input type="checkbox"/> 19"- 21" <input type="checkbox"/>	\$947
	Width: 21"- 23" <input type="checkbox"/>	\$1,112
<input type="checkbox"/>	52083 Acta-Back Standard Stretch Air DEEP (Extra Tall-20")	
	Width: 15"- 17" <input type="checkbox"/> 17"- 19" <input type="checkbox"/> 19"- 21" <input type="checkbox"/>	\$1,043
	Width: 21"- 23" <input type="checkbox"/> 23"- 25" <input type="checkbox"/>	\$1,214

HEADREST OPTIONS - (The Comfort Company)		
<input type="checkbox"/>	50707 Headrest Mounting Plate	\$75
<input type="checkbox"/>	51846 Fully Adjustable Hardware	\$330
<input type="checkbox"/>	51848 Moldable Headrest Pad - Stretch-Air	\$272
	Width: 8" <input type="checkbox"/> 10" <input type="checkbox"/> 14" <input type="checkbox"/>	

TILT / RECLINE OPTIONS		
<input type="checkbox"/>	50706 Reclining Backrest(Manual) Max. Recline - 30°	\$1,200
	<small>Not compatible with head rest extension, adjustable angle pushbar handle, cane handle back, spring back suspension or fold down back. If upholstery is ordered with Recline, it must be tension adjustable.</small>	

ATTENDANT TILT / USER OPERATED SELF - TILT		
<input type="checkbox"/>	50714 Attendant Tilt Only - Right <input type="checkbox"/> Left <input type="checkbox"/>	STD
<input type="checkbox"/>	50716 Self Tilt Only	\$99
	<small>Not available on widths greater than 22. Full Length arm pads only.</small>	
	Right Armrest <input type="checkbox"/> Left Armrest <input type="checkbox"/>	
<input type="checkbox"/>	50718 Dual Activator	\$195
	<small>Not available on widths greater than 22. Full Length arm pads only.</small>	
	Backrest - Mount: Right <input type="checkbox"/> Left <input type="checkbox"/>	
	Armrest - Mount: Right <input type="checkbox"/> Left <input type="checkbox"/>	
<input type="checkbox"/>	50722 Trigger Extension	\$48

FRONT RIGGING		
<input type="checkbox"/>	50723 Omit Front Rigging	NC
<input type="checkbox"/>	50724 Sector Blocks Only (No Front Rigging Supplied)	NC
<input type="checkbox"/>	50725 Plug in 70° Front Rigging	STD
<input type="checkbox"/>	50726 Plug in Adjustable Knee Angle Front Rigging	\$290
<input type="checkbox"/>	50727 Swing Away 60° Front Rigging	NC
<input type="checkbox"/>	50729 Swing Away 70° Front Rigging	NC
<input type="checkbox"/>	50731 Swing Away Elevating Legrests	\$355

FOOT PLATES		
<input type="checkbox"/>	50472 Composite Footplates	STD
<input type="checkbox"/>	50733 Adjustable Angle Footplates	\$140
<input type="checkbox"/>	50735 One Piece FIXED Footplates	\$350
	<small>Available on plug-in front rigging only. 14"-22" wide only.</small>	
<input type="checkbox"/>	50736 One Piece FLIP-UP Footplates	\$395
<input type="checkbox"/>	51832 Interlocking Angle Adjustable Footplate	\$195
	<small>14"-21" wide only.</small>	

FOOTREST/LEGREST ACCESSORIES		
<input type="checkbox"/>	52252 60° & 70° Swing Away Calf Pad	\$135
<input type="checkbox"/>	51680 Swing Away Elevating Legrest Gel Knee Pad	\$99
<input type="checkbox"/>	50740 Adjustable Length Heel Loops	\$55
<input type="checkbox"/>	50741 Bodypoint Padded Calf Strap	\$49
	<small>Sizes greater than 22" are not padded.</small>	
<input type="checkbox"/>	50742 Miller's Dynamic Footrest Mounts For Plug In Front Rigging.	\$190
<input type="checkbox"/>	52253 Miller's Dynamic Footrest Mounts For Swing Away 70° Front Rigging	\$380

EXTENSION TUBE		
<input type="checkbox"/>	50737 Standard Length (13" - 17")	STD
<input type="checkbox"/>	50738 Long Extension Tubes (15" - 19")	\$75
<input type="checkbox"/>	50739 Standard Length (13" - 17") Solid Extension Tubes - Extra Strength	\$99
	<small>Available on Plug-in front Rigging only.</small>	
<input type="checkbox"/>	51454 Long Extension Tubes (15" - 19") Solid Extension Tubes - Extra Strength	\$99
	<small>Available on Plug-in front Rigging only.</small>	

POSITIONING STRAPS		
<input type="checkbox"/>	50743 Bodypoint - Center-Pull Auto Style Hip Belt	\$85
<input type="checkbox"/>	50744 Bodypoint - Padded, Center-Pull Auto Style Hip Belt	\$125
<input type="checkbox"/>	50745 Bodypoint - 4 Point Hip Belt	\$150
<input type="checkbox"/>	50746 Bodypoint - Chest Strap	\$75

MISCELLANEOUS OPTIONS		
<input type="checkbox"/>	51314 One Arm Drive Left <input type="checkbox"/> Right <input type="checkbox"/>	\$3,500
	<small>Max 22" wide, 24" Comp rw required, 18" STF Minimum.</small>	
<input type="checkbox"/>	51281 Protective Lower Padding	\$425
<input type="checkbox"/>	51316 IV Pole	\$295
	<small>Not compatible with TTDS (Transit Tie Down System), Rear Stabilizers, Rear Floor Jacks, Attendant Wheel Locks.</small>	
<input type="checkbox"/>	50747 O2 Holder	\$295
	<small>Not compatible with TTDS (Transit Tie Down System), Rear Stabilizers, Rear Floor Jacks, Attendant Wheel Locks.</small>	
<input type="checkbox"/>	50748 Accessory Tray	call for quote
	<small>May restrict chair specifications.</small>	

CUSTOM MODIFICATIONS - NOTES	
Total Retail Value - This order	\$
Have all desired specs been checked? A complete spec sheet ensures your chair will arrive on time and as you expect.	