



Bentley

PO#	Tag:	Dealer Name:	Dealer Contact:
Mailing Address:		Shipping Address:	
City:	State:	Zip Code:	
Tel:	Fax:		

BENTLEY FRAME TYPE - 20 DEGREES OF TILT			
<input checked="" type="checkbox"/> 50330	Bentley: 20 Degree Manual Tilt Wheelchair	K0009	\$2,895

Confirm Order Via E-mail:	<input type="checkbox"/> YES	<input type="checkbox"/> NO
Send PDG e-News:	<input type="checkbox"/> YES	<input type="checkbox"/> NO
E-Mail:		

CLIENT INFORMATION - (Ensures optimal set-up of chair)		
<input checked="" type="checkbox"/> 50545	Weight (lbs.) _____ Height (ft./in.) _____	NC

FRAME OPTIONS			
<input checked="" type="checkbox"/> 50331	Standard Base & Seat Frame Max Weight Capacity 250 lbs.	NC	
<input type="checkbox"/> 50546	Transit Tie-Down System (TTDS) <small>(TTDS includes 4 factory installed transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks or Extended Base.</small>	K0108	\$250
<input type="checkbox"/> 50547	Heavy Duty Package (Double Gas Strut) <small>Max Weight Capacity 450 lbs. Requires 20" composite or 24" heavy duty spoked wheel. Fixed axles only. Order as "add on" for seat widths up to 22". Heavy duty package included on seat width 23" or wider.</small>	K0108	\$300
<input type="checkbox"/> 51304	Extended Base with Rear Floor Jacks <small>Lifts rear wheels off ground for increased stability. Not compatible with 22" rear wheel, Attendant Wheel Locks, Anti tipper, IV Pole, O2 Holder. Limited ground clearance.</small>		\$600

SEAT WIDTH - (Inches)			
<input type="checkbox"/> 50549	14" Seat Width		NC
<input type="checkbox"/> 50550	15" Seat Width		NC
<input type="checkbox"/> 50551	16" Seat Width		NC
<input type="checkbox"/> 50552	17" Seat Width		NC
<input type="checkbox"/> 50553	18" Seat Width		NC
<input type="checkbox"/> 50554	19" Seat Width		NC
<input type="checkbox"/> 50555	20" Seat Width	E2201	\$300
<input type="checkbox"/> 50556	21" Seat Width	E2201	\$300
<input type="checkbox"/> 50557	22" Seat Width	E2201	\$300

CUSTOM SEAT WIDTH - (Requires 4 Week Delivery)			
<small>Max weight capacity 450 lbs. Requires 20" composite or 24" heavy duty spoked wheel. Fixed axles only.</small>			
<input type="checkbox"/> 50558	23" Seat Width	E2202	\$600
<input type="checkbox"/> 50559	24" Seat Width	E2202	\$600
<input type="checkbox"/> 50560	25" Seat Width	E2202	\$600
<input type="checkbox"/> 50561	26" Seat Width	E2202	\$600
<input type="checkbox"/> 50562	27" Seat Width	E2202	\$600
<input type="checkbox"/> 50563	28" Seat Width		\$900
<input type="checkbox"/> 50564	29" Seat Width		\$900
<input type="checkbox"/> 50565	30" Seat Width		\$900
<input type="checkbox"/> 50566	31" Seat Width		\$900
<input type="checkbox"/> 50567	32" Seat Width		\$900

SEAT DEPTH - (Adjustable)			
<input type="checkbox"/> 50568	16" Seat Depth		NC
<input type="checkbox"/> 50569	17" Seat Depth		NC
<input type="checkbox"/> 50570	18" Seat Depth		NC
<input type="checkbox"/> 50571	19" Seat Depth		NC
<input type="checkbox"/> 50572	20" Seat Depth	E2203	\$300
<input type="checkbox"/> 50573	21" Seat Depth	E2203	\$300
<input type="checkbox"/> 50574	22" Seat Depth	E2203	\$300

FRAME COLOR		
<input type="checkbox"/> 50575	Blue Shimmer	STD
<input type="checkbox"/> 50576	Tracer Red	NC
<input type="checkbox"/> 50577	Black Shimmer	NC
<input type="checkbox"/> 50578	Burnt Copper	NC
<input type="checkbox"/> 50579	Misty Purple	NC
<input type="checkbox"/> 50580	Forest Green	NC
<input type="checkbox"/> 51790	Custom Color: _____ (Contact PDG for Specification)	\$300

ARM TYPE			
<input type="checkbox"/> 50582	No Armrests		STD
<input type="checkbox"/> 50583	Dual Post Fixed Height 12" Full length pad only.		NC
<input type="checkbox"/> 50584	Dual Post Adjustable Height Armrest	E0973	\$245
Option: 8" - 11" <input type="checkbox"/> 10" - 14" <input type="checkbox"/> 12" - 16" <input type="checkbox"/>			
<input type="checkbox"/> 50587	3" Armrest Extension - Forward		\$59
<small>Not compatible with flip-up armrests.</small>			
<input type="checkbox"/> 50588	3" Armrest Extension - Rearward		\$59
<small>Not compatible with flip-up armrests.</small>			
<input type="checkbox"/> 50589	Locking Flip-up Armrests	E0973	\$245
<small>Not recommended for clients with excessive strength or tone. Not compatible with fold down back, clothing guards, back upholstery.</small>			
Set At: 10" <input type="checkbox"/> 11" <input type="checkbox"/> 12" <input type="checkbox"/> 13" <input type="checkbox"/> 14" <input type="checkbox"/>			

ARMPAD OPTIONS			
<input type="checkbox"/> 50594	Desk Length Armpad, 2" X 10" (pair)		NC
<input type="checkbox"/> 50595	Full Length Armpad, 2" X 12" (pair)		NC
<input type="checkbox"/> 50596	3" Wide Arm Pads		\$49
<small>May interfere with wheels in some tilt configuration.</small>			
<input type="checkbox"/> 50596	Full Length 3"x12" (pair)		
<input type="checkbox"/> 50597	Desk Length 3"x12" (pair)		
<input type="checkbox"/> 50598	Gel Armpad		\$115
<small>May interfere with wheels in some tilt configuration.</small>			
<input type="checkbox"/> 50598	Full Length 2"x12" (pair)		
<input type="checkbox"/> 50599	Full Length 3.5"x12" (pair)		
<input type="checkbox"/> 50600	Contour Arm Trough	E2209	\$85
<small>May interfere with wheels in some tilt configuration.</small>			
<input type="checkbox"/> 50600	Right 4.5"x13"		
<input type="checkbox"/> 50601	Left 4.5"x13"		
<input type="checkbox"/> 50602	Flat Hand Pan		\$55
<small>Requires Contour Arm Trough.</small>			
<input type="checkbox"/> 50602	Right 4.5"x5"		
<input type="checkbox"/> 50603	Left 4.5"x5"		
<input type="checkbox"/> 50604	Clothing Guards (Black Plastic)		\$120
<small>Not compatible with flip-up armrests.</small>			
<input type="checkbox"/> 50605	Armrest wideners (pr)		\$375
<small>Not compatible with Clothing Guards.</small>			

SEAT FRAME OPTIONS			
<input type="checkbox"/> 50606	No Seat Pan		STD
<input type="checkbox"/> 50607	Solid Seat Pan	E2231	\$95

PO#	Tag:	Dealer Name:	Dealer Contact:
------------	-------------	---------------------	------------------------

SEAT CUSHION OPTIONS - (The Comfort Company)			
<input type="checkbox"/> 50608	No Seat Cushion		STD
<input type="checkbox"/> 52101	Cushion M2 Zero Elevation Stretch-Air	E2607	\$297
	Width: 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20" <input type="checkbox"/>		
	Depth: 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20" <input type="checkbox"/>		
<input type="checkbox"/> 52116	Spare Cushion Cover		\$60
<input type="checkbox"/> 52110	Cushion M2 Zero Elevation Stretch-Air (Extra Wide)	E2608	\$398
	Width: 22" <input type="checkbox"/> 24" <input type="checkbox"/>		
	Depth: 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20" <input type="checkbox"/>		
<input type="checkbox"/> 52125	Spare Cushion Cover (Extra Wide)		\$60
<input type="checkbox"/>	Incontinent Cushion Liner		\$58

CASTER OPTIONS			
<input type="checkbox"/> 50615	6" x 1-1/4" Urethane		NC
<input type="checkbox"/> 50616	6" x 1-1/4" Pneumatic		\$110
<input type="checkbox"/> 51311	6" x 1-3/4" Urethane		\$75
<input type="checkbox"/> 50617	8" x 1-1/4" Urethane	8" Casters - Compatible only with 16" and 24" rear wheels.	NC
<input type="checkbox"/> 50618	8" x 1-1/4" Pneumatic		\$110
<input type="checkbox"/> 50619	8" x 2" Urethane Caster		\$135
<input type="checkbox"/> 50620	8" x 2" Pneumatic Caster		\$135

REAR WHEEL OPTIONS			
<input type="checkbox"/> 50621	12" Composite Wheel		NC
<input type="checkbox"/> 50622	16" Composite Wheel		NC
<input type="checkbox"/> 50623	20" Composite Wheel	Required for Heavy Duty Package.	NC
<input type="checkbox"/> 51312	20" Spoke Wheel	Required for Heavy Duty Package.	\$75
<input type="checkbox"/> 51461	22" Composite Wheel		NC
<input type="checkbox"/> 51460	22" Spoke Wheel		\$75
<input type="checkbox"/> 50624	24" Composite Wheel		NC
<input type="checkbox"/> 51313	24" Spoke Wheel		\$75
<input type="checkbox"/> 50625	24" Spoked Wheel - Heavy Duty	Required for Heavy Duty Package. Not compatible with Pneumatic rear tires with Flat Free Inserts.	\$350

SEAT HEIGHT OPTIONS						
Rear Wheel	6" Casters					(CIRCLE) desired seat-to-floor height in chart.
12"	13"	14"	15"	16"	nc	
			17"	18"	\$24	
16"	15"	16"	17"	18"	nc	
			19"	20"	\$24	
20"	13"	14"	15"	16"	nc	
			17"	18"	\$24	
22"	14"	15"	16"	17"	nc	
			18"	19"	\$24	
24"	15"	16"	17"	18"	nc	
			19"	20"	\$24	
Rear Wheel	8" Casters					
16"	15"	16"	17"	18"	nc	
			19"	20"	\$24	
24"	15"	16"	17"	18"	nc	
			19"	20"	\$24	

REAR TIRE OPTIONS			
<input type="checkbox"/> 50637	Urethane Tires		STD
<input type="checkbox"/> 50638	Pneumatic Tires		\$60
Pneumatic tires may result in a finished seat height 1/2" higher than stated.			
<input type="checkbox"/> 50639	Pneumatic Tires With Flat-Free Inserts	E2213	\$70
Pneumatic tires may result in a finished seat height 1/2" higher than stated. Flat Free Inserts not compatible with Spoked Wheels.			

AXLE TYPE			
<input type="checkbox"/> 50640	Fixed Rear Axle	Required for Heavy Duty Package.	STD
<input type="checkbox"/> 50641	Quick Release Rear Axles		NC

HAND RIM OPTIONS			
<input type="checkbox"/> 51632	Omit Hand Rim	20", 22" or 24" rear wheels only.	NC
<input type="checkbox"/> 50642	Standard Hand Rims	20", 22" or 24" rear wheels only.	STD
<input type="checkbox"/> 50643	Plastic Coated Hand Rims	Require a 4 to 6 week lead time.	\$95
<input type="checkbox"/> 50644 50645 50649 50650 50651	Plastic Coated Projection Hand Rims	20", 22" or 24" rear wheels only.	\$225
	Number of Projections - 8		
	Projection Length - 1"		
	Select Projection Angle: 0° <input type="checkbox"/> 45° <input type="checkbox"/>		

WHEEL LOCKS			
<input type="checkbox"/> 50653	Push To Lock Wheel Locks		STD
<input type="checkbox"/> 50654	Pull To Lock Wheel Locks		NC
<input type="checkbox"/> 50655	Attendant Wheel Lock (Push To Lock)	Not compatible with Transit Tie Down System.	\$85
<input type="checkbox"/> 50656	User Wheel Locks And Attendant Wheel Locks	Available only on 20" and 24" rear wheels. Not compatible with Transit Tie Down System.	\$125
<input type="checkbox"/> 50657	Wheel Lock Extensions (Pair)	E0971	\$60
Not compatible with Attendant wheel locks.			
<input type="checkbox"/> 50658	Hill Holders		\$350

ANTI-TIP OPTIONS			
<input checked="" type="checkbox"/> 50659	Rear Anti-Tippers (Pair)	Required safety option. E0971	\$135
<input type="checkbox"/> 50660	Front Anti-Tippers (Pair)	Not compatible with 13" Seat to Floor Height.	\$300
<input type="checkbox"/> 50661	Extended Base - with Heavy Duty Rear Stabilizers.	Provides additional rearward stability for clients with excessive movement. Not compatible with Attendant Wheel Locks, Rear Floor Jacks, IV Pole, O2 Holder.	\$600

BACK POST STYLE			
<input type="checkbox"/> 50662	Push Bar Back - 25" High		STD
<input type="checkbox"/> 50663	Push Bar Back - 31" High		\$35
<input type="checkbox"/> 50664	Cane Handle Back - 25" High	Seat widths 14" - 22" only.	\$35
	Not compatible with headrest mounting plate.		
<input type="checkbox"/> 50665	Cane Handle Back - 31" High	Seat widths 14" - 22" only.	\$35
	Not compatible with headrest mounting plate.		

BACK OPTIONS			
<input type="checkbox"/> 50666	Adjustable - Height / Angle Push Bar Handle	Not compatible with headrest mounting plate or cane Handle Back.	\$235
<input type="checkbox"/> 50667	Seat To Back Angle:		STD
	Set At: 90° <input type="checkbox"/> 98° <input type="checkbox"/> 106° <input type="checkbox"/> 114° <input type="checkbox"/> 122° <input type="checkbox"/>		
<input type="checkbox"/> 50673	Backrest Extension 10" High	Not compatible with adjustable height / angle push bar handle.	\$120
<input type="checkbox"/> 50674	Fold Down Back	Not compatible with Reclining Back.	NC
<input checked="" type="checkbox"/> 50675	Spring Back Suspension	Not compatible with Reclining Back.	STD

BACK UPHOLSTERY - (Back height measured from seat pan to top of upholstery)			
<input type="checkbox"/> 50676	No Upholstery		STD
<input type="checkbox"/> 50677	Standard Back Upholstery		\$125
	Set Height At: 21" <input type="checkbox"/> 22" <input type="checkbox"/> 23" <input type="checkbox"/>		
<input type="checkbox"/> 50681	Short Upholstery - (25" frame only)		\$125
	Set Height At: 18" <input type="checkbox"/> 19" <input type="checkbox"/> 20" <input type="checkbox"/>		
<input type="checkbox"/> 50684	Tall Upholstery - (31" frame only)		\$125
	Set Height At: 24" <input type="checkbox"/> 25" <input type="checkbox"/> 26" <input type="checkbox"/>		
<input type="checkbox"/> 50688	Adjustable Tension Upholstery	K0108	\$250
	Set Height At: 21" <input type="checkbox"/> 22" <input type="checkbox"/> 23" <input type="checkbox"/>		

PO#	Tag:	Dealer Name:	Dealer Contact:
------------	-------------	---------------------	------------------------

BACK OPTIONS A - (The Comfort Company)			
<input type="checkbox"/> 52064	Acta-Back Standard Stretch Air		
	Height: 16" <input type="checkbox"/> 18" <input type="checkbox"/>		
	Width: 15"-17" <input type="checkbox"/> 17"-19" <input type="checkbox"/> 19"-21" <input type="checkbox"/>	E2613	\$497
	Width: 21"-23" <input type="checkbox"/>	E2614	\$584
<input type="checkbox"/> 52066	Acta-Back Standard Stretch Air (Extra Tall-20")		
	Width: 15"-17" <input type="checkbox"/> 17"-19" <input type="checkbox"/> 19"-21" <input type="checkbox"/>	E2613	\$565
	Width: 21"-23" <input type="checkbox"/> 23"-25" <input type="checkbox"/>	E2614	\$640

BACK OPTIONS B - (The Comfort Company)			
<input type="checkbox"/> 52081	Acta-Back Standard Stretch Air DEEP		
	Height: 16" <input type="checkbox"/> 18" <input type="checkbox"/>		
	Width: 15"-17" <input type="checkbox"/> 17"-19" <input type="checkbox"/> 19"-21" <input type="checkbox"/>	E2615	\$598
	Width: 21"-23" <input type="checkbox"/>	E2616	\$664
<input type="checkbox"/> 52083	Acta-Back Standard Stretch Air DEEP (Extra Tall-20")		
	Width: 15"-17" <input type="checkbox"/> 17"-19" <input type="checkbox"/> 19"-21" <input type="checkbox"/>	E2615	\$665
	Width: 21"-23" <input type="checkbox"/> 23"-25" <input type="checkbox"/>	E2616	\$739

HEADREST OPTIONS - (The Comfort Company)			
<input type="checkbox"/> 50707	Headrest Mounting Plate		\$75
<input type="checkbox"/> 51846	Fully Adjustable Hardware	E1028	\$211
<input type="checkbox"/> 51848	Moldable Headrest Pad - Stretch-Air	E0955	\$221
	Width: 8" <input type="checkbox"/> 10" <input type="checkbox"/> 14" <input type="checkbox"/>		

TILT / RECLINE OPTIONS			
<input type="checkbox"/> 50706	Reclining Backrest(Manual) Max. Recline - 30" <small>Not compatible with head rest extension, adjustable angle pushbar handle, cane handle back, spring back suspension or fold down back. If upholstery is ordered with Recline, it must be tension adjustable.</small>	E1225	\$830

ATTENDANT TILT / USER OPERATED SELF - TILT			
<input type="checkbox"/> 50714	Attendant Tilt Only - Right <input type="checkbox"/> Left <input type="checkbox"/>		STD
<input type="checkbox"/> 50716	Self Tilt Only <small>Not available on widths greater than 22. Full Length arm pads only.</small>		\$99
	Right Armrest <input type="checkbox"/> Left Armrest <input type="checkbox"/>		
<input type="checkbox"/> 50718	Dual Activator <small>Not available on widths greater than 22. Full Length arm pads only.</small>		\$195
	Backrest - Mount: Right <input type="checkbox"/> Left <input type="checkbox"/>		
	Armrest - Mount: Right <input type="checkbox"/> Left <input type="checkbox"/>		
<input type="checkbox"/> 50722	Trigger Extension		\$48

FRONT RIGGING			
<input type="checkbox"/> 50723	Omit Front Rigging		NC
<input type="checkbox"/> 50724	Sector Blocks Only (No Front Rigging Supplied)		NC
<input type="checkbox"/> 50725	Plug in 70° Front Rigging		STD
<input type="checkbox"/> 50726	Plug in Adjustable Knee Angle Front Rigging		\$240
<input type="checkbox"/> 50727	Swing Away 60° Front Rigging		NC
<input type="checkbox"/> 50729	Swing Away 70° Front Rigging		NC
<input type="checkbox"/> 50731	Swing Away Elevating Legrests	E0990	\$355

FOOT PLATES			
<input type="checkbox"/> 50472	Composite Footplates		STD
<input type="checkbox"/> 50733	Adjustable Angle Footplates	E0040	\$155
<input type="checkbox"/> 50735	One Piece FIXED Footplates <small>Available on plug-in front rigging only. 14"-22" wide only.</small>		\$350
<input type="checkbox"/> 50736	One Piece FLIP-UP Footplates		\$395
<input type="checkbox"/> 51832	Interlocking Angle Adjustable Footplate <small>14"-21" wide only.</small>		\$195

FOOTREST/LEGREST ACCESSORIES			
<input type="checkbox"/> 52252	60° & 70° Swing Away Calf Pad		\$135
<input type="checkbox"/> 51680	Swing Away Elevating Legrest Gel Knee Pad		\$99
<input type="checkbox"/> 50740	Adjustable Length Heel Loops	E0951	\$50
<input type="checkbox"/> 50741	Bodypoint Padded Calf Strap <small>Sizes greater than 22" are not padded.</small>	K0038	\$49
<input type="checkbox"/> 50742	Miller's Dynamic Footrest Mounts For Plug In Front Rigging.		\$190
<input type="checkbox"/> 52253	Miller's Dynamic Footrest Mounts For Swing Away 70° Front Rigging		\$380

EXTENSION TUBE			
<input type="checkbox"/> 50737	Standard Length (13" - 17")		STD
<input type="checkbox"/> 50738	Long Extension Tubes (15" - 19")		\$75
<input type="checkbox"/> 50739	Standard Length (13" - 17") Solid Extension Tubes - Extra Strength <small>Available on Plug-in front Rigging only.</small>		\$99
<input type="checkbox"/> 51454	Long Extension Tubes (15" - 19") Solid Extension Tubes - Extra Strength <small>Available on Plug-in front Rigging only.</small>		\$99

POSITIONING STRAPS			
<input type="checkbox"/> 50743	Bodypoint - Center-Pull Auto Style Hip Belt	E0978	\$40
<input type="checkbox"/> 50744	Bodypoint - Padded, Center-Pull Auto Style Hip Belt	E0978	\$90
<input type="checkbox"/> 50745	Bodypoint - 4 Point Hip Belt	E0978	\$110
<input type="checkbox"/> 50746	Bodypoint - Chest Strap	E0978	\$75

MISCELLANEOUS OPTIONS			
<input type="checkbox"/> 51314	One Arm Drive Left <input type="checkbox"/> Right <input type="checkbox"/> <small>Max 22" wide, 24" Comp rw required, 18" STF Minimum.</small>	E0958	\$3,500
<input type="checkbox"/> 51281	Protective Lower Padding		\$425
<input type="checkbox"/> 51316	IV Pole <small>Not compatible with TTDS (Transit Tie Down System), Rear Stabilizers, Rear Floor Jacks, Attendant Wheel Locks.</small>		\$295
<input type="checkbox"/> 50747	O2 Holder <small>Not compatible with TTDS (Transit Tie Down System), Rear Stabilizers, Rear Floor Jacks, Attendant Wheel Locks.</small>	E2208	\$295
<input type="checkbox"/> 50748	Accessory Tray <small>May restrict chair specifications.</small>	E1029	call for quote

CUSTOM MODIFICATIONS - NOTES	
Total Retail Value - This order	\$
Have all desired specs been checked? A complete spec sheet ensures your chair will arrive on time and as you expect.	
HCPCS codes supplied are for reference only and do not ensure coverage or payment.	