



# Eclipse Express

**Canadian Price List**  
Revised April-2012 Page 1 of 1  
**THE COMFORT COMPANY**

PO#	Tag:	Dealer Name:	Dealer Contact:
Mailing Address:		Shipping Address:	
City:	Province:	Postal Code:	
Tel:	Fax:		

## EXPRESS FORM

**THIS CONFIGURATION INCLUDES THE FOLLOWING  
PRE-SELECTED OPTIONS - NO SUBSTITUTIONS**

Confirm Order Via E-mail:	<input type="checkbox"/> YES	<input type="checkbox"/> NO
Send PDG e-News:	<input type="checkbox"/> YES	<input type="checkbox"/> NO
E-Mail:		

PDG ECLIPSE		
<input checked="" type="checkbox"/> 50908	PDG Eclipse 600, Bariatric Wheelchair	\$3,649
<input type="checkbox"/> 50909	Standard Frame - 600 lbs Weight Limit	STD
<input type="checkbox"/> 50910	Heavy Duty Update - 1000 lbs Weight Limit	\$150

REAR WHEEL (Dealer Adjustable)		
<input checked="" type="checkbox"/> 51047	24" Spoked Wheel - Heavy Duty	STD
<input checked="" type="checkbox"/> 51048	Rear Axle Position	STD
Set At: Rear - 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3 <input type="checkbox"/> 4 <input type="checkbox"/> - Front		

FRAME COLOR		
<input checked="" type="checkbox"/> 50913	Black Shimmer	STD

FRONT CASTER		
<input checked="" type="checkbox"/> 51046	5" Caster With Heavy Duty Steel Fork	STD

SEAT WIDTH - (Inches)		
<input type="checkbox"/> 50920	20" - 26" Seat Adjustable Width	NC
Set Width At: 20" <input type="checkbox"/> 22" <input type="checkbox"/> 24" <input type="checkbox"/> 26" <input type="checkbox"/>		
<input type="checkbox"/> 50928	26" - 32" Seat Adjustable Width	\$200
Set Width At: 26" <input type="checkbox"/> 28" <input type="checkbox"/> 30" <input type="checkbox"/> 32" <input type="checkbox"/>		

REAR TIRE OPTIONS		
<input checked="" type="checkbox"/> 51453	Urethane Tires	STD

SEAT DEPTH - (Dealer Adjustable)		
<input type="checkbox"/> 50974	19" Seat Depth (Adjustable 16" - 20")	NC
<input type="checkbox"/> 50977	22" Seat Depth (Adjustable 20" - 22")	\$300

AXLE TYPE (5/8" Diameter)		
<input checked="" type="checkbox"/> 51051	Quick Release Rear Axles	NC

SEAT TO FLOOR HEIGHT (Adjustable 14.5" - 20")		
<input checked="" type="checkbox"/> 50958	Front Seat to Floor	NC
Set Height At: 15" <input type="checkbox"/> 16" <input type="checkbox"/> 17" <input type="checkbox"/> 18" <input type="checkbox"/>		
<input checked="" type="checkbox"/> 50965	Rear Seat to Floor	NC
Set Height At: 15" <input type="checkbox"/> 16" <input type="checkbox"/> 17" <input type="checkbox"/> 18" <input type="checkbox"/>		

HAND RIMS		
<input checked="" type="checkbox"/> 51052	Aluminum Hand Rims	STD

ARM TYPE		
<input checked="" type="checkbox"/> 51036	Dual Post Adjustable Height 10" - 14"	NC
<input checked="" type="checkbox"/> 51039	Full Length Armpad 2" x 12"	NC

WHEEL LOCKS & Anti-Tippers		
<input checked="" type="checkbox"/> 51054	Push To Lock Wheel Locks	STD
<input type="checkbox"/> 51056	Wheel Lock Extensions (Pair)	\$60
<input checked="" type="checkbox"/> 51058	Rear Anti-Tippers (Pair) <small>Required safety option.</small>	NC

SEATING OPTIONS		
<input checked="" type="checkbox"/> 50980	Solid Seat Pan	STD

FRONT RIGGING		
<input type="checkbox"/> 51062	Swing Away 60° Front Rigging	NC
<input type="checkbox"/> 51063	Swing Away Elevating Legrests <small>Includes standard calf pads.</small>	\$355
<input type="checkbox"/> 51076	Optional Upgrade Oversize 9" x 10" Calf Pad (ELR Only)	\$200

PDG CUSHIONS		
3" Contour Foam / Texture Cover		
<input type="checkbox"/> 50985	up to 24"	\$398
<input type="checkbox"/> 50986	25 To 32"	\$498

FOOT PLATES		
<input checked="" type="checkbox"/> 51067	Adjustable Angle Footplates	NC

BACK POST STYLE		
<input checked="" type="checkbox"/> 51004	Post Style - Push Bar Back - 8° Bend	STD
<input checked="" type="checkbox"/> 51006	Post Height - 25"	STD
<input checked="" type="checkbox"/> 51015	Eclipse Adjustable Back Upholstery 19" Back Height	STD

FOOTREST/LEGREST ACCESSORIES		
<input type="checkbox"/> 51074	Calf Strap	\$60
<input checked="" type="checkbox"/> 51075	Adjustable Length Heel Loops	NC

EXTENSION TUBE		
<input checked="" type="checkbox"/> 51070	Standard Length (13" - 17")	STD

POSITIONING STRAPS (www.bodypoint.com)		
<input type="checkbox"/> 50904	Bodypoint - Center-Pull Auto Style Hip Belt	\$85

CUSTOM MODIFICATIONS - NOTES	
Total Retail Value - This order	\$
Have all desired specs been checked? A complete spec sheet ensures your chair will arrive on time and as you expect.	