



Eclipse

April, 2011
Price List
Canadian \$

Express

International

P.O. #:	Quote No:
Dealer:	
Address:	
City:	State: Zip Code:
Tel:	Fax:
Tag:	

PDG Eclipse

<input checked="" type="checkbox"/>	50908	PDG Eclipse 600, bariatric wheelchair	\$3,949
	50909	Standard Frame (600 lb (272kg) weight limit)	STD
<input type="checkbox"/>	50910	Heavy Duty Upgrade (1000 lb (454kg) weight limit)	\$150

Frame Color

50913	Black Shimmer	STD
-------	---------------	-----

Seat Width

50920	20" - 26" (50.8cm - 66cm) adjustable width, set at:	NC
50926	20" (50.8cm) 22" (55.9cm) 24" (61cm) 26" (66cm)	
50928	26" - 32" (66cm - 81.3cm) adjustable width, set at:	\$200
50934	26" (66cm) 28" (71.1cm) 30" (76.2cm) 32" (81.3cm)	

Seat Depth (Adjustable by Dealer)

<input type="checkbox"/>	50973	19" (48.3cm) seat depth (Adjustable by Dealer 16" - 20" (40.6cm - 50.8cm))	NC
<input type="checkbox"/>	50977	22" (55.9cm) seat depth (Adjustable by Dealer 20" - 22" (50.8cm - 55.9cm))	\$300

Seat to Floor Height - Front (Adjustable by Dealer 14.5" - 20" (36.2cm - 50.8cm))

50957	Set at: X_15" (38.1cm)	NC
-------	------------------------	----

Seat to Floor Height - Rear (Adjustable by Dealer 14.5" - 20" (36.2cm - 50.8cm))

50964	Set at: X_14.5" (36.8cm)	NC
-------	--------------------------	----

Arm Type

<input checked="" type="checkbox"/>	51036	Dual post adjustable height 10" - 14" (25.4cm - 35.6cm)	NC
-------------------------------------	-------	--	----

Armpad Options

51039	Full Length Armpad 2" X 12" (5cm x 30.5cm)	NC
-------	--	----

Seating Options

50980	Solid seat pan	STD
-------	----------------	-----

PDG Cushions

<input type="checkbox"/>	50985	3" (7.5cm) Contour foam / texture cover (up to 24" (61cm) wide)	\$398
<input type="checkbox"/>	50986	3" (7.5cm) Contour foam / texture cover (25" to 32" (63.5cm - 81.3cm) wide)	\$498

MaTRx Cushions

<input type="checkbox"/>	50998	MaTRx-V cushion: 20"w	\$425
<input type="checkbox"/>	51000	MaTRx-V cushion: 21" w - 22" w	\$475
<input type="checkbox"/>	51002	MaTRx-V cushion: 23" w - 30"w	\$625

Back Post Style

	Push Bar Back - 8° bend	
--	-------------------------	--

Back Post Height

	25" (63.5cm)	STD
--	--------------	-----

Back Upholstery - Adjustable Tension

	Eclipse adjustable back upholstery 19" back height	STD
--	--	-----

Rear Wheel (Adjustable by Dealer)

51047	24" spoked wheel - heavy duty	STD
	Rear axle position set at: Rear 1__ 2__ 3_X_ 4__ Front	

Front Caster

51046	5" caster with heavy duty steel fork	STD
-------	--------------------------------------	-----

Rear Tire Options

	Urethane Tires	STD
--	----------------	-----

Axle Type (5/8" diameter)

<input checked="" type="checkbox"/>	51051	Quick release rear axles	NC
-------------------------------------	-------	--------------------------	----

Handrims

51052	Aluminum hand rims	STD
-------	--------------------	-----

Wheel Locks, Rear Anti-tippers

51054	Push to lock wheel locks	STD	
51056	Wheel lock extensions	NC	
<input checked="" type="checkbox"/>	51058	Rear anti-tippers	NC

Front Rigging

<input type="checkbox"/>	51062	Swing away, removable 60° front rigging	NC
<input type="checkbox"/>	51063	Swing away elevating legrests with oversize 9"x10" calf pads & Large adjustable angle footplates	\$710

Footplate

<input checked="" type="checkbox"/>	51067	Large adjustable angle footplates (included with ELR option)	NC
-------------------------------------	-------	--	----

Extension Tubes

51070	Standard length (13" - 17")	STD
-------	-----------------------------	-----

Footrest / Legrest Accessories

<input type="checkbox"/>	51074	Calf strap	\$60
<input checked="" type="checkbox"/>	51075	Adjustable length heel loops	NC

Positioning Straps (www.bodypoint.com)

<input type="checkbox"/>	51118	BodyPoint - centre-pull auto style hip belt	\$85
--------------------------	-------	---	------

Total Retail Value - This order

This is a required option.
Many other options are available. Visit www.PDGmobility.com to order fully-featured products.