



# Eclipse

<b>PO#</b>	<b>Tag:</b>	<b>Dealer Name:</b>	<b>Dealer Contact:</b>
<b>Mailing Address:</b>		<b>Shipping Address:</b>	
<b>City:</b>	<b>State:</b>	<b>Zip Code:</b>	<b>City:</b>
<b>Tel:</b>	<b>Fax:</b>		<b>Tel:</b>
<b>PDG ECLIPSE 600</b>		<input type="checkbox"/> Confirm Order Via E-mail: YES <input type="checkbox"/> NO <input type="checkbox"/> Send PDG e-News: YES <input type="checkbox"/> NO <input type="text"/> E-Mail:	
<input checked="" type="checkbox"/> 50908	PDG Eclipse 600, Bariatric Wheelchair	K0009	\$2,600
<b>CLIENT INFORMATION - (Ensures optimal set-up of chair)</b>			
<input checked="" type="checkbox"/> 50334	Weight (lbs.) _____	Height (ft./in.) _____	NC
<b>FRAME OPTIONS</b>			
<input type="checkbox"/> 50909	<b>Standard Frame</b> Max Weight Capacity 600 lbs.		STD
<input type="checkbox"/> 50910	<b>Heavy Duty Upgrade</b> Max Weight Capacity 1000 lbs, includes reinforced frame.		\$150
<input type="checkbox"/> 50911	<b>Transit Tie-Down System (TTDS)</b> <small>(TTDS includes 4 factory installed transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details.) Available on seat widths 24" or less.</small>		\$250
<b>FRAME COLOR</b>			
<input type="checkbox"/> 50913	<b>Black Shimmer</b>		STD
<input type="checkbox"/> 50914	<b>Tracer Red</b>		NC
<input type="checkbox"/> 50915	<b>Blue Shimmer</b>		NC
<input type="checkbox"/> 51792	<b>Orange Copper</b>		NC
<input type="checkbox"/> 50917	<b>Misty Purple</b>		NC
<input type="checkbox"/> 50918	<b>Forest Green</b>		NC
<input type="checkbox"/> 50919	<b>Custom Color:</b> _____ (Contact PDG for Specification)		\$300
<b>SEAT WIDTH - (Inches)</b>			
<input type="checkbox"/> 50920	<b>20" - 26" Seat Adjustable Width</b>		\$375
	Set Width At: 20" <input type="checkbox"/> 21" <input type="checkbox"/> 22" <input type="checkbox"/> 23" <input type="checkbox"/> 24" <input type="checkbox"/> 25" <input type="checkbox"/> 26" <input type="checkbox"/>		
<input type="checkbox"/> 50928	<b>26" - 32" Seat Adjustable Width</b>	E2201	\$475
	Set Width At: 26" <input type="checkbox"/> 27" <input type="checkbox"/> 28" <input type="checkbox"/> 29" <input type="checkbox"/> 30" <input type="checkbox"/> 31" <input type="checkbox"/> 32" <input type="checkbox"/>		
<input type="checkbox"/> 50935	<b>30" - 36" Seat Adjustable Width</b> <small>Requires extended delivery time.</small>		\$475
	Set Width At: 30" <input type="checkbox"/> 31" <input type="checkbox"/> 32" <input type="checkbox"/> 33" <input type="checkbox"/> 34" <input type="checkbox"/> 35" <input type="checkbox"/> 36" <input type="checkbox"/>		
<input type="checkbox"/> 50942	<b>36" - 42" Seat Adjustable Width</b> <small>Requires extended delivery time.</small>	E2202	\$475
	Set Width At: 36" <input type="checkbox"/> 37" <input type="checkbox"/> 38" <input type="checkbox"/> 39" <input type="checkbox"/> 40" <input type="checkbox"/> 41" <input type="checkbox"/> 42" <input type="checkbox"/>		
<input type="checkbox"/> 50949	<b>42" - 48" Seat Adjustable Width</b> <small>Requires extended delivery time.</small>		\$475
	Set Width At: 42" <input type="checkbox"/> 43" <input type="checkbox"/> 44" <input type="checkbox"/> 45" <input type="checkbox"/> 46" <input type="checkbox"/> 47" <input type="checkbox"/> 48" <input type="checkbox"/>		
<input type="checkbox"/> 50956	<b>48" Seat or Wider</b>		Ask PDG
<b>SEAT DEPTH - (Adjustable)</b>			
<input type="checkbox"/> 50971	<b>16" Seat Depth</b>		NC
<input type="checkbox"/> 50972	<b>17" Seat Depth</b>		NC
<input type="checkbox"/> 50973	<b>18" Seat Depth</b>		NC
<input type="checkbox"/> 50974	<b>19" Seat Depth</b>		NC
<input type="checkbox"/> 50975	<b>20" Seat Depth</b>	E2203	\$400
<input type="checkbox"/> 50976	<b>21" Seat Depth</b>	E2203	\$400
<input type="checkbox"/> 50977	<b>22" Seat Depth</b>	E2204	\$600
<input type="checkbox"/> 50978	<b>23" Seat Depth</b>	E2204	\$600
<input type="checkbox"/> 50979	<b>24" Seat Depth</b>	E2204	\$600
<b>SEAT FRAME OPTIONS</b>			
<input checked="" type="checkbox"/> 50980	<b>Solid Seat Pan</b>		STD
<b>SEAT TO FLOOR HEIGHT</b>			
<input type="checkbox"/> 50957	<b>Front - Height</b> 14.5" Not Available with TTDS. <small>Oversize footplates required at 14.5" and 15" front height.</small>		NC
	Set At: 14.5" <input type="checkbox"/> 15" <input type="checkbox"/> 16" <input type="checkbox"/> 17" <input type="checkbox"/> 18" <input type="checkbox"/>		
<input type="checkbox"/> 50962	Set At: 19" <input type="checkbox"/> 20" <input type="checkbox"/>		\$200
<input type="checkbox"/> 50964	<b>Rear - Height</b> 14.5" Not Available with TTDS. <small>Oversize footplates required at 14.5" and 15" front height.</small>		NC
	Set At: 14.5" <input type="checkbox"/> 15" <input type="checkbox"/> 16" <input type="checkbox"/> 17" <input type="checkbox"/> 18" <input type="checkbox"/>		
<input type="checkbox"/> 50969	Set At: 19" <input type="checkbox"/> 20" <input type="checkbox"/>		\$200
<b>ARM TYPE</b>			
<input type="checkbox"/> 51035	<b>Dual Post Fixed Height 12"</b>		STD
<input type="checkbox"/> 51036	<b>Dual Post Adjustable Height</b>	E0793	\$245
	Option: 10" - 14" <input type="checkbox"/> 12" - 16" <input type="checkbox"/>		
<b>ARMPAD OPTIONS</b>			
<input type="checkbox"/> 51038	<b>Desk Length Armpad, 2" X 10" (pair)</b>		NC
<input type="checkbox"/> 51039	<b>Full Length Armpad, 2" X 12" (pair)</b>		NC
	<b>3" Wide Arm Pads</b> <small>May interfere with wheels in some tilt configuration.</small>		\$49
<input type="checkbox"/> 51040	Full Length 3"x12" (pair)		
<input type="checkbox"/> 51041	Desk Length 3"x12" (pair)		
	<b>Gel Armpad</b> <small>May interfere with wheels in some tilt configuration.</small>		\$115
<input type="checkbox"/> 51042	Full Length 2"x12" (pair)		
<input type="checkbox"/> 51043	Full Length 3"x12" (pair)		
<input type="checkbox"/> 51044	<b>Extend Armrest 3" Forward</b>		\$49
<input type="checkbox"/> 51045	<b>Forward Guards (extend side support)</b>		\$200
<b>PDG CUSHIONS</b>			
<input type="checkbox"/> 50981	<b>2" Basic Foam / Nylon Cover (Up To 24")</b>		\$168
<b>3" Basic Foam / Nylon Cover</b>			
<input type="checkbox"/> 50982	<b>up to 24"</b>		\$198
<input type="checkbox"/> 50983	<b>25 To 32"</b>		\$298
<input type="checkbox"/> 50984	<b>&gt; 32"</b>		\$398
<b>3" Contour Foam / Texture Cover</b>			
<input type="checkbox"/> 50985	<b>up to 24"</b>		\$398
<input type="checkbox"/> 50986	<b>25 To 32"</b>		\$498
<input type="checkbox"/> 50987	<b>&gt; 32"</b>		\$598
<b>3" Contour Pressure Reduction</b>			
<input type="checkbox"/> 50988	<b>up to 24"</b>		\$598
<input type="checkbox"/> 50989	<b>25 To 32"</b>		\$698
<input type="checkbox"/> 50992	<b>&gt; 32"</b>		\$898
<input type="checkbox"/> 50993	<b>Spare Cushion Cover</b>		\$100
	Cushion Depth:: 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20" <input type="checkbox"/> 22" <input type="checkbox"/>		
<input type="checkbox"/> 51730	<b>Custom Cushion Size</b>		Ask PDG



# Eclipse

PO#	Tag:	Dealer Name:	Dealer Contact:
-----	------	--------------	-----------------

BACK POST STYLE		
<input type="checkbox"/> 51004	Push Bar Back - 8° Bend	STD
<input type="checkbox"/> 51007	Cane Handle Back - 8° Bend	NC

BACK POST HEIGHT		
<input type="checkbox"/> 51006	25"	STD
<input type="checkbox"/> 51005	30"	\$75

BACK OPTIONS			
<input type="checkbox"/> 51008	Seat To Back Angle:		NC
	Set At: 85° <input type="checkbox"/> 90° <input type="checkbox"/> 95° <input type="checkbox"/> 100° <input type="checkbox"/>		
<input type="checkbox"/> 51012	Backrest Extension 10" High <small>Extends Sling Upholstery.</small>	E9066	\$125
<input type="checkbox"/> 51013	30° Reclining Backrest <small>Reclining backrest is 1" narrower than seat width. Back frame is not removable.</small>	E1225	\$830

BACK UPHOLSTERY - (Back height measured from seat pan to top of upholstery)			
<input type="checkbox"/> 51014	No Upholstery		NC
<input type="checkbox"/> 51015	Standard Back Upholstery		STD
	Set Height At: 18" <input type="checkbox"/> 19" <input type="checkbox"/> 20" <input type="checkbox"/>		
<input type="checkbox"/> 51019	Short Upholstery		\$75
	Set Height At: 16" <input type="checkbox"/> 17" <input type="checkbox"/>		
<input type="checkbox"/> 51022	Tall Upholstery <small>Requires 30" back posts.</small>		\$75
	Set Height At: 21" <input type="checkbox"/> 22" <input type="checkbox"/> 23" <input type="checkbox"/>		

REAR WHEEL		
<input checked="" type="checkbox"/> 51047	24" Spoked Wheel - Heavy Duty	STD

FRONT CASTER		
<input checked="" type="checkbox"/> 51046	5" Caster With Heavy Duty Steel Fork	STD

REAR TIRE OPTIONS		
<input type="checkbox"/> 51453	Urethane Tires <small>Urethane Tires recommended for heavy loading.</small>	STD
<input type="checkbox"/> 51049	Pneumatic Tires <small>Pneumatic tires may roll off under heavy loading. Pneumatic tires may result in a finished seat height 1/2" higher than stated.</small>	\$125

AXLE TYPE (5/8" Diameter)		
<input type="checkbox"/> 51050	Fixed Rear Axle	STD
<input type="checkbox"/> 51051	Quick Release Rear Axles	NC
<input type="checkbox"/> 51048	Rear Axle Position	NC
	Set At: Rear - 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3 <input type="checkbox"/> 4 <input type="checkbox"/> - Front	

HAND RIMS		
<input type="checkbox"/> 51052	Aluminum Hand Rims	STD
<input type="checkbox"/> 51053	Plastic Coated Hand Rims	\$250

WHEEL LOCKS & Anti-Tippers			
<input type="checkbox"/> 51054	Push To Lock Wheel Locks		STD
<input type="checkbox"/> 51055	Pull To Lock Wheel Locks		NC
<input type="checkbox"/> 51056	Wheel Lock Extensions (Pair)	E0961	\$60
<input type="checkbox"/> 51057	Attendant Wheel Lock (Push To Lock) <small>Not compatible with Transit Tie Down System.</small>		\$125
<input type="checkbox"/> 51058	Rear Anti-Tippers (Pair) <small>Required safety option.</small>	E0971	\$135

FRONT RIGGING			
<input type="checkbox"/> 51059	Omit Front Rigging		STD
<input type="checkbox"/> 51060	Sector Blocks Only (No Front Rigging Supplied)		NC
<input type="checkbox"/> 51062	Swing Away 60° Front Rigging		NC
<input type="checkbox"/> 51063	Swing Away Elevating Legrests		\$355
<input type="checkbox"/> 51065	Center Mount Front Rigging Bracket <small>Compatible with 60° front rigging only.</small>	K0053	\$260

FOOT PLATES			
<input type="checkbox"/> 51066	Composite Footplates <small>Not compatible with 14" &amp; 15" seat to floor.</small>		STD
<input type="checkbox"/> 51067	Oversize Angle Adjustable Footplates <small>Required (no charge) for 14" &amp; 15" Seat to floor</small>	K0040	\$155
<input type="checkbox"/> 51069	One Piece FLIP-UP Footplates (spans footrests) <small>Not compatible with center mounting front rigging.</small>		\$120
<input type="checkbox"/> 52260	Interlocking Angle Adjustable Footplate <small>20"-21" wide only.</small>		\$395

EXTENSION TUBE		
<input type="checkbox"/> 51070	Standard Length (13" - 17")	STD
<input type="checkbox"/> 51071	Long Extension Tubes (15" - 19")	\$75

FOOTREST/LEGREST ACCESSORIES			
<input type="checkbox"/> 51074	Calf Strap	E0038	\$60
<input type="checkbox"/> 51075	Adjustable Length Heel Loops	E0951	\$50
<input type="checkbox"/> 52261	Swing Away Elevating Legrest Gel Knee Pad		\$99
<input type="checkbox"/> 51076	Oversize 9"x10" Calf Pads <small>(elevating legrests only).</small>		\$200

SEAT WIDTH ADJUSTMENT KITS			
<small>Order the kit if you intend on Growing or Shrinking the chair. Kit includes seat pan, back posts, back upholstery, and width adjustment spacers.</small>			
<input type="checkbox"/> 51077	20" <input type="checkbox"/> 21" <input type="checkbox"/> 22" <input type="checkbox"/> 23" <input type="checkbox"/> 24" <input type="checkbox"/> 25" <input type="checkbox"/> 26" <input type="checkbox"/>		\$500
<input type="checkbox"/> 51078	26" <input type="checkbox"/> 27" <input type="checkbox"/> 28" <input type="checkbox"/> 29" <input type="checkbox"/> 30" <input type="checkbox"/> 31" <input type="checkbox"/> 32" <input type="checkbox"/>		\$520
<input type="checkbox"/> 51079	30" <input type="checkbox"/> 31" <input type="checkbox"/> 32" <input type="checkbox"/> 33" <input type="checkbox"/> 34" <input type="checkbox"/> 35" <input type="checkbox"/> 36" <input type="checkbox"/>		\$540
<input type="checkbox"/> 51080	36" <input type="checkbox"/> 37" <input type="checkbox"/> 38" <input type="checkbox"/> 39" <input type="checkbox"/> 40" <input type="checkbox"/> 41" <input type="checkbox"/> 42" <input type="checkbox"/>		\$560
<input type="checkbox"/> 51081	42" <input type="checkbox"/> 43" <input type="checkbox"/> 44" <input type="checkbox"/> 45" <input type="checkbox"/> 46" <input type="checkbox"/> 47" <input type="checkbox"/> 48" <input type="checkbox"/>		\$580

HEADREST OPTIONS - (The Comfort Company)			
<input type="checkbox"/> 51977	Fully Adjustable Hardware	E1028	\$211
<input type="checkbox"/> 51111	Headrest Mounting Plate		\$75
<input type="checkbox"/> 51979	Moldable Headrest Pad - Stretch-Air <small>Width: 6" <input type="checkbox"/> 8" <input type="checkbox"/> 10" <input type="checkbox"/> 14" <input type="checkbox"/></small>	E0955	\$221

POSITIONING STRAPS (www.bodypoint.com)			
<input type="checkbox"/> 50904	Bodypoint - Center-Pull Auto Style Hip Belt	E0978	\$40
<input type="checkbox"/> 50905	Bodypoint - Padded, Center-Pull Auto Style Hip Belt	E0978	\$90

CUSTOM MODIFICATIONS - NOTES	
Total Retail Value - This order	\$
Have all desired specs been checked? A complete spec sheet ensures your chair will arrive on time and as you expect.	
HCPCS codes supplied are for reference only and do not ensure coverage or payment.	