



# FuzeT20

July 9, 10  
Price List

Express

International

P.O. #:	Quote No:	
Dealer:		
Address:		
City:	Province:	Postal Code:
Tel:	Fax:	
Tag:		
Dealer Contact:		

### This configuration includes the following pre-selected options - No substitutions

<b>Fuze Frame Type - T20 Std type frame</b>		
Fuze T20: 20 Degree Manual Tilt Wheelchair W/ capacity = 250 lb (113.5 kg)		\$3,295
Transit Tie Down Eylets (does not include Lap belt)		STD
<b>Select Seat Width</b>		
<input type="checkbox"/> Seat Width 16" (40.6cm) with 600 N Gas Strut		NC
<input type="checkbox"/> Seat Width 18" (45.7cm) with 1000 N Gas Strut		NC
<input type="checkbox"/> Seat Width 20" (50.8cm) with 1200 N Gas Strut		NC
<b>Seat Depth (Adjustable)</b>		
51134 17" (43.2cm) Seat Depth		STD
<b>Frame Color (Seat Frame Grey)</b>		
51140 Black Shimmer		STD
<b>Arm Type</b>		
<input checked="" type="checkbox"/> 51147 Removable Single Post Adjustable Height Armrest 10" - 14" (25.4cm - 35.6cm)		NC
<b>Armpad Options</b>		
51155 Full Length Armpad, 2" X 12" (pair)		STD
<b>Seat Pan Options</b>		
<input checked="" type="checkbox"/> 51165 Solid Seat Pan		NC
<b>Seat Cushion Options</b>		
51166 No Seat Cushion		STD
<b>Seat Height (Dealer can adjust)</b>		
13" (32.5cm) or 16" (40.0cm) (see rear wheel options)		STD
<b>Caster Options</b> Note: If changing seat height, some settings require different components.		
51174 5" x 1-1/4" Urethane with 20" rear wheels		STD
51179 6" x 1-1/4" Urethane with 22" rear wheels		STD
<b>Rear Wheel Options</b>		
<input type="checkbox"/> 51185 20" rear wheels - Seat height = 13" (33cm) - Adjusts to 17" (43.2cm)		STD
<input type="checkbox"/> 51187 22" rear wheels - Seat height = 16" (40.6cm) - Adjusts 15" to 18" (38.1cm - 45.7cm)		STD
<b>Rear Tire Options</b>		
51189 Urethane Tires		STD
<b>Axle Type</b>		
51192 Quick Release Rear Axles		STD
<b>Hand Rim Options</b>		
51194 Steel Hand Rims		STD
<b>Wheel Locks, Anti-Tippers</b>		
51206 Push To Lock Wheel Locks		STD
<input type="checkbox"/> 51209 Wheel Lock Extensions		\$60
<input checked="" type="checkbox"/> 51210 Rear Anti-Tippers (Pair)		NC
<b>Back Post Style (7/8" Diameter tubing)</b>		
51213 Push Bar Back, 8" bend		STD
<b>Back Post Height (Seat Pan To Handle)</b>		
51216 25" (63.5cm)		STD
51217 Angle/Height Adjustable Push Bar		NC
<b>Back Upholstery</b>		
51226 No Upholstery		STD
<b>Tilt / Recline Options</b>		
252 Attendant Tilt		STD
<b>Front Rigging</b>		
<input type="checkbox"/> Plug-in 70 degree Front Rigging		STD
<input type="checkbox"/> 51257 Swing Away Removable 70° Front Rigging		STD
<input type="checkbox"/> 51259 Swing Away Elevating Legrests		\$355
<b>Foot Plates</b>		
<input type="checkbox"/> 51260 Composite Footplates		STD
<input type="checkbox"/> 51261 Adjustable Angle Footplates		\$140
<b>Footrest/Legrest Accessories</b>		
<input checked="" type="checkbox"/> 51263 Adjustable Length Heel Loops		NC
<b>Extension Tube</b>		
51265 Standard Length (13" - 17")		STD
<b>Positioning Straps</b>		
<input type="checkbox"/> 51267 BodyPoint - center-pull auto style hip belt		\$85

Total Retail Value - This order \$ \_\_\_\_\_

This is a required option.

Many other options are available. Visit [www.pdgmobility.com](http://www.pdgmobility.com) to view and order full-featured products.