

<b>PO#</b>	<b>Tag:</b>	<b>Dealer Name:</b>	<b>Dealer Contact:</b>
<b>Mailing Address:</b>		<b>Shipping Address:</b>	
<b>City:</b>	<b>Province:</b>	<b>Postal Code:</b>	<b>City:</b>
<b>Tel:</b>	<b>Fax:</b>		<b>Fax:</b>

  

<b>Confirm Order Via E-mail:</b> <input type="checkbox"/> YES <input type="checkbox"/> NO		
<b>Send PDG e-News:</b> <input type="checkbox"/> YES <input type="checkbox"/> NO		
<b>E-Mail:</b>		

FUZE FRAME TYPE		
<input checked="" type="checkbox"/> 51120	<b>Fuze T20: 20 Degree Manual Tilt Wheelchair</b> <small>Max Weight Capacity 250 lbs.</small>	\$3,995
<input type="checkbox"/> 51121	<b>Client Information</b> <small>(Ensures Optimal Set-Up Of Chair)</small>	NC
	Weight (Lbs.) _____ Height (Ft./In.) _____	

FUZE FRAME OPTIONS		
<input checked="" type="checkbox"/> 51122	<b>Standard Growing Base &amp; Seat Frame</b> <small>Max Weight Capacity 250 lbs.</small>	STD
<input type="checkbox"/> 51123	<b>Transit Tie-Down System (TTDS)</b> <small>(TTDS includes 4 factory installed transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details.)</small>	NC

GAS STRUT OPTIONS		
<input type="checkbox"/> 51306	<b>600N Gas Strut</b> <small>Recommended for: 15" and 16" Seat Depths. Up to 150 pound client weight.</small>	NC
<input type="checkbox"/> 51307	<b>1000N Gas Strut</b> <small>Recommended for: 17" and 18" Seat Depths. Up to 190 pound client weight.</small>	NC
<input type="checkbox"/> 51308	<b>1200N Gas Strut</b> <small>Recommended for: 19" and 20" Seat Depths. Up to 250 pound client weight.</small>	NC

SEAT WIDTH - (Adjustable) (1" Diameter Tubing)		
<input type="checkbox"/> 51124	<b>13" - 15" Seat Width</b> Set At: 13" <input type="checkbox"/> 14" <input type="checkbox"/> 15" <input type="checkbox"/>	NC
<input type="checkbox"/> 51126	<b>15" - 19" Seat Width</b> Set At: 15" <input type="checkbox"/> 16" <input type="checkbox"/> 17" <input type="checkbox"/> 18" <input type="checkbox"/> 19" <input type="checkbox"/>	NC
<input type="checkbox"/> 51131	<b>20" - 22" Seat Width</b> Set At: 20" <input type="checkbox"/> 21" <input type="checkbox"/> 22" <input type="checkbox"/> <small>Max Weight Capacity 300 lbs. Includes double gas strut (600N)</small>	\$300

SEAT DEPTH - (Dealer Adjustable 15"-20")		
<input type="checkbox"/> 51132	<b>15" Seat Depth</b>	NC
<input type="checkbox"/> 51133	<b>16" Seat Depth</b>	NC
<input type="checkbox"/> 51134	<b>17" Seat Depth</b>	NC
<input type="checkbox"/> 51135	<b>18" Seat Depth</b>	NC
<input type="checkbox"/> 51136	<b>19" Seat Depth</b>	NC
<input type="checkbox"/> 51137	<b>20" Seat Depth</b>	\$300

FRAME COLOR - (Seat Frame Grey)		
<input type="checkbox"/> 51138	<b>Blue Shimmer</b>	STD
<input type="checkbox"/> 51139	<b>Tracer Red</b>	NC
<input type="checkbox"/> 51140	<b>Black Shimmer</b>	NC
<input type="checkbox"/> 51795	<b>Orange Copper</b>	NC
<input type="checkbox"/> 51142	<b>Misty Purple</b>	NC
<input type="checkbox"/> 51143	<b>Forest Green</b>	NC
<input type="checkbox"/> 51144	<b>Custom Color:</b> _____ (Contact PDG for Specification)	\$300

ARM TYPE		
<input type="checkbox"/> 51145	<b>No Armrests</b>	STD
<input type="checkbox"/> 51146	<b>Single Post Fixed Height Armrest - 12"</b>	NC
<input type="checkbox"/> 51147	<b>Single Post Adjustable Height Armrest - 10" - 14"</b>	\$295
<input type="checkbox"/> 51148	<b>Locking Flip-up Armrests</b> <small>Not recommended for clients with excessive strength or tone. Not compatible with fold down back, curved armpad, back upholstery.</small> Set At: 10" <input type="checkbox"/> 11" <input type="checkbox"/> 12" <input type="checkbox"/> 13" <input type="checkbox"/> 14" <input type="checkbox"/>	\$275
	<b>Dual Post Adjustable Height Armrest</b>	\$295
<input type="checkbox"/> 50157	Option: 8" - 11"	
<input type="checkbox"/> 51449	Option: 10" - 14"	
<input type="checkbox"/> 51450	Option: 12" - 16"	

ARMPAD OPTIONS		
<input type="checkbox"/> 51153	<b>Curved Armpad (pair)</b> <small>Not compatible with Lap Tray, Locking Flip-up Armrests or Dual Post Armrests.</small>	STD
<input type="checkbox"/> 51154	<b>Desk Length Armpad, 2" X 10" (pair)</b>	NC
<input type="checkbox"/> 51155	<b>Full Length Armpad, 2" X 12" (pair)</b>	NC
	<b>3" Wide Arm Pads</b> <small>May interfere with wheels in some tilt configuration.</small>	\$49
<input type="checkbox"/> 51156	Full Length 3"x12" (pair)	
<input type="checkbox"/> 51157	Desk Length 3"x12" (pair)	
	<b>Gel Armpad</b> <small>May interfere with wheels in some tilt configuration.</small>	\$115
<input type="checkbox"/> 51158	Full Length 2"x12" (pair)	
<input type="checkbox"/> 51159	Full Length 3.5"x12" (pair)	
	<b>Contour Arm Trough</b> <small>Not compatible with Lap Tray, Locking Flip-up Armrests or Dual Post Armrests. May interfere with wheels in some tilt configuration.</small>	\$85
<input type="checkbox"/> 51161	Left (4.5" X 13")	
<input type="checkbox"/> 51160	Right (4.5" X 13")	
	<b>Flat Hand Pad</b> <small>Must be ordered with contour arm trough.</small>	\$55
<input type="checkbox"/> 51163	Left (4.5" X 5")	
<input type="checkbox"/> 51162	Right (4.5" X 5")	

SEAT FRAME OPTIONS		
<input type="checkbox"/> 51164	<b>No Seat Pan</b>	STD
<input type="checkbox"/> 51165	<b>Solid Seat Pan</b>	\$140

<b>PO#</b>	<b>Tag:</b>	<b>Dealer Name:</b>	<b>Dealer Contact:</b>
------------	-------------	---------------------	------------------------

SEAT CUSHION OPTIONS - (The Comfort Company)		
<input type="checkbox"/> 51166	No Seat Cushion	STD
<input type="checkbox"/> 52227	Cushion M2 Zero Elevation Stretch-Air	\$577
Width: 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20" <input type="checkbox"/>		
Depth: 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20" <input type="checkbox"/>		
<input type="checkbox"/> 52239	Spare Cushion Cover	\$109
<input type="checkbox"/> 52236	Cushion M2 Zero Elevation Stretch-Air (Extra Wide)	\$677
Width: 22" <input type="checkbox"/>		
Depth: 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20" <input type="checkbox"/>		
<input type="checkbox"/> 52248	Spare Cushion Cover (Extra Wide)	\$109
<input type="checkbox"/> 5	Incontinent Cushion Liner	\$85

CASTER OPTIONS		
<input type="checkbox"/> 51174	5" x 1-1/4" Urethane	NC
<input type="checkbox"/> 51175	6" x 1-1/4" Urethane	NC
<input type="checkbox"/> 51176	6" x 1-1/4" Pneumatic	\$120
<input type="checkbox"/> 51451	6" x 1-3/4" Semi-Pneumatic	\$85
<input type="checkbox"/> 51179	8" x 1-1/4" Urethane	NC
<input type="checkbox"/> 51180	8" x 1-1/4" Pneumatic	\$120
<input type="checkbox"/> 51181	8" x 2" Urethane Caster	\$135
<input type="checkbox"/> 51182	8" x 2" Pneumatic Caster	\$135

REAR WHEEL OPTIONS		
<input type="checkbox"/> 51183	12" Composite Wheel	NC
<input type="checkbox"/> 51184	16" Composite Wheel	NC
<input type="checkbox"/> 51185	20" Composite Wheel	NC
<input type="checkbox"/> 51468	20" Spoke Wheel	\$75
<input type="checkbox"/> 51186	22" Composite Wheel	NC
<input type="checkbox"/> 51467	22" Spoke Wheel	\$75
<input type="checkbox"/> 51187	24" Composite Wheel	NC
<input type="checkbox"/> 52274	24" Spoke Wheel	\$75
<input type="checkbox"/> 51188	24" Spoked Wheel - Heavy Duty	\$350
<input type="checkbox"/> 52275	Rear Wheel Camber 5°	NC

REAR TIRE OPTIONS		
<input type="checkbox"/> 51189	Urethane Tires	STD
<input type="checkbox"/> 51190	Pneumatic Tires <small>Pneumatic tires may result in a finished seat height 1/2" higher than stated.</small>	NC
<input type="checkbox"/> 51191	Pneumatic Tires With Flat-Free Inserts <small>Pneumatic tires may result in a finished seat height 1/2" higher than stated. Flat Free Inserts not compatible with Spoked Wheels.</small>	\$200

AXLE TYPE		
<input type="checkbox"/> 51192	Quick Release Rear Axles	STD
<input type="checkbox"/> 51193	Fixed Rear Axles	NC

HAND RIM OPTIONS		
<input type="checkbox"/> 51634	Omit Hand Rims <small>20", 22" or 24" rear wheels only. Not recommended for weights over 195 lbs.</small>	NC
<input type="checkbox"/> 51194	Standard Hand Rims <small>20", 22" or 24" rear wheels only.</small>	STD
<input type="checkbox"/> 51195	Plastic Coated Hand Rims <small>20", 22" or 24" rear wheels only.</small>	\$110
<input type="checkbox"/> 51196	Plastic Coated Projection Hand Rims <small>20", 22" or 24" rear wheels only.</small>	\$250
51198	Number of Projections - 8 Projection Length - 1" Select Projection Angle: 0° <input type="checkbox"/> 45° <input type="checkbox"/>	
51205		
51201		
51202		

SEAT HEIGHT OPTIONS (Standard T20)											
(CIRCLE) desired seat-to-floor height in chart below.											
Rear Wheel	5" Casters					Rear Wheel	6" Casters / 6" Fork				
	12"	13"	-	-	-		12"	-	-	-	-
16"	13"	14"	15"	-	-	16"	14"	15"	-	-	
20"	13"	14"	15"	16"	17"	20"	14"	15"	16"	17"	
22"	-	14"	15"	16"	17"	22"	14"	15"	16"	17"	
24"	-	-	15"	16"	17"	24"	-	15"	16"	17"	
Rear Wheel	6" Casters / 8" Fork					Rear Wheel	8" Casters / 8" Fork				
	16"	15"	-	-	-		16"	-	-	-	-
20"	15"	16"	17"	-	-	20"	16"	17"	-	-	
22"	15"	16"	17"	18"	-	22"	16"	17"	18"	-	
24"	15"	16"	17"	18"	19"	24"	16"	17"	18"	19"	

SEAT HEIGHT OPTIONS (Narrow T20)									
Note: The Narrow T20 is recommended only when overall width is a concern. The overall width is 2" narrower overall than the Stadar T20. Some wheel and seat combinations offer less than 20° tilt, but all offer a minimum of 12° tilt.									
(CIRCLE) desired seat-to-floor height in chart below.									
Rear Wheel	5" Casters				Rear Wheel	6" Casters / 6" Fork			
	12"	-	13"	14"		12"	14"	-	-
16"	14"	15"	16"	-	16"	14"	15"	16"	
20"	16"	17"	-	-	20"	16"	17"	18"	
22"	17"	-	-	-	22"	17"	18"	-	
24"	-	-	-	-	24"	18"	-	-	
Rear Wheel	6" Casters / 8" Fork				Rear Wheel	8" Casters / 8" Fork			
	12"	-	-	-		12"	-	-	-
16"	-	15"	16"	-	16"	16"	-	-	
20"	16"	17"	18"	-	20"	16"	17"	18"	
22"	17"	18"	19"	-	22"	17"	18"	19"	
24"	18"	19"	-	-	24"	18"	19"	20"	

WHEEL LOCKS, ANTI-TIPPERS		
<input type="checkbox"/> 51206	Push To Lock Wheel Locks	STD
<input type="checkbox"/> 51207	Pull To Lock Wheel Locks	NC
<input type="checkbox"/> 51208	Attendant Wheel Lock	No Extention. \$99
<input type="checkbox"/> 51336	User & Attendant Wheel Lock	\$99
<input type="checkbox"/> 51209	Wheel Lock Extensions	\$60
<input checked="" type="checkbox"/> 51210	Rear Anti-Tippers (Pair)	Required safety option. \$99

HEADREST OPTIONS - (The Comfort Company)		
<input type="checkbox"/> 51443	Headrest Mounting Plate	\$75
<input type="checkbox"/> 51967	Fully Adjustable Hardware	\$330
<input type="checkbox"/> 51969	Moldable Headrest Pad - Stretch-Air	\$272
Width: 8" <input type="checkbox"/> 10" <input type="checkbox"/> 14" <input type="checkbox"/>		

BACK POST STYLE (7/8" Diameter Tubing)		
<input type="checkbox"/> 51211	Push Bar Back, Straight	STD
<input type="checkbox"/> 51212	Cane Handle Back, Straight	NC
<input type="checkbox"/> 51213	Push Bar Back, 8° Bend	\$75
<input type="checkbox"/> 51214	Cane Handle Back, 8° Bend	\$75

BACK POST HEIGHT (Seat Pan to Handle)		
<input type="checkbox"/> 51215	20"	STD
<input type="checkbox"/> 51216	25"	STD

