



# Fuze T20

<b>PO#</b>	<b>Tag:</b>	<b>Dealer Name:</b>	<b>Dealer Contact:</b>
<b>Mailing Address:</b>		<b>Shipping Address:</b>	
<b>City:</b>	<b>State:</b>	<b>Zip Code:</b>	<b>City:</b>
<b>Tel:</b>	<b>Fax:</b>	<b>State:</b>	<b>Zip Code:</b>
<b>Weight (Lbs.)</b> _____ <b>Height (Ft./In.)</b> _____		<b>E-Mail:</b> _____	

FUZE FRAME TYPE			
<input checked="" type="checkbox"/> 51120	Fuze T20: 20 Degree Manual Tilt Wheelchair <small>Max Weight Capacity 250 lbs.</small>	K0009	\$2,895
<input type="checkbox"/> 51121	Client Information <small>(Ensures Optimal Set-Up Of Chair)</small>		NC
Weight (Lbs.) _____ Height (Ft./In.) _____			

<b>Confirm Order Via E-mail:</b>	<input type="checkbox"/> YES	<input type="checkbox"/> NO
<b>Send PDG e-News:</b>	<input type="checkbox"/> YES	<input type="checkbox"/> NO

FUZE FRAME OPTIONS			
<input checked="" type="checkbox"/> 51122	Standard Growing Base & Seat Frame <small>Max Weight Capacity 250 lbs.</small>		STD
<input type="checkbox"/> 51123	Transit Tie-Down System (TTDS) <small>(TTDS includes 4 factory installed transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details.)</small>	K0108	\$250

ARM TYPE			
<input type="checkbox"/> 51145	No Armrests		STD
<input type="checkbox"/> 51146	Single Post Fixed Height Armrest - 12"		NC
<input type="checkbox"/> 51147	Single Post Adjustable Height Armrest - 10" - 14"	E0973	\$230
<input type="checkbox"/> 51148	Locking Flip-up Armrests <small>Not recommended for clients with excessive strength or tone. Not compatible with fold down back, curved armpad, back upholstery.</small>	E0973	\$230
Set At: 10" <input type="checkbox"/> 11" <input type="checkbox"/> 12" <input type="checkbox"/> 13" <input type="checkbox"/> 14" <input type="checkbox"/>			
	Dual Post Adjustable Height Armrest	E0973	\$350
<input type="checkbox"/> 50157	Option: 8" - 11"		
<input type="checkbox"/> 51449	Option: 10" - 14"		
<input type="checkbox"/> 51450	Option: 12" - 16"		

GAS STRUT OPTIONS			
<input type="checkbox"/> 51306	600N Gas Strut <small>Recommended for: 15" and 16" Seat Depths. Up to 150 pound client weight.</small>		NC
<input type="checkbox"/> 51307	1000N Gas Strut <small>Recommended for: 17" and 18" Seat Depths. Up to 190 pound client weight.</small>		NC
<input type="checkbox"/> 51308	1200N Gas Strut <small>Recommended for: 19" and 20" Seat Depths. Up to 250 pound client weight.</small>		NC

SEAT WIDTH - (Adjustable) (1" Diameter Tubing)			
<input type="checkbox"/> 51124	13" - 15" Seat Width		NC
Set At: 13" <input type="checkbox"/> 14" <input type="checkbox"/> 15" <input type="checkbox"/>			
<input type="checkbox"/> 51126	15" - 19" Seat Width		NC
Set At: 15" <input type="checkbox"/> 16" <input type="checkbox"/> 17" <input type="checkbox"/> 18" <input type="checkbox"/> 19" <input type="checkbox"/>			
<input type="checkbox"/> 51131	20" - 22" Seat Width	E2201	\$400
Set At: 20" <input type="checkbox"/> 21" <input type="checkbox"/> 22" <input type="checkbox"/> <small>Max Weight Capacity 300 lbs. Includes double gas strut (600N)</small>			

ARMPAD OPTIONS			
<input type="checkbox"/> 51153	Curved Armpad (pair) <small>Not compatible with Lap Tray, Locking Flip-up Armrests or Dual Post Armrests.</small>		STD
<input type="checkbox"/> 51154	Desk Length Armpad, 2" X 10" (pair)		NC
<input type="checkbox"/> 51155	Full Length Armpad, 2" X 12" (pair)		NC
	3" Wide Arm Pads <small>May interfere with wheels in some tilt configuration.</small>		\$49
<input type="checkbox"/> 51156	Full Length 3"x12" (pair)		
<input type="checkbox"/> 51157	Desk Length 3"x12" (pair)		
	Gel Armpad <small>May interfere with wheels in some tilt configuration.</small>		\$115
<input type="checkbox"/> 51158	Full Length 2"x12" (pair)		
<input type="checkbox"/> 51159	Full Length 3.5"x12" (pair)		
	Contour Arm Trough <small>Not compatible with Lap Tray, Locking Flip-up Armrests or Dual Post Armrests. May interfere with wheels in some tilt configuration.</small>	E2290	\$85
<input type="checkbox"/> 51161	Left (4.5" X 13")		
<input type="checkbox"/> 51160	Right (4.5" X 13")		
	Flat Hand Pad <small>Must be ordered with contour arm trough.</small>		\$55
<input type="checkbox"/> 51163	Left (4.5" X 5")		
<input type="checkbox"/> 51162	Right (4.5" X 5")		

SEAT DEPTH - (Dealer Adjustable 15"-20")			
<input type="checkbox"/> 51132	15" Seat Depth		NC
<input type="checkbox"/> 51133	16" Seat Depth		NC
<input type="checkbox"/> 51134	17" Seat Depth		NC
<input type="checkbox"/> 51135	18" Seat Depth		NC
<input type="checkbox"/> 51136	19" Seat Depth		NC
<input type="checkbox"/> 51137	20" Seat Depth	E2203	\$300

FRAME COLOR - (Seat Frame Grey)			
<input type="checkbox"/> 51138	Blue Shimmer		STD
<input type="checkbox"/> 51139	Tracer Red		NC
<input type="checkbox"/> 51140	Black Shimmer		NC
<input type="checkbox"/> 51795	Orange Copper		NC
<input type="checkbox"/> 51142	Misty Purple		NC
<input type="checkbox"/> 51143	Forest Green		NC
<input type="checkbox"/> 51144	Custom Color: _____ (Contact PDG for Specification)		\$300

SEAT FRAME OPTIONS			
<input type="checkbox"/> 51164	No Seat Pan		STD
<input type="checkbox"/> 51165	Solid Seat Pan		\$95

<b>PO#</b>	<b>Tag:</b>	<b>Dealer Name:</b>	<b>Dealer Contact:</b>
------------	-------------	---------------------	------------------------

SEAT CUSHION OPTIONS - (The Comfort Company)			
<input type="checkbox"/> 51166	No Seat Cushion		STD
<input type="checkbox"/> 52227	Cushion M2 Zero Elevation Stretch-Air	E2607	\$297
	Width: 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20" <input type="checkbox"/>		
	Depth: 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20" <input type="checkbox"/>		
<input type="checkbox"/> 52239	Spare Cushion Cover		\$60
<input type="checkbox"/> 52236	Cushion M2 Zero Elevation Stretch-Air (Extra Wide)	E2608	\$398
	Width: 22" <input type="checkbox"/>		
	Depth: 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20" <input type="checkbox"/>		
<input type="checkbox"/> 52248	Spare Cushion Cover (Extra Wide)		\$60
<input type="checkbox"/> 5	Incontinent Cushion Liner		\$58

CASTER OPTIONS		
<input type="checkbox"/> 51174	5" x 1-1/4" Urethane	NC
<input type="checkbox"/> 51175	6" x 1-1/4" Urethane	NC
<input type="checkbox"/> 51176	6" x 1-1/4" Pneumatic	\$75
<input type="checkbox"/> 51451	6" x 1-3/4" Semi-Pneumatic	\$85
<input type="checkbox"/> 51179	8" x 1-1/4" Urethane	NC
<input type="checkbox"/> 51180	8" x 1-1/4" Pneumatic	\$75
<input type="checkbox"/> 51181	8" x 2" Urethane Caster	\$145
<input type="checkbox"/> 51182	8" x 2" Pneumatic Caster	\$145

REAR WHEEL OPTIONS		
<input type="checkbox"/> 51183	12" Composite Wheel	NC
<input type="checkbox"/> 51184	16" Composite Wheel	NC
<input type="checkbox"/> 51185	20" Composite Wheel	NC
<input type="checkbox"/> 51468	20" Spoke Wheel	\$75
<input type="checkbox"/> 51186	22" Composite Wheel	NC
<input type="checkbox"/> 51467	22" Spoke Wheel	\$75
<input type="checkbox"/> 51187	24" Composite Wheel	NC
<input type="checkbox"/> 52274	24" Spoke Wheel	\$75
<input type="checkbox"/> 51188	24" Spoked Wheel - Heavy Duty	\$315
<input type="checkbox"/> 52275	Rear Wheel Camber 5°	NC

REAR TIRE OPTIONS		
<input type="checkbox"/> 51189	Urethane Tires	STD
<input type="checkbox"/> 51190	Pneumatic Tires <small>Pneumatic tires may result in a finished seat height 1/2" higher than stated.</small>	\$60
<input type="checkbox"/> 51191	Pneumatic Tires With Flat-Free Inserts <small>Pneumatic tires may result in a finished seat height 1/2" higher than stated. Flat Free Inserts not compatible with Spoked Wheels.</small>	E2213 \$70

AXLE TYPE		
<input type="checkbox"/> 51192	Quick Release Rear Axles	STD
<input type="checkbox"/> 51193	Fixed Rear Axles	NC

HAND RIM OPTIONS		
<input type="checkbox"/> 51634	Omit Hand Rims <small>20", 22" or 24" rear wheels only. Not recommended for weights over 195 lbs.</small>	NC
<input type="checkbox"/> 51194	Standard Hand Rims <small>20", 22" or 24" rear wheels only.</small>	STD
<input type="checkbox"/> 51195	Plastic Coated Hand Rims <small>20", 22" or 24" rear wheels only.</small>	\$98
<input type="checkbox"/> 51196 51198 51205 51201 51202	Plastic Coated Projection Hand Rims <small>20", 22" or 24" rear wheels only.</small> Number of Projections - 8 Projection Length - 1" Select Projection Angle: 0° <input type="checkbox"/> 45° <input type="checkbox"/>	\$250

SEAT HEIGHT OPTIONS (Standard T20)										
(CIRCLE) desired seat-to-floor height in chart below.										
Rear Wheel	5" Casters					Rear Wheel	6" Casters / 6" Fork			
	12"	13"	-	-	-		12"	-	-	-
16"	13"	14"	15"	-	-	16"	14"	15"	-	-
20"	13"	14"	15"	16"	17"	20"	14"	15"	16"	17"
22"	-	14"	15"	16"	17"	22"	14"	15"	16"	17"
24"	-	-	15"	16"	17"	24"	-	15"	16"	17"
Rear Wheel	6" Casters / 8" Fork					Rear Wheel	8" Casters / 8" Fork			
	16"	15"	-	-	-		16"	-	-	-
20"	15"	16"	17"	-	-	20"	16"	17"	-	-
22"	15"	16"	17"	18"	-	22"	16"	17"	18"	-
24"	15"	16"	17"	18"	19"	24"	16"	17"	18"	19"

SEAT HEIGHT OPTIONS (Narrow T20)									
Note: The Narrow T20 is recommended only when overall width is a concern. The overall width is 2" narrower overall than the Stadard T20. Some wheel and seat combinations offer less than 20° tilt, but all offer a minimum of 12° tilt.									
(CIRCLE) desired seat-to-floor height in chart below.									
Rear Wheel	5" Casters				Rear Wheel	6" Casters / 6" Fork			
	12"	-	13"	14"		12"	14"	-	-
16"	14"	15"	16"	16"	14"	15"	16"		
20"	16"	17"	-	20"	16"	17"	18"		
22"	17"	-	-	22"	17"	18"	-		
24"	-	-	-	24"	18"	-	-		
Rear Wheel	6" Casters / 8" Fork				Rear Wheel	8" Casters / 8" Fork			
	12"	-	-	-		12"	-	-	-
16"	-	15"	16"	16"	16"	-	-		
20"	16"	17"	18"	20"	16"	17"	18"		
22"	17"	18"	19"	22"	17"	18"	19"		
24"	18"	19"	-	24"	18"	19"	20"		

WHEEL LOCKS, ANTI-TIPPERS		
<input type="checkbox"/> 51206	Push To Lock Wheel Locks	STD
<input type="checkbox"/> 51207	Pull To Lock Wheel Locks	NC
<input type="checkbox"/> 51208	Attendant Wheel Lock	No Extention. \$85
<input type="checkbox"/> 51336	User & Attendant Wheel Lock	\$99
<input type="checkbox"/> 51209	Wheel Lock Extensions	E0961 \$60
<input checked="" type="checkbox"/> 51210	Rear Anti-Tippers (Pair)	Required safety option. E0971 \$135

HEADREST OPTIONS - (The Comfort Company)		
<input type="checkbox"/> 51443	Headrest Mounting Plate	\$75
<input type="checkbox"/> 51967	Fully Adjustable Hardware	E1028 \$211
<input type="checkbox"/> 51969	Moldable Headrest Pad - Stretch-Air	E0955 \$221
	Width: 8" <input type="checkbox"/> 10" <input type="checkbox"/> 14" <input type="checkbox"/>	

BACK POST STYLE (7/8" Diameter Tubing)		
<input type="checkbox"/> 51211	Push Bar Back, Straight	STD
<input type="checkbox"/> 51212	Cane Handle Back, Straight	NC
<input type="checkbox"/> 51213	Push Bar Back, 8° Bend	\$75
<input type="checkbox"/> 51214	Cane Handle Back, 8° Bend	\$75

BACK POST HEIGHT (Seat Pan to Handle)		
<input type="checkbox"/> 51215	20"	STD
<input type="checkbox"/> 51216	25"	STD

<b>PO#</b>	<b>Tag:</b>	<b>Dealer Name:</b>	<b>Dealer Contact:</b>
------------	-------------	---------------------	------------------------

BACK POST OPTIONS		
<input type="checkbox"/> 51217	<b>Angle/Height Adjustable Push Bar</b> <small>Not compatible with Cane Handle back.</small>	\$235
<input type="checkbox"/> 51218	<b>Removable Cane Handle Push Bar</b> <small>Not compatible with Push Bar back</small>	\$75
<input type="checkbox"/> 51219	<b>Seat To Back Angle:</b> Set At: 90° <input type="checkbox"/> 96° <input type="checkbox"/> 102° <input type="checkbox"/> 108° <input type="checkbox"/> 114° <input type="checkbox"/> 120° <input type="checkbox"/>	STD
<input type="checkbox"/> 51225	<b>Fold Down Back</b> <small>No 21" or 22" Depth.</small>	NC
<input type="checkbox"/> 52276	<b>Fixed Back</b>	\$75

BACK UPHOLSTERY - (Not compatible with Cane Handle Back)		
<input type="checkbox"/> 51226	<b>No Upholstery</b>	STD
<input type="checkbox"/> 51227	<b>Standard Back Upholstery</b> <small>Back height measured from seat pan to top of upholstery.</small> Set At: 18" <input type="checkbox"/> 19" <input type="checkbox"/> 20" <input type="checkbox"/> 21" <input type="checkbox"/> 22" <input type="checkbox"/> 23" <input type="checkbox"/>	\$125
<input type="checkbox"/> 51234	<b>Adjustable Tension Upholstery</b> <small>K0108</small> Set At: 21" <input type="checkbox"/> 22" <input type="checkbox"/> 23" <input type="checkbox"/>	\$250

BACK OPTIONS A - (The Comfort Company)		
<input type="checkbox"/> 52192	<b>Acta-Back Standard Stretch Air</b> <small>E2613</small> Height: 16" <input type="checkbox"/> 18" <input type="checkbox"/> Width: 15"-17" <input type="checkbox"/> 17"-19" <input type="checkbox"/> 19"-21" <input type="checkbox"/>	\$497
<input type="checkbox"/> 52194	<b>Acta-Back Standard Stretch Air (Extra Tall-20")</b> <small>E2613</small> Width: 15"-17" <input type="checkbox"/> 17"-19" <input type="checkbox"/> 19"-21" <input type="checkbox"/>	\$565

BACK OPTIONS B - (The Comfort Company)		
<input type="checkbox"/> 52208	<b>Acta-Back Stretch Air DEEP</b> <small>E2615</small> Height: 16" <input type="checkbox"/> 18" <input type="checkbox"/> Width: 15"-17" <input type="checkbox"/> 17"-19" <input type="checkbox"/> 19"-21" <input type="checkbox"/>	\$598
<input type="checkbox"/> 52110	<b>Acta-Back Standard Stretch Air DEEP (Extra Tall-20")</b> <small>E2615</small> Width: 15"-17" <input type="checkbox"/> 17"-19" <input type="checkbox"/> 19"-21" <input type="checkbox"/>	\$665

TILT / RECLINE OPTIONS		
<input checked="" type="checkbox"/> 51252	<b>Attendant Tilt Only</b>	STD

FRONT RIGGING		
<input type="checkbox"/> 51254	<b>Omit Front Rigging</b>	NC
<input type="checkbox"/> 51255	<b>Sector Blocks Only (No Front Rigging Supplied)</b>	NC
<input type="checkbox"/> 51257	<b>Swing Away 70° Front Rigging</b>	STD
<input type="checkbox"/> 51452	<b>Swing Away 90° Front Rigging (3"-9")</b> <small>Not compatible with 8" casters.</small>	NC
<input type="checkbox"/> 51259	<b>Swing Away Elevating Legrests</b> <small>E0990</small>	\$355
<input type="checkbox"/> 51481	<b>Plug in 70° Front Rigging</b> <small>Must have frame for plug in front rigging.</small>	NC
<input type="checkbox"/> 52277	<b>Plug in Adjustable Knee Angle Front Rigging</b> <small>Must have frame for plug in front rigging.</small>	\$240

FOOT PLATES		
<input type="checkbox"/> 51260	<b>Composite Footplates</b>	STD
<input type="checkbox"/> 51261	<b>Adjustable Angle Footplates</b> <small>K0040</small>	\$155
<input type="checkbox"/> 51834	<b>Interlocking Angle Adjustable Footplate</b> <small>14"-21" wide only</small>	\$195

FOOTREST/LEGREST ACCESSORIES		
<input type="checkbox"/> 51263	<b>Adjustable Length Heel Loops</b> <small>E0951</small>	\$50
<input type="checkbox"/> 51264	<b>Calf Strap</b> <small>K0038</small>	\$60
<input type="checkbox"/> 52278	<b>Swing Away 70° Calf Pad</b>	\$135
<input type="checkbox"/> 51338	<b>Miller's Dynamic footrest Mounts for Plug in Front Rigging</b>	\$190
<input type="checkbox"/> 52279	<b>Miller's Dynamic footrest Mounts for Swing Away 70° Front Rigging</b>	\$380
<input type="checkbox"/> 52280	<b>Swing Away 70° Gel Padding</b>	\$99
<input type="checkbox"/> 51799	<b>Swing Away Elevating Legrest Gel Knee Pad</b>	\$99

EXTENSION TUBE		
<input type="checkbox"/> 51265	<b>Standard Length (12" - 16")</b> <small>Not compatible with 90°</small>	STD
<input type="checkbox"/> 51266	<b>Long Extension Tubes (15"-19")</b> <small>Not compatible with 90°</small> <small>K0043</small>	\$75
<input type="checkbox"/> 52281	<b>Standard Length (13" - 17") Solid Extension Tubes - Extra Strangth</b> <small>Available on Plug-in front Rigging only.</small>	\$99
<input type="checkbox"/> 52282	<b>Long Extension Tubes (15" - 19") Solid Extension Tubes - Extra Strangth</b> <small>Available on Plug-in front Rigging only.</small>	\$99
<input type="checkbox"/> 51785	<b>90° Extension Tubes (11" - 16")</b> <small>Not compatible with 60°, 70°, ELR.</small>	STD

POSITIONING STRAPS		
<input type="checkbox"/> 51267	<b>BodyPoint - Center-pull auto style hip belt</b> <small>E0978</small>	\$40
<input type="checkbox"/> 51268	<b>BodyPoint - Padded, center-pull auto style hip belt</b> <small>E0978</small>	\$90
<input type="checkbox"/> 51269	<b>BodyPoint - 4 point hip belt</b> <small>E0978</small>	\$110
<input type="checkbox"/> 51270	<b>BodyPoint - Chest Strap</b> <small>E0960</small>	\$75

MISCELLANEOUS OPTIONS		
<input type="checkbox"/> 51489	<b>Clear Lap Tray, slide-on</b>	\$490

CUSTOM MODIFICATIONS - NOTES	
<b>Total Retail Value - This order</b>	\$
Have all desired specs been checked? A complete spec sheet ensures your chair will arrive on time and as you expect. HCPCS codes supplied are for reference only and do not ensure coverage or payment.	