



Fuze T50

PO#	Tag:	Dealer Name:	Dealer Contact:
------------	-------------	---------------------	------------------------

Mailing Address:				Shipping Address:			
City:	State:	Zip Code:		City:	State:	Zip Code:	
Tel:		Fax:		Tel:		Fax:	

FUZE FRAME TYPE			
<input checked="" type="checkbox"/> 50497	Fuze T50: 50 Degree Manual Tilt Wheelchair <small>Max Weight Capacity 250 lbs.</small>	E1161	\$2,895

FUZE FRAME OPTIONS			
<input checked="" type="checkbox"/> 51281	Standard Growing Base & Seat Frame <small>Max Weight Capacity 250 lbs.</small>	STD	
<input type="checkbox"/> 51982	Frame for Plug in Front Rigging		\$300
<input type="checkbox"/> 50728	Transit Tie-Down System (TTDS) <small>(TTDS includes 4 factory installed transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details.)</small>	K0108	\$250

SEAT WIDTH - (Adjustable) (1" Diameter Tubing)			
<input type="checkbox"/> 51769	13" - 15" Seat Width <small>Set At: 13" <input type="checkbox"/> 14" <input type="checkbox"/> 15" <input type="checkbox"/></small>	NC	
<input type="checkbox"/> 51772	15" - 20" Seat Width <small>Set At: 15" <input type="checkbox"/> 16" <input type="checkbox"/> 17" <input type="checkbox"/> 18" <input type="checkbox"/> 19" <input type="checkbox"/></small>	NC	
	<small>Set At: 20" <input type="checkbox"/> 21" <input type="checkbox"/> 22" <input type="checkbox"/></small>	E2201	\$300
<small>21" & 22" Max Weight capacity 350 lbs. Fixed crossbar, non adjustable.</small>			

SEAT DEPTH - (Adjustable)			
<input type="checkbox"/> 50757	15" Seat Depth	NC	
<input type="checkbox"/> 50758	16" Seat Depth	NC	
<input type="checkbox"/> 50759	17" Seat Depth	NC	
<input type="checkbox"/> 50760	18" Seat Depth	NC	
<input type="checkbox"/> 50761	19" Seat Depth	NC	
<input type="checkbox"/> 50762	20" Seat Depth	E2203	\$300
<input type="checkbox"/> 51476	21" Seat Depth <small>Requires 25" Push bar Back Straight, angle/height adjustable push bar handle. Not compatible with Reclining Back.</small>	E2203	\$300
<input type="checkbox"/> 51477	22" Seat Depth <small>Requires 25" Push bar Back Straight, angle/height adjustable push bar handle. Not compatible with Reclining Back.</small>	E2203	\$300

FRAME COLOR - (Seat Frame Grey)			
<input type="checkbox"/> 50763	Blue Shimmer	STD	
<input type="checkbox"/> 50764	Tracer Red	NC	
<input type="checkbox"/> 50765	Black Shimmer	NC	
<input type="checkbox"/> 51794	Orange Copper	NC	
<input type="checkbox"/> 50767	Misty Purple	NC	
<input type="checkbox"/> 50768	Forest Green	NC	
<input type="checkbox"/> 50769	Custom Color: _____ (Contact PDG for Specification)		\$300

ARM TYPE			
<input type="checkbox"/> 50770	No Armrests	STD	
<input type="checkbox"/> 50771	Single Post Fixed Height Armrest - 12"	NC	
<input type="checkbox"/> 50772	Single Post Adjustable Height Armrest - 10" - 14"	E0973	\$230
<input type="checkbox"/> 50773	Locking Flip-up Armrests <small>Not recommended for clients with excessive strength or tone. Not compatible with fold down back, curved armpad, back upholstery.</small>	E0973	\$230
	<small>Set At: 10" <input type="checkbox"/> 11" <input type="checkbox"/> 12" <input type="checkbox"/> 13" <input type="checkbox"/> 14" <input type="checkbox"/></small>		
	Dual Post Adjustable Height Armrest	E0973	\$295
<input type="checkbox"/> 51457	Option: 10" - 14"		
<input type="checkbox"/> 51458	Option: 12" - 16"		

Confirm Order Via E-mail:	<input type="checkbox"/> YES	<input type="checkbox"/> NO
Send PDG e-News:	<input type="checkbox"/> YES	<input type="checkbox"/> NO
E-Mail: _____		

ARMPAD OPTIONS			
<input type="checkbox"/> 50778	Curved Armpad (pair) <small>Not compatible with Lap Tray, Locking Flip-up Armrests or Dual Post Armrests.</small>	STD	
<input type="checkbox"/> 50779	Desk Length Armpad, 2" X 10" (pair)	NC	
<input type="checkbox"/> 50780	Full Length Armpad, 2" X 12" (pair)	NC	
	3" Wide Arm Pads <small>May interfere with wheels in some tilt configuration.</small>		\$49
<input type="checkbox"/> 50781	Full Length 3"x12" (pair)		
<input type="checkbox"/> 50782	Desk Length 3"x12" (pair)		
	Gel Armpad <small>May interfere with wheels in some tilt configuration.</small>		\$115
<input type="checkbox"/> 50783	Full Length 2"x12" (pair)		
<input type="checkbox"/> 50784	Full Length 3.5"x12" (pair)		
	Contour Arm Trough <small>Not compatible with Lap Tray, Locking Flip-up Armrests or Dual Post Armrests. May interfere with wheels in some tilt configuration.</small>		\$85
<input type="checkbox"/> 50786	Left (4.5" X 13")		
<input type="checkbox"/> 50785	Right (4.5" X 13")		
	Flat Hand Pad <small>Must be ordered with contour arm trough.</small>		\$55
<input type="checkbox"/> 50787	Left (4.5" X 5")		
<input type="checkbox"/> 50788	Right (4.5" X 5")		

SEAT FRAME OPTIONS			
<input type="checkbox"/> 50789	No Seat Pan	STD	
<input type="checkbox"/> 50790	Solid Seat Pan	E2231	\$95

SEAT CUSHION OPTIONS - (The Comfort Company)			
<input type="checkbox"/> 50791	No Seat Cushion	STD	
<input type="checkbox"/> 52167	Cushion M2 Zero Elevation Stretch-Air <small>Width: 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20" <input type="checkbox"/></small> <small>Depth: 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20" <input type="checkbox"/></small>	E2607	\$297
<input type="checkbox"/> 52179	Spare Cushion Cover		\$60
<input type="checkbox"/> 52176	Cushion M2 Zero Elevation Stretch-Air (Extra Wide) <small>Width: 22" <input type="checkbox"/></small> <small>Depth: 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20" <input type="checkbox"/></small>	E2608	\$398
<input type="checkbox"/> 52188	Spare Cushion Cover (Extra Wide)		\$60
<input type="checkbox"/> 5	Incontinent Cushion Liner		\$58

CASTER OPTIONS			
<input type="checkbox"/> 50804	5" x 1-1/4" Urethane	NC	
<input type="checkbox"/> 50805	6" x 1-1/4" Urethane	NC	
<input type="checkbox"/> 50806	6" x 1-1/4" Pneumatic		\$75
<input type="checkbox"/> 51334	6" x 1-3/4" Urethane Wider	K0075	\$85
<input type="checkbox"/> 50807	8" x 1-1/4" Urethane	NC	
<input type="checkbox"/> 50808	8" x 1-1/4" Pneumatic		\$75
<input type="checkbox"/> 50809	8" x 2" Urethane Caster		\$145
<input type="checkbox"/> 50810	8" x 2" Pneumatic Caster	K0074	\$145

PO#	Tag:	Dealer Name:	Dealer Contact:
-----	------	--------------	-----------------

REAR WHEEL OPTIONS		
<input type="checkbox"/> 50811	12" Composite Wheel	NC
<input type="checkbox"/> 50812	16" Composite Wheel	NC
<input type="checkbox"/> 50813	20" Composite Wheel	NC
<input type="checkbox"/> 51339	20" Spoke Wheel <small>Not compatible with Pneumatic with Flat Free Inserts.</small>	\$75
<input type="checkbox"/> 50814	22" Composite Wheel	NC
<input type="checkbox"/> 51340	22" Spoke Wheel <small>Not compatible with Pneumatic with Flat Free Inserts.</small>	\$75
<input type="checkbox"/> 50815	24" Composite Wheel	NC
<input type="checkbox"/> 51341	24" Spoke Wheel <small>Not compatible with Pneumatic with Flat Free Inserts.</small>	\$75
<input type="checkbox"/> 50816	24" Spoked Wheel - Heavy Duty <small>Not compatible with Pneumatic with Flat Free Inserts.</small>	\$315
<input type="checkbox"/> 52267	Rear Wheel Camber 5°	NC

SEAT HEIGHT OPTIONS			
(CIRCLE) desired seat-to-floor height in chart below.			
<small>All dimensions are approximate and depend on rear wheel type and pre-tilt setting (shipped standard with 3° pre-tilt). Actual seat height can vary by as much as 1" but can be adjusted down by removing the pre-tilt collar/bumper.</small>			
Rear Wheel	5" Casters		
12"	-	13"	14"
16"	14"	15"	16"
20"	16"	17"	-
22"	17"	-	-
Rear Wheel	6" Casters / 6" Fork		
12"	14"	15"	16"
16"	14"	15"	16"
20"	16"	17"	18"
22"	17"	18"	-
24"	18"	-	-
Rear Wheel	6" Casters / 8" Fork		
12"	18" Inverted		
	19" Inverted		
16"	-	15"	16"
20"	16"	17"	18"
22"	17"	18"	19"
24"	18"	19"	-
Rear Wheel	8" Casters / 8" Fork		
12"	18" Inverted		
	19" Inverted		
16"	20" Inverted		
20"	16"	17"	18"
22"	17"	18"	19"
24"	18"	19"	20"

REAR TIRE OPTIONS		
<input type="checkbox"/> 50817	Urethane Tires	STD
<input type="checkbox"/> 50818	Pneumatic Tires <small>Pneumatic tires may result in a finished seat height 1/2" higher than stated.</small>	\$60
<input type="checkbox"/> 50819	Pneumatic Tires With Flat-Free Inserts <small>Pneumatic tires may result in a finished seat height 1/2" higher than stated. Flat Free Inserts not compatible with Spoked Wheels.</small>	E2213 \$70

AXLE TYPE		
<input type="checkbox"/> 50820	Quick Release Rear Axles	STD
<input type="checkbox"/> 50821	Fixed Rear Axles	NC

HAND RIM OPTIONS		
<input type="checkbox"/> 51633	Omit Hand Rims	20", 22" or 24" rear wheels only. NC
<input type="checkbox"/> 50822	Standard Hand Rims	20", 22" or 24" rear wheels only. NC
<input type="checkbox"/> 50823	Plastic Coated Hand Rims	20", 22" or 24" rear wheels only. \$110
<input type="checkbox"/> 50824 50825 50829 50830 50831	Plastic Coated Projection Hand Rims Number of Projections - 8 Projection Length - 1" Select Projection Angle: 0° <input type="checkbox"/> 45° <input type="checkbox"/>	\$250

WHEEL LOCKS, ANTI-TIPPERS		
<input type="checkbox"/> 51976	Foot Operated Wheel Locks <small>Not compatible with 1 arm drive.</small>	\$200
<input type="checkbox"/> 50833	Push To Lock Wheel Locks	STD
<input type="checkbox"/> 50834	Pull To Lock Wheel Locks	NC
<input type="checkbox"/> 50835	Attendant Wheel Locks <small>Not compatible with IV Pole and O2 Holder, No Extension.</small>	\$99
<input type="checkbox"/> 51335	Attendant & User Wheel Locks <small>Not compatible with IV Pole and O2 Holder, No Extension.</small>	\$99
<input type="checkbox"/> 50836	Wheel Lock Extensions	E0961 \$60
<input checked="" type="checkbox"/> 50837	Rear Anti-Tippers (Pair) <small>Required safety option.</small>	E0971 \$135

HEADREST OPTIONS - (The Comfort Company)		
<input type="checkbox"/> 51303	Headrest Mounting Plate	\$75
<input type="checkbox"/> 51934	Fully Adjustable Hardware	E1028 \$211
<input type="checkbox"/> 51936	Moldable Headrest Pad - Stretch-Air Width: 8" <input type="checkbox"/> 10" <input type="checkbox"/> 14" <input type="checkbox"/>	E0955 \$221

BACK UPHOLSTERY		
<input type="checkbox"/> 50853	No Upholstery	STD
<input type="checkbox"/> 50854	Standard Back Upholstery <small>Back height measured from seat pan to top of upholstery.</small> Set At: 18" <input type="checkbox"/> 19" <input type="checkbox"/> 20" <input type="checkbox"/> 21" <input type="checkbox"/> 22" <input type="checkbox"/> 23" <input type="checkbox"/>	\$125
<input type="checkbox"/> 50861	Adjustable Tension Upholstery Set At: 21" <input type="checkbox"/> 22" <input type="checkbox"/> 23" <input type="checkbox"/>	K0108 \$250

BACK OPTIONS A - (The Comfort Company)		
<input type="checkbox"/> 52132	Acta-Back Standard Stretch Air Height: 16" <input type="checkbox"/> 18" <input type="checkbox"/> Width: 15"-17" <input type="checkbox"/> 17"-19" <input type="checkbox"/> 19"-21" <input type="checkbox"/> Width: 21"-23" <input type="checkbox"/>	E2613 \$497 E2614 \$584
<input type="checkbox"/> 52134	Acta-Back Standard Stretch Air (Extra Tall-20") Width: 15"-17" <input type="checkbox"/> 17"-19" <input type="checkbox"/> 19"-21" <input type="checkbox"/> Width: 21"-23" <input type="checkbox"/>	E2613 \$565 E2614 \$640

BACK OPTIONS B - (The Comfort Company)		
<input type="checkbox"/> 52184	Acta-Back Standard Stretch Air DEEP Height: 16" <input type="checkbox"/> 18" <input type="checkbox"/> Width: 15"-17" <input type="checkbox"/> 17"-19" <input type="checkbox"/> 19"-21" <input type="checkbox"/> Width: 21"-23" <input type="checkbox"/>	E2615 \$598 E2616 \$664
<input type="checkbox"/> 52150	Acta-Back Standard Stretch Air DEEP (Extra Tall-20") Width: 15"-17" <input type="checkbox"/> 17"-19" <input type="checkbox"/> 19"-21" <input type="checkbox"/> Width: 21"-23" <input type="checkbox"/>	E2615 \$665 E2616 \$739

BACK POST STYLE (7/8" Diameter Tubing)		
<input type="checkbox"/> 50838	Push Bar Back, Straight	STD
<input type="checkbox"/> 50839	Cane Handle Back, Straight	NC
<input type="checkbox"/> 50840	Push Bar Back, 8° Bend	\$75
<input type="checkbox"/> 50841	Cane Handle Back, 8° Bend	\$75

PO#	Tag:	Dealer Name:	Dealer Contact:
-----	------	--------------	-----------------

BACK POST HEIGHT (Seat Pan to Handle)		
<input type="checkbox"/> 50842	20"	STD
<input type="checkbox"/> 50843	25"	STD

BACK POST OPTIONS		
<input type="checkbox"/> 50844	Angle/Height Adjustable Push Bar <small>Not compatible with Cane Handle back.</small>	\$235
<input type="checkbox"/> 50845	Removable Cane Handle Push Bar <small>Not compatible with Push Bar back.</small>	\$75
<input type="checkbox"/> 50846	Seat To Back Angle: Set At: 90° <input type="checkbox"/> 96° <input type="checkbox"/> 102° <input type="checkbox"/> 108° <input type="checkbox"/> 114° <input type="checkbox"/> 120° <input type="checkbox"/>	STD
<input type="checkbox"/> 50852	Fold Down Back <small>No 21" or 22" Depth.</small>	NC
<input type="checkbox"/> 52251	Fixed Back	\$75

TILT / RECLINE OPTIONS		
<input checked="" type="checkbox"/> 50879 50884	Tilt Limiting Device - Limits Tilt At < 50° Set Tilt Limit At: 20° <input type="checkbox"/> 25° <input type="checkbox"/> 30° <input type="checkbox"/> 35° <input type="checkbox"/> 40° <input type="checkbox"/> 45° <input type="checkbox"/> 50° <input type="checkbox"/>	NC
<input type="checkbox"/> 51309 51324 51323	Power Tilt <small>Tilt range 4° to 45°</small> Toggle Switch Left <input type="checkbox"/> Right <input type="checkbox"/> Button Switch Left <input type="checkbox"/> Right <input type="checkbox"/>	E1002 \$3,400
<input type="checkbox"/> 51342	Reclining Backrest(Manual) Max. Recline - 30° <small>Not compatible with fold down back or cane handle back. Requires 25" Angle / Height Adjustable Push Bar with 8" bend. Not Available with 15", 16", 21", or 22" seat depth. If upholstery is ordered with recline it must be adjustable tension.</small>	E1225 \$830

FRONT RIGGING		
<input type="checkbox"/> 50891	Omit Front Rigging	NC
<input type="checkbox"/> 50892	Sector Blocks Only (No Front Rigging Supplied)	NC
<input type="checkbox"/> 50894	Swing Away 70° Front Rigging	STD
<input type="checkbox"/> 51343	Swing Away 90° Front Rigging (3"-9") <small>Not compatible with 8" casters.</small>	NC
<input type="checkbox"/> 50896	Swing Away Elevating Legrests	E0990 \$355
<input type="checkbox"/> 51983	Plug in 70° Front Rigging <small>Must have frame for plug in front rigging.</small>	NC
<input type="checkbox"/> 51984	Plug in Adjustable Knee Angle Front Rigging <small>Must have frame for plug in front rigging.</small>	\$240

FOOT PLATES		
<input type="checkbox"/> 50898	Adjustable Angle Footplates	K0040 \$155
<input type="checkbox"/> 51474	One Piece Flip-up Footplates	\$395
<input type="checkbox"/> 50897	Compostite Footplates	STD
<input type="checkbox"/> 51765	Interlocking Angle Adjustable Footplate <small>14"-21" wide only.</small>	\$195

FOOTREST/LEGREST ACCESSORIES			
<input type="checkbox"/> 50900	Adjustable Length Heel Loops	E0951	\$50
<input type="checkbox"/> 50901	Calf Strap	K0038	\$60
<input type="checkbox"/> 52270	Swing Away 70° Calf Pad		\$135
<input type="checkbox"/> 51337	Miller's Dynamic footrest Mounts for Swing Away 70° Front Rigging		\$380
<input type="checkbox"/> 52271	Swing Away 70° Gel Padding		\$99
<input type="checkbox"/> 51791	Swing Away Elevating Legrest Gel Knee Pad		\$99

EXTENSION TUBE			
<input type="checkbox"/> 50902	Standard Length (12" - 16") <small>Not compatible with 90°</small>		STD
<input type="checkbox"/> 50903	Long Extension Tubes (15" - 19") <small>Not compatible with 90°</small>	K0043	\$75
<input type="checkbox"/> 52268	Standard Length (13" - 17") Solid Extension Tubes - Extra Strangth <small>Available on Plug-in front Rigging only.</small>		\$99
<input type="checkbox"/> 52269	Long Extension Tubes (15" - 19") Solid Extension Tubes - Extra Strangth <small>Available on Plug-in front Rigging only.</small>		\$99
<input type="checkbox"/> 51783	90° Extension Tubes (11" - 16") <small>Not compatible with 60°, 70°, ELR.</small>		STD

POSITIONING STRAPS			
<input type="checkbox"/> 51118	BodyPoint - Center-pull auto style hip belt	E0978	\$40
<input type="checkbox"/> 51119	BodyPoint - Padded, center-pull auto style hip belt	E0978	\$90
<input type="checkbox"/> 50906	BodyPoint - 4 point hip belt	E0978	\$110
<input type="checkbox"/> 50907	BodyPoint - Chest Strap	E0978	\$75

MISCELLANEOUS OPTIONS			
<input type="checkbox"/> 51470	Stealth Flip-Down Abductor	E0957	\$489
<input type="checkbox"/> 51319	Quick Release Bracket - Accessory Mount <small>Not Compatible with Transit Tie Down System (TTDS).</small>	E1028	\$200
<input type="checkbox"/> 51986	Hip Gel Guides <small>Not Compatible with Dual Post Arms.</small>		\$325
<input type="checkbox"/> 51344	One Arm Drive Left <input type="checkbox"/> Right <input type="checkbox"/> <small>Requires 24" composite rear wheel. Forward Wheel Position. Available with 18" and 20" seat-to-floor only. Max 22" wide.</small>	E0958	\$3,500
<input type="checkbox"/> 51345	IV Pole	E1029	\$350
<input type="checkbox"/> 51346	Ventilator Tray <small>Minimum seat width is 16". Reduces ground clearance. Restricted to 20", 22", and 24" rear wheels only. May become unstable when no occupant in the chair. Not compatatible with O2 holder.</small>	E1029	\$295
<input type="checkbox"/> 51347	Ventilator Bracket <small>Compatible with LTV1 (lap topce) ventilator.</small>	E1029	\$95
<input type="checkbox"/> 51348	Ventilator Battery Tray <small>Minimum seat width is 16". Reduces ground clearance. Fits U1 Battery Only. Restricted to 20", 22", and 24" rear wheels, with seat height restrictions.</small>	E1029	\$95
<input type="checkbox"/> 51428	O2 Holder <small>May add to overall width.</small>	E2208	\$350
<input type="checkbox"/> 51475	Clear Lap Tray, slide-on		\$490

CUSTOM MODIFICATIONS - NOTES	
Total Retail Value - This order	\$
Have all desired specs been checked? A complete spec sheet ensures your chair will arrive on time and as you expect. HCPCS codes supplied are for reference only and do not ensure coverage or payment.	