



Fuze T50 Jr

PO#	Tag:	Dealer Name:	Dealer Contact:
Mailing Address:		Shipping Address:	
City:	State:	Zip Code:	
Tel:	Fax:		
FUZE FRAME TYPE <input checked="" type="checkbox"/> 51491 Fuze T50JR: 50 Degree Manual Tilt Wheelchair E1233 \$2,850		Confirm Order Via E-mail: <input type="checkbox"/> YES <input type="checkbox"/> NO Send PDG e-News: <input type="checkbox"/> YES <input type="checkbox"/> NO E-Mail: _____	

FUZE FRAME OPTIONS			
<input checked="" type="checkbox"/> 51492	Standard Growing Base & Seat Frame <small>Max Weight Capacity 180 lbs.</small>		STD
<input type="checkbox"/> 51493	Transit Tie-Down System (TTDS) <small>(TTDS includes 4 factory installed transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details.)</small>	K0108	\$250

SEAT WIDTH - (Adjustable) (1" Diameter Tubing)			
<input type="checkbox"/> 51494	12" - 15" Seat Width		NC
Set At: 12" <input type="checkbox"/> 13" <input type="checkbox"/> 14" <input type="checkbox"/> 15" <input type="checkbox"/>			
<input type="checkbox"/> 51498	15" - 18" Seat Width		NC
Set At: 15" <input type="checkbox"/> 16" <input type="checkbox"/> 17" <input type="checkbox"/> 18" <input type="checkbox"/>			

SEAT DEPTH - (Adjustable)			
<input type="checkbox"/> 51502	13" Seat Depth		NC
<input type="checkbox"/> 51503	14" Seat Depth		NC
<input type="checkbox"/> 51504	15" Seat Depth		NC
<input type="checkbox"/> 51505	16" Seat Depth		NC
<input type="checkbox"/> 51506	17" Seat Depth		NC
<input type="checkbox"/> 51507	18" Seat Depth	E2203	NC

FRAME COLOR - (Seat Frame Grey)			
<input type="checkbox"/> 51508	Blue Shimmer		NC
<input type="checkbox"/> 51509	Tracer Red		NC
<input type="checkbox"/> 51510	Black Shimmer		NC
<input type="checkbox"/> 51796	Orange Copper		NC
<input type="checkbox"/> 51512	Misty Purple		NC
<input type="checkbox"/> 51619	Forest Green		NC
<input type="checkbox"/> 51513	Custom Color: _____ (Contact PDG for Specification)		\$300

ARM TYPE			
<input type="checkbox"/> 51514	No Armrests		STD
<input type="checkbox"/> 51515	Single Post Fixed Height Armrest - 12"		NC
<input type="checkbox"/> 51516	Single Post Adjustable Height Armrest - 10" - 14"	E0973	\$250
<input type="checkbox"/> 51519	Locking Flip-up Armrests <small>Not recommended for clients with excessive strength or tone. Not compatible with fold down back, curved armpad, back upholstery.</small>	E0973	\$250
Set At: 8" <input type="checkbox"/> 9" <input type="checkbox"/> 10" <input type="checkbox"/> 11" <input type="checkbox"/> 12" <input type="checkbox"/> 13" <input type="checkbox"/> 14" <input type="checkbox"/>			
	Dual Post Adjustable Height Armrest	E0973	\$295
<input type="checkbox"/> 51527	Option: 8" - 11"		
<input type="checkbox"/> 51528	Option: 10" - 14"		
<input type="checkbox"/> 51529	Option: 12" - 16"		

ARMPAD OPTIONS			
<input type="checkbox"/> 51530	Curved Armpad (pair) <small>Not compatible with Lap Tray, Locking Flip-up Armrests or Dual Post Armrests.</small>		STD
<input type="checkbox"/> 51531	Desk Length Armpad, 2" X 10" (pair)		NC
<input type="checkbox"/> 51532	Full Length Armpad, 2" X 12" (pair)		NC
	3" Wide Arm Pads <small>May interfere with wheels in some tilt configuration.</small>		\$49
<input type="checkbox"/> 51533	Full Length 3"x12" (pair)		
<input type="checkbox"/> 51534	Desk Length 3"x10" (pair)		
	Gel Armpad <small>May interfere with wheels in some tilt configuration.</small>		\$115
<input type="checkbox"/> 51535	Full Length 2"x12" (pair)		
<input type="checkbox"/> 51536	Full Length 3.5"x12" (pair)		
	Contour Arm Trough <small>Not compatible with Lap Tray, Locking Flip-up Armrests or Dual Post Armrests. May interfere with wheels in some tilt configuration.</small>	E2209	\$85
<input type="checkbox"/> 51537	Left (4.5" X 13")		
<input type="checkbox"/> 51538	Right (4.5" X 13")		
	Flat Hand Pad <small>Must be ordered with contour arm trough.</small>		\$55
<input type="checkbox"/> 51539	Left (4.5" X 5")		
<input type="checkbox"/> 51540	Right (4.5" X 5")		

SEAT FRAME OPTIONS			
<input type="checkbox"/> 51541	No Seat Pan		STD
<input type="checkbox"/> 51542	Solid Seat Pan	E2231	\$95

CASTER OPTIONS			
<input type="checkbox"/> 51544	5" x 1-1/4" Urethane		NC
<input type="checkbox"/> 51545	6" x 1-1/4" Urethane		NC
<input type="checkbox"/> 51546	6" x 1-1/4" Pneumatic		\$75
<input type="checkbox"/> 51547	6" x 1-3/4" Semi-Pneumatic	K0075	\$85
<input type="checkbox"/> 51548	8" x 1-1/4" Urethane		NC
<input type="checkbox"/> 51549	8" x 1-1/4" Pneumatic		\$75
<input type="checkbox"/> 51550	8" x 2" Urethane Caster	K0074	\$145
<input type="checkbox"/> 51551	8" x 2" Pneumatic Caster		\$145

PO#	Tag:	Dealer Name:	Dealer Contact:
------------	-------------	---------------------	------------------------

REAR WHEEL OPTIONS		
<input type="checkbox"/> 51552	12" Composite Wheel	NC
<input type="checkbox"/> 51553	16" Composite Wheel	NC
<input type="checkbox"/> 51554	20" Composite Wheel	NC
<input type="checkbox"/> 51555	20" Spoke Wheel	\$75
<input type="checkbox"/> 51556	22" Composite Wheel	NC
<input type="checkbox"/> 51557	22" Spoke Wheel	\$75
<input type="checkbox"/> 51558	24" Composite Wheel	NC
<input type="checkbox"/> 51559	24" Spoke Wheel	\$75
<input type="checkbox"/> 5	Rear Wheel Camber 5°	NC

SEAT HEIGHT OPTIONS

(CIRCLE) desired seat-to-floor height in chart below.

All dimensions are approximate and depend on rear wheel type and pre-tilt setting (shipped standard with 3° pre-tilt). Actual seat height can vary by as much as 1" but can be adjusted down by removing the pre-tilt collar/bumper.

Rear Wheel	5" Casters				Rear Wheel	6" Casters / 6" Fork			
12"	-	13"	14"		12"	14"	18" Inverted		
16"	14"	15"	16"		16"	14"	15"	16"	-
20"	16"	17"	-		20"	16"	17"	18"	-
22"	17"	-	-		22"	17"	18"	-	-
					24"	18"	-	-	-

Rear Wheel	6" Casters / 8" Fork				Rear Wheel	8" Casters / 8" Fork			
12"	18" Inverted				12"	18" Inverted			
	19" Inverted					19" Inverted			
16"	-	15"	16"		16"	20" Inverted			
20"	16"	17"	18"		20"	16"	17"	18"	
22"	17"	18"	19"		22"	17"	18"	19"	
24"	18"	19"	-		24"	18"	19"	20"	

REAR TIRE OPTIONS		
<input type="checkbox"/> 51560	Urethane Tires	STD
<input type="checkbox"/> 51561	Pneumatic Tires <small>Pneumatic tires may result in a finished seat height 1/2" higher than stated.</small>	\$60
<input type="checkbox"/> 51562	Pneumatic Tires With Flat-Free Inserts <small>Pneumatic tires may result in a finished seat height 1/2" higher than stated. Flat Free Inserts not compatible with Spoked Wheels.</small>	E2213 \$95

AXLE TYPE		
<input type="checkbox"/> 51563	Quick Release Rear Axles	STD
<input type="checkbox"/> 51564	Fixed Rear Axles	NC

HAND RIM OPTIONS		
<input type="checkbox"/> 51635	Omit Hand Rims	NC
<input type="checkbox"/> 51565	Standard Hand Rims	STD
<input type="checkbox"/> 52283	Plastic Coated Hand Rims	\$95
<input type="checkbox"/> 52284 52285 52286 52288	Plastic Coated Projection Hand Rims <small>20", 22" or 24" rear wheels only.</small>	\$250
	Number of Projections - 8	
	Projection Length - 1"	
	Select Projection Angle: 0° <input type="checkbox"/> 45° <input type="checkbox"/>	

WHEEL LOCKS, ANTI-TIPPERS		
<input type="checkbox"/> 51798	Foot Operated Wheel Locks <small>Not compatible with 1 arm drive.</small>	\$200
<input type="checkbox"/> 51567	Push To Lock Wheel Locks	STD
<input type="checkbox"/> 51568	Pull To Lock Wheel Locks	\$85
<input type="checkbox"/> 51569	Attendant Wheel Lock <small>No Extension.</small>	\$99
<input type="checkbox"/> 51570	Wheel Lock Extensions <small>E0961</small>	\$70
<input checked="" type="checkbox"/> 51571	Rear Anti-Tippers (Pair) <small>Required safety option. E0971</small>	\$135

BACK POST STYLE (7/8" Diameter Tubing)		
<input type="checkbox"/> 51572	Push Bar Back, Straight	STD
<input type="checkbox"/> 51573	Cane Handle Back, Straight <small>Not compatible with back upholstery.</small>	NC
<input type="checkbox"/> 51574	Push Bar Back, 8° Bend	\$75
<input type="checkbox"/> 51575	Cane Handle Back, 8° Bend <small>Not compatible with back upholstery.</small>	\$75

BACK POST HEIGHT (Seat Pan to Handle)		
<input type="checkbox"/> 51576	20"	STD
<input type="checkbox"/> 51577	25"	NC

BACK POST OPTIONS		
<input type="checkbox"/> 51578	Angle/Height Adjustable Push Bar <small>Not compatible with Cane Handle back.</small>	\$235
<input type="checkbox"/> 51579	Removable Cane Handle Push Bar <small>Not compatible with Push Bar back.</small>	\$75
<input type="checkbox"/> 51580	Seat To Back Angle: <small>Set At: 90° <input type="checkbox"/> 96° <input type="checkbox"/> 102° <input type="checkbox"/> 108° <input type="checkbox"/> 114° <input type="checkbox"/> 120° <input type="checkbox"/></small>	STD
<input type="checkbox"/> 51586	Fold Down Back	STD
<input type="checkbox"/> 52289	Fixed Back	\$75

HEADREST OPTIONS - (The Comfort Company)		
<input type="checkbox"/> 51588	Headrest Mounting Plate	\$75
<input type="checkbox"/> 52254	Fully Adjustable Hardware <small>E1028</small>	\$211
<input type="checkbox"/> 52255	Moldable Headrest Pad - Stretch-Air <small>E0955</small> <small>Width: 8" <input type="checkbox"/> 10" <input type="checkbox"/> 14" <input type="checkbox"/></small>	\$221

TILT / RECLINE OPTIONS		
<input checked="" type="checkbox"/> 51620	Tilt Limiting Device - Limits Tilt At < 50° <small>Set Tilt Limit At: 20° <input type="checkbox"/> 25° <input type="checkbox"/> 30° <input type="checkbox"/> 35° <input type="checkbox"/> 35° <input type="checkbox"/> 40° <input type="checkbox"/> 45° <input type="checkbox"/> 50° <input type="checkbox"/></small>	STD

FRONT RIGGING		
<input type="checkbox"/> 51602	Omit Front Rigging	NC
<input type="checkbox"/> 51603	Sector Blocks Only (No Front Rigging Supplied)	NC
<input type="checkbox"/> 52290	External Front Rigging Mount	\$75
<input type="checkbox"/> 51604	Swing Away 70° Front Rigging	NC
<input type="checkbox"/> 51605	Swing Away 90° Front Rigging (3"-9") <small>Not compatible with 8" casters.</small>	STD
<input type="checkbox"/> 51606	Swing Away Elevating Legrests <small>E0990</small> <small>Not available on 12" or 13" wide.</small>	\$295



Fuze T50 Jr

PO#	Tag:	Dealer Name:	Dealer Contact:
-----	------	--------------	-----------------

FOOT PLATES		
<input type="checkbox"/> 51607	Composite Footplates <small>Not compatible with 12" wide chair.</small>	STD
<input type="checkbox"/> 51608	Adjustable Angle Footplates K0040	\$175
<input type="checkbox"/> 51609	One Piece Flip-up Footplates	\$395
<input type="checkbox"/> 51839	Interlocking Angle Adjustable Footplate 14"-21" wide only	\$200

FOOTREST/LEGREST ACCESSORIES		
<input type="checkbox"/> 51610	Adjustable Length Heel Loops E0951	\$50
<input type="checkbox"/> 51611	Calf Strap K0038	\$60
<input type="checkbox"/> 52291	Swing Away 70° Calf Pad <small>Not available on 12" or 13" wide.</small>	\$135
<input type="checkbox"/> 52292	Miller's Dynamic footrest Mounts for Swing Away 70° Front Rigging	\$380
<input type="checkbox"/> 52293	Swing Away 70° Gel Padding	\$99
<input type="checkbox"/> 51806	Swing Away Elevating Legrest Gel Knee Pad	\$99

EXTENSION TUBE		
<input type="checkbox"/> 51613	Standard Length (13" - 17") <small>Not compatible with 90°</small>	STD
<input type="checkbox"/> 51614	Long Extension Tubes (15" - 19") <small>Not compatible with 90°</small> K0043	\$75
<input type="checkbox"/> 51784	90° Extension Tubes (11" - 16") <small>Not compatible with 60°, 70°, ELR.</small>	NC

POSITIONING STRAPS		
<input type="checkbox"/> 51615	BodyPoint - Center-pull auto style hip belt E0978	\$49
<input type="checkbox"/> 51616	BodyPoint - Padded, center-pull auto style hip belt E0978	\$90
<input type="checkbox"/> 51617	BodyPoint - 4 point hip belt E0978	\$110
<input type="checkbox"/> 51618	BodyPoint - Chest Strap E0960	\$75

MISCELLANEOUS OPTIONS		
<input type="checkbox"/> 51797	Clear lap tray, slide on.	\$490

CUSTOM MODIFICATIONS - NOTES	
Total Retail Value - This order	\$
Have all desired specs been checked? A complete spec sheet ensures your chair will arrive on time and as you expect.	
HCPCS codes supplied are for reference only and do not ensure coverage or payment.	