

ACCESSORIES



Fuze 21/22" Seat Width

- Weight capacity – 350 lbs
- Fixed cross bar



Fuze Bodypoint Lapbelts

- Made for PDG
- Allows for on-chair, easy adjustment



Fuze Fixed Back on T50 / T20



Fuze Foot-operated Wheel Lock

- Compatible with all rear wheels
- Step-on function
- Limits forward wheel position



Fuze Vent Tray with Separate Battery Tray

- Minimum seat width – 16"
- Accommodates U1 Gel battery



Fuze Inverted Axle Plate

- Allows for higher seat heights c/w small rear wheels
- 12" rear wheel up to 20" seat height



Fuze Hip Guides c/w mounting hardware

- Mounts with removeable hardware
- Includes adjustable gel pad



Fuze Camber

- Allows for 5° camber
- Available on all wheel sizes



Fuze JR External Front Rigging

- Increased foot positioning space
- Accommodates lower leg deformities



WHEELCHAIRS FOR FUNCTION AND INDEPENDENCE

www.pdgmobility.com



Fuze JR 90° Front Rigging

- Adjusts 5½" to 10½"
- Compatible with adjustable angle footplates



Bentley Heavy Duty Anti-tips



70° Front Rigging with Calf Pad

- Ideal for foot positioning in tilt
- Prevents need for ELR's



Fuze Jr. 90° Front Rigging c/w extensions

- adjusts 11½" to 15½"
- compatible with adjustable angle footplates



Bentley Rear Floor Jacks

- Raises rear wheels
- Stabilizes chair



**Gel Knee Protectors
(optional ELR Accessory)**

- Prevents skin breakdown
- Improves lower leg positioning



Locking Adjustable Angle Footplates

- Prevents footplates from separating in middle
- Compatible with Adjustable Angle footplates



Stellar 'Stentley' Strut Bracket

- shortened strut bracket
- improved heel clearance for foot propelling
- 45° tilt



Stellar HD

- 450 lb weight capacity
- 600 lbs capacity, with tilt limitations
- 45° tilt, with tilt limiter
- 22" – 30" seat width
- Includes 3 gas struts for functional tilt action
- Forward front casters
- 24" Heavy Duty Spoked Wheels



USA
 PDG Product Design Group
 PMB # 173 - 1124 Fir Avenue
 Blaine, WA
 98230, USA

CANADA
 PDG Product Design Group
 Unit 103 - 318 East Kent Avenue South
 Vancouver BC
 V5X 4N6

T: +1 604 323 9220
 F: +1 604 323 9097
 Toll Free (in N. America): 888-858-4422
 Toll Free Fax (in N. America): 877-323-9093
 www.pdgmobility.com