

# eNews

May 2008



Fuze T50



## The Fuze T50 Chosen by MDEA As an Award Winner

### Great News!

The Fuze T50 has been selected as a winner by the 2008 Medical Design Excellence Awards (MDEA) judging panel. MDEA is the premier awards program for the medical community. MDEA exclusively recognizes contributions and advances in medical product design.

The Fuze T50 took top honors in the Rehabilitation and Assistive Technology Products category.

Here are some features of the Fuze that contributed to PDG winning this award:

- We placed a simple tilt pivot point at the occupant's center of gravity. This makes for easy, effortless tilting regardless of client weight.
- The well-balanced tilt mechanism benefits both patient and care-giver.
- The wheelchair is easily adjusted. Most changes can be made using one Allen key.
- The Fuze can be easily converted from a 20° tilt to a 50° tilt. This allows the user to keep the same wheelchair as their needs change over time.

For more information about this great wheelchair and other PDG products, call us or visit PDG online at [www.pdgmobility.com](http://www.pdgmobility.com).

### PDG On the Road



June 9 - 12 ..... Heartland ..... Waterloo IA ..... Booth 407



June 11 - 13 ..... Autonomic ..... Paris, France ..... Booth F100

