



A NEW Math Equation



PDG + 70 DEGREE
+ NRC
(NO ROTATION CHANNLED)
X FOOTRESTS
= BEAUTIFUL !!!!!

NEW 70° NRC's (NO ROTATION CHANNLED) FRONT RIGGING!

PDG is proud to announce it is **NOW** shipping our new swing-away 70° removeable front rigging with composite or angle-adjustable footplates for the Fuze T20 and Fuze T50. The forged and cast aluminium channel design allows the extension tube to fit perfectly inside the hanger tube creating a non-rotational system.

The front rigging are easily adjusted. The footplates will not sag. Screw heads and nuts are recessed in the channel design to prevent potential injury to the wheelchair user. Most adjustments can be done with one allen key. PDG's new composite footplates come standard with forward and rear mounting choices.

* Watch for the Bentley and Stellar version coming soon!

pdg This spring, see us at the following events:

...On the Road

- Posture and Mobility Conference in Warwick, UK from April 14th to 16th
- NSM 2009 in Nashville, TN on April 26 in Booth 12
- Naidex in Birmingham, UK from April 28th to 30th in Stand C96
- CSMC in Toronto, ON from May 12th to 14th in Booth 801



Click on the HME logo to register!

Come join us at HME Expo in Baltimore! As a valued PDG customer we are pleased to provide you with free admission to the expo hall at HME EXPO.

Register today at <http://www.hmeexpo.com> and use promotional code: HME10075.

