

## FUZE T50 JR

# ...Now Available for Kids



The FUZE T50 JR

labeled FUZE T50 JR. This version of the product is built around the same architecture as the well known FUZE T50.

Some nice features for this product:

- the same effortless tilt as the FUZE T50 adult version
- shorter wheelbase, making navigation in tight spaces easier
- very fast to set-up and fit end-user needs

Prescribers using the Fuze T50 JR like the rigid frame design combined with the ability for the chair to grow with the child. It is easy to adjust and provides excellent service with its simple pivot-point design and single-trigger pivot activator.

Responding to demand from prescribers and dealers familiar with the FUZE family of products, we recently completed a smaller size FUZE chair for younger clients. Smaller clients are more likely to need lower armrests, shorter seat depths and 90 degree front rigging. To help prescribers pick options commonly needed for kids and young adults, we set up a specialized order form

## New Childrens' Lap Belt



In collaboration with Bodypoint ©, PDG now has a specialized lap-belt that makes belt length adjustment easy with the child seated in the wheelchair. Available in either a two-point or four-point style.



Feb 23:	Medichair Dealers National show	Toronto
March 3 – 5:	ISS	Nashville
May 22:	NSM show	Nashville
April 21:	Medichair Halton Show	Halton
May 17 -19:	CSMC	Toronto

**Have questions? Call us.**  
**We'll be happy to tell you more.**

T: Toll Free in N. America: 1.888.858.4422  
email: [info@PDGmobility.com](mailto:info@PDGmobility.com)