

Stellar and Stellar HD Manual Tilt Wheelchair

Ideal for individuals requiring dynamic tilt and/or with limited stamina.



- Tilts 45 degrees, rotating about a point near the front of the seat
- Improves posture and prevents sliding forward
- Excellent access to rear wheels for effortless propulsion
- Adjusts to fit individual needs
- Many options available



Basic Components



Ideal Access to Push-rims



- Rear wheels are infinitely adjustable to ensure easy access to push-rims
- Adjust rear wheels forward to reduce load on casters – makes the chair very easy to propel



Sit at a Table While in Tilt



- Seat tilts about a point located near the front of the seat
- Knees rise only $\frac{3}{4}$ " (1.9cm) when seat is tilted 20° , or 1-1/2" (3.8cm) when seat is tilted to 45°



Basic Options

Very adjustable frame to ensure a perfect fit

- Seat depth 16" (40.6cm) - 22" (55.9cm)
- Seat height 14" (35.6cm) - 20" (50.8cm)
- Seat width 14" (35.6cm) - 22" (55.9cm)
- HD seat width 23" (58.4cm) – 32" (81.3cm)
- Weight capacity 250 lb (113kg)
- HD weight capacity up to 600 lb (272kg)
- Rear wheel options 12" to 24" (HD rear wheels 24")



Stellar HD – Weight limit 600 lb (272kg)



- In-stock sizes up to 32" (81.3cm) seat width
- Attendant operated tilt with ease



Front Rigging Options



Choose swing-away front rigging or other types appropriate for user-specific needs



Plug-in front rigging eliminates potential obstructions when foot propelling.



Knee angle adjustable front rigging.



Back Canes Options



Choose from:

- Standard push-bar
- Adjustable angle
- Cane style



Transit Tie-Down Option

Rear tie-down mount



Forward tie-down mount



Thanks for Taking a Look at the Stellar



To learn more about PDG product applications, please visit www.PDGmobility.com and look for Case Histories

