



Stellar Express

Canadian Price List
Revised April-2012 Page 1 of 1
THE COMFORT COMPANY

PO#	Tag:	Dealer Name:	Dealer Contact:
Mailing Address:		Shipping Address:	
City:	Province:	Postal Code:	
Tel:	Fax:		

EXPRESS FORM

**THIS CONFIGURATION INCLUDES THE FOLLOWING
PRE-SELECTED OPTIONS - NO SUBSTITUTIONS**

Confirm Order Via E-mail:	<input type="checkbox"/> YES	<input type="checkbox"/> NO
Send PDG e-News:	<input type="checkbox"/> YES	<input type="checkbox"/> NO
E-Mail:		

STELLAR FRAME TYPE - 45 DEGREES OF TILT		
<input checked="" type="checkbox"/> 50331	Stellar 45 Degree Manual Tilt Wheelchair Max Weight Capacity 250 lbs.	\$2,750
<input checked="" type="checkbox"/> 50332	Transit Tie-Down Eyelets Does not include WC 19 Lap Belt. (TTDS includes 4 factory installed transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks or Extended Base.	STD

ANTI-TIP OPTIONS		
<input checked="" type="checkbox"/> 50445	Rear Anti-Tippers (Required Safety option).	STD

SEAT WIDTH		
<input type="checkbox"/> 50337	Seat Width 16"	NC
<input type="checkbox"/> 50339	Seat Width 18"	NC
<input type="checkbox"/> 50341	Seat Width 20"	NC

BACK POST STYLE		
<input checked="" type="checkbox"/> 50448	Push Bar Back - 25" High	STD
<input checked="" type="checkbox"/> 50452	Angle/Height Adjustable Push Bar Handle	STD

SEAT DEPTH - (Dealer Adjustable 16"-21")		
<input checked="" type="checkbox"/> 50355	17" Seat Depth	STD

BACK OPTIONS (Adjustable)		
<input checked="" type="checkbox"/> 50453	Seat To Back Angle: Standard: 90° Dealer Adjustable: 96°/ 102°/ 108°/ 114°/ 120°	STD

FRAME COLOR		
<input checked="" type="checkbox"/> 50363	Black Shimmer	STD

BACK UPHOLSTERY		
<input checked="" type="checkbox"/> 50462	No Upholstery	STD

ARM TYPE - (Removable & Adjustable Height)		
<input checked="" type="checkbox"/> 50372	Dual Post Adjustable Height, 10" - 14"	STD
<input checked="" type="checkbox"/> 50383	Full Length Armpad, 2" X 12" (pair)	STD

ATTENDANT TILT		
<input checked="" type="checkbox"/> 50504	Attendant Tilt - Right Hand Only	STD

SEAT PAN OPTIONS		
<input checked="" type="checkbox"/> 50393	Solid Seat Pan	STD
<input checked="" type="checkbox"/> 50394	No Seat Cushion	STD

FRONT RIGGING		
<input type="checkbox"/> 50524	Swing Away 70° Front Rigging	STD
<input type="checkbox"/> 50526	Swing Away Elevating Legrest	\$355

REAR WHEEL / CASTER / SEAT TO FLOOR		
<input type="checkbox"/> 50409	20" Composite Wheel, 6" x 1-1/4" Urethane 15" STF (Adjustable STF 14" or 16")	STD
<input type="checkbox"/> 51463	22" Composite Wheel, 6" x 1-1/4" Urethane 16" STF (Adjustable STF 15" or 17")	NC

FOOT PLATES		
<input type="checkbox"/> 50528	Compostite Footplates	STD
<input type="checkbox"/> 50529	Adjustable Angle Footplates	\$140

REAR TIRE OPTIONS		
<input checked="" type="checkbox"/> 50419	Urethane Tires	STD

FOOTREST/LEGREST ACCESSORIES		
<input checked="" type="checkbox"/> 50536	Adjustable Length Heel Loops	STD

AXLE TYPE		
<input checked="" type="checkbox"/> 50426	Fixed Rear Axles	STD

EXTENSION TUBE (Adjustable)		
<input checked="" type="checkbox"/> 50533	Standard Length (13" - 17")	STD

HAND RIM OPTIONS		
<input checked="" type="checkbox"/> 50428	Steel Hand Rims	STD

MISCELLANEOUS OPTIONS		
<input type="checkbox"/> 50539	Bodypoint - Center Pull Auto Style Hip Belt	\$85
<input type="checkbox"/> 50427 50460	QR Rear Axles & QR Fold Down Back - Accommodates Folding For Auto Transport	\$80

WHEEL LOCKS		
<input checked="" type="checkbox"/> 50439	Push To Lock Wheel Locks	STD
<input type="checkbox"/> 50443	Wheel Lock Extensions (pair)	\$60

NOTES	
Total Retail Value - This order	\$
Have all desired specs been checked? A complete spec sheet ensures your chair will arrive on time and as you expect.	
Many other options are available. Visit www.pdgmobility.com to view and order full-featured products.	