



Stellar

Express

Ph: +1 888 858 4422 July 9, 10
 Ph: +1 604 323 9220 Price List
 Fax: +1 877 323 9093

International
 In Cdn \$

P.O.#:	Quote #:
Dealer:	
Address:	
City:	State/Prov:
Tel:	Fax:
Tag:	
Dealer Contact:	

This configuration includes the following pre-selected options - No substitutions

Stellar Frame Type - 45 Degrees of Tilt		
<input checked="" type="checkbox"/>	50331 Stellar manual tilt wheelchair Wt Capacity = 250 lb (113.5 kg)	\$2,950
<input type="checkbox"/>	50332 Transit Tie Down Eylets (does not include WC 19 Lap belt)	NC

Select Seat Width (Inches)		
<input type="checkbox"/>	50337 16" seat width (40.6 cm)	NC
<input type="checkbox"/>	50339 18" seat width (45.7 cm)	NC
<input type="checkbox"/>	50341 20" seat width (50.8 cm)	NC

Seat Depth - Dealer can adjust 16" - 21" (40.6cm - 53.3cm)		
<input type="checkbox"/>	50355 17" seat depth (43 cm)	NC

Frame Color		
<input type="checkbox"/>	50361 Black Shimmer	STD

Arm Type - Removable & Adjustable Height		
<input checked="" type="checkbox"/>	50371 Dual post adjustable height 10" - 14" (25.4cm - 35.6cm)	NC

Armrest Options		
<input type="checkbox"/>	50383 Full length armapd, 2" x 12" (pair)	NC

Seat Pan Options		
<input checked="" type="checkbox"/>	50393 Solid seat pan	NC

Seat Height Options (Dealer can adjust)		
<input type="checkbox"/>	15" (38.1cm) or 16" (40.6cm) (see rear wheel options)	STD

Caster Options		
<input type="checkbox"/>	50401 6" x 1-1/4" urethane	NC

Rear Wheel Options		
<input type="checkbox"/>	50409 20" rear wheels - Seat height = 15" (38.1 cm) - Adjustable to 14" (35.6cm) or 16" (40.6cm)	STD
<input type="checkbox"/>	50410 22" rear wheels - Seat height = 16" (40.6cm) - Adjustable to 15" (38.1cm) - 17" (43.2cm)	NC

Rear Tire Options		
<input type="checkbox"/>	50419 Urethane tires	STD

Axle Type		
<input type="checkbox"/>	50426 Fixed rear axles	STD

Hand Rim Options		
<input type="checkbox"/>	50428 Steel hand rims	STD

Wheel Locks		
<input type="checkbox"/>	50539 Push to lock wheel locks	STD
<input type="checkbox"/>	50443 Wheel lock extensions (pair)	\$60

Anti-Tip Options		
<input checked="" type="checkbox"/>	50445 Rear anti-tippers	NC

Back Post Style		
<input type="checkbox"/>	50448 Push Bar Back - 25" high (63.5 cm)	STD
<input type="checkbox"/>	50452 Adjustable height / angle push bar handle	NC

Back Options (Adjustable)		
<input type="checkbox"/>	50453 Seat to back angle - set at: 90° <u>X</u> Adjustable to 96°, 102°, 108°, 114°, or 120°	STD

Back Upholstery		
<input type="checkbox"/>	50462 No back upholstery	STD

Attendant Tilt / User Operated Self - Tilt		
<input type="checkbox"/>	50504 Attendant tilt only: right <u>X</u>	STD

Front Rigging 38		
<input type="checkbox"/>	50524 Swing away removable 70° front rigging	NC
<input type="checkbox"/>	50526 Swing away elevating legrest	\$355

Foot Plates		
<input type="checkbox"/>	50528 Composite footplates	STD
<input type="checkbox"/>	50529 Adjustable angle footplates	\$140

Extension Tube (Adjustable)		
<input type="checkbox"/>	50533 Standard length 13" - 17" (33cm - 43.2cm)	STD

Foot Rest/Legrest Accessories		
<input checked="" type="checkbox"/>	50536 Adjustable length heel loops	NC

Other Options		
<input type="checkbox"/>	50743 BodyPoint - center-pull auto style hip belt	\$85
<input type="checkbox"/>	QR rear axles & QR fold-down back - accomodates folding for auto transport	\$80

Total Retail Value - This order \$ _____

This is a required option
 Many other options are available. Visit www.pdgmobility.com to view and order full-featured products.