

<b>PO#</b>	<b>Tag:</b>	<b>Dealer Name:</b>	<b>Dealer Contact:</b>
<b>Mailing Address:</b>		<b>Shipping Address:</b>	
<b>City:</b>	<b>Province:</b>	<b>Postal Code:</b>	<b>City:</b>
<b>Tel:</b>	<b>Fax:</b>		<b>Tel:</b>
<b>CLIENT INFORMATION - (Ensures optimal set-up of chair)</b>		<b>Confirm Order Via E-mail:</b> <input type="checkbox"/> YES <input type="checkbox"/> NO <b>Send PDG e-News:</b> <input type="checkbox"/> YES <input type="checkbox"/> NO <b>E-Mail:</b>	
<input checked="" type="checkbox"/> 50334	<b>Weight (lbs.)</b> _____ <b>Height (ft./in.)</b> _____	NC	

STELLAR FRAME TYPE - 45 DEGREES OF TILT		
<input checked="" type="checkbox"/> 50331	<b>Stellar Manual Tilt Wheelchair</b> <small>Max Weight Capacity 250 lbs.</small>	\$3,995

FRAME OPTIONS		
<input type="checkbox"/> 704	<b>Heavy Duty Package - Max 22" Wide</b> <small>Max Weight Capacity 300 lbs, includes secondary gas strut.</small>	\$200
<input type="checkbox"/> 51985	<b>Frame Upgrade/Heavy Duty Package</b> <small>Available on 23" &amp; 24" seat width only. Max Weight Capacity 450 lbs, includes frame modification for dual gas strut.</small>	\$300
<b>STELLAR HD (If selected - Skip to Page 4, for HD Specifications)</b>		\$1200
<input type="checkbox"/> 51768	<b>Stellar HD</b> <small>Max Weight Capacity 450 lbs. Requires 20" composite or 24" heavy duty spoked wheel. Longer frame length for optional performance.</small>	
<input type="checkbox"/> 51789	<b>Stellar with Shortened Strut Bracket</b> <small>Increases heel clearance for foot propeling.</small>	\$75
<input type="checkbox"/> 51305	<b>Extended Base with Rear Floor Jacks</b> <small>Lifts rear wheels off ground for increased stability. Not compatible with Attendant Wheel Locks, 22" rw, Anti tipper, IV Pole, O2 Holder. Limited Ground clearance.</small>	\$600
<input type="checkbox"/> 50332	<b>Transit Tie-Down System (TTDS)</b> <small>(TTDS includes 4 factory installed transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details.) Not compatible with Attendant Wheel Locks or Extended Base.</small>	\$95

SEAT WIDTH - (Inches)		
<input type="checkbox"/> 50335	14" Seat Width	NC
<input type="checkbox"/> 50336	15" Seat Width	NC
<input type="checkbox"/> 50337	16" Seat Width	NC
<input type="checkbox"/> 50338	17" Seat Width	NC
<input type="checkbox"/> 50339	18" Seat Width	NC
<input type="checkbox"/> 50340	19" Seat Width	NC
<input type="checkbox"/> 50341	20" Seat Width	NC
<input type="checkbox"/> 50342	21" Seat Width	\$195
<input type="checkbox"/> 50343	22" Seat Width	\$195
<input type="checkbox"/> 50344	<b>23" Seat Width</b> <small>Max weight cap. 300lbs. Requires 20" Composite or 24" Heavy duty spoked wheel. Fixed axle only. Requires frame upgrade / heavy duty package above.</small>	\$600
<input type="checkbox"/> 50345	<b>24" Seat Width</b> <small>Max weight cap. 300lbs. Requires 20" Composite or 24" Heavy duty spoked wheel. Fixed axle only. Requires frame upgrade / heavy duty package above.</small>	\$600

SEAT DEPTH - (Adjustable)		
<input type="checkbox"/> 50354	16" Seat Depth	NC
<input type="checkbox"/> 50355	17" Seat Depth	NC
<input type="checkbox"/> 50356	18" Seat Depth	NC
<input type="checkbox"/> 50357	19" Seat Depth	NC
<input type="checkbox"/> 50458	20" Seat Depth	NC
<input type="checkbox"/> 50359	21" Seat Depth	NC
<input type="checkbox"/> 50360	22" Seat Depth	NC

FRAME COLOR		
<input type="checkbox"/> 50361	Blue Shimmer	NC
<input type="checkbox"/> 50362	Tracer Red	NC
<input type="checkbox"/> 50363	Black Shimmer	NC
<input type="checkbox"/> 51793	Orange Copper	NC
<input type="checkbox"/> 50365	Misty Purple	NC
<input type="checkbox"/> 50366	Forest Green	NC
<input type="checkbox"/> 50367	Custom Color: _____ (Contact PDG for Specification)	\$300

ARM TYPE		
<input type="checkbox"/> 50368	No Armrests	STD
<input type="checkbox"/> 50369	Dual Post Fixed Height 12" <small>Full length pad only.</small>	NC
<input type="checkbox"/> 50371	Dual Post Adjustable Height Armrest <small>Option: 8" - 11" <input type="checkbox"/> 10" - 14" <input type="checkbox"/> 12" - 16" <input type="checkbox"/></small>	\$295
<input type="checkbox"/> 50374	3" Armrest Extension - Forward <small>Not compatible with flip-up armrests.</small>	\$59
<input type="checkbox"/> 50375	3" Armrest Extension - Rearward <small>Not compatible with flip-up armrests.</small>	\$59
<input type="checkbox"/> 50377	Locking Flip-up Armrests <small>Not recommended for clients with excessive strength or tone. Not compatible with fold down back, clothing guards, back upholstery or Stellar HD.</small> <small>Set At: 10" <input type="checkbox"/> 11" <input type="checkbox"/> 12" <input type="checkbox"/> 13" <input type="checkbox"/> 14" <input type="checkbox"/></small>	\$295

ARMPAD OPTIONS		
<input type="checkbox"/> 50382	Desk Length Armpad, 2" X 10" (pair)	NC
<input type="checkbox"/> 50383	Full Length Armpad, 2" X 12" (pair)	NC
<input type="checkbox"/> 50384	3" Wide Arm Pads <small>May interfere with wheels in some tilt configuration.</small>	\$49
<input type="checkbox"/> 50384	Full Length 3"x12" (pair)	
<input type="checkbox"/> 50385	Desk Length 3"x12" (pair)	
<input type="checkbox"/> 50386	Gel Armpad <small>May interfere with wheels in some tilt configuration.</small>	\$115
<input type="checkbox"/> 50386	Full Length 2"x12" (pair)	
<input type="checkbox"/> 50387	Full Length 3.5"x12" (pair)	
<input type="checkbox"/> 50370	Contour Arm Trough <small>May interfere with wheels in some tilt configuration.</small>	\$85
<input type="checkbox"/> 50370	Right 4.5"x13"	
<input type="checkbox"/> 50388	Left 4.5"x13"	
<input type="checkbox"/> 50376	Flat Hand Pan <small>Requires Contour Arm Trough.</small>	\$55
<input type="checkbox"/> 50376	Right 4.5"x5"	
<input type="checkbox"/> 50389	Left 4.5"x5"	
<input type="checkbox"/> 50390	Clothing Guards (Black Plastic) <small>Not compatible with flip-up armrests.</small>	\$205
<input type="checkbox"/> 50391	Armrest wideners (pr) <small>Limits Tilt, Not available with 20" and 24" rear wheels.</small>	\$375

SEAT FRAME OPTIONS		
<input type="checkbox"/> 50392	No Seat Pan	STD
<input type="checkbox"/> 50393	Solid Seat Pan	\$140

<b>PO#</b>	<b>Tag:</b>	<b>Dealer Name:</b>	<b>Dealer Contact:</b>
------------	-------------	---------------------	------------------------

SEAT CUSHION OPTIONS - (The Comfort Company)		
<input type="checkbox"/> 50394	No Seat Cushion	NC
<input type="checkbox"/> 52029	Cushion M2 Zero Elevation Stretch-Air	\$577
	Width: 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20" <input type="checkbox"/>	
	Depth: 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20" <input type="checkbox"/>	
<input type="checkbox"/> 52044	Spare Cushion Cover	\$109
<input type="checkbox"/> 52038	Cushion M2 Zero Elevation Stretch-Air (Extra Wide)	\$677
	Width: 22" <input type="checkbox"/> 24" <input type="checkbox"/>	
	Depth: 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20" <input type="checkbox"/>	
<input type="checkbox"/> 52053	Spare Cushion Cover (Extra Wide)	\$109
<input type="checkbox"/> 5	Incontinent Cushion Liner	\$85

CASTER OPTIONS			
<input type="checkbox"/> 50401	6" x 1-1/4" Urethane		NC
<input type="checkbox"/> 50402	6" x 1-1/4" Pneumatic		\$110
<input type="checkbox"/> 51459	6" x 1-3/4" Urethane Caster		\$75
<input type="checkbox"/> 50403	8" x 1-1/4" Urethane	8" Casters - Compatible only with 16" and 24" rear wheels.	NC
<input type="checkbox"/> 50404	8" x 1-1/4" Pneumatic		\$110
<input type="checkbox"/> 50405	8" x 2" Urethane Caster		\$135
<input type="checkbox"/> 50406	8" x 2" Pneumatic Caster		\$135

REAR WHEEL OPTIONS		
<input type="checkbox"/> 50407	12" Composite Wheel	NC
<input type="checkbox"/> 50408	16" Composite Wheel	NC
<input type="checkbox"/> 50409	20" Composite Wheel	Required for Heavy Duty Package. NC
<input type="checkbox"/> 51325	20" Spoke Wheel	Required for Heavy Duty Package. \$75
<input type="checkbox"/> 51463	22" Composite Wheel	NC
<input type="checkbox"/> 51462	22" Spoke Wheel	\$75
<input type="checkbox"/> 50410	24" Composite Wheel	NC
<input type="checkbox"/> 51326	24" Spoke Wheel	\$75
<input type="checkbox"/> 50411	24" Spoked Wheel - Heavy Duty	Required for Heavy Duty Package. Not compatible with Pneumatic rear tires with Flat Free Inserts. \$350

SEAT HEIGHT OPTIONS				
(CIRCLE) desired seat-to-floor height in chart below.				
Rear Wheel	6" Casters			
12"	14"	15"	16"	nc
16"	16"	17"	18"	\$60
20"	14"	15"	16"	nc
22"	15"	16"	17"	nc
24"	16"	17"	18"	nc
	17"	18"	19"	\$60
	18"	19"	20"	\$60
	19"	20"		\$60

REAR TIRE OPTIONS		
<input type="checkbox"/> 50419	Urethane Tires	STD
<input type="checkbox"/> 50424	Pneumatic Tires	NC
	Pneumatic tires may result in a finished seat height 1/2" higher than stated.	
<input type="checkbox"/> 50425	Pneumatic Tires With Flat-Free Inserts	\$200
	Pneumatic tires may result in a finished seat height 1/2" higher than stated. Not compatible with Spoked Wheels.	

AXLE TYPE		
<input type="checkbox"/> 50426	Fixed Rear Axle	STD
<input type="checkbox"/> 50427	Quick Release Rear Axles	NC

HAND RIM OPTIONS		
<input type="checkbox"/> 51631	Omit Hand Rim	20", 22" or 24" rear wheels only. Not recommended for weights over 195 lbs. NC
<input type="checkbox"/> 50428	Standard Hand Rims	20", 22" or 24" rear wheels only. NC
<input type="checkbox"/> 50429	Plastic Coated Hand Rims	20" or 24" rear wheels only. \$95
<input type="checkbox"/> 50430	Plastic Coated Projection Hand Rims	20", 22" or 24" rear wheels only. \$225
<input type="checkbox"/> 50431		
<input type="checkbox"/> 50435		
<input type="checkbox"/> 50436		
<input type="checkbox"/> 50437		
	Number of Projections - 8	
	Projection Length - 1"	
	Select Projection Angle: 0° <input type="checkbox"/> 45° <input type="checkbox"/>	

WHEEL LOCKS		
<input type="checkbox"/> 50439	Push To Lock Wheel Locks	STD
<input type="checkbox"/> 50440	Pull To Lock Wheel Locks	NC
<input type="checkbox"/> 50441	Attendant Wheel Lock (Push To Lock)	Not compatible with Transit Tie Down System. \$85
<input type="checkbox"/> 50442	User Wheel Locks And Attendant Wheel Locks	Available only on 20" and 24" rear wheels. Not compatible with Transit Tie Down System. \$150
<input type="checkbox"/> 50443	Wheel Lock Extensions (Pair)	Not compatible with Attendant wheel locks. \$60
<input type="checkbox"/> 50444	Hill Holders	\$350

ANTI-TIP OPTIONS		
<input checked="" type="checkbox"/> 50445	Rear Anti-Tippers (Pair)	Required safety option. \$99
<input type="checkbox"/> 50446	Front Anti-Tippers	\$300
<input type="checkbox"/> 50447	Heavy Duty Rear Anti-Tipper	Provides additional rearward stability for clients with excessive movement. Not compatible with Attendant Wheel Locks. \$600

BACK POST STYLE		
<input type="checkbox"/> 50448	Push Bar Back - 25" High	STD
<input type="checkbox"/> 50449	Push Bar Back - 31" High	\$35
<input type="checkbox"/> 50450	Cane Handle Back - 25" High	Seat widths 14" - 22" only. Not compatible with headrest mounting plate. \$35
<input type="checkbox"/> 50451	Cane Handle Back - 31" High	Seat widths 14" - 22" only. Not compatible with headrest mounting plate. \$35

BACK OPTIONS		
<input type="checkbox"/> 50452	Adjustable - Height / Angle Push Bar Handle	Not compatible with headrest mounting plate or cane handle back. \$395
<input type="checkbox"/> 50453	Seat To Back Angle:	STD
	Set At: 90° <input type="checkbox"/> 98° <input type="checkbox"/> 106° <input type="checkbox"/> 114° <input type="checkbox"/> 122° <input type="checkbox"/>	
<input type="checkbox"/> 50459	Backrest Extension 10" High	Not compatible with adjustable height / angle push bar handle. \$150
<input type="checkbox"/> 50460	Fold Down Back	Not compatible with Reclining Back. NC
<input type="checkbox"/> 50461	Spring Back Suspension	Not compatible with Reclining Back. \$105

BACK UPHOLSTERY - (Back height measured from seat pan to top of upholstery)		
<input type="checkbox"/> 50462	No Upholstery	STD
<input type="checkbox"/> 50463	Standard Back Upholstery	\$125
	Set Height At: 21" <input type="checkbox"/> 22" <input type="checkbox"/> 23" <input type="checkbox"/>	
<input type="checkbox"/> 50470	Short Upholstery - (25" frame only)	\$125
	Set Height At: 18" <input type="checkbox"/> 19" <input type="checkbox"/> 20" <input type="checkbox"/>	
<input type="checkbox"/> 50471	Tall Upholstery - (31" frame only)	\$125
	Set Height At: 24" <input type="checkbox"/> 25" <input type="checkbox"/> 26" <input type="checkbox"/>	
<input type="checkbox"/> 50476	Adjustable Tension Upholstery	\$250
	Set Height At: 21" <input type="checkbox"/> 22" <input type="checkbox"/> 23" <input type="checkbox"/>	
	Not available in widths greater than 22".	

<b>PO#</b>	<b>Tag:</b>	<b>Dealer Name:</b>	<b>Dealer Contact:</b>
------------	-------------	---------------------	------------------------

BACK OPTIONS A - (The Comfort Company)			
<input type="checkbox"/> 51988	<b>Acta-Back Standard Stretch Air</b>		
	Height: 16" <input type="checkbox"/> 18" <input type="checkbox"/>		
	Width: 15"-17" <input type="checkbox"/> 17"-19" <input type="checkbox"/> 19"-21" <input type="checkbox"/>		\$815
	Width: 21"-23" <input type="checkbox"/> 23"-25" <input type="checkbox"/>		\$898
<input type="checkbox"/> 51990	<b>Acta-Back Standard Stretch Air (Extra Tall-20")</b>		
	Width: 15"-17" <input type="checkbox"/> 17"-19" <input type="checkbox"/> 19"-21" <input type="checkbox"/>		\$898
	Width: 21"-23" <input type="checkbox"/> 23"-25" <input type="checkbox"/>		\$980

BACK OPTIONS B - (The Comfort Company)			
<input type="checkbox"/> 52007	<b>Acta-Back Standard Stretch Air DEEP</b>		
	Height: 16" <input type="checkbox"/> 18" <input type="checkbox"/>		
	Width: 15"-17" <input type="checkbox"/> 17"-19" <input type="checkbox"/> 19"-21" <input type="checkbox"/>		\$947
	Width: 21"-23" <input type="checkbox"/> 23"-25" <input type="checkbox"/>		\$1,112
<input type="checkbox"/> 52009	<b>Acta-Back Standard Stretch Air DEEP (Extra Tall-20")</b>		
	Width: 15"-17" <input type="checkbox"/> 17"-19" <input type="checkbox"/> 19"-21" <input type="checkbox"/>		\$1,043
	Width: 21"-23" <input type="checkbox"/> 23"-25" <input type="checkbox"/>		\$1,214

HEADREST OPTIONS - (The Comfort Company)			
<input type="checkbox"/> 50496	Headrest Mounting Plate		\$75
<input type="checkbox"/> 51890	Fully Adjustable Hardware		\$330
<input type="checkbox"/> 51892	<b>Moldable Headrest Pad - Stretch-Air</b>		\$272
	Width: 8" <input type="checkbox"/> 10" <input type="checkbox"/> 14" <input type="checkbox"/>		

TILT / RECLINE OPTIONS			
<input type="checkbox"/> 50494	<b>Reclining Backrest(Manual) Max. Recline - 30°</b> <small>Not compatible with cane handle back, spring back suspension or fold down back.If upholstery is ordered with Recline, it must be tension adjustable.</small>		\$1,200

ATTENDANT TILT / USER OPERATED SELF - TILT			
<input type="checkbox"/> 50504	<b>Attendant Tilt Only</b>		STD
	Right <input type="checkbox"/> Left <input type="checkbox"/>		
<input type="checkbox"/> 50506	<b>Self Tilt Only</b>		\$99
	<small>Not available on widths greater than 22". Full Length arm pads only.</small>		
	Right Armrest <input type="checkbox"/> Left Armrest <input type="checkbox"/>		
<input type="checkbox"/> 50508	<b>Dual Activator: Backrest/Armrest</b>		\$195
	<small>Not available on widths greater than 22. Full Length arm pads only.</small>		
	Backrest - Mount: Right <input type="checkbox"/> Left <input type="checkbox"/> Armrest - Mount: Right <input type="checkbox"/> Left <input type="checkbox"/>		
<input type="checkbox"/> 50512	<b>Trigger Extension</b>		\$48
<input type="checkbox"/> 50513	<b>Tilt Limiting Device</b>		STD
	Set At: 20° <input type="checkbox"/> 25° <input type="checkbox"/> 30° <input type="checkbox"/> 35° <input type="checkbox"/> 40° <input type="checkbox"/> 45° <input type="checkbox"/>		

FRONT RIGGING		
<input type="checkbox"/> 50519	Omit Front Rigging	NC
<input type="checkbox"/> 50520	Sector Blocks Only (No Front Rigging Supplied)	NC
<input type="checkbox"/> 50521	<b>Plug in 70° Front Rigging</b> <small>Must have frame for plug in front rigging.</small>	STD
<input type="checkbox"/> 50522	<b>Plug in Adjustable Knee Angle Front Rigging</b> <small>Must have frame for plug in front rigging.</small>	\$290
<input type="checkbox"/> 50523	Swing Away 60° Front Rigging	NC
<input type="checkbox"/> 50524	Swing Away 70° Front Rigging	NC
<input type="checkbox"/> 50526	Swing Away Elevating Legrests	\$355

FOOT PLATES		
<input type="checkbox"/> 50528	Composite Footplates	STD
<input type="checkbox"/> 50529	Adjustable Angle Footplates	\$140
<input type="checkbox"/> 50531	<b>One Piece FIXED Footplates</b> <small>Available on plug-in front rigging only. 14"-22" wide only.</small>	\$350
<input type="checkbox"/> 50532	<b>One Piece FLIP-UP Footplates</b>	\$395
<input type="checkbox"/> 51841	<b>Interlocking Angle Adjustable Footplate</b> <small>14"-21" wide only</small>	\$195

FOOTREST/LEGREST ACCESSORIES		
<input type="checkbox"/> 52059	60° & 70° Swing Away Calf Pad	\$135
<input type="checkbox"/> 51761	Swing Away Elevating Legrest Gel Knee Pad	\$99
<input type="checkbox"/> 52060	70° Swing Away Gel Padding	\$99
<input type="checkbox"/> 50536	Adjustable Length Heel Loops	\$55
<input type="checkbox"/> 50537	<b>Bodypoint Padded Calf Strap</b> <small>Sizes greater than 22" are not padded.</small>	\$49
<input type="checkbox"/> 50538	<b>Miller's Dynamic Footrest Mounts For Plug In Front Rigging.</b>	\$190
<input type="checkbox"/> 52061	<b>Miller's Dynamic Footrest Mounts For Swing Away 70° Front Rigging</b>	\$380

EXTENSION TUBE		
<input type="checkbox"/> 50533	Standard Length (13" - 17")	STD
<input type="checkbox"/> 50534	Long Extension Tubes (15" - 19")	\$75
<input type="checkbox"/> 50535	<b>Standard Length (13" - 17") Solid Extension Tubes - Extra Strength</b> <small>Available on Plug-in front Rigging only.</small>	\$99
<input type="checkbox"/> 51455	<b>Long Extension Tubes (15" - 19") Solid Extension Tubes - Extra Strength</b> <small>Available on Plug-in front Rigging only.</small>	\$99

POSITIONING STRAPS (www.bodypoint.com)		
<input type="checkbox"/> 50539	Bodypoint - Center-Pull Auto Style Hip Belt	\$85
<input type="checkbox"/> 50540	Bodypoint - Padded, Center-Pull Auto Style Hip Belt	\$125
<input type="checkbox"/> 50541	Bodypoint - 4 Point Hip Belt	\$150
<input type="checkbox"/> 50542	Bodypoint - Chest Strap	\$75

MISCELLANEOUS OPTIONS		
<input type="checkbox"/> 51471	Protective Lower Padding	\$425
<input type="checkbox"/> 51284	<b>IV Pole</b> <small>Not compatible with TTDS (Transit Tie Down System), Rear Stabilizers, Rear Floor Jacks, Attendant Wheel Locks.</small>	\$295
<input type="checkbox"/> 50543	<b>O2 Holder</b> <small>Not compatible with TTDS (Transit Tie Down System), Rear Stabilizers, Rear Floor Jacks, Attendant Wheel Locks.</small>	\$295
<input type="checkbox"/> 50544	<b>Accessory Tray</b> <small>Model and dimension must be supplied.</small>	\$350
	Width: ____ x Height: ____ x Length: ____	

PO#	Tag:	Dealer Name:	Dealer Contact:
-----	------	--------------	-----------------

**STELLAR HD ONLY**  
**THIS CONFIGURATION INCLUDES THE FOLLOWING PRE-SELECTED OPTIONS - NO SUBSTITUTIONS**

CLIENT INFORMATION - (Ensures optimal set-up of chair)		
<input checked="" type="checkbox"/> 50334	Weight (lbs.) _____ Height (ft./in.) _____	NC
STELLAR FRAME TYPE - 45 DEGREES OF TILT		
<input checked="" type="checkbox"/> 50331	<b>Stellar Manual Tilt Wheelchair</b> <small>Max Weight Capacity 250 lbs.</small>	\$3,995

STANDARD OPTIONS INCLUDE		
<input checked="" type="checkbox"/> 51768	<b>Stellar HD</b> <small>Max Weight Capacity 450 lbs. Longer frame length for optional performance.</small>	\$1200
<input checked="" type="checkbox"/> 52258	<b>450 lbs Weight Capacity</b> Up to 600 lbs with 25° tilt limiter.	STD
<input checked="" type="checkbox"/> 52259	<b>45 Degrees of Tilt</b> Limited to 25° over 450 lbs.	STD
<input checked="" type="checkbox"/> 51472 51973 51974 51975	<b>Default Gas Struts</b> <input type="checkbox"/> up to 250 lbs = 2 Gas Struts - 1200{N} <input type="checkbox"/> 250 lbs - 350 lbs = 3 Gas Struts - 800{N} <input type="checkbox"/> 350 lbs - 425 lbs = 3 Gas Struts - 1000{N} <input type="checkbox"/> 425 lbs - 500 lbs = 3 Gas Struts - 1200{N}	STD
<input type="checkbox"/> 51800	<b>Seat Width</b> Set At: 22" <input type="checkbox"/> 24" <input type="checkbox"/> 26" <input type="checkbox"/> 28" <input type="checkbox"/> 30" <input type="checkbox"/> 32" <input type="checkbox"/>	STD
<input type="checkbox"/> 50356	<b>Seat Depth - (Dealer Adjustable)</b> Set At: 18" <input type="checkbox"/> 19" <input type="checkbox"/> 20" <input type="checkbox"/> 21" <input type="checkbox"/> 22" <input type="checkbox"/>	STD
<input type="checkbox"/> 50361	<b>Frame Color</b> <input type="checkbox"/> Blue Shimmer <input type="checkbox"/> Black Shimmer <input type="checkbox"/> Tracer Red <input type="checkbox"/> Orange Copper <input type="checkbox"/> Misty Purple <input type="checkbox"/> Forest Green <input type="checkbox"/> Custom Color: _____ (Contact PDG for Specification)	STD \$300
<input checked="" type="checkbox"/> 50372	<b>Dual Post Adjustable Height Arm Rests, 10" - 14"</b>	\$295
<input checked="" type="checkbox"/> 50393	<b>Solid Seat Pan</b>	\$140
<input checked="" type="checkbox"/> 50394	<b>No Seat Cushion</b>	STD
<input type="checkbox"/> 50414	<b>Seat to Floor Height</b> Set At: 16" <input type="checkbox"/> 17" <input type="checkbox"/> 18" <input type="checkbox"/>	STD
<input checked="" type="checkbox"/> 50411	<b>24" Spoked Heavy Duty Rear Wheel</b>	\$350
<input checked="" type="checkbox"/> 50419	<b>Urethane Tire</b>	STD
<input checked="" type="checkbox"/> 50456	<b>Forward Set, 5" Heavy Duty Steel Fork with Urethane Caster</b>	STD
<input checked="" type="checkbox"/> 50428	<b>Aluminum Hand Rim</b>	STD
<input checked="" type="checkbox"/> 50426	<b>Fixed Rear Axles</b>	STD

ARMPAD OPTIONS		
<input type="checkbox"/> 50382	Desk Length Armpad, 2" X 10" (pair)	NC
<input type="checkbox"/> 50383	Full Length Armpad, 2" X 12" (pair)	NC
	<b>3" Wide Arm Pads</b> <small>May interfere with wheels in some tilt configuration.</small>	\$49
<input type="checkbox"/> 50384	Full Length 2"x12" (pair)	
<input type="checkbox"/> 50385	Desk Length 2"x12" (pair)	
	<b>Gel Armpad</b> <small>May interfere with wheels in some tilt configuration.</small>	\$115
<input type="checkbox"/> 50386	Full Length 2"x12" (pair)	
<input type="checkbox"/> 50387	Full Length 3.5"x12" (pair)	
	<b>Contour Arm Trough</b> <small>May interfere with wheels in some tilt configuration.</small>	\$85
<input type="checkbox"/> 50370	Right 4.5"x13"	
<input type="checkbox"/> 50388	Left 4.5"x13"	
	<b>Flat Hand Pan</b> <small>Requires Contour Arm Trough.</small>	\$55
<input type="checkbox"/> 50376	Right 4.5"x5"	
<input type="checkbox"/> 50389	Left 4.5"x5"	
<input type="checkbox"/> 50390	<b>Clothing Guards (Black Plastic)</b> <small>Not compatible with flip-up armrests.</small>	\$205

STANDARD OPTIONS INCLUDE		
<input checked="" type="checkbox"/> 50439	<b>Push To Lock Wheel Locks</b>	STD
<input checked="" type="checkbox"/> 50445	<b>Rear Anti-Tippers</b>	\$99
<input checked="" type="checkbox"/> 50448	<b>25" Push Back Bar</b>	STD
<input checked="" type="checkbox"/> 50453	<b>Seat To Back Angle 90°</b>	STD
<input type="checkbox"/> 50476	<b>Adjustable Tension Back Upholstery</b> Set Height At: 21" <input type="checkbox"/> 22" <input type="checkbox"/> 23" <input type="checkbox"/>	\$250
<input checked="" type="checkbox"/> 50504	<b>Attendant Tilt Only</b>	STD
<input type="checkbox"/> 50523	<b>Swing Away 60° Front Rigging</b>	STD
<input type="checkbox"/> 50526	<b>Swing Away Elevating Legrests</b>	\$355
<input checked="" type="checkbox"/> 50533	<b>Standard Length Extension Tube (13" - 17")</b>	STD
<input checked="" type="checkbox"/> 50529	<b>Oversized Angle Adjustable Footplates</b>	STD

CUSTOM MODIFICATIONS - NOTES	
Total Retail Value - This order	\$
Have all desired specs been checked? A complete spec sheet ensures your chair will arrive on time and as you expect.	