

Manual **TILT** Wheelchairs by



FUZE T50 



FUZE T20 



FUZE T50 JR. 



STELLAR 



BENTLEY 

For pressure relief and positioning

- Tilts to 50°
- Center of gravity tilt
- Easy to adjust axle position
- Effortless tilt
- Seat width: 12" (30.4 cm) - 20" (50.8) adjustable
- Seat depth: 15" (38.1 cm) - 22" (55.9 cm) adjustable
- Seat height: 13" (33 cm) - 20" (50.8) adjustable
- Backrest angle: 90° to 120°
- Weight limit: 250 lbs (113.4 kg)

For positioning and foot propulsion

- Tilts to 20°
- Excellent for foot and hand propulsion in and out of tilt
- Easy to adjust axle position
- Pivot point at the front so knees remain level while in tilt
- Seat width: 13" (33 cm) - 20" (50.8) adjustable
- Seat depth: 15" (38.1 cm) - 20" (50.8) adjustable
- Seat height: 13" (33 cm) - 20" (50.8) adjustable
- Backrest angle: 90° to 120°
- Weight limit: 250 lbs (113.4 kg)

For pressure relief and positioning in smaller sizes

- Tilts to 50°
- Center of gravity tilt
- Easy to adjust axle position
- Effortless tilt
- Seat width: 12" (30.4 cm) - 18" (45.7 cm) adjustable
- Seat depth: 13" (33 cm) - 20" (50.8 cm) adjustable
- Seat height: 13" (33 cm) - 20" (50.8 cm) adjustable
- Backrest angle: 90° to 120°
- Weight limit: 180 lbs (81.6 kg)

For pressure relief and positioning

- Tilts to 45°
- Excellent for hand propulsion
- Easy to adjust axle position
- Pivot point near the front so knees stay low
- Seat width: 14" (35.6 cm) - 32" (81.3 cm) non-adjustable
- Seat depth: 16" (40.6 cm) - 22" (55.9 cm) adjustable
- Seat height: 14" (35.6 cm) - 20" (50.8 cm) adjustable
- Backrest angle: 90° to 122°
- Weight limit: 250 lbs (113.4 kg); option to 450 lbs (203.1 kg)



For positioning and postural support For clients with agitation

- Tilts to 20°
- Excellent for foot and hand propulsion in and out of tilt
- Easy to adjust axle position
- Pivot point at the front so knee and foot height doesn't change when tilted
- Seat width: 14" (35.6 cm) - 32" (81.3 cm) non-adjustable
- Seat depth: 16" (40.6 cm) - 22" (55.9 cm) adjustable
- Seat height: 13" (33 cm) - 20" (50.8 cm) adjustable
- Backrest angle: 90° to 122°
- Weight limit: 250 lbs (113.4 kg); option to 450 lbs (203.1 kg)
- High agitation option available



T: +1 604.323.9220
Toll Free (in N. America): 888.858.4422
info@pdgmobility.com

www.PDGmobility.com

 All wheelchairs are WC 19 approved.
 Denotes pivot point for tilt on all chairs
Custom Wheelchair sizes and options are available