



STELLAR GL

AADL Cat.# 321

Order Form - AADL
Prices in Canadian Dollars. October 1, 2021

Customer and Order Information

Purchase Order *	Tag	v.1007
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in
		*required

Base Chair

<input checked="" type="checkbox"/> 1901000	Stellar GL 45° Manual Tilt Wheelchair	\$ 3,735
	Also available: the lightweight Stellar GLT with aluminum frame. See separate order form. pdgmobility.com/glt	

Weight Capacity

Chair Weight Capacity		
<input type="checkbox"/> 1902000	250 lbs Capacity Available for 14" to 22" Seat Width.	
<input type="checkbox"/> 1902050	300 lbs Capacity Available for 14" to 22" Seat Width.	\$ 200
<input type="checkbox"/> 1905006	450 lbs Capacity Available for 20" to 32" Seat Width. Not compatible with TTDS. The Stellar HD offers superior strength and performance. See separate order form.	\$ 500
User Weight Range Required for setup of tilt uplift force.		
<input type="checkbox"/> 1902430	Up to 75 lbs	
<input type="checkbox"/> 1902431	75 lbs to 100 lbs	
<input type="checkbox"/> 1902432	100 lbs to 115 lbs	
<input type="checkbox"/> 1902433	115 lbs to 150 lbs	
<input type="checkbox"/> 1902434	150 lbs to 200 lbs	
<input type="checkbox"/> 1902435	200 lbs to 250 lbs	
<input type="checkbox"/> 1902436	250 lbs to 300 lbs	
<input type="checkbox"/> 1902437	300 lbs to 375 lbs	
<input type="checkbox"/> 1902438	375 lbs to 450 lbs	

Frame Options

<input type="checkbox"/> 1902200	Extended Base with Rear Floor Jacks	\$ 900
	Lifts Rear Wheels off ground for increased stability. Limited ground clearance. Not compatible with Attendant Wheel Lock, HD Anti-Tipper, IV Pole, O2 Holder.	
<input type="checkbox"/> 1902300	Transit Tie Down System (TTDS)	
	TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Not compatible with Seat Widths greater than 22". Tested to max. Weight Capacity of 300 lbs.	

Seat Width

Overall chair width is +9.5" with 12/16" Rear Wheels and +10.5" with 20/22/24" Rear Wheels.

Standard Seat Widths		
<input type="checkbox"/> 1903014	14" Seat Width	
<input type="checkbox"/> 1903015	15" Seat Width	
<input type="checkbox"/> 1903016	16" Seat Width	
<input type="checkbox"/> 1903017	17" Seat Width	
<input type="checkbox"/> 1903018	18" Seat Width	
<input type="checkbox"/> 1903019	19" Seat Width	
<input type="checkbox"/> 1903020	20" Seat Width	
<input type="checkbox"/> 1903021	21" Seat Width	
<input type="checkbox"/> 1903022	22" Seat Width	
Custom Seat Widths Requires 450 lbs Weight Capacity Frame. Extended delivery lead times may apply.		
<input type="checkbox"/> 1903023	23" Seat Width	
<input type="checkbox"/> 1903024	24" Seat Width	
<input type="checkbox"/> 1903025	25" Seat Width	\$ 375
<input type="checkbox"/> 1903026	26" Seat Width	\$ 375
<input type="checkbox"/> 1903027	27" Seat Width	
<input type="checkbox"/> 1903028	28" Seat Width	\$ 700
<input type="checkbox"/> 1903029	29" Seat Width	
<input type="checkbox"/> 1903030	30" Seat Width	\$ 700
<input type="checkbox"/> 1903031	31" Seat Width	\$ 700
<input type="checkbox"/> 1903032	32" Seat Width	\$ 700

Seat Depth

Standard Seat Depths (Adjustable)		
<input type="checkbox"/> 1904016	16" Seat Depth	
<input type="checkbox"/> 1904017	17" Seat Depth	
<input type="checkbox"/> 1904018	18" Seat Depth	
<input type="checkbox"/> 1904019	19" Seat Depth	
<input type="checkbox"/> 1904020	20" Seat Depth	
Custom Seat Depths		
<input type="checkbox"/> 1904021	21" Seat Depth	
<input type="checkbox"/> 1904022	22" Seat Depth	\$ 400

Frame Color			
PDG Premium Colors			
<input type="checkbox"/> 1905000	Tracer Red	<input type="checkbox"/> 1905001	Gloss Black Shimmer
<input type="checkbox"/> 1905002	Silver	<input type="checkbox"/> 1905003	Fine Texture Black
<input type="checkbox"/> 1905004	Orange Fusion	<input type="checkbox"/> 1905005	Karl's Imagination
<input type="checkbox"/> 1905006	Blue Shimmer	<input type="checkbox"/> 1905007	Gloss White
<input type="checkbox"/> 1905008	Dark Titanium	<input type="checkbox"/> 1905009	Geyser Green
<input type="checkbox"/> 1905010	Lollypop Blue	<input type="checkbox"/> 1905011	Lollypop Red
<input type="checkbox"/> 1905012	Bubble Gum Pink	<input type="checkbox"/> 1905013	Berry Purple
<input type="checkbox"/> 1905014	Daffodil Yellow	<input type="checkbox"/> 1905200	Custom Color Contact PDG for specifications.
			\$ 300

Armrests				
<input type="checkbox"/> 1906000	No Armrests			
Single Post Height Adjustable Armrests			Not compatible with the 450 lbs Weight Capacity Frame.	\$ 295
<input type="checkbox"/> 1906111	8"-11" Range	<input type="checkbox"/> 1906114	10"-14" Range	
Dual Post Height Adjustable Armrests			Required for the 450 lbs Weight Capacity Frame.	\$ 295
<input type="checkbox"/> 1906211	8"-11" Range	<input type="checkbox"/> 1906214	10"-14" Range	<input type="checkbox"/> 1906216
<input type="checkbox"/> 1906216	12"-16" Range			
Fixed Height Armrests (12")			Full Length Armpad only.	
<input type="checkbox"/> 1906100	Single Post	<input type="checkbox"/> 1906200	Dual Post	
	Armrest Height If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.			
	Seat Pan to top of Armrest: _____		Does not include height of Armpad	
	Armrest Position Indicate preference for Armrest position. If left blank, the default setting will be used.			
In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.				
<input type="checkbox"/> 1906401	Rear	<input type="checkbox"/> 1906402	Mid	<input type="checkbox"/> 1906403
<input type="checkbox"/> 1906403	Front			
Locking Flip-up			Not recommended for users with excessive strength or tone.	\$ 295
Not compatible with Fold Down Back, Clothing Guards, Back Upholstery. Available with Full Length and Desk Length Armpads.				
<input type="checkbox"/> 1906310	Set at 10"	<input type="checkbox"/> 1906311	Set at 11"	<input type="checkbox"/> 1906312
<input type="checkbox"/> 1906312	Set at 12"			
<input type="checkbox"/> 1906313	Set at 13"	<input type="checkbox"/> 1906314	Set at 14"	

Armpads			
Regular Armpads			
<input type="checkbox"/> 1907100	Desk Length - 2" x 10"	<input type="checkbox"/> 1907101	Full Length - 2" x 14"
Wide Armpads			May interfere with wheels in some tilt configurations.
<input type="checkbox"/> 1907110	Desk Length - 3" x 10"	<input type="checkbox"/> 1907111	Full Length - 3" x 12"
Gel Armpads			May interfere with wheels in some tilt configurations.
<input type="checkbox"/> 1907120	Full Length - 2" x 12"	<input type="checkbox"/> 1907121	Full Length - 3.5" x 12"
Contour Arm Trough with Hand Pad			May interfere with wheels in some tilt configurations.
<input type="checkbox"/> 1907130	Right (4.5" x 18")	<input type="checkbox"/> 1907131	Left (4.5" x 18")

Armrest Options			
Armrest Wideners Move Armrests outward by 1" on each side. Do not increase seat size.			
<input type="checkbox"/> 1907210	For Single Post Armrests	May cause tilt interference at lower Seat Height and/or Seat Depth settings.	
<input type="checkbox"/> 1907220	For Dual Post Armrests	Limits tilt. Not compatible with 20", 22", or 24" Rear Wheels, Clothing Guards.	
Clothing Guards Not compatible with Flip-up Armrests, Armrest Wideners.			
<input type="checkbox"/> 1907310	For Single Post Armrests	Available for 10"-14" Height range only.	
<input type="checkbox"/> 1907320	For Dual Post Armrests		
Armrest Extensions			Not compatible with Flip-up or Single Post Armrests.
<input type="checkbox"/> 1907400	3" Extension - Forward	<input type="checkbox"/> 1907401	3" Extension - Rearward

Seat Pan

- 1908000 **No Seat Pan**
- 1908001 **Solid Aluminum Seat Pan** \$ 135

Seat Height

Don't forget **Select desired Seat-to-Floor Height** STF Height is adjustable within the range of the chosen Subframe Plate (Low/High).

	Low							High						
	RW	13"	14"	15"	16"	17"	18"	RW	15"	16"	17"	18"	19"	20"
6" Casters	12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
8" Casters	RW	13"	14"	15"	16"	17"	18"	RW	15"	16"	17"	18"	19"	20"
	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Casters

- Urethane Casters** 8" Casters are compatible only with 16" and 24" Rear Wheels.
- 1911061 6" x 1.25" Caster size 1911081 8" x 1.25" Caster size
 - 1911063 6" x 1.50" Caster size \$ 85
 - 1911084 8" x 2.00" Caster size \$ 135
- Pneumatic Casters** 8" Casters are compatible only with 16" and 24" Rear Wheels.
- 1911161 6" x 1.25" Caster size \$ 110
 - 1911181 8" x 1.25" Caster size \$ 120
 - 1911184 8" x 2.00" Caster size \$ 150
 - 1911260 **Frog Legs 6" x 1.40" Soft Roll Casters** \$ 95
 - 1911400 **Frog Legs Suspension Forks** Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster. \$ 395

Rear Wheels

- Composite Wheels**
- 1912012 12" Wheel size 1912016 16" Wheel size 1912020 20" Wheel size
 - 1912022 22" Wheel size 1912024 24" Wheel size
- Spoke Wheels** \$ 75
- 1912120 20" Wheel size 1912122 22" Wheel size 1912124 24" Wheel size
 - 1912224 **24" Heavy Duty Spoke Wheels** \$ 350
- Rear Wheel Position** Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.
In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.
- 1912401 Rear 1912402 Mid 1912403 Front

Axle Type

- 1913000 **Quick Release Rear Axles**
- 1913010 **Fixed Rear Axles**

Rear Tires

- 1914000 **Urethane Tires**
- 1914010 **Pneumatic Tires** May result in a finished Seat Height 1/2" greater than stated. \$ 60
- 1914020 **Pneumatic Tires with Flat-Free Inserts** May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels. \$ 95

Handrims		Handrims are only available with 20", 22", and 24" Rear Wheels.
<input type="checkbox"/>	1915000 No Handrims	
<input type="checkbox"/>	1915010 Standard Aluminum Handrims	
<input type="checkbox"/>	1915020 Plastic Coated Handrims	\$ 150

Wheel Locks & Brakes				
User Wheel Locks				
<input type="checkbox"/>	1916010 Push to Lock	<input type="checkbox"/>	1916020 Pull to Lock	
<input type="checkbox"/>	1916030 Attendant Wheel Locks (Push to Lock)	Not compatible with TTDS or the 300 & 450 lbs Weight Capacity Frames.		\$ 85
<input type="checkbox"/>	1916040 User & Attendant Wheel Locks	Available with 20", 22", or 24" Rear Wheels only. Not compatible with Drum Brakes.		\$ 125
<input type="checkbox"/>	2416300 Drum Brakes	Dual trigger activation (right/left side activated independently). Available with 20" Rear Wheels only.		\$ 795

Wheel Lock Options			
<input type="checkbox"/>	1916100 Wheel Lock Extensions	Not compatible with Attendant Wheel Lock.	\$ 60
<input type="checkbox"/>	1916200 Hill Holders	Requires Pneumatic Rear Tires with Flat-Free Inserts.	\$ 450

Anti-Tippers			
<input checked="" type="checkbox"/>	1917010 Rear Anti-Tippers	Required Safety Option.	\$ 75
<input type="checkbox"/>	1917020 Front Anti-Tippers		\$ 795
<input type="checkbox"/>	1917030 Heavy Duty Rear Anti-Tipper	Provides additional rearward stability for users with excessive movement. Not compatible with Attendant Wheel Locks, IV Pole, O2 Holder.	\$ 895

Back Posts					
Back Post Style					
<input type="checkbox"/>	1918010 Push Bar Back	<input type="checkbox"/>	1918020 Cane Handle Back	Seat Widths of 14"-22" only. Not compatible with Headrest Mounting Plate.	
Back Post Height		Measured from the Seat Pan to the Handle.			
<input type="checkbox"/>	1918425 Short (22")	<input type="checkbox"/>	1918431 Tall (25")		
Seat to Back Angle					
<input type="checkbox"/>	1919090 Set at 90°	<input type="checkbox"/>	1919098 Set at 98°	<input type="checkbox"/>	1919106 Set at 106°
<input type="checkbox"/>	1919114 Set at 114°	<input type="checkbox"/>	1919122 Set at 122°		

Back Options				
Angle/Height Adjustable Push Bar Handle		Not compatible with Headrest Mounting Plate or Cane Handle Back.	\$ 150	
<input type="checkbox"/>	1919010 Button on Outside	<input type="checkbox"/>	1919011 Button on Inside	
<input type="checkbox"/>	1919210 10" Backrest Extension	Requires the Angle/Height Adjustable Push Bar Handle.		\$ 150
<input type="checkbox"/>	1919212 Headrest Mounting Plate	Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.		\$ 95
<input type="checkbox"/>	1919220 Fold Down Back	Not compatible with Reclining Back.		
<input type="checkbox"/>	1919230 Spring Back Suspension	Not compatible with Reclining Back.		\$ 105
<input type="checkbox"/>	1919240 Seating Dynamics Dynamic Rocker Back Interface	Not compatible with Fold Down Back, Spring Back Suspension, 30° Recline, Cane Handle. Requires Adjustable Tension Upholstery or rigid Backrest. Finished Back Post Height 3" taller than selected Back Post Height.		\$ 2,995

Back Upholstery					
<input type="checkbox"/>	1920000 No Back Upholstery				
Standard Upholstery		Not available with Cane Handle Back.	\$ 195		
<input type="checkbox"/>	1920418 Short (18" Height)	<input type="checkbox"/>	1920421 Regular (21" Height)	<input type="checkbox"/>	1020424 Tall (24" Height)
<input type="checkbox"/>	1920450 Tension Adjustable Upholstery	Height of Upholstery is 16". Not available with Cane Handle Back.		\$ 295	

See next page for more options >

Position of Upholstery Select the desired position of the top of the Upholstery relative to the Seat Pan.

Style	22" Back Post								
	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

Style	25" Back Post								
	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tall							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Tilt Activation

Attendant Tilt only Select trigger mounting location below.

- 1923011 Right Back Post 1923012 Left Back Post

User Tilt only Not available on widths greater than 22". Full Length Arm Pads only. Select trigger mounting location below. \$ 100

- 1923021 Right Armrest 1923022 Left Armrest

Attendant & User Tilt Not available on widths greater than 22". Full Length Arm Pads only. Select trigger mounting locations below. \$ 195

- 1923031 Right Back Post 1923032 Left Back Post
 1923041 Right Armrest 1923042 Left Armrest

1923110 **Tilt Trigger Extension** \$ 50

- 1923400 **Power Tilt** 250 lbs Weight Capacity. 45° tilt range. Includes 2 batteries for off-chair charging. pdgmobility.com/power-tilt \$ 1,795

Minimum Seat Width is 16". Power Tilt is compatible with either O2 Holder or IV Pole on the same chair, not both.

- Toggle Switch 1923421 Installed on Right side 1923422 Installed on Left side
 Button Switch 1923431 Installed on Right side 1923432 Installed on Left side
 Stealth Egg Switch 1923440 Not installed

Optional for Egg Switch: 1923445 Swing Away Mounting Hardware Not compatible with all Headrest makes and models. \$ 625

Specify Headrest to be used: _____

Tilt Limit Optional. If no tilt limit is selected, the full tilt range will be available.

- 1923220 Set at 20° 1923230 Set at 30° 1923240 Set at 40°

Recline

- 1923310 **Reclining Backrest (Manual)** \$ 895

30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery.

Front Rigging

- 1925000 **No Front Rigging** 1925001 **Sector Blocks only**

Easy-Out™ Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame. pdgmobility.com/easy-out

- 1925410 Fixed 70° \$ 95

- 1925420 Universal Angle Adjustable (0°-110°) Set at angle: _____ (use increments of 10°) \$ 295

Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings.

Swing Away Front Rigging Legacy Design. Not compatible with Easy-Out.

- 1925130 Swing Away 60° 1925140 Swing Away 70°

- 1925200 Swing Away Elevating Legrests Includes Calf Pads. \$ 375

- 1925300 1" Offset Hangers for Swing Away \$ 75

- 1925110 **Standard Plug In 70° Front Rigging** Legacy Design. Not compatible with Easy-Out.

Footplates

- 1926010 **Composite Footplates**

- 1926020 **Aluminum Angle Adjustable Footplates** \$ 165

- 1926030 **One Piece Fixed Footplate** Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable. \$ 350


- 1926040 **One Piece Flip-up Footplate** Requires selection of Composite or Aluminum Angle Adjustable. \$ 395

- 1926050 **Interlocking Angle Adjustable Footplates** Available for Seat Widths up to 21" only. \$ 195

Footrest/Legrest Accessories

Calf Pads Not compatible with Gel Padding.		
<input type="checkbox"/>	1927012 For Easy-Out™ Front Rigging Includes reinforced premium Calf Pad mount.	\$ 200
<input type="checkbox"/>	1927010 For Swing Away Front Rigging	\$ 135
<input type="checkbox"/>	1927020 Gel Padding Not compatible with Calf Pad or Easy-Out Front Rigging.	\$ 125
Gel Knee Pads		
<input type="checkbox"/>	1927032 For Easy-Out™ Angle Adjustable Front Rigging	\$ 200
<input type="checkbox"/>	1927030 For Swing Away ELR Front Rigging	\$ 100
<input type="checkbox"/>	1927040 Length Adjustable Heel Loops	\$ 50
<input type="checkbox"/>	1927050 Bodypoint Padded Calf Strap	\$ 195
Miller's Dynamic Footrest Mounts Not compatible with Easy-Out Front Rigging.		
<input type="checkbox"/>	1927060 For Standard Plug In Front Rigging	\$ 445
<input type="checkbox"/>	1927065 For Swing Away 70° Front Rigging	\$ 575

Front Rigging Extension Tubes

Regular Extension Tubes		
<input type="checkbox"/>	1928010 Short 14"-18" for Easy-Out Front Rigging 13"-17" for Swing Away and Standard Plug In Front Rigging	
<input type="checkbox"/>	1928020 Long 18"-24" for Easy-Out Front Rigging 15"-19" for Swing Away and Standard Plug In Front Rigging	\$ 75
Extra Strength Extension Tubes (Solid) Available on Plug In Front Rigging only. Not compatible with Easy-Out Front Rigging.		
<input type="checkbox"/>	1928110 Standard (13"-17")	<input type="checkbox"/> 1928120 Long (15"-19")
	Lower Leg Length If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.	
	Seat Pan to Foot Plates: _____	Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.

Miscellaneous Options

<input type="checkbox"/>	1930080 Gel Hip Guides Reduces available Seat Width by 2".	\$ 395
<input type="checkbox"/>	1930110 IV Pole Not compatible with HD Anti Tipper, Attendant Wheel Locks.	\$ 295
<input type="checkbox"/>	1930120 O2 Holder Not compatible with HD Anti Tipper, Attendant Wheel Locks.	\$ 395
<input type="checkbox"/>	1930170 Adhesive Back Hook Strip, 2" x 30" For Cushion mounting.	\$ 10

On-Chair Seating & Positioning

 Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. pdgmobility.com/seating

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.