

Customer and Order Information

Purchase Order *	Tag
Dealer/Provider Name *	
Purchaser Name and Contact Information *	
Billing Address *	
Shipping Address * *required	

Frame Type

Elevation Ultralight Manual Wheelchair Extended delivery lead times may apply. This wheelchair is not recommended for children under 12 years.

<input checked="" type="checkbox"/> 2001000	Ultralight Manual Rigid Wheelchair	\$ 3,900
<input checked="" type="checkbox"/> 2001001	User Controlled Anterior/Posterior Back Angle Recline System	\$ 850
<input checked="" type="checkbox"/> 2001002	User Controlled Anterior Tilt & Elevation System	\$ 1,800

Frame Options

2002000 **Standard Base and Seat Frame** 250 lbs Weight Capacity.

User Weight Range Required for selection of Gas Struts (Force measured in Newton - N)

<input type="checkbox"/> 2002150	65 lbs to 90 lbs (150 N)	<input type="checkbox"/> 2002200	90 lbs to 140 lbs (200 N)
<input type="checkbox"/> 2002300	140 lbs to 200 lbs (300 N)	<input type="checkbox"/> 2002400	200 lbs to 250 lbs (400 N)

Seat Width

Standard Seat Widths

<input type="checkbox"/> 2003012	12" Seat Width	<input type="checkbox"/> 2003013	13" Seat Width	<input type="checkbox"/> 2003014	14" Seat Width
<input type="checkbox"/> 2003015	15" Seat Width	<input type="checkbox"/> 2003016	16" Seat Width	<input type="checkbox"/> 2003017	17" Seat Width
<input type="checkbox"/> 2003018	18" Seat Width	<input type="checkbox"/> 2003019	19" Seat Width		

Seat Depth

Standard Seat Depths

<input type="checkbox"/> 2004014	14" Seat Depth	<input type="checkbox"/> 2004015	15" Seat Depth	<input type="checkbox"/> 2004016	16" Seat Depth
<input type="checkbox"/> 2004017	17" Seat Depth	<input type="checkbox"/> 2004018	18" Seat Depth	<input type="checkbox"/> 2004019	19" Seat Depth
<input type="checkbox"/> 2004020	20" Seat Depth				

Frame Colour

PDG Premium Colours

<input type="checkbox"/> 2005000	Tracer Red	<input type="checkbox"/> 2005001	Gloss Black	<input type="checkbox"/> 2005002	Silver
<input type="checkbox"/> 2005003	Black Shimmer	<input type="checkbox"/> 2005004	Orange Fusion	<input type="checkbox"/> 2005005	Karl's Imagination
<input type="checkbox"/> 2005006	Blue Shimmer	<input type="checkbox"/> 2005007	Gloss White	<input type="checkbox"/> 2005008	Dark Titanium
<input type="checkbox"/> 2005009	Geyser Green	<input type="checkbox"/> 2005010	Lollypop Blue	<input type="checkbox"/> 2005011	Lollypop Red
<input type="checkbox"/> 2005012	Bubble Gum Pink	<input type="checkbox"/> 2005013	Berry Purple	<input type="checkbox"/> 2005014	Safety Yellow
<input type="checkbox"/> 2005200	Custom Colour	Contact PDG for specifications.			\$ 300

Arm Type

<input type="checkbox"/> 2006000	No Armrest	
<input type="checkbox"/> 2006001	Rigid Armrest	Not compatible with Push to Lock, Fold Down Back. \$ 150

Seat Option

<input type="checkbox"/> 2008110	Standard Seat Sling
<input type="checkbox"/> 2008100	No Seat Sling Third-party Seat Sling must be installed prior to wheelchair use.

Cushion

<input checked="" type="checkbox"/> 2009000	No Cushion
---	-------------------

Seat Height

24" Wheel Size					
<input type="checkbox"/> 2010018	18" STF	<input type="checkbox"/> 2010019	19" STF	<input type="checkbox"/> 2010020	20" STF
25" Wheel Size Not compatible with 4" Caster with Frog Leg Uni-tine Fork.					
<input type="checkbox"/> 2010118	18.5" STF	<input type="checkbox"/> 2010119	19.5" STF	<input type="checkbox"/> 2010120	20.5" STF
26" Wheel Size Not compatible with Frog Leg Uni-tine Fork.					
<input type="checkbox"/> 2010219	19" STF	<input type="checkbox"/> 1610220	20" STF	<input type="checkbox"/> 1610221	21" STF

Casters

Composite					
<input type="checkbox"/> 2011040	4" x 1.0" Composite Sport	<input type="checkbox"/> 2011042	4" x 1.5" Composite Softroll	<input type="checkbox"/> 2011052	5" x 1.5" Composite Softroll
					\$ 125
Aluminium					
<input type="checkbox"/> 2011140	4" x 1.0" Aluminium Sport	<input type="checkbox"/> 2011142	4" x 1.5" Aluminium Softroll	<input type="checkbox"/> 2011152	5" x 1.5" Aluminium Softroll
					\$ 195
Caster Fork & Fork Suspension					
<input type="checkbox"/> 2011410	Standard Aluminium Fork				
<input type="checkbox"/> 2011420	Frog Legs Uni-tine Fork				
<input type="checkbox"/> 2011421	Frog Legs Ultra-Sport Suspension Forks Compatible with 4" Caster only.				
					\$ 450

Rear Wheels

Spinergy Wheels					
<input type="checkbox"/> 2012310	Spinergy Wire				
<input type="checkbox"/> 2012311	Spinergy SPOX Everyday				\$ 745
<input type="checkbox"/> 2012312	Spinergy SPOX Sport				\$ 745
<input type="checkbox"/> 2012313	Spinergy LX Everyday				\$ 895
<input type="checkbox"/> 2012314	Spinergy LX Sport				\$ 895
<input type="checkbox"/> 2012315	Spinergy Lite Extreme LX Flexrims LX Wheels only. Not compatible with 26" Wheel size.				\$ 1,245
Select Wheel Camber					
<input type="checkbox"/> 2012400	0° Camber	<input type="checkbox"/> 2012402	2° Camber		
<input type="checkbox"/> 2012404	4° Camber	<input type="checkbox"/> 2012406	6° Camber		



Axle Type

<input checked="" type="checkbox"/> 2013000	Quick Release Rear Axles
---	---------------------------------

Rear Tires

<input type="checkbox"/> 2014110	Primo Pneumatic Tires	
<input type="checkbox"/> 2014120	Schwalbe Marathon Plus Evolution	\$ 125

Hand Rims

<input type="checkbox"/> 2015000	No Hand Rims	
<input type="checkbox"/> 2015100	Aluminum Clear Anodized	
<input type="checkbox"/> 2015110	Plastic Coated	\$ 135
<input type="checkbox"/> 2015120	Natural Fit with Standard Grip	\$ 425

Wheel Locks

<input type="checkbox"/> 2016010	Push to Lock Not compatible with Rigid Armrest.
<input type="checkbox"/> 2016090	Composite Scissor Locks

Anti-Tippers

<input type="checkbox"/>	2017010	Rear Anti-Tippers	\$ 145
--------------------------	---------	--------------------------	--------

Back Options

<input type="checkbox"/>	2019020	Cane Style Back Post	\$ 195
<input type="checkbox"/>	2019221	Rigid Back	
<input type="checkbox"/>	2019220	Fold Down Back Not compatible with Rigid Armrest.	\$ 150
Seat Back Height			
<input type="checkbox"/>	2019300	Short (8.5"-12") - Rigidizer Bar Height = 5.5"	
<input type="checkbox"/>	2019301	Medium (10.5"-16") - Rigidizer Bar Height = 7.5"	
<input type="checkbox"/>	2019302	Tall (16.5"-20") - Rigidizer Bar Height = 9.5"	
<input type="checkbox"/>	2019310	Custom Rigidizer Bar Height _____ Must be at least 4.5' above Seat Sling.	\$ 150
<input type="checkbox"/>	2019311	Deep Rigidizer Bar For installing third-party Backrest.	\$ 150

Back Upholstery

<input type="checkbox"/>	2020010	Standard Back Upholstery	
<input type="checkbox"/>	2020000	No Back Upholstery Third-party Back Upholstery must be installed prior to wheelchair use.	

Footrest

<input type="checkbox"/>	2025010	Standard Footrest	
<input type="checkbox"/>	2025011	High Mount Footrest	

Foot Plates

<input type="checkbox"/>	2026020	Aluminum Angle Adjustable	
<input type="checkbox"/>	2026060	Flip-up	

Calf Straps

<input type="checkbox"/>	2027050	Bodypoint - Padded Calf Strap	\$ 105
--------------------------	---------	--------------------------------------	--------

Positioning Straps

<input checked="" type="checkbox"/>	2029005	Standard Lap Belt	
Bodypoint - Hip Belts			
Push button with flat-mount end fitting. Recommended for mounting through slots in side guards, to allow "anti-shear" feature to be used.			
<input type="checkbox"/>	2029010	Center-Pull Auto Style	\$ 105
<input type="checkbox"/>	2029015	Center-Pull Auto Style, Padded	\$ 155
<input type="checkbox"/>	2029200	Bodypoint - EvoFlex Pelvic Stabilizer	\$ 275
Push button recommended for when anti-shear feature is not used. Mounted directly to the seat tubes using the included seat tube mounting bracket.			

Side Guards

<input type="checkbox"/>	2030400	No Side Guards	
Non-Lifting Side Guards			
<input type="checkbox"/>	2030410	Rigid - Small	
<input type="checkbox"/>	2030411	Rigid - Regular	\$ 295
<input type="checkbox"/>	2030415	Folding - Regular	\$ 350
Lifting (Anti-Shear) Side Guards			
<input type="checkbox"/>	2030420	Rigid - Regular	\$ 600
<input type="checkbox"/>	2030425	Folding - Regular	
<input type="checkbox"/>	2030426	Folding - Large Not compatible with Rigid Armrest.	\$ 650

Miscellaneous Options			
<input type="checkbox"/>	2030450	Alber Wheel Mount Adapter	\$ 600
Quad Triggers			
<input type="checkbox"/>	2030461	Inboard Mounting	\$ 115
<input type="checkbox"/>	2030462	Outboard Mounting	\$ 115
<input type="checkbox"/>	2030463	Thumb Loops	\$ 130

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.



Scan this code or go to pdgmobility.com/products/elevation to learn more about the **Elevation Ultralight Manual Wheelchair**



  PDG Mobility