

Customer and Order Information

Purchase Order *	Tag	v.0101
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in

Base Chair

Fuze T50n 50° Manual Tilt Wheelchair

- | | | |
|---|--|----------|
| <input checked="" type="checkbox"/> 2601000 | Fuze T50 50° Manual Tilt Wheelchair | \$ 4,300 |
| <input checked="" type="checkbox"/> 2601008 | Narrow Frame Upgrade Reduces overall chair width to seat width plus only 4". | \$ 200 |

Frame Type & Weight Capacity

Fixed Width Frame Select a Fixed Seat Width below.

- | | | |
|----------------------------------|---|--------|
| <input type="checkbox"/> 2602010 | 300 lbs Weight Capacity | |
| <input type="checkbox"/> 2602011 | 400 lbs Weight Capacity Not compatible with Transit Tie Down System. | |
| <input type="checkbox"/> 2602025 | Adjustable Width Frame 250 lbs Weight Capacity. Select an Adjustable Seat Width below. | \$ 100 |

Frame Options

- | | | |
|----------------------------------|---|--------|
| <input type="checkbox"/> 2602222 | Take-Apart Frame The Take-Apart frame allows attendants to take the Fuze apart in seconds without tools. pdgmobility.com/take-apart | \$ 500 |
|----------------------------------|---|--------|
- This feature reduces the weight of the heaviest piece of the wheelchair to 29 lbs or less. Perfect for transporting it in the trunk of a car or SUV. Not compatible with Power Tilt, Ventilator Tray, Ventilator Battery Tray, DRB, User Wheel Locks. Take-Apart rear crossbar is not width adjustable. PDG will provide one no-charge crossbar in a different width over the lifetime of each chair ordered with adjustable width frame (can be orderd with the chair or later). Note: Crossbar for the T50n is 4" shorter than the Seat Width.



- | | | |
|----------------------------------|---------------------------------------|--|
| <input type="checkbox"/> 2602300 | Transit Tie Down System (TTDS) | |
|----------------------------------|---------------------------------------|--|
- TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Tested to max. Weight Capacity of 300 lbs.

Seat Width

Overall chair width is +4.0" for 12" RW

Fixed Seat Widths Order with Fixed Frame.

- | | | | | | |
|----------------------------------|----------------|----------------------------------|----------------|----------------------------------|----------------|
| <input type="checkbox"/> 2603016 | 16" Seat Width | <input type="checkbox"/> 2603017 | 17" Seat Width | <input type="checkbox"/> 2603018 | 18" Seat Width |
| <input type="checkbox"/> 2603019 | 19" Seat Width | <input type="checkbox"/> 2603020 | 20" Seat Width | <input type="checkbox"/> 2603021 | 21" Seat Width |
| <input type="checkbox"/> 2603022 | 22" Seat Width | | | | |
| <input type="checkbox"/> 2603023 | 23" Seat Width | <input type="checkbox"/> 2603024 | 24" Seat Width | | \$ 375 |

Adjustable Seat Widths Order with Adjustable Frame.

- | | | | |
|----------------------------------|----------------------------|----------------------------------|----------------------------|
| <input type="checkbox"/> 2603216 | 16"-19" Range - Set at 16" | <input type="checkbox"/> 2603217 | 16"-19" Range - Set at 17" |
| <input type="checkbox"/> 2603218 | 16"-19" Range - Set at 18" | <input type="checkbox"/> 2603219 | 16"-19" Range - Set at 19" |
| <input type="checkbox"/> 2603319 | 19"-20" Range - Set at 19" | <input type="checkbox"/> 2603320 | 19"-20" Range - Set at 20" |



FUZE T50n

Order Form - Canada

Prices in Canadian Dollars. January 1, 2024

Seat Depth

Standard Seat Depths (Adjustable)

- | | | |
|---|---|---|
| <input type="checkbox"/> 2604015 15" Seat Depth | <input type="checkbox"/> 2604016 16" Seat Depth | <input type="checkbox"/> 2604017 17" Seat Depth |
| <input type="checkbox"/> 2604018 18" Seat Depth | <input type="checkbox"/> 2604019 19" Seat Depth | <input type="checkbox"/> 2604020 20" Seat Depth |

Custom Seat Depths

Requires Stroller Handle Back - Straight. Not compatible with Reclining or Fold-down Back.

\$ 400

- | | |
|---|---|
| <input type="checkbox"/> 2604021 21" Seat Depth | <input type="checkbox"/> 2604022 22" Seat Depth |
|---|---|

Frame Color

PDG Premium Colors

View color swatches on PDG's website. pdgmobility.com/colors

- | | | |
|---|---|---|
| <input type="checkbox"/> 2605001 Black Shimmer | <input type="checkbox"/> 2605015 Sky Blue | <input type="checkbox"/> 2605016 Poppy Red |
| <input type="checkbox"/> 2605003 Fine Texture Black | <input type="checkbox"/> 2605010 Lollipop Blue | <input type="checkbox"/> 2605011 Lollipop Red |
| <input type="checkbox"/> 2605017 Arctic Grey | <input type="checkbox"/> 2605018 Tropical Teal | <input type="checkbox"/> 2605019 Tangerine Orange |
| <input type="checkbox"/> 2605020 Olive Green | <input type="checkbox"/> 2605021 Lush Lilac | <input type="checkbox"/> 2605012 Bubble Gum Pink |
| <input type="checkbox"/> 2605022 Leaf Green | <input type="checkbox"/> 2605023 Boysenberry Purple | <input type="checkbox"/> 2605014 Daffodil Yellow |

- | | |
|---|---------------------------------|
| <input type="checkbox"/> 2605200 Custom Color | Contact PDG for specifications. |
|---|---------------------------------|

\$ 300

Armrests

- | |
|--|
| <input type="checkbox"/> 2606000 No Armrests |
|--|

Height Adjustable Armrests (Single Post)

\$ 295

- | | |
|--|--|
| <input type="checkbox"/> 2606114 10"-14" Range | <input type="checkbox"/> 2606116 12"-16" Range |
|--|--|

- | | |
|--|--------------------------|
| <input type="checkbox"/> 2606100 Fixed Height Armrests (12") | Full Length Armpad only. |
|--|--------------------------|



Armrest Height

If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.

Seat Pan to top of Armrest: _____ Does not include height of Armpad



Armrest Position

Indicate preference for Armrest position. If left blank, the default setting will be used.

In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

- | | |
|--------------------------------------|--|
| <input type="checkbox"/> 2606402 Mid | <input type="checkbox"/> 2606403 Front |
|--------------------------------------|--|

Locking Flip-up Armrests

Not recommended for users with excessive strength or tone. Not compatible with Fold-down Back.

\$ 295

- | | | |
|---|---|---|
| <input type="checkbox"/> 2606310 Set at 10" | <input type="checkbox"/> 2606311 Set at 11" | <input type="checkbox"/> 2606312 Set at 12" |
| <input type="checkbox"/> 2606313 Set at 13" | <input type="checkbox"/> 2606314 Set at 14" | |

Armpads

Regular Armpads

- | | |
|---|---|
| <input type="checkbox"/> 2607100 Desk Length - 2" x 10" | <input type="checkbox"/> 2607101 Full Length - 2" x 14" |
|---|---|

Waterfall Armpads

\$ 115

- | | |
|---|---|
| <input type="checkbox"/> 2607105 Desk Length - 2" x 10" | <input type="checkbox"/> 2607106 Full Length - 2" x 14" |
|---|---|

Wide Armpads

May interfere with wheels in some tilt configurations.

\$ 200

- | | |
|---|---|
| <input type="checkbox"/> 2607110 Desk Length - 3" x 10" | <input type="checkbox"/> 2607111 Full Length - 3" x 12" |
|---|---|

Gel Armpads

May interfere with wheels in some tilt configurations.

\$ 135

- | | |
|---|---|
| <input type="checkbox"/> 2607120 Full Length - 2" x 12" | <input type="checkbox"/> 2607121 Full Length - 3.5" x 12" |
|---|---|

Contour Arm Trough with Hand Pad

\$ 125

Not compatible with Lap Tray, Locking Flip-up Armrests. May interfere with wheels in some tilt configurations.

- | | |
|---|--|
| <input type="checkbox"/> 2607130 Right (4.5" x 18") | <input type="checkbox"/> 2607131 Left (4.5" x 18") |
|---|--|

Armrest Options

Armrest Wideners

Move Armrests outward by up to 1" on each side. Seat size does not change.

- | | |
|--|--------|
| <input type="checkbox"/> 2607210 1" Total Additional Width | \$ 95 |
| <input type="checkbox"/> 2607211 2" Total Additional Width | \$ 195 |



FUZE T50n


Order Form - Canada

Prices in Canadian Dollars. January 1, 2024

Seat Pan

- ☐ 2608000 **No Seat Pan**
- ☐ 2608001 **Solid Aluminum Seat Pan** \$ 135

Seat Height

 **Don't forget** **Select desired Seat-to-Floor Height** STF Height is adjustable.

12" Rear Wheels					
Caster size	16"	17"	18"	19"	20"
5"	<input type="checkbox"/>	<input type="checkbox"/>			
6"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
8"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
6" Soft Roll	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

Casters

Urethane Casters


- ☐ 2611051 5" x 1.25" Caster size ☐ 2611061 6" x 1.25" Caster size ☐ 2611081 8" x 1.25" Caster size
- ☐ 2611063 6" x 1.50" Caster size \$ 85
- ☐ 2611084 8" x 2.00" Caster size \$ 135

Pneumatic Casters

- ☐ 2611161 6" x 1.25" Caster size \$ 110
- ☐ 2611181 8" x 1.25" Caster size \$ 120
- ☐ 2611184 8" x 2.00" Caster size \$ 150
- ☐ 2611260 **Frog Legs 6" x 1.50" Aluminum Soft-Roll Casters** \$ 275

Rear Wheels

☒ 2612012 **12" Composite Wheels**

 **Rear Wheel Position** Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.
In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

☐ 2612401 Rear More stable ☐ 2612402 Mid ☐ 2612403 Front Less stable

Axle Type

- ☐ 2613000 **Quick Release Rear Axles**
- ☐ 2613010 **Fixed Rear Axles**

Rear Tires

- ☐ 2614000 **Urethane Tires**
- ☐ 2614010 **Pneumatic Tires** May result in a finished Seat Height 1/2" greater than stated. \$ 60
- ☐ 2614020 **Pneumatic Tires with Flat-Free Inserts** May result in a finished Seat Height 1/2" greater than stated. \$ 95

Handrims

- ☒ 2615000 **No Handrims**

Wheel Locks & Brakes

- User Wheel Locks** Not compatible with Take-Apart Frame. \$ 150
- ☐ 2616010 Push to Lock ☐ 2616020 Pull to Lock
- ☐ 2616030 **Attendant Wheel Locks (Push to Lock)**
- ☐ 2616005 **Foot Operated Wheel Lock** Rear-mounted unilateral single-activation lock. Not compatible with Ventilator Tray. \$ 595
- ☐ 2616300 **Drum Brakes** Available with Urethane Tires only. Compatible with 17.5", 18.5", 19.5" STF Height only. **COMING SOON** \$ 795
- User-activated Unilateral Wheel Lock** Pull to Lock. Available up to 22" Seat Width only. \$ 795
- ☐ 2616310 Left Activation ☐ 2616311 Right Activation



FUZE T50n

Order Form - Canada

Prices in Canadian Dollars. January 1, 2024

Wheel Lock Options

☐ 2616100 **Wheel Lock Extensions** \$ 60

Anti-Tippers

☒ 2617010 **Rear Anti-Tippers** Required Safety Option. \$ 75

Back Posts

Back Post Style

☐ 2618100 Stroller Handle Back ☐ 2618105 Cane Handle Back Not compatible with Headrest Mounting Plate.

Back Post Type

☐ 2618200 Straight ☐ 2618205 8° Bend

Back Post Height

 Measured from the Seat Pan to the Handle.

☐ 2618420 Short (17") ☐ 2618425 Tall (23")

Back Post Folding

 Fold-down Back is not compatible with Dynamic Rocker Back Interface, Reclining Backrest.

☐ 2619220 Fold-down ☐ 2619221 Non-folding

Seat to Back Angle

☐ 2619090 Set at 90° ☐ 2619096 Set at 98° ☐ 2619102 Set at 106°

☐ 2619108 Set at 114° ☐ 2619114 Set at 122°

Back Options

Angle/Height Adjustable Stroller Handle

 Not compatible with Headrest Mounting Plate or Cane Handle Back.

☐ 2619010 Button on Outside ☐ 2619011 Button on Inside \$ 150

Removable Cane Handle Push Bar

 Not compatible with Stroller Handle Back.

☐ 2619050 \$ 140

Endomorph Back

 Backrest will be narrower than Seat by total of 2x the Offset selected below. Minimum width of Backrest is 12".

☐ 2618102 1" Offset per side ☐ 2618103 1.5" Offset per side ☐ 2618104 2" Offset per side \$ 395

Headrest Mounting Plate

 Not compatible with Angle/Height Adjustable Stroller Handle, Cane Handle Back. Not required if mounting to a solid Back.

☐ 2619212 \$ 95

Seating Dynamics Dynamic Rocker Back Interface

☐ 2619240 \$ 2,995

Not compatible with Take-Apart Frame, Fold-down Back, 30° Recline, Standard Upholstery. Only available with Stroller Handle - 8° Bend.

Finished Back Post Height is 7" taller than selected Back Post Height. Maximum compatible Seat Depth is 20". 90° and 98° Seat Angle only.

Back Upholstery

☐ 2620000 **No Back Upholstery**

Standard Upholstery

 Not available with Cane Handle Back.

☐ 2620418 Short (18" Height) ☐ 2620421 Regular (21" Height) \$ 195

Tension Adjustable Upholstery

 Height of Upholstery is 16". Not available with Cane Handle Back.

☐ 2620450 \$ 295

Position of Upholstery

 Select the desired position of the top of the Upholstery relative to the Seat Pan.

17" Back Post									
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>								

23" Back Post									
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
Regular					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
Tension Adj				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

Tilt Activation

Attendant Tilt only

 Select trigger mounting location below.

☐ 2623011 Right Back Post ☐ 2623012 Left Back Post

See next page for more options >



FUZE T50n

Order Form - Canada

Prices in Canadian Dollars. January 1, 2024

Power Tilt 40° tilt range. Includes 2 batteries for off-chair charging. pdgmobility.com/power-tilt

<input type="checkbox"/>	2623400	Standard Power Tilt (250 lbs Weight Capacity)			\$ 1,795
<input type="checkbox"/>	2623401	Standard Power Tilt (350 lbs Weight Capacity)	Requires Fixed Frame.		\$ 2,095
<input type="checkbox"/>	2623410	Power Tilt with Manual Tilt Override (250 lbs Weight Capacity)	Allows for instant manual release of tilt in case of power tilt failure.		\$ 2,295
<input type="checkbox"/>	2623411	Power Tilt with Manual Tilt Override (350 lbs Weight Capacity)	Allows for instant manual release of tilt in case of power tilt failure. Requires Fixed Frame.		\$ 2,595
		Toggle Switch	<input type="checkbox"/> 2623421 Installed on Right side	<input type="checkbox"/> 2623422 Installed on Left side	
		Button Switch	<input type="checkbox"/> 2623431 Installed on Right side	<input type="checkbox"/> 2623432 Installed on Left side	
		Stealth Egg Switch	<input type="checkbox"/> 2623440 Not installed		
		Optional for Egg Switch:	<input type="checkbox"/> 2623445 Swing Away Mounting Hardware	Not compatible with all Headrest makes and models.	\$ 625
			Specify Headrest to be used: _____		

Tilt Limit Optional. If no tilt limit is selected, the full tilt range will be available.

<input type="checkbox"/>	2623220	Set at 20°	<input type="checkbox"/>	2623230	Set at 30°	<input type="checkbox"/>	2623240	Set at 40°
--------------------------	---------	------------	--------------------------	---------	------------	--------------------------	---------	------------

Recline

<input type="checkbox"/>	2623310	Reclining Backrest (Manual)		\$ 895
30° maximum recline. Not compatible with Cane Handle Back, Fold-down Back. Not available with 22" Seat Depth. Finished Back Post Height is 3" taller than selected Back Post Height.				

Front Rigging

<input type="checkbox"/>	2625000	No Front Rigging		<input type="checkbox"/>	2625001	Sector Blocks only	
Easy-Out™ Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame. pdgmobility.com/easy-out							
<input type="checkbox"/>	2625410	Fixed 70°				\$	95
<input type="checkbox"/>	2625420	Universal Angle Adjustable (0°-110°)		Set at angle: _____ (use increments of 10°)		\$	295
Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings.							
Swing Away Front Rigging Legacy Design. Not compatible with Easy-Out. Requires pre-tilt of 5°.							
<input type="checkbox"/>	2625130	Swing Away 60°		<input type="checkbox"/>	2625140	Swing Away 70°	
<input type="checkbox"/>	2625160	Swing Away 90°		Includes 11"-16" Extension Tubes. Not compatible with 8" Casters.			\$ 75
<input type="checkbox"/>	2625200	Swing Away Elevating Legrests		Includes Calf Pads.			\$ 375
<input type="checkbox"/>	2625300	1" Offset Hangers for Swing Away					\$ 75
<input type="checkbox"/>	2625110	Standard Plug In 70° Front Rigging		Legacy Design. Not compatible with Easy-Out.			

Footplates

<input type="checkbox"/>	2626010	Composite Footplates	
<input type="checkbox"/>	2626020	Aluminum Angle Adjustable Footplates	\$ 165
<input type="checkbox"/>	2626030	One Piece Fixed Footplate	Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable. \$ 350
<input type="checkbox"/>	2626040	One Piece Flip-up Footplate	Requires selection of Composite or Aluminum Angle Adjustable. \$ 395
<input type="checkbox"/>	2626050	Interlocking Angle Adjustable Footplates	Available for Seat Widths up to 21" only. \$ 195

Footrest/Legrest Accessories

Calf Pad Not compatible with Gel Padding.				
<input type="checkbox"/>	2627012	For Easy-Out™ Front Rigging	Includes reinforced premium Calf Pad mount.	\$ 200
<input type="checkbox"/>	2627010	For Swing Away Front Rigging		\$ 135
<input type="checkbox"/>	2627020	Gel Padding	Not compatible with Calf Pad or Easy-Out Front Rigging.	\$ 125
Gel Knee Pad				
<input type="checkbox"/>	2627032	For Easy-Out™ Angle Adjustable Front Rigging		\$ 200
<input type="checkbox"/>	2627030	For Swing Away ELR Front Rigging		\$ 100
<input type="checkbox"/>	2627040	Length Adjustable Heel Loops		\$ 50
<input type="checkbox"/>	2627050	Bodypoint Padded Calf Strap	Sizes greater than 22" are not padded.	\$ 195



FUZE T50n

Order Form - Canada

Prices in Canadian Dollars. January 1, 2024

Front Rigging Extension Tubes

Regular Extension Tubes

No selection required for 90° Front Rigging

- ☐ 2628010 Short 14"-18" for Easy-Out Front Rigging 13"-17" for Swing Away and Standard Plug In Front Rigging
- ☐ 2628020 Long 18"-24" for Easy-Out Front Rigging 15"-19" for Swing Away and Standard Plug In Front Rigging

\$ 75



Lower Leg Length

If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.

Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.

Miscellaneous Options

Stealth Flip Down Abductor

Requires Solid Seat Pan.

\$ 650

- ☐ 2630060 Small Pad 3"x4"x3"x2" size ☐ 2630062 Large Pad 5"x5"x4"x3" size

QR Bracket - Accessory Mount

\$ 200

Gel Hip Guides

Reduces available Seat Width by 2".

\$ 495

Hip Belt Mount

Provides improved mounting options for Bodypoint Hip Belt or other third-party belts.

\$ 55

Ventilator Tray

Reduces ground clearance. Chair may be unstable when unoccupied. Minimum Seat Width is 20".
Not compatible with Power Tilt, Foot Operated Wheel Locks.

\$ 295

Ventilator Bracket

Compatible with LTV1 (laptop) ventilator.

\$ 295

Clear Lap Tray (Slide On)

Available for Seat Widths of 17"-22" only. Not compatible with Contour Arm Trough.

\$ 490

Q'Straint Lap Belt

Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" Belt Length.

\$ 245

Wheelchair Backpack

Black with PDG Logo. 17"x14" size.

\$ 145

Adhesive Back Hook Strip, 2" x 30"

For Cushion mounting.

\$ 10

On-Chair Seating & Positioning

Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. pdgmobility.com/seating

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.