

Customer and Order Information	
Purchase Order *	Tag
Dealer/Provider Name *	
Purchaser Name and Contact Information *	
Billing Address *	
Shipping Address * <span style="float: right;">*required</span>	

Frame Type	
<b>Elevation Ultralight Manual Wheelchair</b> Extended delivery lead times may apply. This wheelchair is not recommended for children under 12 years.	
<input checked="" type="checkbox"/> 2001000 Ultralight Manual Rigid Wheelchair	<i>K0005 / E1235 (Pediatric sizes)</i> \$ 3,900
<input checked="" type="checkbox"/> 2001001 User Controlled Anterior/Posterior Back Angle Recline System	<i>E1225</i> \$ 850
<input checked="" type="checkbox"/> 2001002 User Controlled Anterior Tilt & Elevation System	<i>K0108</i> \$ 1,800

Frame Options	
<input checked="" type="checkbox"/> 2002000 <b>Standard Base and Seat Frame</b> 250 lbs Weight Capacity.	
<b>User Weight Range</b> Required for selection of Gas Struts (Force measured in Newton - N)	
<input type="checkbox"/> 2002150 65 lbs to 90 lbs (150 N)	<input type="checkbox"/> 2002200 90 lbs to 140 lbs (200 N)
<input type="checkbox"/> 2002300 140 lbs to 200 lbs (300 N)	<input type="checkbox"/> 2002400 200 lbs to 250 lbs (400 N)

Seat Width		
<b>Standard Seat Widths</b>		
<input type="checkbox"/> 2003012 12" Seat Width	<input type="checkbox"/> 2003013 13" Seat Width	<input type="checkbox"/> 2003014 14" Seat Width
<input type="checkbox"/> 2003015 15" Seat Width	<input type="checkbox"/> 2003016 16" Seat Width	<input type="checkbox"/> 2003017 17" Seat Width
<input type="checkbox"/> 2003018 18" Seat Width	<input type="checkbox"/> 2003019 19" Seat Width	

Seat Depth		
<b>Standard Seat Depths</b>		
<input type="checkbox"/> 2004014 14" Seat Depth	<input type="checkbox"/> 2004015 15" Seat Depth	<input type="checkbox"/> 2004016 16" Seat Depth
<input type="checkbox"/> 2004017 17" Seat Depth	<input type="checkbox"/> 2004018 18" Seat Depth	<input type="checkbox"/> 2004019 19" Seat Depth
<input type="checkbox"/> 2004020 20" Seat Depth		<i>E2203</i> \$ 135

Frame Color		
<b>PDG Premium Colors</b>		
<input type="checkbox"/> 2005000 Tracer Red	<input type="checkbox"/> 2005001 Gloss Black	<input type="checkbox"/> 2005002 Silver
<input type="checkbox"/> 2005003 Black Shimmer	<input type="checkbox"/> 2005004 Orange Fusion	<input type="checkbox"/> 2005005 Karl's Imagination
<input type="checkbox"/> 2005006 Blue Shimmer	<input type="checkbox"/> 2005007 Gloss White	<input type="checkbox"/> 2005008 Dark Titanium
<input type="checkbox"/> 2005009 Geyser Green	<input type="checkbox"/> 2005010 Lollypop Blue	<input type="checkbox"/> 2005011 Lollypop Red
<input type="checkbox"/> 2005012 Bubble Gum Pink	<input type="checkbox"/> 2005013 Berry Purple	<input type="checkbox"/> 2005014 Safety Yellow
<input type="checkbox"/> 2005200 <b>Custom Color</b> Contact PDG for specifications.		\$ 300

Arm Type	
<input type="checkbox"/> 2006000 <b>No Armrest</b>	
<input type="checkbox"/> 2006001 <b>Rigid Armrest</b> Not compatible with Push to Lock, Fold Down Back.	\$ 150

Seat Option	
<input type="checkbox"/> 2008110 <b>Standard Seat Sling</b>	
<input type="checkbox"/> 2008100 <b>No Seat Sling</b> Third-party Seat Sling must be installed prior to wheelchair use.	

## Cushion

2009000 **No Cushion**

## Seat Height

### 24" Wheel Size

2010018 18" STF  2010019 19" STF  2010020 20" STF

### 25" Wheel Size

Not compatible with 4" Caster with Frog Leg Uni-tine Fork.

2010118 18.5" STF  2010119 19.5" STF  2010120 20.5" STF

### 26" Wheel Size

Not compatible with Frog Leg Uni-tine Fork.

2010219 19" STF  1610220 20" STF  1610221 21" STF

## Casters

### Composite

2011040 4" x 1.0" Composite Sport  2011042 4" x 1.5" Composite Softroll  2011052 5" x 1.5" Composite Softroll

### Aluminium

2011140 4" x 1.0" Aluminum Sport  2011142 4" x 1.5" Aluminum Softroll  2011152 5" x 1.5" Aluminum Softroll

\$ 125

### Caster Fork & Fork Suspension

2011410 Standard Aluminum Fork

2011420 Frog Legs Uni-tine Fork

\$ 195

2011421 Frog Legs Ultra-Sport Suspension Forks Compatible with 4" Caster only.

E1015 \$ 450

## Rear Wheels

### Spinergy Wheels

2012310 Spinergy Wire

2012311 Spinergy SPOX Everyday

\$ 745

2012312 Spinergy SPOX Sport

\$ 745

2012313 Spinergy LX Everyday

\$ 895

2012314 Spinergy LX Sport

\$ 895

2012315 **Spinergy Lite Extreme LX Flexrims** LX Wheels only. Not compatible with 26" Wheel size.

\$ 1,245



### Select Wheel Camber

2012400 0° Camber

2012402 2° Camber

2012404 4° Camber

2012406 6° Camber

## Axle Type

2013000 **Quick Release Rear Axles**

## Rear Tires

2014110 **Primo Pneumatic Tires**

E2211 & E2212

2014120 **Schwalbe Marathon Plus Evolution**



E2211 & E2212 \$ 125

## Hand Rims

2015000 **No Hand Rims**

2015100 **Aluminum Clear Anodized**

2015110 **Plastic Coated**

\$ 135

2015120 **Natural Fit with Standard Grip**

\$ 425

## Wheel Locks

2016010 **Push to Lock** Not compatible with Rigid Armrest.

2016090 **Composite Scissor Locks**

## Anti-Tippers

<input type="checkbox"/>	2017010	<b>Rear Anti-Tippers</b>	E0971	\$	145
--------------------------	---------	--------------------------	-------	----	-----

## Back Options

<input type="checkbox"/>	2019020	<b>Cane Style Back Post</b>		\$	195
<input type="checkbox"/>	2019221	<b>Rigid Back</b>			
<input type="checkbox"/>	2019220	<b>Fold Down Back</b>	Not compatible with Rigid Armrest.	\$	150
<b>Seat Back Height</b>					
<input type="checkbox"/>	2019300	Short (8.5"-12") - Rigidizer Bar Height = 5.5"			
<input type="checkbox"/>	2019301	Medium (10.5"-16") - Rigidizer Bar Height = 7.5"			
<input type="checkbox"/>	2019302	Tall (16.5"-20") - Rigidizer Bar Height = 9.5"			
<input type="checkbox"/>	2019310	<b>Custom Rigidizer Bar Height</b> _____ Must be at least 4.5' above Seat Sling.		\$	150
<input type="checkbox"/>	2019311	<b>Deep Rigidizer Bar</b>	For installing third-party Backrest.	\$	150

## Back Upholstery

<input type="checkbox"/>	2020010	<b>Standard Back Upholstery</b>	E2611		
<input type="checkbox"/>	2020000	<b>No Back Upholstery</b>	Third-party Back Upholstery must be installed prior to wheelchair use.		

## Footrest

<input type="checkbox"/>	2025010	<b>Standard Footrest</b>			
<input type="checkbox"/>	2025011	<b>High Mount Footrest</b>			


## Foot Plates

<input type="checkbox"/>	2026020	<b>Aluminum Angle Adjustable</b>	K0040		
<input type="checkbox"/>	2026060	<b>Flip-up</b>			

## Calf Straps

<input type="checkbox"/>	2027050	<b>Bodypoint - Padded Calf Strap</b>	K0038	\$	105
--------------------------	---------	--------------------------------------	-------	----	-----

## Positioning Straps

<input checked="" type="checkbox"/>	2029005	<b>Standard Lap Belt</b>			
<b>Bodypoint - Hip Belts</b>					
Push button with flat-mount end fitting. Recommended for mounting through slots in side guards, to allow "anti-shear" feature to be used.					
<input type="checkbox"/>	2029010	Center-Pull Auto Style		E0978	\$ 105
<input type="checkbox"/>	2029015	Center-Pull Auto Style, Padded		E0978	\$ 155
<input type="checkbox"/>	2029200	<b>Bodypoint - EvoFlex Pelvic Stabilizer</b>		E0978	\$ 275
Push button recommended for when anti-shear feature is not used. Mounted directly to the seat tubes using the included seat tube mounting bracket.					

## Side Guards

<input type="checkbox"/>	2030400	<b>No Side Guards</b>			
<b>Non-Lifting Side Guards</b>					
<input type="checkbox"/>	2030410	Rigid - Small	<input type="checkbox"/>	2030411	Rigid - Regular
					\$ 295
<input type="checkbox"/>	2030415	Folding - Regular			\$ 350
<b>Lifting (Anti-Shear) Side Guards</b>					
<input type="checkbox"/>	2030420	Rigid - Regular			\$ 600
<input type="checkbox"/>	2030425	Folding - Regular	<input type="checkbox"/>	2030426	Folding - Large
					Not compatible with Rigid Armrest.
					\$ 650

Miscellaneous Options			
<input type="checkbox"/>	2030450	<b>Alber Wheel Mount Adapter</b>	A9900 \$ 600
<b>Quad Triggers</b>			
<input type="checkbox"/>	2030461	Inboard Mounting	\$ 115
<input type="checkbox"/>	2030462	Outboard Mounting	\$ 115
<input type="checkbox"/>	2030463	Thumb Loops	\$ 130

**Order Notes**

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions.  
Use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.



Scan this code or go to  
[pdgmobility.com/products/elevation](http://pdgmobility.com/products/elevation)  
to learn more about the  
**Elevation Ultralight Manual Wheelchair**

