

Customer and Order Information		
Purchase Order *	Tag	v.070118
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Billing Address *		
Shipping Address *		*required

Frame Type		
<input checked="" type="checkbox"/> 1601000	Fuze T20 20° Manual Tilt Wheelchair	E1161 \$ 2,995

Frame Options		
<input type="checkbox"/> 1602010	Fixed Base and Seat Frame 250 lbs Weight Capacity. Select Fixed Seat Width below.	
<input type="checkbox"/> 1602025	Adjustable Base and Seat Frame 250 lbs Weight Capacity. Select Adjustable Seat Width below.	\$ 100
<input type="checkbox"/> 1602300	Transit Tie Down System (TTDS) TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details.	
User Weight Range Required for selection of Gas Struts (Force measured in Newton - N)		
<input type="checkbox"/> 1602406	Up to 150 lbs (600 N)	
<input type="checkbox"/> 1602410	150 lbs to 200 lbs (1000 N)	
<input type="checkbox"/> 1602412	200 lbs to 250 lbs (1200 N)	

Seat Width		
		Overall chair width is +9.5" for 12/16" RW and +10.5" for 20/22/24" RW
Fixed Seat Widths Order with Fixed Base and Seat Frame.		
<input type="checkbox"/> 1603015	15" Seat Width	
<input type="checkbox"/> 1603016	16" Seat Width	
<input type="checkbox"/> 1603017	17" Seat Width	
<input type="checkbox"/> 1603018	18" Seat Width	
<input type="checkbox"/> 1603019	19" Seat Width	
<input type="checkbox"/> 1603020	20" Seat Width	
<input type="checkbox"/> 1603021	21" Seat Width	
<input type="checkbox"/> 1603022	22" Seat Width	E2201 \$ 325
Adjustable Seat Widths Order with Adjustable Base and Seat Frame.		
<input type="checkbox"/> 1603112	12"-15" Range - Set at 12"	
<input type="checkbox"/> 1603113	12"-15" Range - Set at 13"	
<input type="checkbox"/> 1603114	12"-15" Range - Set at 14"	
<input type="checkbox"/> 1603115	12"-15" Range - Set at 15"	
<input type="checkbox"/> 1603215	15"-20" Range - Set at 15"	
<input type="checkbox"/> 1603216	15"-20" Range - Set at 16"	
<input type="checkbox"/> 1603217	15"-20" Range - Set at 17"	
<input type="checkbox"/> 1603218	15"-20" Range - Set at 18"	
<input type="checkbox"/> 1603219	15"-20" Range - Set at 19"	
<input type="checkbox"/> 1603220	15"-20" Range - Set at 20"	E2201 \$ 325
<input type="checkbox"/> 1603321	21"-22" Range - Set at 21"	
<input type="checkbox"/> 1603322	21"-22" Range - Set at 22"	E2201 \$ 325

Seat Depth		
Standard Seat Depths (Adjustable)		
<input type="checkbox"/> 1604015	15" Seat Depth	
<input type="checkbox"/> 1604016	16" Seat Depth	
<input type="checkbox"/> 1604017	17" Seat Depth	
<input type="checkbox"/> 1604018	18" Seat Depth	
<input type="checkbox"/> 1604019	19" Seat Depth	
<input type="checkbox"/> 1604020	20" Seat Depth	E2203 \$ 475

Frame Color		
PDG Premium Colors		
<input type="checkbox"/> 1605000	Tracer Red	
<input type="checkbox"/> 1605001	Gloss Black	
<input type="checkbox"/> 1605002	Silver	
<input type="checkbox"/> 1605003	Black Shimmer	
<input type="checkbox"/> 1605004	Orange Fusion	
<input type="checkbox"/> 1605005	Karl's Imagination	
<input type="checkbox"/> 1605006	Blue Shimmer	
<input type="checkbox"/> 1605007	Gloss White	
<input type="checkbox"/> 1605008	Dark Titanium	
<input type="checkbox"/> 1605009	Geyser Green	
<input type="checkbox"/> 1605010	Lollypop Blue	
<input type="checkbox"/> 1605011	Lollypop Red	
<input type="checkbox"/> 1605012	Bubble Gum Pink	
<input type="checkbox"/> 1605013	Berry Purple	
<input type="checkbox"/> 1605014	Safety Yellow	
<input type="checkbox"/> 1605200	Custom Color Contact PDG for specifications.	\$ 300

Arm Type			
<input type="checkbox"/>	1606000	No Armrest	
<input type="checkbox"/>	1606100	Single Post Fixed Height (12")	
		Single Post Adjustable Height	8"-11" not compatible with 22" and 24" Rear Wheels. E0973 \$ 255
<input type="checkbox"/>	1606111	8"-11"	<input type="checkbox"/> 1606114 10"-14"
<input type="checkbox"/>	1606200	Dual Post Fixed Height (12")	Full Length Armpad only.
		Dual Post Adjustable Height	8"-11" not compatible with 22" and 24" Rear Wheels. E0973 \$ 255
<input type="checkbox"/>	1606211	8"-11"	<input type="checkbox"/> 1606214 10"-14" <input type="checkbox"/> 1606216 12"-16"
		Locking Flip-up	Not recommended for users with excessive strength or tone. Not compatible with Fold Down Back, Back Upholstery. E0973 \$ 255
<input type="checkbox"/>	1606310	Set at 10"	<input type="checkbox"/> 1606311 Set at 11" <input type="checkbox"/> 1606312 Set at 12"
<input type="checkbox"/>	1606313	Set at 13"	<input type="checkbox"/> 1606314 Set at 14"

Armrests			
<input type="checkbox"/>	1607095	Curved Armpads	Not compatible with Lap Tray, Locking Flip-up Armrests, Dual Post Armrests, 22" and 24" Rear Wheels.
		Regular Armpads	
<input type="checkbox"/>	1607100	Desk Length - 2" x 10"	<input type="checkbox"/> 1607101 Full Length - 2" x 14"
<input type="checkbox"/>	1607102	Full Length Premium Grade - 2" x 14"	\$ 35
		Wide Armpads	May interfere with wheels in some tilt configurations. \$ 115
<input type="checkbox"/>	1607110	Desk Length - 3" x 10"	<input type="checkbox"/> 1607111 Full Length - 3" x 12"
		Gel Armpads	May interfere with wheels in some tilt configurations. \$ 125
<input type="checkbox"/>	1607120	Full Length - 2" x 12"	<input type="checkbox"/> 1607121 Full Length - 3.5" x 12"
		Contour Arm Trough with Hand Pad	May interfere with wheels in some tilt configurations. E2209 \$ 125
<input type="checkbox"/>	1607130	Right (4.5" x 18")	<input type="checkbox"/> 1607131 Left (4.5" x 18")

Seat Pan			
<input type="checkbox"/>	1608000	No Seat Pan	
<input type="checkbox"/>	1608001	Solid Aluminum Seat Pan	E2231 \$ 120

Cushion			
<input type="checkbox"/>	1609000	No Cushion	
		nxt BioFit	
<input type="checkbox"/>	1609101	14"-20" width, 16"-20" depth	E2607 \$ 595
<input type="checkbox"/>	1609102	22" width, 16"-20" depth	E2608 \$ 625
		Matrx VI	
<input type="checkbox"/>	1609111	14"-20" width, 16"-20" depth	E2607 \$ 550
<input type="checkbox"/>	1609112	22" width, 18"-20" depth	E2608 \$ 600
<input type="checkbox"/>	1609115	HD - 18"-22" width, 18"-20" depth	E2608 \$ 1,150
		ROHO Quadro Select High Profile	
<input type="checkbox"/>	1609121	14"-20" width, 16"-20" depth	E2624 \$ 650
<input type="checkbox"/>	1609122	22" width, 16"-20" depth	E2625 \$ 750
		ROHO High Profile	
<input type="checkbox"/>	1609125	14"-20" width, 16"-20" depth	E2622 \$ 650
<input type="checkbox"/>	1609126	22" width, 16"-20" depth	E2623 \$ 750
		Varilite Reflex	
<input type="checkbox"/>	1609131	14"-20" width, 16"-20" depth	E2603 \$ 375
		Cushion Size	
<input type="checkbox"/>	1609400	Specify Size (Width" x Depth")	

Seat Height

Select desired Seat-to-Floor Height

Standard T20 Options

RW	5" Casters						
12"	13"						
16"	13"	14"	15"				
20"	13"	14"	15"	16"	17"		
22"		14"	15"	16"	17"		
24"			15"	16"	17"		

RW	6" Casters & 6" Forks						
16"		14"	15"				
20"		14"	15"	16"	17"		
22"		14"	15"	16"	17"	18"	
24"			15"	16"	17"	18"	

RW	6" Casters & 8" Forks						
16"			15"				
20"			15"	16"	17"		
22"			15"	16"	17"	18"	
24"			15"	16"	17"	18"	19"

RW	8" Casters & 8" Forks						
20"				16"	17"		
22"				16"	17"	18"	
24"				16"	17"	18"	19"

Narrow T20 Options

Reduces the overall width of the Chair by 2".
Recommended only if the width of the chair is a concern. Some STF combinations offer less than 20° of tilt, but all offer a minimum of 12°.

RW	5" Casters						
12"	13"	14"					
16"		14"	15"	16"			
20"				16"	17"		
22"					17"		

RW	6" Casters & 6" Forks						
12"		14"					
16"		14"	15"	16"			
20"				16"	17"	18"	
22"					17"	18"	
24"						18"	

RW	6" Casters & 8" Forks						
16"			15"	16"			
20"				16"	17"	18"	
22"					17"	18"	19"
24"						18"	19"

RW	8" Casters & 8" Forks						
16"				16"			
20"				16"	17"	18"	
22"					17"	18"	19"
24"						18"	19"

- Standard
 Narrow

Wheel size _____

Caster size _____

Fork size _____

STF height _____

Casters

Urethane

- 1611051 5" x 1.25" Caster size
- 1611061 6" x 1.25" Caster size
- 1611063 6" x 1.50" Caster size \$ 85
- 1611081 8" x 1.25" Caster size
- 1611084 8" x 2.00" Caster size \$ 135

Pneumatic

- 1611161 6" x 1.25" Caster size E2214 & E2215 \$ 110
- 1611181 8" x 1.25" Caster size E2214 & E2215 \$ 120
- 1611184 8" x 2.00" Caster size E2214 & E2215 \$ 150

- 1611260 **6" Soft Roll** \$ 95

Rear Wheels

Composite

- 1612012 12" Wheel size
- 1612016 16" Wheel size
- 1612020 20" Wheel size
- 1612022 22" Wheel size
- 1612024 24" Wheel size

Spoke

- 1612120 20" Wheel size \$ 75
- 1612122 22" Wheel size
- 1612124 24" Wheel size
- 1612224 **24" Heavy Duty Spoke** \$ 350

Axle Type

- 1613000 **Quick Release Rear Axles**
- 1613010 **Fixed Rear Axles**

Rear Tires			
<input type="checkbox"/>	1614000	Urethane	
<input type="checkbox"/>	1614010	Pneumatic Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated.	E2211 & E2212 \$ 60
<input type="checkbox"/>	1614020	Pneumatic with Flat-Free Inserts Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoke Wheels.	E2211 & E2213 \$ 95

Hand Rims			
<input type="checkbox"/>	1615000	No Hand Rims	
<input type="checkbox"/>	1615010	Standard 20", 22", and 24" Rear Wheels only.	
<input type="checkbox"/>	1615020	Plastic Coated 20", 22", and 24" Rear Wheels only.	\$ 135

Wheel Locks			
<input type="checkbox"/>	1616010	Push to Lock	
<input type="checkbox"/>	1616020	Pull to Lock	
<input type="checkbox"/>	1616030	Attendant Wheel Lock (Push to Lock) Not compatible with Wheel Lock Extensions.	\$ 85
<input type="checkbox"/>	1616040	User & Attendant Wheel Lock	\$ 125
<input type="checkbox"/>	1616100	Wheel Lock Extensions Not compatible with Attendant Wheel Lock.	E0961 \$ 60

Anti-Tippers			
<input checked="" type="checkbox"/>	1617010	Rear Anti-Tippers Required Safety Option.	E0971 \$ 60

Back Post Style			
Push Bar Back			
<input type="checkbox"/>	1618010	Straight	
<input type="checkbox"/>	1618030	8° Bend	\$ 75
Cane Handle Back Not compatible with Headrest Mounting Plate.			
<input type="checkbox"/>	1618020	Straight	
<input type="checkbox"/>	1618040	8° Bend	\$ 75
Back Post Height Seat Pan to Handle.			
<input type="checkbox"/>	1618420	20" Height	
<input type="checkbox"/>	1618425	25" Height	

Back Options			
Angle/Height Adjustable Push Bar Handle Not compatible with Headrest Mounting Plate or Cane Handle Back. \$ 150			
<input type="checkbox"/>	1619010	Button on Outside	
<input type="checkbox"/>	1619011	Button on Inside	
<input type="checkbox"/>	1619050	Removable Cane Handle Push Bar Not compatible with Push Bar Back.	\$ 140
Seat to Back Angle			
<input type="checkbox"/>	1619090	Set at 90°	
<input type="checkbox"/>	1619096	Set at 96°	
<input type="checkbox"/>	1619102	Set at 102°	
<input type="checkbox"/>	1619108	Set at 108°	
<input type="checkbox"/>	1619114	Set at 114°	
<input type="checkbox"/>	1619120	Set at 120°	
<input type="checkbox"/>	1619220	Fold Down Back	
<input type="checkbox"/>	1619221	Non Fold Down Back	\$ 75

Back Upholstery			
<input type="checkbox"/>	1620000	No Back Upholstery	
Standard Back height measured from Seat Pan to top of Upholstery. Not available with Cane Handle Back. E2611 \$ 155			
<input type="checkbox"/>	1620018	Set at 18"	
<input type="checkbox"/>	1620019	Set at 19"	
<input type="checkbox"/>	1620020	Set at 20"	
<input type="checkbox"/>	1620021	Set at 21"	
<input type="checkbox"/>	1620022	Set at 22"	
<input type="checkbox"/>	1620023	Set at 23"	
Adjustable Tension Back height measured from Seat Pan to top of Upholstery. Not available with Cane Handle Back. E2611 \$ 295			
<input type="checkbox"/>	1620321	Set at 21"	
<input type="checkbox"/>	1620322	Set at 22"	
<input type="checkbox"/>	1620323	Set at 23"	

Backrest

<input type="checkbox"/>	1621000	No Backrest			
nxt Xtend					
<input type="checkbox"/>	1621100	16"-20" width, 16"-20" height		E2620	\$ 690
<input type="checkbox"/>	1621105	22" width, 16"-20" height		E2621	\$ 720
nxt Xtend Deep					
<input type="checkbox"/>	1621101	16"-20" width, 16"-20" height		E2620	\$ 690
<input type="checkbox"/>	1621106	22" width, 16"-20" height		E2621	\$ 720
Matrx Elite					
<input type="checkbox"/>	1621110	14"-20" width, 10"-20" height		E2620	\$ 825
<input type="checkbox"/>	1621112	HD - 22" width, 16"-20" height		E2621	\$ 1,095
<input type="checkbox"/>	1621116	Deep - 14"-20" width, 12"-20" height		E2620	\$ 950
ROHO Agility					
<input type="checkbox"/>	1621120	Min Contour 14"-20" width, 10"-18" height		E2620	\$ 825
<input type="checkbox"/>	1621122	Mid Contour 14"-20" width, 10"-18" height		E2620	\$ 850
<input type="checkbox"/>	1621124	Max Contour 14"-20" width, 10"-18" height		E2620	\$ 950
Varilite Icon					
<input type="checkbox"/>	1621130	Tall - 14"-20" width, 14"-24" height		E2615	\$ 850
<input type="checkbox"/>	1621132	Deep - 14"-20" width, 17"-24" height		E2620	\$ 895
Backrest Size					
	1621400	Specify Size (Width" x Height")	_____		

Tilt/Recline Options

Attendant Tilt only					
<input type="checkbox"/>	1623011	Trigger on Right Back Post	<input type="checkbox"/>	1623012	Trigger on Left Back Post

Headrest

<input type="checkbox"/>	1624000	No Headrest			
<input type="checkbox"/>	1624010	Headrest Mounting Plate			\$ 75
Not compatible with Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.					
<input type="checkbox"/>	1624100	nxt MultiFit Head Support Hardware		E1028	\$ 235
nxt MultiFit Head Support Pad					
<input type="checkbox"/>	1624105	Regular (9" x 5")	<input type="checkbox"/>	1624106	Large (10" x 5")
<input type="checkbox"/>	1624107	X-Large (14" x 5")		E0955	\$ 245
Stealth Headrest Hardware					
<input type="checkbox"/>	1624110	Multi-Axis HMO475		E1028	\$ 150
<input type="checkbox"/>	1624111	Link TWB480		E1028	\$ 325
Stealth Headrest Pad					
<input type="checkbox"/>	1624115	10" width	<input type="checkbox"/>	1624116	14" width

Front Rigging

<input type="checkbox"/>	1625000	No Front Rigging	<input type="checkbox"/>	1625001	Sector Blocks only		
PDG Easy-Out™ Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame.							
<input type="checkbox"/>	1625410	Fixed 70°					\$ 95
<input type="checkbox"/>	1625420	Swing In/Away Universal Angle Adjustable (0°-110°)	Set at (use increments of 10°): _____				\$ 195
<input type="checkbox"/>	1625430	Swing In/Away Elevating Legrests	High knee pivot. Includes reinforced premium Calf Pad mount.			E0990	\$ 495
Standard Front Rigging Legacy Design. Not compatible with Easy-Out™.							
<input type="checkbox"/>	1625140	Swing Away 70°					
<input type="checkbox"/>	1625160	Swing Away 90°	Not compatible with 8" Casters.				\$ 55
<input type="checkbox"/>	1625200	Swing Away Elevating Legrests				E0990	\$ 375
<input type="checkbox"/>	1625300	1" Offset Hangers for Swing Away					\$ 75

Foot Plates

<input type="checkbox"/>	1626010	Composite					
<input type="checkbox"/>	1626020	Aluminum Angle Adjustable				K0040	\$ 155
<input type="checkbox"/>	1626050	Interlocking Angle Adjustable	Seat Widths of 14"-21" only.			K0040	\$ 195

Footrest/Legrest Accessories

<input type="checkbox"/>	1627010	Swing Away 70° Calf Pad	Not compatible with 70° Gel Padding.				\$ 135
<input type="checkbox"/>	1627020	Swing Away 70° Gel Padding	Not compatible with 70° Swing Away Calf Pad, Easy-Out Front Rigging.				\$ 100
<input type="checkbox"/>	1627030	Swing Away ELR Gel Knee Pad					\$ 100
<input type="checkbox"/>	1627040	Adjustable Length Heel Loops				E0951	\$ 50
<input type="checkbox"/>	1627050	Bodypoint - Padded Calf Strap				K0038	\$ 105

Extension Tubes

<input type="checkbox"/>	1528010	Short	14"-18" for Easy-Out Fixed & SA Angle Adjustable	11"-15" for Easy-Out SA ELR	13"-17" for Standard Front Rigging		
<input type="checkbox"/>	1528020	Long	18"-24" for Easy-Out Fixed & SA Angle Adjustable	15"-21" for Easy-Out SA ELR	15"-19" for Standard Front Rigging		\$ 75
<input type="checkbox"/>	1628200	For 90° FR (11"-16")	Not compatible with Easy-Out, Standard 70° or ELR				\$ 75

Positioning Straps

<input type="checkbox"/>	1629005	Lap Belt Mount				E1028	\$ 55
Bodypoint - Hip Belts							
<input type="checkbox"/>	1629010	Center-Pull Auto Style				E0978	\$ 105
<input type="checkbox"/>	1629015	Center-Pull Auto Style, Padded				E0978	\$ 155
<input type="checkbox"/>	1629020	4-Point, Padded				E0978	\$ 175
<input type="checkbox"/>	1629030	Bodypoint - Chest Strap				E0960	\$ 85



Miscellaneous Options

<input type="checkbox"/>	1630080	Gel Hip Guides	Not compatible with Dual Post Arms. Reduces available Seat Width by 2".			E0956 & E1028	\$ 325
<input type="checkbox"/>	1630110	IV Pole	15" STF and higher.			E1029	\$ 295
<input type="checkbox"/>	1630120	O2 Holder	15" STF and higher.			E2208	\$ 395
<input type="checkbox"/>	1630160	Clear Lap Tray (Slide On)	Not compatible with 8"- 11" Dual Post, Curved Armpads, Contour Arm Trough. Seat Widths of 15"-20" only. Premium Grade or Gel Armpads recommended to ensure snug fit of the Tray.			E0950	\$ 490
<input type="checkbox"/>	1630170	Adhesive Back Hook Strip, 2" x 30"	For Cushion mounting.				\$ 10
<input type="checkbox"/>	1630301	Fuze T20 to T50 Conversion Kit				K0108	\$ 995

Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions.
Use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.



Scan this code or go to
pdgmobility.com/products/fuze-T20/
to learn more about the
Fuze T20 20° Manual Tilt Wheelchair



  PDG Mobility