



STELLAR GLT

Order Form - Canada

Prices in Canadian Dollars. January 1, 2024

Customer and Order Information

Purchase Order *	Tag	v.0227
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in
		*required

Base Chair

Stellar GLT 45° Manual Tilt Wheelchair

- | | | |
|---|--|----------|
| <input checked="" type="checkbox"/> 2501000 | Stellar GL 45° Manual Tilt Wheelchair | \$ 4,300 |
| <input checked="" type="checkbox"/> 2501009 | Lightweight Aluminum Frame Upgrade Aluminum Frame is not recommended for agitated users. | \$ 300 |
| <input type="checkbox"/> Choose the standard Stellar GL with a steel frame for agitated or heavier users. See separate order form. | | |

Weight Capacity

- | | | | |
|---|--------------------------------------|----------------------------------|--------------------|
| <input checked="" type="checkbox"/> 2502000 | 250 lbs Chair Weight Capacity | | |
| User Weight Range Required for setup of tilt uplift force. | | | |
| <input type="checkbox"/> 1902430 | Up to 75 lbs | <input type="checkbox"/> 2502431 | 75 lbs to 100 lbs |
| <input type="checkbox"/> 1902433 | 115 lbs to 150 lbs | <input type="checkbox"/> 2502434 | 150 lbs to 200 lbs |
| | | <input type="checkbox"/> 2502432 | 100 lbs to 115 lbs |
| | | <input type="checkbox"/> 2502435 | 200 lbs to 250 lbs |

Frame Options

- | | |
|---|---------------------------------------|
| <input type="checkbox"/> 2502300 | Transit Tie Down System (TTDS) |
| TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. | |

Seat Width

Overall chair width is +9.5" with 12/16" Rear Wheels and +10.5" with 20/22/24" Rear Wheels.

Standard Seat Widths

- | | | | | | |
|----------------------------------|----------------|----------------------------------|----------------|----------------------------------|----------------|
| <input type="checkbox"/> 2503014 | 14" Seat Width | <input type="checkbox"/> 2503015 | 15" Seat Width | <input type="checkbox"/> 2503016 | 16" Seat Width |
| <input type="checkbox"/> 2503017 | 17" Seat Width | <input type="checkbox"/> 2503018 | 18" Seat Width | <input type="checkbox"/> 2503019 | 19" Seat Width |
| <input type="checkbox"/> 2503020 | 20" Seat Width | <input type="checkbox"/> 2503021 | 21" Seat Width | <input type="checkbox"/> 2503022 | 22" Seat Width |

Seat Depth

Standard Seat Depths (Adjustable)

- | | | | | | |
|----------------------------------|----------------|----------------------------------|----------------|----------------------------------|----------------|
| <input type="checkbox"/> 2504016 | 16" Seat Depth | <input type="checkbox"/> 2504017 | 17" Seat Depth | <input type="checkbox"/> 2504018 | 18" Seat Depth |
| <input type="checkbox"/> 2504019 | 19" Seat Depth | <input type="checkbox"/> 2504020 | 20" Seat Depth | | |

Custom Seat Depths

- | | | | | |
|----------------------------------|----------------|----------------------------------|----------------|--------|
| <input type="checkbox"/> 2504021 | 21" Seat Depth | <input type="checkbox"/> 2504022 | 22" Seat Depth | \$ 400 |
|----------------------------------|----------------|----------------------------------|----------------|--------|

Frame Color

PDG Premium Colors



View color swatches on PDG's website.

pdgmobility.com/colors

- | | | | | | |
|----------------------------------|--------------------|----------------------------------|--------------------|----------------------------------|------------------|
| <input type="checkbox"/> 2505001 | Black Shimmer | <input type="checkbox"/> 2505015 | Sky Blue | <input type="checkbox"/> 2505016 | Poppy Red |
| <input type="checkbox"/> 2505003 | Fine Texture Black | <input type="checkbox"/> 2505010 | Lollypop Blue | <input type="checkbox"/> 2505011 | Lollypop Red |
| <input type="checkbox"/> 2505017 | Arctic Grey | <input type="checkbox"/> 2505018 | Tropical Teal | <input type="checkbox"/> 2505019 | Tangerine Orange |
| <input type="checkbox"/> 2505020 | Olive Green | <input type="checkbox"/> 2505021 | Lush Lilac | <input type="checkbox"/> 2505012 | Bubble Gum Pink |
| <input type="checkbox"/> 2505022 | Leaf Green | <input type="checkbox"/> 2505023 | Boysenberry Purple | <input type="checkbox"/> 2505014 | Daffodil Yellow |

- | | | | |
|----------------------------------|---------------------|---------------------------------|--------|
| <input type="checkbox"/> 2505200 | Custom Color | Contact PDG for specifications. | \$ 300 |
|----------------------------------|---------------------|---------------------------------|--------|

Armrests

<input type="checkbox"/> 2506000	No Armrests				
Height Adjustable Armrests (Single Post)				\$ 295	
<input type="checkbox"/> 2506111	8"-11" Range	<input type="checkbox"/> 2506114	10"-14" Range	<input type="checkbox"/> 2506116	12"-16" Range
<input type="checkbox"/> 2506100	Fixed Height Armrests (12") Full Length Armpad only.				
	Armrest Height If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.				
	Seat Pan to top of Armrest: _____		Does not include height of Armpad		
	Armrest Position Indicate preference for Armrest position. If left blank, the default setting will be used.				
In the interest of maintaining chair stably and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.					
<input type="checkbox"/> 2506401	Rear	<input type="checkbox"/> 2506402	Mid	<input type="checkbox"/> 2506403	Front

Armpads

Regular Armpads					
<input type="checkbox"/> 2507100	Desk Length - 2" x 10"	<input type="checkbox"/> 2507101	Full Length - 2" x 14"		
Waterfall Armpads					
<input type="checkbox"/> 2507105	Desk Length - 2" x 10"	<input type="checkbox"/> 2507106	Full Length - 2" x 14"		\$ 115
Wide Armpads May interfere with wheels in some tilt configurations.					
<input type="checkbox"/> 2507110	Desk Length - 3" x 10"	<input type="checkbox"/> 2507111	Full Length - 3" x 12"		\$ 200
Gel Armpads May interfere with wheels in some tilt configurations.					
<input type="checkbox"/> 2507120	Full Length - 2" x 12"	<input type="checkbox"/> 2507121	Full Length - 3.5" x 12"		\$ 135
Contour Arm Trough with Hand Pad May interfere with wheels in some tilt configurations.					
<input type="checkbox"/> 2507130	Right (4.5" x 18")	<input type="checkbox"/> 2507131	Left (4.5" x 18")		\$ 125


Armrest Options

Armrest Wideners Move Armrests outward by up to 1" on each side. Seat size does not change. May cause tilt interference.					
<input type="checkbox"/> 2507210	1" Total Additional Width				\$ 95
<input type="checkbox"/> 2507211	2" Total Additional Width				\$ 195
<input type="checkbox"/> 2507310	Clothing Guards Not compatible with 8"-11" Height range.				\$ 110
Armrest Extensions Not compatible with Desk Length Armpads.					
<input type="checkbox"/> 2507400	3" Extension - Forward	<input type="checkbox"/> 2507401	3" Extension - Rearward		\$ 75

Seat Pan

<input type="checkbox"/> 2508000	No Seat Pan	
<input type="checkbox"/> 2508001	Solid Aluminum Seat Pan	\$ 135

Seat Height

 **Don't forget**

Select desired Seat-to-Floor Height STF Height is adjustable within the range of the chosen Subframe Plate (Low/High).

6" Casters	Low								High							
	RW	13"	14"	15"	16"	17"	18"	RW	15"	16"	17"	18"	19"	20"		
	12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
8" Casters	Low								High							
	RW	13"	14"	15"	16"	17"	18"	RW	15"	16"	17"	18"	19"	20"		
	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

Casters

Urethane Casters 8" Casters are compatible only with 16" and 24" Rear Wheels.

<input type="checkbox"/> 2511061	6" x 1.25" Caster size	<input type="checkbox"/> 2511081	8" x 1.25" Caster size		
<input type="checkbox"/> 2511063	6" x 1.50" Caster size			\$	85
<input type="checkbox"/> 2511084	8" x 2.00" Caster size			\$	135

Pneumatic Casters 8" Casters are compatible only with 16" and 24" Rear Wheels.

<input type="checkbox"/> 2511161	6" x 1.25" Caster size			\$	110
<input type="checkbox"/> 2511181	8" x 1.25" Caster size			\$	120
<input type="checkbox"/> 2511184	8" x 2.00" Caster size			\$	150

<input type="checkbox"/> 2511260	Frog Legs 6" x 1.50" Aluminum Soft-Roll Casters			\$	275
----------------------------------	--	--	--	----	-----

<input type="checkbox"/> 2511400	Frog Legs Suspension Forks	Only compatible with 16 and 24" Rear Wheels. Requires 6" Soft-Roll Caster.		\$	595
----------------------------------	-----------------------------------	--	--	----	-----

Rear Wheels

Composite Wheels

<input type="checkbox"/> 2512012	12" Wheel size	<input type="checkbox"/> 2512016	16" Wheel size	<input type="checkbox"/> 2512020	20" Wheel size
<input type="checkbox"/> 2512022	22" Wheel size	<input type="checkbox"/> 2512024	24" Wheel size		

Spoke Wheels

<input type="checkbox"/> 2512120	20" Wheel size	<input type="checkbox"/> 2512122	22" Wheel size	<input type="checkbox"/> 2512124	24" Wheel size	\$	75
----------------------------------	----------------	----------------------------------	----------------	----------------------------------	----------------	----	----



Rear Wheel Position Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.

In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

<input type="checkbox"/> 2512401	Rear	<input type="checkbox"/> 2512402	Mid	<input type="checkbox"/> 2512403	Front
----------------------------------	------	----------------------------------	-----	----------------------------------	-------

Axle Type

<input type="checkbox"/> 2513000	Quick Release Rear Axles
----------------------------------	---------------------------------

<input type="checkbox"/> 2513010	Fixed Rear Axles
----------------------------------	-------------------------

Rear Tires

<input type="checkbox"/> 2514000	Urethane Tires		
----------------------------------	-----------------------	--	--

<input type="checkbox"/> 2514010	Pneumatic Tires	May result in a finished Seat Height 1/2" greater than stated.	\$	60
----------------------------------	------------------------	--	----	----

<input type="checkbox"/> 2514020	Pneumatic Tires with Flat-Free Inserts	May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels.	\$	95
----------------------------------	---	---	----	----

Handrims

Handrims are only available with 20", 22", and 24" Rear Wheels.

<input type="checkbox"/> 2515000	No Handrims		
----------------------------------	--------------------	--	--

<input type="checkbox"/> 2515010	Standard Aluminum Handrims		
----------------------------------	-----------------------------------	--	--

<input type="checkbox"/> 2515020	Plastic Coated Handrims	Not compatible with Drum Brakes.	\$	150
----------------------------------	--------------------------------	----------------------------------	----	-----

<input type="checkbox"/> 2515030	Projection Handrims	Composite Wheels only. 10 Projections, 30°, 1-1.8" Length.	\$	495
----------------------------------	----------------------------	--	----	-----

Wheel Locks & Brakes

User Wheel Locks

<input type="checkbox"/> 2516010	Push to Lock	<input type="checkbox"/> 2516020	Pull to Lock
----------------------------------	--------------	----------------------------------	--------------

<input type="checkbox"/> 2516030	Attendant Wheel Locks (Push to Lock)		\$	85
----------------------------------	---	--	----	----

<input type="checkbox"/> 2516040	User & Attendant Wheel Locks	Not compatible with Drum Brakes.	\$	125
----------------------------------	---	----------------------------------	----	-----

<input type="checkbox"/> 2516300	Drum Brakes	Available with 12" Composite or 20" Spoke Rear Wheels. Urethane Tires only.	12" WHEELS COMING SOON	\$	795
----------------------------------	--------------------	---	-------------------------------	----	-----

User-activated Unilateral Wheel Lock Pull to Lock.

<input type="checkbox"/> 2516310	Left Activation	<input type="checkbox"/> 2516311	Right Activation	\$	795
----------------------------------	-----------------	----------------------------------	------------------	----	-----



<input type="checkbox"/> 2523110	Tilt Trigger Extension		\$ 50
Tilt Limit Optional. If no tilt limit is selected, the full tilt range will be available.			
<input type="checkbox"/> 2523220	Set at 20°	<input type="checkbox"/> 2523230 Set at 30°	<input type="checkbox"/> 2523240 Set at 40°
Power Tilt 250 lbs Weight Capacity. 45° tilt range. Includes 2 batteries for off-chair charging. pdgmobility.com/power-tilt Minimum Seat Width is 16". Power Tilt is compatible with either O2 Holder or IV Pole on the same chair, not both.			
<input type="checkbox"/> 2523400	Standard Power Tilt		\$ 1,795
<input type="checkbox"/> 2523410	Power Tilt with Manual Tilt Override	Allows for instant manual release of tilt in case of power tilt failure.	\$ 2,295
	Toggle Switch	<input type="checkbox"/> 2523421 Installed on Right side	<input type="checkbox"/> 2523422 Installed on Left side
	Button Switch	<input type="checkbox"/> 2523431 Installed on Right side	<input type="checkbox"/> 2523432 Installed on Left side
	Stealth Egg Switch	<input type="checkbox"/> 2523440 Not installed	
	Optional for Egg Switch:	<input type="checkbox"/> 2523445 Swing Away Mounting Hardware	Not compatible with all Headrest makes and models.
Specify Headrest to be used: _____			

Recline

<input type="checkbox"/> 2523310	Reclining Backrest (Manual)	\$ 895
30° maximum recline. Not compatible with Backrest Extension, Fold-down Back, Tension Adjustable Upholstery. Finished Back Post Height is 3" taller than selected Back Post Height. Limits Tilt by up to 20°.		

Front Rigging

<input type="checkbox"/> 2525000	No Front Rigging	<input type="checkbox"/> 2525001 Sector Blocks only	
Easy-Out™ Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame. pdgmobility.com/easy-out			
<input type="checkbox"/> 2525410	Fixed 70°		\$ 95
<input type="checkbox"/> 2525420	Universal Angle Adjustable (0°-110°)	Set at angle: _____ (use increments of 10°)	\$ 295
Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings.			
Swing Away Front Rigging Legacy Design. Not compatible with Easy-Out.			
<input type="checkbox"/> 2525130	Swing Away 60°	<input type="checkbox"/> 2525140 Swing Away 70°	
<input type="checkbox"/> 2525200	Swing Away Elevating Legrests	Includes Calf Pads.	\$ 375
<input type="checkbox"/> 2525300	1" Offset Hangers for Swing Away		\$ 75
<input type="checkbox"/> 2525110	Standard Plug In 70° Front Rigging	Legacy Design. Not compatible with Easy-Out.	

Footplates

<input type="checkbox"/> 2526010	Composite Footplates	
<input type="checkbox"/> 2526020	Aluminum Angle Adjustable Footplates	\$ 165
<input type="checkbox"/> 2526030	One Piece Fixed Footplate	Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.
<input type="checkbox"/> 2526040	One Piece Flip-up Footplate	Requires selection of Composite or Aluminum Angle Adjustable.
<input type="checkbox"/> 2526050	Interlocking Angle Adjustable Footplates	Available for Seat Widths up to 21" only.

Footrest/Legrest Accessories

Calf Pads Not compatible with Gel Padding.		
<input type="checkbox"/> 2527012	For Easy-Out™ Front Rigging	Includes reinforced premium Calf Pad mount.
<input type="checkbox"/> 2527010	For Swing Away Front Rigging	
<input type="checkbox"/> 2527020	Gel Padding	Not compatible with Calf Pad or Easy-Out Front Rigging.
Gel Knee Pads		
<input type="checkbox"/> 2527032	For Easy-Out™ Angle Adjustable Front Rigging	
<input type="checkbox"/> 2527030	For Swing Away ELR Front Rigging	
<input type="checkbox"/> 2527040	Length Adjustable Heel Loops	
<input type="checkbox"/> 2527050	Bodypoint Padded Calf Strap	



Front Rigging Extension Tubes

Regular Extension Tubes

- | | | | | | |
|----------------------------------|-------|------------------------------------|---|----|----|
| <input type="checkbox"/> 2528010 | Short | 14"-18" for Easy-Out Front Rigging | 13"-17" for Swing Away and Standard Plug In Front Rigging | | |
| <input type="checkbox"/> 2528020 | Long | 18"-24" for Easy-Out Front Rigging | 15"-19" for Swing Away and Standard Plug In Front Rigging | \$ | 75 |



Lower Leg Length If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.

Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.

Miscellaneous Options

Spoke Guards Available for 20"-24" Composite Rear Wheels. Made from durable and lightweight aluminum.

\$ 195

- | | | | |
|----------------------------------|-------|----------------------------------|-------------|
| <input type="checkbox"/> 2530020 | Black | <input type="checkbox"/> 2530021 | Frame Color |
|----------------------------------|-------|----------------------------------|-------------|

Stealth Flip Down Abductor Requires Solid Seat Pan.

\$ 650

- | | | | | | |
|----------------------------------|-----------|------------------|----------------------------------|-----------|------------------|
| <input type="checkbox"/> 2530060 | Small Pad | 3"x4"x3"x2" size | <input type="checkbox"/> 2530062 | Large Pad | 5"x5"x4"x3" size |
|----------------------------------|-----------|------------------|----------------------------------|-----------|------------------|

- | | | | | |
|----------------------------------|-----------------------|-------------------------------------|----|-----|
| <input type="checkbox"/> 2530080 | Gel Hip Guides | Reduces available Seat Width by 2". | \$ | 495 |
|----------------------------------|-----------------------|-------------------------------------|----|-----|

- | | | | | |
|----------------------------------|----------------|--|----|-----|
| <input type="checkbox"/> 2530110 | IV Pole | Not compatible with Attendant Wheel Locks. Requires frame modification and cannot be added to chair later. | \$ | 295 |
|----------------------------------|----------------|--|----|-----|

- | | | | | |
|----------------------------------|------------------|--|----|-----|
| <input type="checkbox"/> 2530120 | O2 Holder | Not compatible with Attendant Wheel Locks. Requires frame modification and cannot be added to chair later. | \$ | 395 |
|----------------------------------|------------------|--|----|-----|

- | | | | | |
|----------------------------------|---------------------------|--|----|-----|
| <input type="checkbox"/> 2530200 | Q'Straint Lap Belt | Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" Belt Length. | \$ | 245 |
|----------------------------------|---------------------------|--|----|-----|

- | | | | | |
|----------------------------------|----------------------------|------------------------------------|----|-----|
| <input type="checkbox"/> 2530222 | Wheelchair Backpack | Black with PDG Logo. 17"x14" size. | \$ | 145 |
|----------------------------------|----------------------------|------------------------------------|----|-----|

- | | | | | |
|----------------------------------|---|-----------------------|----|----|
| <input type="checkbox"/> 2530170 | Adhesive Back Hook Strip, 2" x 30" | For Cushion mounting. | \$ | 10 |
|----------------------------------|---|-----------------------|----|----|

On-Chair Seating & Positioning

- Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. pdgmobility.com/seating

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.